# ARTROMOT®-SP3

SW 1.X

- DE · Gebrauchsanleitung
- GB · Operating Instructions
- FR · Mode d'emploi
- ES · Instrucciones de uso
- IT · Istruzioni per l'uso
- NL · Gebruiksaanwijzing

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<https://www.enovis-medtech.de/unternehmen/qualitaet-und-verantwortung>

# Contents

<b>Description of the ARTROMOT®-SP3</b>	<b>III</b>
<b>ARTROMOT®-SP3 setup illustrations</b>	<b>VIII</b>
<b>Symbol overview ARTROMOT®-SP3</b>	<b>V</b>
<b>1. How to use the CPM device</b>	<b>39</b>
1.1 Fields of application	39
1.2 Therapy objectives	39
1.3 Indications	39
1.4 Contraindications	39
<b>2. Description of the ARTROMOT®-SP3</b>	<b>40</b>
2.1 Description of the device components	41
2.2 Description of the control pendant	42
2.3 Explanation of symbols	45
2.4 Explanation of symbols (connections and nameplate)	46
<b>3. Safety information</b>	<b>48</b>
<b>4. Device setup</b>	<b>52</b>
4.1 Connecting the device, performance check	52
4.2 Adjusting the device to the patient	53
<b>5. Setting the treatment values</b>	<b>54</b>
5.1 General information on programming the ARTROMOT®-SP3	54
5.2 Programming the ARTROMOT®-SP3 Standard models	55
5.3 Treatment value details – Standard models	56
5.4 Programming the ARTROMOT®-SP3 Comfort models	59
5.5 Protocol details – Comfort model	60
<b>6. Care, Maintenance, Transport, Conversion</b>	<b>64</b>
6.1 Care / re-use	64
6.2 Maintenance (fuse replacement)	65
6.3 Transport	67
6.4 Conversion	68
<b>7. Environmental protection statement</b>	<b>69</b>
<b>8. Specifications</b>	<b>69</b>
<b>9. Contact</b>	<b>71</b>
<b>10. Technical service</b>	<b>72</b>
10.1 Technical hotline	72
10.2 Shipment	72
10.3 Spare parts	72
<b>11. CE-Declaration of Conformity</b>	<b>73</b>
<b>12. IEC 60601-1-2 and IEC 60601-1-11</b>	<b>226</b>
12.1 Electromagnetic emissions	226
12.2 Electromagnetic immunity	227
12.3 RF (Radio Frequency) Interference	229

# 1. How to use the CPM device

## 1.1 Fields of application

The **ARTROMOT®-SP3** is a motor-operated Continuous Passive Motion (CPM) device providing motion to the ankle joint.

Suitable for use in hospitals, clinics, general practices and rental services as well as in the patient's home, it is an important supplement to medical and therapeutic treatment.

## 1.2 Therapy objectives

CPM therapy with the **ARTROMOT®-SP3** is mainly used to prevent the negative effects of immobilization, to allow patients to regain painless mobility of joints at an early stage and to promote healing and achieve a positive functional result.

Other objectives of therapy include:

- improvement of joint metabolism
- prevention of joint stiffness
- promotion of the regeneration and healing of cartilage and damaged ligaments
- faster hematoma/fluid resorption
- improved lymph and blood circulation
- thrombosis and embolism prophylaxis

## 1.3 Indications

The CPM device is indicated in the treatment of most injuries and diseases of the ankle joint as well as in the postoperative treatment after joint surgery and in the treatment of joint diseases. Examples:

- joint distortion and contusion
- arthrotomy and arthroscopy procedures in combination with synovectomy, arthrolysis, or other intra-articular interventions
- mobilization of joints in anesthetized patients
- exercise-stable, surgically restored fractures, pseudoarthrosis and osteotomy
- myoplasty
- arthroplasties including implantation of endoprostheses

## 1.4 Contraindications

**Do NOT use** the **ARTROMOT®-SP3** on patients with:

- acute inflammatory processes in the joints, unless on the order of a physician
- spastic paralysis
- unstable osteosynthesis

## 2. Description of the ARTROMOT®-SP3

The motorized CPM device provides passive motion to the ankle joint as follows:

plantar flexion / dorsal extension  
50° / 0° / 40°

inversion / eversion  
40° / 0° / 20°

It can be reconfigured for use on either side and allows patients to exercise in bed and while sitting on a chair.

### Note!

To unambiguously represent the current position of the CPM device, the values for plantar flexion and inversion are marked with the symbol "–" both on the display and in this document.

These are some of the outstanding features of the **ARTROMOT®-SP3**:

- anatomically correct adjustment
- physiological movements
- control pendant for precise adjustment of patient-specific therapy values
- symbols for easy operation of the control pendant
- chip card for storage of the programmed therapy parameters

### Biocompatibility

Those parts of the **ARTROMOT®-SP3** that come into contact with the patient when the device is used as intended, are designed to fulfill the biocompatibility requirements of the applicable standards.

### Essential Performance Features

- The programmed angles are accurately maintained with a tolerance of +/- 2°.
- The programmed speeds are accurately maintained with a tolerance of +/- 5%.
- The selected mode and mechanical settings do not change during operation.

### Frequent functions:

- a) Unpacking (device and accessories)
- b) Assemble and block footplate on footplate bracket
- c) Assemble and block motor B and footplate on moving bracket
- d) Device connections (mains power connection, connecting the control pendant and the motion element)
- e) Mechanical adjustments of the motion element on the rotational axes
- f) Programming the control pendant (range of motion, speed, functions)
- g) Storage


## 2.1 Description of the device components

**Note: Please fold out page III!**

1. Clamping levers to adjust the height of the calf support assembly
2. Calf support assembly
3. Motor A (center of rotation for dorsal extension/plantar flexion)
4. Bracket for control pendant
5. Clamping lever to adjust the inclination of the carriage
6. Aiming pin for adjustment of the rotational axis of the upper ankle joint
7. Footplate
8. Knurled screw for adjustment of the footplate height
9. Knurled screw to adjust the footplate to the leg length
10. Motor B
11. Footplate bracket
12. Thumbscrew to secure motor B on the moving bracket for dorsal extension / plantar flexion
13. Bore hole to mount motor B to the moving bracket for dorsal extension / plantar flexion
14. CPM base
15. Moving bracket for dorsal extension / plantar flexion
16. Power switch (ON/OFF)
17. Fuse cap
18. Connection for power cord
19. Locking device to adjust the height of the carriage
20. Perforation to set the height of the carriage
21. Patient chip card<sup>1</sup>
22. Control pendant
23. Rating plate
24. Connection for control pendant


## 2.2 Description of the control pendant

### 2.2.1 Control pendant in normal mode:


<sup>1</sup> Devices of the ARTROMOT®-SP3 product family with patient chip card only.

## 2.2.2 Control pendant in MENU selection mode

## 2.2.3 Control pendant in ROM programming mode

## 2.2.4 Control pendant in general programming mode


## 2.3 Explanation of symbols

Also refer to symbol overview on page V.

### Standard protocols:

	Dorsal extension
	Plantar flexion
	Inversion
	Eversion
	Pause
	Timer
	Speed
	New patient
	Reverse on load - motor A
	Reverse on load - motor B
	Motor A ON/OFF
	Motor B ON/OFF
	Transport setting
A+B	Synchronized/ non-synchronized mode
	Total therapy time
	Service menu

### Comfort protocols:

	Warm-up protocol
AIII B	Isolation protocol
	Dorsal extension / plantar flexion therapy documentation
	Inversion / eversion therapy documentation
	Stretching in dorsal extension
	Stretching in plantar flexion
	Stretching in inversion
	Stretching in eversion

## 2.4 Explanation of symbols (connections and nameplate)

Alternating current

for serial numbers < 20,000:  
protection class I equipment. The medical device must be connected to a system with protective earth conductor!

for serial numbers > 20,000:  
protection class II equipment. The medical device has a double or reinforced insulation.

Type B applied part

Power switch OFF

Power switch ON

The year and month of manufacture are indicated next to this factory symbol.

The name next to this factory symbol is the manufacturer.

The number next to this symbol is the article reference number

Device complies with Council Directive 93/42/EEC about medical devices, tested and approved by DQS Medizinprodukte GmbH.

The number next to this symbol is the serial number.

Caution! Observe warnings set forth in instruction manual!

Refer to instruction manual! ISO 7010-M002

Refer to instruction manual! ISO 7000-1641

Do not dispose with unsorted municipal waste.

Protect from moisture.

## IP21

The IP rating indicates the level of protection and thus the suitability of the device for use under different ambient conditions.

The rating IP21 means:

2 is the level of protection against contact and solid objects.

The digit 2 means:

- Protection against contact: protected against contact with a finger.
- Protection against solid objects: protected against solid objects (diameter of 12.5 mm and greater).

1 is the level of protection against moisture.

The digit 1 means: protection against vertically falling water drops.

Observe temperature limits (storage)! ISO 7000-0632

The overall weight of the device is indicated next to this symbol.

### 3. Safety information

#### Definitions

Be sure to read the safety information before using the CPM device. The safety statements are classified as follows:

#### **Danger!**

indicates an imminent hazard. If not avoided, this hazard will result in death or serious injury.

#### **Warning!**

indicates a hazard. If not avoided, the hazard can result in death or serious injury.

#### **Caution!**

indicates a potential hazard. If not avoided, the hazard may result in minor injury and/or product/property damage.

#### Safety information

#### **Danger!**

Explosion hazard –

The **ARTROMOT®-SP3** is not designed for **use in areas of rooms used for medical purposes** where an explosion hazard may occur. An explosion hazard may result from the use of flammable anesthetics, skin cleansing agents, and disinfectants.

#### **Warning!**

Patient hazard –

- Only **Authorized individuals** are allowed to operate the **ARTROMOT®-SP3**. Individuals are authorized after receiving training in the operation of the device and reading these operating instructions.
- Before using the CPM device, the operator must ascertain that it is in correct working order and operating condition. The cables and connectors, in particular, must be checked for signs of damage. Damaged parts must be replaced immediately, before use.
- **Before therapy**, a **test run** consisting of several exercise cycles must be completed, first without and then with the patient. Check that all setting screws are tightened.
- Stop therapy immediately when you have doubts about the device settings and/or the therapy protocol.

## Warning!

Patient hazard –

- It is important that the patient's **position is anatomically correct**. Therefore, carefully verify the following settings/positions:
  1. Carriage angle
  2. Ankle joint axis/carriage axis motor A
  3. Footplate adjustment
  4. Calf support assembly
- Movements must **not** cause **any pain or irritation**.
- The patient must be **fully conscious** while being instructed in the use of the CPM device and when using the device.
- The **choice** of the **therapy parameters** to program and of the **therapy protocols** to use is restricted to the responsible **physician or therapist**. The doctor or therapist must decide on a case-by-case basis whether the device can be used with each particular patient.
- The patient must be familiar with the functions of the **ARTROMOT®-SP3 control pendant** which must be **within easy reach** of the patient, allowing him or her to stop therapy if needed. **Patients unable to operate the control pendant**, e.g. paralytic patients, must always be supervised by specially trained staff during therapy.
- After data storage, write the patient's name on the **patient chip card**<sup>1</sup>. The card should only be used for this particular patient. If the **patient chip card**<sup>1</sup> is used for another patient, be sure to **delete the previous patient's data** from the card first (see: section 5.3, "New patient"). **Use original chip cards only.**<sup>1</sup>
- The **ARTROMOT®-SP3** must only be used with **accessories** approved by DJO.
- Modifications to the medical device described in this document without the manufacturer's written consent are prohibited.

- Stability of the CPM device must always be ensured while it is in use. The ARTROMOT-SP3 must only be set up on surfaces that guarantee its stability. Very soft or instable surfaces (such as waterbeds) are NOT suitable.
- Do not allow **parts of the body or any objects** (such as blankets, cushions or cables) to get caught **in the moving parts** of the CPM device.

## Warning!

- Special care must be taken when small children and babies are present while the device is in use! Sufficient distance to the device is mandatory for their safety!
- Never leave the device unattended when it is switched on! Switch the device off and disconnect the power line from the wall outlet!
- After use, store the device in a safe place! Ensure device stability also during storage!

<sup>1</sup> Devices of the **ARTROMOT®-SP3** product family **with patient chip card only**.

## Warning!

Shock hazard –

Strictly observe the following warnings. Failure to do so endangers the lives of the patient, the user, and other persons involved.

- Allow the **ARTROMOT®-SP3** to reach room temperature **before use**. If the device was transported at **temperatures below 0°C (32°F)**, allow it to stand at room temperature for about 2 hours, until all condensation has disappeared.
- The **ARTROMOT®-SP3** device must only be operated in dry rooms.
- When disconnecting the device from the power line, remove the plug from the wall outlet first, before disconnecting the cable from the device.
- When connecting the device to other equipment or when creating a medical system, check that the sum of leakage currents will not cause any hazard. Please contact DJO, if you have questions in this matter.
- Do not use multiple portable socket outlets (MPSO) to connect the device to the power line. The **ARTROMOT®-SP3** must be connected to a **properly installed wall outlet with a non-fused earthed wire**. Before connecting the power cord, it must be completely unrolled and placed such that it will not get caught in the moving parts of the device.
- Before cleaning and service interventions, **disconnect the device from the power line** by removing the power cord from the wall outlet.
- Liquids must not be allowed to enter the CPM device or the control pendant. If liquids have entered into the devices, the **ARTROMOT®-SP3** must be immediately checked by a service technician, before it can be reused.

## Warning!

Equipment malfunction –

- Magnetic and electrical fields are capable of interfering with the proper performance of the device. For this reason make sure that all external devices operated in the vicinity of the device comply with the relevant EMC requirements. X-ray equipment, MRI devices, radio systems, cell phones, etc. are possible sources of interference as they may emit higher levels of electromagnetic radiation. Keep the device away from such equipment and verify its performance before use.
- Refer **repair and maintenance** to authorized persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.
- **Route all cables below the device frame to either side**, ensuring that they cannot get caught by the moving parts during operation.
- Check the **ARTROMOT®-SP3** for damage and loose connections at least **once a year**. Damaged and worn parts must immediately be replaced with original spare parts by authorized staff.

## Caution!

Preventing chafing and pressure sores –

When your patient is **adipose**, very **tall** or **very short**, be sure to prevent chafing and pressure sores. Place the leg concerned in a moderate abductive position if deemed appropriate.

## **Caution!**

Equipment damage –

- Check that the voltage and frequency ratings of your local **power line** are those indicated on the nameplate.
- The leg support element withstands a **maximum continuous load of 20 kg**.
- Do not allow **any objects** (such as blankets, cushions, or cables) to get caught in the **moving parts** of the device.
- Do not expose the ARTROMOT®-SP3 to direct sunlight, because some of the components may warm up excessively.

Furthermore, the device must be set up at a safe distance from radiators to prevent excessive temperature rises.

- The presence of children, pets and rodents does not normally impair the functioning of the device. However, avoid contamination of the device by children or animals, from dust and lint, and keep them at a safe distance from the device. The safety statements set forth apply.

## 4. Device setup

**Note:** For a better understanding of each step, please fold out pages III and VIII.

### 4.1 Connecting the device, performance check

The equipment supplied includes these items:

base unit, footplate bracket (11), footplate (7), control pendant (22), patient chip card (for ARTROMOT®-SP3 with chip card only), power cord (not shown), operating instructions

1. Set up the device and connect the required cables by reversing the transport instructions given in section 6.3.
2. Turn the **power switch** (16) on.

#### **Caution!**

Equipment damage –

Be sure to lift the moving bracket for dorsal extension / plantar flexion (15) a little, relieving the strain on it, before you either press the locking pins for the height adjustment together or loosen the clamping lever for adjustment of the inclination (5).

3. Release the locking device (19). To do so, first relieve the strain on the carriage by lifting the moving bracket for dorsal extension / plantar flexion (15) a little. Then press the two black pins of the locking device together simultaneously (Fig. A).
4. Set the carriage to an angle of approx. 60°, then release the pins of the locking device. Check that the pins engage in the perforation (20).
5. Again, hold the moving bracket for dorsal extension / plantar flexion (15) until it forms an angle of 45° with the base. Close the clamping lever to lock the carriage in that position (Fig. B).

6. Follow these steps to set the carriage to the **home position**:

#### **ARTROMOT®-SP3 without patient chip card**

Press the **MENU** button on the control pendant twice to access programming level 2 (you advance one level each time you press the button).

Press the "New Patient" parameter button **→ 0 ←** and select this function (a check mark appears in the circle next to the function).

Press the **START** button. The carriage automatically enters the home position.

#### **ARTROMOT®-SP3 with patient chip card**

##### – **Initial adjustment for new patients**

Insert the original patient chip card (21) into the control pendant (22).

Press the **MENU** button on the control pendant three times to access programming level 3 (you advance one level each time you press the button).

Press the "New patient" parameter button **→ 0 ←** and select this function (a check mark appears in the circle next to the function)

Press the **START** button.

The CPM device automatically enters the **home position**.

##### – **Adjustment with programmed chip card**

Insert the original patient chip card (21) into the control pendant (22).

Press the **START** button.

The CPM device automatically enters the **home position**.

#### **Performance check**

If the control pendant can be operated as described above and the **ARTROMOT®-SP3** enters the home position (for home position

values, refer to sections 5.3 and 5.5), the device has passed the performance check.

The device also runs performance checks regularly during operation. This is what happens, if a problem is identified:

- An audio signal sounds.
- the device switches off immediately
- The message "ERR", accompanied by a code number (e.g. ERR 5), appears on the display.

In this situation, you may attempt to restart the device by turning it briefly off and on again with the power switch. With the device switched off, check that all plugs are correctly connected. Switch the CPM device on again: if the error message persists, have the device inspected by a Service technician, before using it again.

#### 4.2 Adjusting the device to the patient

Before any treatment is possible, adjust the device as follows:

1. Before you adjust the ARTROMOT®-SP3 to the patient, you may have to convert the device for use on the left or right foot (see section 6.4 Conversion).
2. First lift the moving bracket for dorsal extension / plantar flexion (15) a little to hold it securely and to remove the strain. Simultaneously press the two black pins of the locking device (19) together. Now you can adjust the carriage to the correct height. Release the pins of the locking device (19). Check that the pins engage in the perforation (20) (Fig. A).
3. Continue to hold the moving bracket for dorsal extension / plantar flexion (15) and release the clamping lever for adjustment of the inclination (5). Adjust the correct inclination for the carriage. Close the clamping lever to lock the carriage in that position (III. B).
4. Place the patient's leg on the carriage and adjust the footplate (7) to the length of the patient's leg. To do so, loosen the knurled screw (9) and adjust the footplate accordingly. It is the purpose of this setting to adjust the rotational axis of the CPM device to the compromise center of rotation of the upper ankle joint. Check that the foot is firmly positioned against the footplate. Tighten the knurled screw (9) (III. C)

5. Loosen the clamping lever adjusting the height of the calf support assembly (1) (Fig. D) and the knurled screw adjusting the height of the footplate (8) (Fig. E). Set the calf support assembly and the footplate to a level where the rotational axis of the carriage and the compromise center of rotation of the upper ankle joint coincide. Close the clamping lever (1) and tighten the knurled screw (8)
6. Check the settings and make sure that the patient is comfortably positioned. If necessary, repeat steps 2 to 5 to achieve a perfect setup.

#### Note!

- The aiming pin (6) is a valuable aid to adapt the carriage's rotational axis to the compromise center of rotation of the upper ankle joint. It identifies the correct compromise axis (Fig. F)
- If the CPM device is to be used on a supine patient, it is recommended to select a hip angle of approx. 40° at the ARTROMOT®-SP3.

#### ⚠ Caution!

Equipment damage –

- Cover the **leg support assemblies with disposable tissues**, when using the ARTROMOT®-SP3 after surgery. This helps prevent discoloration.
- Before you change the carriage inclination and/or height, always lift the moving bracket for dorsal extension / plantar flexion a little.

#### ⚠ Caution!

Patient hazard –

- Ensure that the rotational axes of the carriage and of the ankle joint coincide both in the vertical and in the horizontal plane (Fig. F).
- Before you change the carriage inclination and/or height, always lift the moving bracket for dorsal extension / plantar flexion a little.

## 5. Setting the treatment values

### Warning!

Patient hazard –

Before therapy, a test run consisting of several exercise cycles must be completed without the patient. Then, repeat the test run with the patient and check that the movements do not cause any pain.

### Note!

See also 2.2 and 2.3 as well as page V!

### 5.1 General information on programming the ARTROMOT®-SP3

1. You activate the programming mode by briefly pressing the **MENU** button on the control pendant.
2. The treatment parameters and functions are allocated to four (Standard model) or six (Comfort model) programming levels (four per level).

To be able to program a parameter you will have to access the corresponding programming level. This is also done with the **MENU** button. You advance one level each time you press the button. The code M1, M2, etc. that appears in the middle of the display indicates the programming level.

3. You activate the treatment parameters and functions with the four parameter buttons below the display. The symbols above the **4 parameter buttons** indicate the assigned parameters and functions.

This is what happens when you press one of the parameter buttons to select a parameter:

- The corresponding symbol appears on the display in a larger format.
- The set value is displayed.
- The symbol above the parameter button is highlighted.

4. With the + / - buttons (plus/minus) you change the displayed value. When you press and hold the button, the value will change at a higher rate.

Some of the (special) functions can only be enabled and disabled. This is done by pressing the corresponding parameter button or with the + / - buttons. Active parameters are identified with a check mark in the circle next to the symbol.

5. Having programmed all parameters, press the **STOP** button to save the values.
6. Next, press the **START** button: The ARTROMOT®-SP3 will check the set values, moves to the middle position and stops.
7. Press the **START** button again to start therapy.

In the **synchronized mode**, the carriage will then move from the middle position to the maximum eversion value. Afterwards, the carriage will simultaneously move to the maximum inversion value and to the maximum plantar flexion value, and then simultaneously to the maximum dorsal extension value and the value halfway between inversion and eversion. After reaching this position, the cycle restarts: maximum eversion at the same time as halfway between dorsal extension / plantar flexion, etc.

After activation of the **START** button in the **non-synchronized mode**, the motors will perform random movements, each motor reversing after reaching the maximum values.

## Note!

- Refer to sections 5,3 and 5,5 for a description of the parameters.
- To **view the set parameter values**, press the corresponding parameter button. Before, however, you have to press the STOP button and access the correct menu level.
- To prevent accidental changes of the parameter settings, **lock the controls** by simultaneously pressing the + (**plus**) and – (**minus**) buttons.

Press both buttons again to unlock.

- Selecting the "New Patient" function will automatically delete the data on the patient chip card. When you have finished programming the device and press the STOP button, the settings will automatically also be saved to the patient chip card.
- **Emergency stop function:** The ARTROMOT®-SP3 will stop immediately, when any of the controls is pressed during therapy. When the START button is pressed in the **synchronized mode**, the carriage will move to the position halfway between the set values, and stops. Patient treatment can be resumed by pressing the START button. When the START button is pressed again in the **non-synchronized mode**, the carriage will reverse its direction.
- The manufacturer recommends an application period of up to 1 hour per therapy session.
- For a uniform, smooth and gentle transition between the two directions of motion, the speed is automatically reduced before reaching the reversal point and, on reaching the reversal point, it is continuously increased until it reaches the set speed.

## Patients with a programmed chip card

- First, complete the mechanical adjustments.
- Then, insert the chip card (the patient is not yet positioned on the CPM device).
- Press the **START** button: the device will move to the middle position as specified by the parameters stored on the chip card and stops.
- Subsequently, place the patient on the CPM device and press the **START** button to initiate therapy.

## 5.2 Programming the ARTROMOT®-SP3 Standard models

Different programming levels are provided to program the **ARTROMOT®-SP3-Standard** models.

You change between levels by pressing the **MENU** button.

The display always indicates the currently selected level.

The following **treatment values, settings** and **information** can be entered/viewed on the control pendant (22):

### LEVEL 1:

- Dorsal extension 
- Plantar flexion 
- Inversion 
- Eversion 

### LEVEL 2:

- Pause 
- Timer (therapy time) 
- Speed 
- New patient 

### LEVEL 3:

- Reverse on load - motor A
- Reverse on load - motor B
- Motor A ON/OFF
- Motor B ON/OFF

### LEVEL 4:

- Transport setting
- Synchronized/  
non-synchronized mode
- Total therapy time
- Service menu

### Note!

- While you adjust the extension/  
flexion values, the CPM device will  
move to the set range. This allows  
you to easily and quickly determine  
the ROM where the patient does not  
experience pain.
- Special functions can be programmed  
and retrieved with **ARTROMOT®-SP3  
Comfort** devices (see sections 5.4  
and 5.5).
- Only **ARTROMOT®-SP3 chip card**  
models allow therapy protocols to be  
**saved** to an inserted chip card.
- The angle entered last for the respec-  
tive directions of motion will be saved  
in each case.

## 5.3 Treatment value details – Standard models

- To access the appropriate programming  
level, press the **MENU button** repeatedly.
- The treatment parameters are selected  
with the different **parameter buttons**.
- To change the treatment values, press  
the **+ / - buttons**.
- To enable/disable a function, press the  
respective **parameter button** again.
- To save the entered values, press the  
**STOP button**.

### LEVEL 1:

- **Dorsal extension** ↑  
maximum value: **40 degrees**
- **Plantar flexion** ↓  
maximum value: **- 50 degrees**
- **Inversion** ← **40 degrees**
- **Eversion** → **20 degrees**

### Note!

- The programmed values and the  
values measured at the patient's foot  
may deviate slightly.
- To ensure a physiological motion, the  
carriage moves to the following posi-  
tions in the synchronized mode in the  
following sequence:  
maximum dorsal value at the same  
time as halfway between inversion /  
eversion  
maximum eversion value at the same  
time as halfway between dorsal  
extension / plantar flexion  
maximum plantar flexion at the same  
time as maximum inversion
- To ensure that the movements you  
program are physiological in nature,  
we recommend to program the car-  
riage in the following order:  
1. Programming the maximum dorsal  
extension value

- Resting the foot by programming a plantar flexion angle of  $0^\circ$

- Programming the maximum eversion value

- Resting the foot by programming an inversion angle of  $0^\circ$

- Programming the maximum plantar flexion value

- Programming the maximum inversion value

- When the **START** button is pressed in the synchronized mode, the device will first move from the middle position to the maximum eversion value. Afterwards, the carriage will simultaneously move to the maximum inversion value and to the maximum plantar flexion value, and then simultaneously to the maximum dorsal extension value and to the value half way between inversion and eversion. After reaching this position, the cycle restarts: maximum eversion at the same time as halfway between dorsal extension / plantar flexion, etc.

## LEVEL 2:

### ■ Pauses

Pauses occur at each programmed maximum value.

These are the **three pause points**:

maximum dorsal value at the same time as halfway between inversion / eversion

maximum eversion value at the same time as halfway between dorsal extension / plantar flexion

maximum plantar flexion at the same time as maximum inversion

Pauses are adjustable in steps of 1 second between 0 and 59 seconds.

**Default setting: no pause**

### ■ Timer (therapy time)

**Default setting is continuous operation.**

A clock symbol in the upper right-hand corner of the display identifies the continuous mode of operation. The clock indicates the elapsed therapy time.

In the **continuous mode**, the device must be stopped with the **STOP** button.

However, you can also select therapy durations of **1 to 59 minutes in steps of 1 minute** and of **1 to 24 hours in steps of 30 minutes**. When the time has elapsed, the device switches **automatically off** and stops in the middle between the set limit values.

In this case, a circle replaces the clock symbol. The circle fills as the therapy time progresses.

### ■ Speed

The speed can be adjusted in steps of 5 % between 5 % and 100 %.

**Default setting: 50 %**

### ■ New patient

With this function, the CPM device will move to the home position. Select the function and press the **START** button. The device enters the home position and existing therapy parameters will be deleted.

On **ARTROMOT®-SP3 devices with a patient chip card** the factory defaults are restored. All values stored on the chip card will be deleted.

The carriage stops in the home position, dorsal extension / plantar flexion =  $0^\circ$  and inversion / eversion =  $5^\circ$ .

The "New Patient" function (home position) selects the following settings:

- Dorsal extension:  $20^\circ$
- Plantar flexion:  $-20^\circ$
- Inversion:  $-5^\circ$
- Eversion:  $15^\circ$
- Pauses: 0
- Timer: continuous operation
- Speed: 50 %
- Reverse on load - motor A: 25
- Reverse on load - motor B: 25
- Motor A: ON

- Motor B: ON
- Synchronized mode: enabled
- Total therapy time: 0
- Special functions disabled

### LEVEL 3:

#### ■ Reverse on load - motor A (feature for patient safety)

In the **synchronized mode**, the carriage automatically moves to the middle position when the patient's resistance (load) exceeds the set value.

In the **non-synchronized mode**, the carriage automatically moves in the opposite direction when the patient's resistance (load) exceeds the set value.

Adjustable levels for reverse on load feature: 1 - 25. At level 1, very low resistance will cause the device to reverse; at level 25, a high resistance is required to initiate the reversal.

**Default setting: Level 25**

#### ■ Reverse on load - motor B (feature for patient safety)

In the **synchronized mode**, the carriage automatically moves to the middle position when the patient's resistance (load) exceeds the set value.

In the **non-synchronized mode**, the carriage automatically moves in the opposite direction when the patient's resistance (load) exceeds the set value.

Adjustable levels for reverse on load feature: 1 - 25. At level 1, very low resistance will cause the device to reverse; at level 25, a high resistance is required to initiate the reversal.

**Default setting: Level 25**

### Caution!

Patient hazard –

The reverse-on-load feature is only a safety measure in the event of cramps, spasms, locked joints and similar conditions. The manufacturer cannot be held liable for misuse of this feature.

#### ■ Motor A ON/OFF A

To permit fully isolated movements, the motors can individually be turned on and off. Motor A performs dorsal extension and plantar flexion movements motor B performs eversion and inversion.

For an isolated eversion / inversion movement, set motor A to the appropriate position (dorsal extension / plantar flexion), then turn the motor off.


In the normal operating mode, the display will show the symbol for the disabled  dorsal / plantar movement .

**Default setting: motor A ON**

#### ■ Motor B ON/OFF B

To permit fully isolated movements, the motors can individually be turned on and off. Motor A performs dorsal extension and plantar flexion movements motor B performs eversion and inversion.

For an isolated dorsal extension / plantar flexion movement, set motor B to the appropriate position (inversion / eversion), then turn the motor off.

In the normal operating mode, the display will show the symbol for the inversion / eversion movement  disabled .

**Default setting: motor B ON**

### Note!

Please note that one motor (A or B) must always be turned on. Otherwise the display will show the message:  OFF 

### LEVEL 4:

#### ■ Transport setting

With this function, the device will move to a position optimally suited for packing. Select the function and press the START button. The carriage moves to the transport position. In the top left-hand corner of the display you will see  (also refer to chapter 6 Transport)

#### ■ Synchronized/non-synchronized mode A+B

Motors A and B can run in synchronized or non-synchronized mode.

**Synchronized:**

Motors A and B perform a synchronized movement according to the ankle joint's physiological motion pattern:

From the middle position the carriage first goes to the maximum eversion value. Afterwards, the carriage will simultaneously move to the maximum inversion value and to the maximum plantar flexion value, and then simultaneously to the maximum dorsal extension value and to the value halfway between inversion and eversion. After reaching this position, the cycle restarts: maximum eversion at the same time as halfway between dorsal extension / plantar flexion, etc.

During operation the synchronized mode is indicated in the top left-hand corner by the symbol .

**non-synchronized:**

Both motors run independently of each other within the selected ranges of motion.

To enable the non-synchronized mode, disable the synchronized mode.

During operation the non-synchronized mode is indicated in the top left-hand corner by the symbol .

**Default setting: synchronized mode enabled**

**⚠ Caution!**

Patient hazard –

The synchronized operation is recommended as a general rule. However, medical and/or therapeutic indications may call for the non-synchronized mode. In the non-synchronized mode, particular caution and attention on the part of the physician/therapist is required to rule out any risk for the patient.

**■ Total therapy time  ****ARTROMOT®-SP3 models without chip card**

The total therapy time is the added sum of operating hours. If the device is used by only one patient, this time is equivalent to the duration of all the patient's therapy sessions.

Under menu item "total therapy time" of **ARTROMOT®-SP3 models with chip card**, you can view each patient's total

therapy time (duration of all the patient's therapy sessions).

**Deleting the stored therapy time**

Press and hold the parameter button for 5 seconds or select the New Patient function.

**■ Service MENU **

For service purposes only, refer to Service Manual.

**Reminder:**

You save the set parameter values by pressing the **STOP** button.

**5.4 Programming the ARTROMOT®-SP3 Comfort models**

**ARTROMOT®-SP3 Comfort** models offer two more programming levels for additional functions.

The programming levels are selected in the same way as on Standard models.

Programming levels 1, 2, 3 and 6 are identical with programming levels 1 to 4 of the Standard models.

All special functions are disabled in the default settings.

The following **treatment values, settings** and **information** can be entered/viewed on the control pendant (22):

**LEVEL 1:**


- Dorsal extension 
- Plantar flexion 
- Inversion 
- Eversion 

**MENU****LEVEL 2:**

- Pause 
- Timer (therapy time) 
- Speed 
- New patient 

### LEVEL 3:

- Reverse on load - motor A
- Reverse on load - motor B
- Motor A ON/OFF
- Motor B ON/OFF

### LEVEL 4:

- Warm-up protocol
- Isolation protocol
- Dorsal extension / plantar flexion therapy documentation
- Inversion / eversion therapy documentation


### LEVEL 5:

- Stretching in dorsal extension
- Stretching in plantar flexion
- Stretching in inversion
- Stretching in eversion

### LEVEL 6:

- Transport setting
- Synchronized/ non-synchronized mode
- Total therapy time
- Service menu

## 5.5 Protocol details – Comfort model

- To access the appropriate programming level, press the **MENU button** repeatedly.
- The treatment parameters are selected with the different **parameter buttons**.
- To change the treatment values, press the **+ / – buttons**.
- To enable/disable a function, press the respective **parameter button** again.
- To save the entered values, press the **STOP button**.

All **special functions** are disabled in the **default settings**.

**Note!**

**LEVEL 1:** equivalent to level 1 of the Standard model (see: 5.3)

**LEVEL 2:** equivalent to level 2 of the Standard model (see: 5.3)

**LEVEL 3:** equivalent to level 3 of the Standard model (see: 5.3)

**LEVEL 6:** equivalent to level 4 of the Standard model (see: 5.3)

### LEVEL 4:

#### ■ Warm-up protocol


During warm-up, the patient will slowly become used to the set maximum dorsal extension / plantar flexion and inversion / eversion values, starting from the middle position.

The warm-up protocol starts in the middle between the two maximum values set for dorsal extension / plantar flexion and inversion / eversion.

In the beginning, the range of inversion / eversion motion increases with each cycle, until the programmed maximum values are reached after a total of 15 cycles. In this phase, motor A (dorsal extension / plantar flexion) is inactive.

Then the carriage stops again in the middle position and starts increasing the range of dorsal extension / plantar flexion motion with each cycle, until the programmed maximum values are reached after a total of 15 cycles. In this phase, motor B (inversion / eversion) is inactive.

After warm-up the device enters the normal mode.

During warm-up, the display will show the symbol  in the upper, left-hand corner.

**Default setting: disabled**

## ■ Isolation protocol AIIIB

For this special function, both motors are turned on, but they will never move simultaneously.

These are the steps of the special function:

- At first motor B moves to the programmed maximum values for dorsal extension and plantar flexion for 10 cycles, then it stops. In this phase, motor A is disabled (symbol for motor A: )
- Then motor A move to the programmed maximum values for dorsal extension and plantar flexion for 10 cycles and stops. The full cycle will start over with motor A. Motor B disabled while motor A is operating (symbol for motor B: )
- The stop position at the end of each tenth cycle (for both dorsal extension / plantar flexion and inversion / eversion) can be selected in steps of 25% between 0% and 100% of the programmed maximum range of motion.
- Steps 1 and 2 can be repeated as often as needed. You can stop treatment with the STOP button; after the programmed therapy duration, it will stop automatically.

The isolation protocol of the CPM device is identified with the  symbol in the upper, left-hand corner of the display.


**Default setting: disabled**

Follow these steps to program the special function:

- Access menu level 4, then press the parameter button to select the special function AIIIB .
- The display will show the symbol for the special function, its status (circle with or without check mark) as well as two boxes with the selected percentages for motor A and motor B stop positions. The circle appears in a box.
- Activate the function by pressing the "+" button (plus). Deactivate the

function by pressing the "-" button (minus). In this case it is not possible to activate/deactivate the function by pressing the parameter button again.

- A check mark appears in the circle which means that the function is active.
- Then press the parameter button for the isolation protocol to select the upper numeric field where you enter the percentage value for the stop position of motor A.
- To modify the value, press the "+" (plus) or "-" (minus) button. Select 75%, for example.
- Press the parameter button again to select the lower numeric field where you enter the stop position for motor B.
- To modify the value, press the "+" (plus) or "-" (minus) button. Select 25%, for example.
- Save the programmed values by pressing the STOP button, then press the START button to begin the treatment.

### This is what the settings mean:

Dorsal extension / plantar flexion is treated ten times first (motor A).

Motor B is positioned at 25% of the programmed range of motion for inversion and eversion.

Subsequently, inversion / eversion is treated ten times (motor B). Motor A is positioned at 75% of the programmed range of motion for dorsal extension and plantar flexion.

### Note!

The percentage values can be modified only when the special function is enabled (check mark).

### ■ Dorsal extension / plantar flexion therapy documentation

**ARTROMOT®-SP3 Comfort** devices with patient chip card have a special documentation function which provides a log of all therapy sessions.

The device run times as well as the range of motion of the sessions are recorded.

The collected data are presented graphically in the form of two curves in a coordinate system (X-axis = range of motion / Y-axis = time) where the upper curve illustrates the trend of the dorsal extension movement and the lower curve the trend of the plantar flexion movement.

### ■ Inversion / eversion therapy documentation

**ARTROMOT®-SP3 Comfort** devices with chip card have a special documentation function which provides a log of all therapy sessions.

The device run times as well as the range of motion of the sessions are recorded.

The collected data are presented graphically in the form of a coordinate system (X-axis = range of motion / Y-axis = time) where the upper curve illustrates the trend of the eversion movement and the lower curve the trend of the inversion movement.

## LEVEL 5:

### ■ Stretching in dorsal extension

With the special "stretching in dorsal extension" function the carriage gently stretches the joint by bending the foot upward.

Starting at the middle position the carriage will first move to the programmed plantar flexion value and then to the programmed dorsal extension value.

Subsequently the carriage reverses 5° towards plantar flexion and then moves very slowly back again to the programmed dorsal extension value (display ). After that it attempts to attain another 5°, moving even slower than before (display .

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the car-

riage moves in the opposite direction.

This stretch cycle is repeated ten times.

After that the carriage moves to the programmed maximum plantar flexion value and restarts the cycle "stretching in dorsal extension".

**Default setting: disabled**

### Note!

- If a pause has been programmed, the carriage will stop for the pause each time the maximum stretching value is attained.
- Enabling the stretching in dorsal extension function automatically disables motor B. The only movement performed is dorsal extension / plantar flexion.
- It is not possible to combine a second stretching function with the special "stretching in dorsal extension" function.
- During operation, the display will show the symbol  in the upper, left-hand corner.

### ■ Stretching in plantar flexion

With the special "stretching in plantar flexion" function the carriage gently stretches the joint by bending the foot downward.

Starting at the middle position the carriage will first move to the programmed dorsal extension value and then to the programmed plantar flexion value.

Subsequently the carriage reverses 5° towards dorsal extension and then moves very slowly back again to the programmed plantar flexion value (display ). After that it attempts to attain another 5°, moving even slower than before (display .

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the carriage moves in the opposite direction.

This stretch cycle is repeated ten times.

After that the carriage moves to the programmed maximum dorsal extension value and restarts the cycle "stretching in plantar flexion".

**Default setting: disabled**




## Note!

- If a pause has been programmed, the carriage will stop for the pause each time the maximum stretching value is attained.
- Enabling the "stretching in plantar flexion" function automatically disables motor B. The only movement performed is dorsal extension/plantar flexion.
- It is not possible to combine a second stretching function with the special "stretching in plantar flexion" function.
- During operation, the display will show the symbol  in the upper, left-hand corner.

### ■ Stretching in inversion

With the special "stretching in inversion" function the carriage gently stretches the joint towards the inside.

Starting at the middle position the carriage will first move to the programmed eversion value and then to the programmed inversion value.

Subsequently the carriage reverses 5° towards eversion and then moves very slowly back again to the programmed inversion value (display:  for the right foot and  for the left foot). After that it attempts to attain another 5°, moving even slower than before (display:  for the right foot and  for the left foot).

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the carriage moves in the opposite direction.

This stretch cycle is repeated ten times.

After that the carriage moves to the programmed maximum eversion value and restarts the cycle "stretching in inversion".

**Default setting: disabled**

## Note!




- If a pause has been programmed, the carriage will stop for the pause each time the maximum stretching value is attained.
- Enabling the "stretching in inversion" function automatically disables motor A. The only movement performed is inversion / eversion.

- It is not possible to combine a second stretching function with the special "stretching in inversion" function.
- During operation, the display will show the symbol  in the upper, left-hand corner.

### ■ Stretching in eversion

With the special "stretching in eversion" function the carriage gently stretches the joint towards the outside.

Starting at the middle position the carriage will first move to the programmed inversion value and then to the programmed eversion value.

Subsequently the carriage reverses 5° towards inversion and then moves very slowly back again to the programmed eversion value (display:  for the right foot and  for the left foot). After that it attempts to attain another 5°, moving even slower than before (display:  for the right foot and  for the left foot).

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the carriage moves in the opposite direction.

This stretch cycle is repeated ten times.

After that the carriage moves to the programmed maximum inversion value and restarts the cycle "stretching in eversion".

**Default setting: disabled**

## Note!

- If a pause has been programmed, the carriage will stop for the pause each time the maximum stretching value is attained.
- Enabling the "stretching in eversion" function automatically disables motor A. The only movement performed is inversion / eversion.
- It is not possible to combine a second stretching function with the special "stretching in eversion" function.
- During operation, the display will show the symbol  in the upper, left-hand corner.

### Reminder:

You save the selected parameter values by pressing the **STOP** button.

## 6. Care, Maintenance, Transport, Conversion

### 6.1 Care / re-use

The **ARTROMOT®-SP3** is suitable for re-use if the following points are observed.

#### **Warning!**

Shock hazard –

Unplug the device from the power line before cleaning.

Shock hazard, equipment damage –

Liquids must not enter the device or the control pendant.

- The **ARTROMOT®-SP3** must be **disinfected by wiping down** with a disinfectant. Thus, it complies with the special hygiene standards for medical technical equipment.
- The **enclosure** and removable **leg support assemblies** can be cleaned with commonly used **disinfectants** and mild **household detergents**.
- Use only a **damp cloth** to wipe the CPM device down.

#### **Warning!**

Patient hazard –

Patient contamination

- Before using the device on another patient, be sure to clean and disinfect it according to the instructions given here.

#### **Note!**

The manufacturer recommends using only special disinfectants approved for medical devices and with the characteristics specified below under "Caution".

Disinfection with the following disinfectants and in compliance with the instructions for use has been evaluated and approved by the manufacturer in the risk assessment:

- DESCOSEPT AF lemon  
(Part no.: 00-311L-xxx)  
BAuA Reg.No.: N-55153, C€-0482  
Manufacturer: Dr. Schumacher GmbH  
Contact time: 2 minutes minimum

#### **Caution!**

Equipment damage –

- The plastic material used is not resistant to mineral acids, formic acid, phenols, cresols, oxidants and strong organic or inorganic acids with a pH value below 4.
- Use only clear disinfectants to prevent discoloration of the material.
- Do not expose the CPM device to strong ultraviolet radiation (sunlight) or fire.

## 6.2 Maintenance (fuse replacement)

### Check before each use

Visually inspect the device for signs of mechanical damage before each use.

If you detect damage or malfunctions that may impair the safety of the patient or of the operator, have the device repaired before using it.

### Technical inspections

Due to their design, the bearings and joints of the **ARTROMOT®-SP3** do not require maintenance and all materials are protected from corrosion, BUT: For safety, the devices require regular maintenance. To maintain the functional and operational safety, check all components for damage and loose connections at least **once a year**.

These checks should be performed by persons with adequate training and experience. Damaged and worn parts must immediately be replaced with original spare parts by authorized staff.

### **Warning!**

Patient hazard, equipment malfunction and damage –

- Refer repair and maintenance to authorized persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.

ORMED GmbH will make all documents required for servicing, such as circuit diagrams, parts lists, descriptions or calibration instructions, available to authorized experts.

The inspections can be carried out by DJO Technical Service within the framework of a service agreement. Contact DJO for details.

**Other than that, the manufacturer does not require any other regular maintenance.**

## Note!

Regarding further technical or other inspections or inspection intervals, observe applicable local requirements, such as IEC 62353, DGUV3 or similar regulations to be observed by users of medical equipment or electric devices.

## Fuse replacement

### **⚠ Warning!**

Patient hazard, equipment malfunction and damage –

The replacement of fuses must be referred to specialists as defined in IEC 60364 or other applicable standards (e.g. biomedical technicians, electricians, electronics installers).

#### **Fuses used must be**

**S# < 20,000: T1A L250Vac**

**S# >20,000: T2A H250Vac.**

Before replacing fuses, turn off the ARTRO-MOT®-SP3 and disconnect the device from the power line.

Use an appropriate tool to remove the fuse holder situated between the power switch and the power connector (Fig. 1). Replace the fuses and reinsert the fuse holder (Fig. 2). Ensure that the fuse holder properly locks into place.


Fig. 1


Fig. 2

## 6.3 Transport

The following adjustments need to be carried out to prepare the ARTROMOT®-SP3 for transport:

1. Activate the "transport setting" function →  in the menu (also refer to 5.3) and start the ARTROMOT®-SP3.
2. Remove the footplate (7) as well as its bracket (11). To do so, simply pull it from the square-end guide at motor B. (Fig. G)
3. For transport of the device, motor B must be mounted on the right side of the moving bracket (15) (marked R). Refer to section 6.4 Conversion, if required.
4. Turn off the ARTROMOT®-SP3 power switch.
5. Disconnect the power cord and the control pendant.
6. Relieve the strain on the carriage by lifting it a little at the moving bracket for dorsal extension / plantar flexion (15) and open the clamping lever which adjusts the inclination of the carriage (5). Select an angle where the carriage is parallel with the adjusted height. Tighten the clamping lever to lock the carriage in this position.

7. Again lift the moving bracket for dorsal extension / plantar flexion (15) a little to remove the strain. Then loosen the locking device for the height adjustment by simultaneously pressing the two black pins of the locking device together (19). Select a height where the locking pins engage in the top hole (marked: Transport).

8. Only use the original shipping box to transport the device. Ormed GmbH cannot be held liable for transport damage, if the device is not shipped in its original shipping box.
9. Put the control pendant (22) in the supplied small box.
10. Place the footplate bracket (11) with the footplate (7) and the power cord in the extra boxes provided.
11. Place the ARTROMOT®-SP3 on the bottom of the shipping box. Then pack the two extra boxes for the control pendant and the footplate + bracket in the shipping box along with the padding material, as shown below.

Extra box with control pendant

Extra box with footplate/bracket and power cord


## 6.4 Conversion

The ARTROMOT®-SP3 can be used for the left and for the right ankle joint. However, the device must be converted first. This is done very quickly.

1. Activate the "transport setting" function →  in the menu (also refer to 5.3) and start the ARTROMOT®-SP3.
2. Remove the footplate (7) as well as its bracket (11). To do so, pull it from the square-end guide at motor B. (Fig. G)
3. Loosen thumbscrew (12) used to mount motor B to the moving bracket (15). Then pull the motor upward and remove (Fig. H).
4. Mount motor B on the opposite side of the carriage, using the thumbscrew.

5. Push the footplate bracket including the footplate into the square-end aperture on motor B. Observe the marks: red dot on square rod of the footplate bracket to match red dot on square-end aperture on motor B, mark L or R on the footplate bracket to match mark L or R on the motor housing.


### Note!

- To assist in the conversion, motor B is equipped with a side identification function.

The following symbols may appear on the display during conversion: (example)

conversion from left to right:

Motor on the left-hand side: 

Motor removed from bracket: 

Conversion to the right-hand side completed: 

- If, after the conversion, the footplate bracket, i.e., the footplate is at an inadmissible angle, the display shows:

 **START**

Please press the START button and the carriage will automatically move to the permitted range and stop.

## 7. Environmental protection statement

The product described in this operation manual must not be disposed of with unsorted household or municipal waste. It requires separate disposal. Please contact DJO for information about the possible recycling of the product.

The expected service life of the device and all supplied parts and accessories is 6 years. The application beyond this period is the responsibility of the user.

## 8. Specifications

<b>Model:</b>	ARTROMOT®-SP3
<b>Device name:</b>	ARTROMOT®-SP3 Standard, Part no.: 80.00.035 ARTROMOT®-SP3 Standard Chip, Part no.: 80.00.036 ARTROMOT®-SP3 Comfort, Part no.: 80.00.037 ARTROMOT®-SP3 Comfort Chip, Part no.: 80.00.038
<b>Input ratings:</b>	100 – 240 Vac / 50 – 60 Hz tolerance -15% to +10%
<b>Current consumption:</b>	100V      240V
Standby (ON):	5VA      5VA
Operation (maximum):	33VA      33VA 330mA    140mA
<b>Fuses:</b>	serial numbers < 20,000: T1A L250 Vac serial numbers >20,000: T2A H250 Vac according to IEC 60127-2/3, UL 248-14, CSA C22.2 no. 248.14
<b>Protection class:</b>	serial numbers < 20,000: I serial numbers > 20,000: II
<b>IP class:</b>	serial numbers < 20,000: IPX0 serial numbers > 20,000: IP21
<b>Applied part:</b>	Type B
<b>Max. load on carriage:</b>	20 kg
<b>Dimensions (transport)</b>	
length:	78.0 cm
width:	42.0 cm
height:	39.5 cm

<b>Accuracy of measured values</b>	
Goniometer in the measuring range:	from -50° to +40° in plantar flexion/ dorsal extension
Accuracy:	+/- 2°
<b>Weight:</b>	11 kg
<b>Materials used:</b>	ABS, POM, PUR, PA, FR 4, aluminum, stainless steel, brass
<b>MDD:</b>	class IIa
<b>Standards compliance:</b>	93/42/EEC IEC 60601-1 IEC 60601-1-6 IEC 60601-1-9 IEC 60601-1-11 (S# >20,000) IEC 62366 IEC 62304 EN ISO 14971 ANSI/UL 60601-1 (S# <20,000) CAN/CSA C22.2 No. 601.1-M90 (S# <20,000) ANSI AAMI ES 60601-1 1st Edition (S# >20,000) CAN CSA 22.2 No. 60601-1-08 (S# >20,000)
<b>EMC</b> (electromagnetic compatibility)	IEC 60601-1-2
<b>Manufactured in compliance with:</b>	EN ISO 13485
<b>Ambient conditions (storage, transport)</b>	
Temperature:	-25 ° to +70 °
Relative humidity:	at 70 °C up to 93%, no condensation
Atmospheric pressure:	500 hPa to 1060 hPa
<b>Ambient conditions (operation)</b>	
Temperature:	+5 °C to +40 °C
Relative humidity:	15 % to 93 %
Atmospheric pressure:	700 hPa to 1060 hPa

## 9. Contact

We would be happy to answer any questions you may have about our products and services.

### **ARTROMOT® International**

Please contact your local dealer, the DJO Headquarters in the USA, the DJO International Headquarters in England, or DJO Germany.

### **DJO International**

Headquarters: DJO UK Ltd.  
1a Guildford Business Park  
Guildford  
Surrey, GU2 8XG, United Kingdom  
Telephone: +44 (0) 1483 459659  
Fax: +44 (0) 1483 459470  
Email: info@DJOglobal.eu  
Website: www.DJOglobal.eu

### **DJO Headquarters**

DJO, LLC  
2900 Lake Vista Drive  
Lewisville, TX 75067, USA  
Telephone: +1 800 321 9549  
Fax: +1 800 936 6569  
Email: webmaster@DJOglobal.com  
Website: www.DJOglobal.com

### **Manufacturer / Headquarters Germany**

ORMED GmbH  
A company of Enovis  
Bötzingen Straße 90  
79111 Freiburg, Germany

Telephone: +49 (0) 1805 167 633  
Fax: +49 (0) 1805 367 633  
Email: kundenservice@enovis.com  
Website: www.enovis-medtech.de

### **Quarters Austria**

DJO  
ORTHOMED Medizintechnik GmbH  
Annagasse 5/1/Top 4  
1010 Vienna, Austria  
Telephone: +43 1 5320834  
Fax: +43 1 532083431  
Email: orthomed@orthomed.cc  
Website: www.orthomed.cc

### **Warranty:**

2 years (mechanical parts)  
2 years (electronics)

## 10. Technical service

### 10.1 Technical hotline

Do you have any technical questions?  
Do you need technical service?

Telephone: +49-180-5-1 ormed de  
+49-180-5-1 67 63 33  
Facsimile: +49-180-5-3 ormed de  
+49-180-5-3 67 63 33

### 10.2 Shipment

To prevent damage during transport, only use the original shipping box. These boxes can be obtained from DJO.

Before packing the therapy unit, set it to the transport position (see chapter 6.3).

### 10.3 Spare parts

Refer to the Service Manual for the most recent list of spare parts.

When ordering spare parts, always specify:

- Item
- Description
- Part number
- Quantity
- Serial number of the device

#### Note!

Refer repairs to authorized, specially trained staff.

DJO offers service training for your staff.

#### Note!

Country-specific power cords are available from DJO. Please contact DJO or your DJO distributor.

Surcharges may apply in certain cases to spare parts ordered in low quantities.

Item	Description	Part no.	Qty
1.	Patient chip card	0.0034.035	<input type="checkbox"/>
2.	Marker pen for patient chip card	0.0031.006	<input type="checkbox"/>
3.	Power cord EU version H05VV-F3G 3 × 1 mm <sup>2</sup> length: 4 m	0.0034.118	<input type="checkbox"/>

## 11. CE-Declaration of Conformity

The certificate can be download from:

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<https://www.enovis-medtech.de/unternehmen/qualitaet-und-verantwortung>

## 12. IEC 60601-1-2 and IEC 60601-1-11

For **ARTROMOT®-SP3** devices with serial numbers > #25.000.

The **ARTROMOT®-SP3** device is subject to particular precautions regarding electromagnetic compatibility (EMC). The device must be installed and put into service strictly in compliance with the EMC directives set forth in the accompanying documents.

Portable and mobile RF communication systems may affect the **ARTROMOT®-SP3**.

The **ARTROMOT®-SP3** should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, **ARTROMOT®-SP3** should be observed to verify normal operation in the configuration in which it will be used.

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If it is necessary to replace assemblies or cables, only the manufacturer's original parts must be used to ensure continued compliance with EMC requirements after repair. This requirement applies to the power supply unit, cables and cable lengths, drive unit consisting of the motor and the control system, the programming unit incl. the coiled cable and the connector.

Regarding the electromagnetic compatibility, there is no maintenance work required during the expected service life of the product.

The product designation **ARTROMOT®-SP3** used in the text below includes all product variants.

### **WARNING!**

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

## 12.1 Electromagnetic emissions

### Guidance and Manufacturer's Declaration – Electromagnetic Emissions

**ARTROMOT®-SP3** is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-SP3** unit is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions to CISPR 11	Group 1	<b>ARTROMOT®-SP3</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions to CISPR 11	Class B	<b>ARTROMOT®-SP3</b> is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions to IEC 61000-3-2	complies with class A requirements	
Voltage fluctuations/flicker emissions to IEC 61 000-3-3	complies	

## 12.2 Electromagnetic immunity

### Guidance and Manufacturer's Declaration – Electromagnetic Immunity

**ARTROMOT®-SP3** is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-SP3** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment – guidance
Electrostatic discharge (ESD) to IEC 61 000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/bursts to IEC 61 000-4-4	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power should be that of a typical commercial or hospital environment.
Surges to IEC 61 000-4-5	± 1 kV line-to-line voltage	± 1 kV line-to-line voltage	Mains power should be that of a typical commercial or hospital environment.
Voltage dips To IEC 61000-4-11	0% $U_T$ ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°  0% $U_T$ ; 1 cycle and 70% $U_T$ ; 25/30 cycles  Single phase at 0°	0% $U_T$ ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°  0% $U_T$ ; 1 cycle and 70% $U_T$ ; 25/30 cycles  Single phase at 0°	Mains power should be that of a typical commercial or hospital environment. If the user of the <b>ARTROMOT®-SP3</b> unit requires continued operation during power mains interruptions, it is recommended that the <b>ARTROMOT®-SP3</b> unit be powered from an uninterruptible power supply or a battery.
Voltage interruptions to IEC 61000-4-11	0% $U_T$ ; 250/300 cycles	0% $U_T$ ; 250/300 cycles	
Power frequency (50/60 Hz) magnetic field to IEC 61 000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.

**Note:**  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and Manufacturer's Declaration – Electromagnetic Immunity

**ARTROMOT®-SP3** is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-SP3** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment – guidance
Conducted RF to IEC 61000-4-6	6 V 0,15 MHz-80 MHz 80% AM at 1 kHz	6 V 0,15 MHz-80 MHz 80% AM at 1 kHz	
Radiated RF to IEC 61000-4-3	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communications equipment to IEC 61000-4-3	27 V/m 380-390 MHz PM 18 Hz  28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine  9 V/m 704-787 MHz PM 217 Hz  28 V/m 800-960 MHz PM 18 Hz  28 V/m 1,7-1,99 GHz PM 217 Hz  28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz  9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	27 V/m 380-390 MHz PM 18 Hz  28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine  9 V/m 704-787 MHz PM 217 Hz  28 V/m 800-960 MHz PM 18 Hz  28 V/m 1,7-1,99 GHz PM 217 Hz  28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz  9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ARTROMOT®-SP3** exceeds the applicable **ARTROMOT®-SP3** unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **ARTROMOT®-SP3** unit.

## 12.3 RF (Radio Frequency) Interference

EMC disturbance can affect the performance of the **ARTROMOT®-SP3**. To prevent the user from unacceptable risks the device runs performance checks regularly during operation. This is what happens, if a problem is identified:

- an audio signal sounds
- the unit stops immediately
- the message ERROR and an error code appear on the display

In this case try to quit the error with the STOP key. If the error persists you may attempt to restart the unit by turning it briefly off and on with the power switches.

### **WARNING!**

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the **ARTROMOT®-SP3**, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

### **Notice!**

(for home use)

Portable and mobile RF communication systems may affect the **ARTROMOT®-SP3** device.

For this reason, make sure that wireless communication equipment, such as wireless home network devices, cell phones, cordless telephones and their base stations and walkie-talkies operate at a minimum distance of 3,3 m from the device. (Calculated on the basis of the maximum power output of a typical cell phone of 2 W).

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### Symbol Overview standard models ARTROMOT®-SP3

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## Symbol Overview standard models ARTROMOT®-SP3

### EBENE 1/LEVEL 1:

Dorsalextension  
dorsal extension

Plantarflexion  
plantar flexion

Inversion  
inversion

Eversion  
eversion

### EBENE 2/LEVEL 2:

Pause  
pause

Timer  
timer

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Neuer Patient  
new patient

### EBENE 3/LEVEL 3:

Lastumkehr Motor A  
load reversal motor A

Lastumkehr Motor B  
load reversal motor B

Motor A Ein/Aus  
motor A ON/OFF

Motor B Ein/Aus  
motor B ON/OFF

### EBENE 4/LEVEL 4:

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Isolationsprogramm  
isolation protocol

Therapieverlaufsdokumentation  
Dorsal-/  
Plantarflexion  
dorsal extension /  
plantar flexion  
therapy  
documentation

Therapieverlaufsdokumentation  
Inversion/Eversion  
inversion / eversion  
therapy  
documentation

### EBENE 5/LEVEL 5:

Dehnung  
Dorsalextension  
stretching in  
dorsalextension

Dehnung  
Plantarflexion  
stretching in  
plantar flexion

Dehnung Inversion  
stretching in inversion

Dehnung Eversion  
stretching in eversion

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





## AMENDMENT TO THE INSTRUCTION FOR USE (IFU)

This document outlines the amendment to the Instruction for Use (IFU) for necessary updates and changes.

### List of the affected devices:

Reference Number	Article Name	GTIN
80.00.040	<b>ARTROMOT®</b> -K1 Standard	00888912403191
80.00.041	<b>ARTROMOT®</b> -K1 Standard Chip	00888912403207
80.00.042	<b>ARTROMOT®</b> -K1 Comfort	00888912403214
80.00.043	<b>ARTROMOT®</b> -K1 Comfort Chip	00888912403221
80.00.044	<b>ARTROMOT®</b> -K1 Standard Basic	00190446751823
80.00.045	<b>ARTROMOT®</b> -K1 Classic	00888912403238
80.00.047	<b>ARTROMOT®</b> -K1 Classic Basic	00190446751830
80.00.023	<b>ARTROMOT®</b> -S3 Standard	00888912403108
80.00.024	<b>ARTROMOT®</b> -S3 Comfort	00888912403115
80.00.033	<b>ARTROMOT®</b> -E2 Compact	00888912403146
80.00.031	<b>ARTROMOT®</b> -E2	00888912403139
80.00.035	<b>ARTROMOT®</b> -SP3 Standard	00888912403153
80.00.036	<b>ARTROMOT®</b> -SP3 Standard Chip	00888912403160
80.00.037	<b>ARTROMOT®</b> -SP3 Comfort	00888912403177
80.00.038	<b>ARTROMOT®</b> -SP3 Comfort Chip	00888912403184
80.00.070	<b>ARTROMOT®</b> ACTIV E-K	00888912403283
80.00.084	<b>ARTROMOT®</b> -S4 Comfort	00190446278764
80.00.022	OptiFlex®-S	00190446794783
80.00.046	OptiFlex®-K1 Classic	00190446796015
80.00.048	OptiFlex®-K1 Standard	00190446796039
80.00.049	OptiFlex®-K1 Comfort Chip	00190446796046

### First change: Warning

New warnings have been implemented due to labeling changes that add or strengthen a warning, precaution, or include additional information.

Safety warnings:

#### **Warning!**

"The treatment of minors or persons with limited physical, sensory or mental abilities (e.g. movement disorders, epilepsy, seizures, mental disorders, behavioural disorders or neurodevelopmental disorders) may only be carried out under the constant observation of attendants (e.g. specialist staff, legal guardians or custodians)

– After data storage, write the patient's name on the patient chip card. The card should only be used for this particular patient. If the patient chip card is used for another patient, be sure to delete the previous patient data from the card first. Use original chip cards only."

## Second change: Label

### Symbols

Table 1: Symbols changes

Symbol change	Old version	New version	Description of symbol acc. ISO 15223 + 2014/
UDI	New requirement		Indicates a carrier that contains unique device identifier information. This symbol identifies the UDI carrier, including the AIDC and human readable information. (01) GTIN code of device (11) Date of production (21) serial number of device
Manufacturer + Date of production + country of manufacture	 Name Address   2015	 DE YYYY-MM	Indicates the medical device manufacturer Indicates the date when the medical device was manufactured DE is the two or three letter country code as defined in ISO 3166-1. Indicates the country of manufacture
Manufacturer's name	<b>ORMED GmbH</b> A DJO Company	<b>ORMED GmbH</b> a company of Enovis	Name of the manufacturer
Manufacturer's address spelling	Böszinger Straße 90 79111 Freiburg / Germany	Böszinger Strasse 90 79111 Freiburg / Germany	Address of the manufacturer
Norm/Standard ETL classification	Conforms to ANSI/ AAMI ES STD 60601-1 Certified CAN/CSA C22.2 No. 60601-1	 InterTek 5029964 MEE	Medical electrical equipment relevant standards according ETL classification – control number above

### Other changes

Manufacturers contact was changed to: Tel. +49 761 4566-01 · Fax +49 761 4566-5501, medizintechnik@enovis.com, enovis-medtech.de

### Third change: Electromagnetic emissions, Electromagnetic immunity and RF (radio frequency) interference

The EMC table provides a detailed comparison of electromagnetic compatibility as per the latest standards.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
ARTROMOT® CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the ARTROMOT® CPM device unit is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions to CISPR 11: 2015+A1:2016+A2:2019	Group 1	<b>ARTROMOT®</b> CPM device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions to CISPR 11: 2015+A1:2016+A2:2019	Class B	<b>ARTROMOT®</b> CPM device is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions to IEC 61000-3-2: 2005+A1:2008 +A2:2009	complies with class A requirements	
Voltage fluctuations/flicker emissions to IEC 61000-3-3: 2013	complies	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

**ARTROMOT®** CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®** CPM device unit is used in such an environment.

Immunity test	IEC 60601-1-2:2014 +A1:2020		Electromagnetic environment – guidance
	test level Home healthcare	Compliance level Home healthcare	
Electrostatic discharge (ESD) to IEC 61000-4-2:2008	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/bursts to IEC 61000-4-4:2012	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power should be that of a typical commercial or hospital environment.
Surges to IEC 61000-4-5:2014	± 1 kV line-to-line voltage	± 1 kV line-to-line voltage	Mains power should be that of a typical commercial or hospital environment.
Voltage dips To IEC 61000-4-11:2017	0% Ut ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% Ut ; 1 cycle and 70% Ut ; 25/30 cycles Single phase at 0°	0% Ut ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% Ut ; 1 cycle and 70% Ut ; 25/30 cycles Single phase at 0°	Mains power should be that of a typical commercial or hospital environment. If the user of the <b>ARTROMOT®</b> CPM device unit requires continued operation during power mains interruptions, it is recommended that the <b>ARTROMOT®</b> CPM device unit be powered from an uninterruptible power supply or a battery.
Voltage interruptions to IEC 61000-4-11:2017	0% Ut; 250/300 cycles	0% Ut; 250/300 cycles	
Power frequency (50/60 Hz) magnetic field to IEC 61000-4-8:2009	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical commercial or hospital environment.

Note: Ut is the a.c. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

**ARTROMOT®** CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®** CPM device unit is used in such an environment.

Immunity test	IEC 60601-1-2:2014 +A1:2020		Electromagnetic environment – guidance
	test level Home healthcare	Compliance level Home healthcare	

Conducted RF to IEC 61000-4-6:2013	3 V 0.15 MHz-80 MHz 80% AM at 1kHz 6V 0.15 MHz-80 MHz 80% AM at 1kHz at ISM frequency bands and Amateur radio bands	3 V 0.15 MHz-80 MHz 80% AM at 1kHz 6V 0.15 MHz-80 MHz 80% AM at 1kHz at ISM frequency bands and Amateur radio bands	
Radiated RF to IEC 61000-4-3: 2006+A1:2007 +A2:2010	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communica- tions equipment to IEC 61000-4-3: 2006+A1:2007 +A2:2010	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	
Radiated fields in close proximity IEC 61000-4-39:2017	8 A/m 30 kHz No modulation (CW) 65 A/m 134.2 kHz PM 2.1 kHz 7,5 A/m 13.56 MHz PM 50 kHz	8 A/m 30 kHz No modulation (CW) 65 A/m 134.2 kHz PM 2.1 kHz 7,5 A/m 13.56 MHz PM 50 kHz	

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ARTROMOT®** CPM device exceeds the applicable **ARTROMOT®** CPM device unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **ARTROMOT®** CPM device unit.

The **ARTROMOT®** CPM movementsplint is intended for operation in the specified electromagnetic environment. The customer or user of the **ARTROMOT®** CPM device should ensure that it is used in such an environment. The required interference levels of the specified standards are met.

Hersteller/Manufacturer:  
ORMED GmbH  
A company of Enovis  
Bötzingen Straße 90 · 79111 Freiburg · Germany  
Tel. +49 761 4566 -01 · Fax +49 761 4566 -5501  
medizintechnik@enovis.com  
[www.enovis-medtech.de](http://www.enovis-medtech.de)