# ARTROMOT®-S3

SW:V3.X

- DE · Gebrauchsanleitung
- GB · Operating Instructions
- FR · Mode d'emploi
- ES · Instrucciones de uso
- IT · Istruzioni per l'uso
- NL · Gebruiksaanwijzing

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**Anhang:**  
ARTROMOT®-S3  
ARTROMOT®-S4  
ARTROMOT® ACTIVE-K  
ARTROMOT®-K1  
ARTROMOT®-SP3  
ARTROMOT®-E2

# Contents

<b>Device description ARTROMOT®-S3</b>	<b>III</b>
<b>Figures ARTROMOT®-S3</b>	<b>VIII</b>
<b>Symbol overview ARTROMOT®-S3</b>	<b>V</b>
<b>1. How to use the CPM device</b>	<b>43</b>
1.1 Fields of application	43
1.2 Therapy objectives	43
1.3 Indications	43
1.4 Contraindications	43
<b>2. ARTROMOT®-S3 description</b>	<b>44</b>
2.1 Description of the device components	45
2.2 Description of the programming unit	46
2.3 Explanation of symbols	49
2.4 Explanation of symbols (connections and nameplate)	50
<b>3. Safety information</b>	<b>52</b>
<b>4. Device setup</b>	<b>55</b>
4.1 Connecting the device, performance check	55
4.2 Adjusting the device to the patient	56
<b>5. Setting the treatment values</b>	<b>58</b>
5.1 General information on programming ARTROMOT®-S3	58
5.2 Programming ARTROMOT®-S3 Standard	60
5.3 Treatment value details – Standard models	60
5.4 Programming ARTROMOT®-S3 Comfort	64
5.5 Protocol details – Comfort model	65
5.6 Examples of use/sample protocols	69
<b>6. Care, Maintenance, Transport, Conversion</b>	<b>71</b>
6.1 Care/re-use	71
6.2 Maintenance (fuse replacement)	71
6.3 Transport	73
6.4 Conversion	74
<b>7. Environmental Protection Statement</b>	<b>75</b>
<b>8. Specifications</b>	<b>75</b>
<b>9. Contact</b>	<b>77</b>
<b>10. Technical Service</b>	<b>78</b>
10.1 Technical hotline	78
10.2 Shipment	78
10.3 Spare parts	78
<b>11. CE-Declaration of Conformity</b>	<b>79</b>
<b>12. IEC 60601-1-2 and IEC 60601-1-11</b>	<b>248</b>
12.1 Electromagnetic emissions	248
12.2 Electromagnetic immunity	249
12.3 Recommended separation distances	251

# 1. How to use the CPM device

## 1.1 Fields of application

**ARTROMOT®-S3** is a motor-operated Continuous Passive Motion (CPM) device providing motion to the shoulder joint.

Suitable for use in hospitals, clinics, general practices and rental services as well as in the patient's home, it is an important supplement to medical and therapeutic treatment.

## 1.2 Therapy objectives

CPM therapy with **ARTROMOT®-S3** is mainly used to prevent the negative effects of immobilization, to allow patients to regain painless mobility of joints at an early stage and to promote healing and achieve a positive functional result.

Other objectives of therapy include:

- improvement of joint metabolism
- prevention of joint stiffness (arthrofibrosis)
- promotion of the regeneration and healing of cartilage and damaged ligaments
- faster hematoma/fluid resorption
- improved lymph and blood circulation
- thrombosis and embolism prophylaxis

## 1.3 Indications

The CPM device is indicated in the treatment of most injuries and postoperative conditions and diseases of the shoulder joint.

Examples:

- joint distortion and contusion
- arthrotomy and arthroscopy procedures in combination with synovectomy, arthrolysis or other intra-articular interventions
- all types of arthroplasty
- mobilization of joints in anesthetized patients
- operative treatment of fractures, pseudoarthrosis, if exercise-stable
- decompression surgery (acromioplasty)
- endoprosthetic implants
- soft tissue surgery in the armpit and the shoulder girdle
- tumor surgery in the shoulder region

## 1.4 Contraindications

**Do NOT use ARTROMOT®-S3** on patients with:

- acute inflammatory processes in the joints, unless on the order of a physician
- spastic paralysis
- unstable osteosynthesis

## 2. ARTROMOT®-S3 description

The motorized CPM device supports the following movements of the shoulder joint:

**Adduction/Abduction** 0°-30°-175°

**Internal/External rotation** 90°-0°-90°

**Extension/Flexion (Elevation)** 0°-30°-175°  
with elbow bent between 60° and 90°

### Horizontal

**Adduction/ Abduction** 0°-120°  
(manual adjustment only)

It can be reconfigured for use on either side.

### Note!

To allow the current position of the CPM device to be clearly indicated, the internal rotation values on the display and throughout this document are identified with "-".

### Note!

Elevation denotes the upward movement of the arm beyond 90° abduction.

Extension and flexion are set by entering the adduction and abduction values

These are some of the outstanding **ARTROMOT®-S3** features:

- anatomically correct setup
- physiological movements
- maximum possible ranges of motion
- programming unit for precise adjustment of patient-specific therapy values
- chip card for storage of the programmed therapy parameters
- easy transport

### Biocompatibility

Those parts of the **ARTROMOT®-S3** device that come into contact with the patient when the device is used as intended, are designed to fulfil the biocompatibility requirements of the applicable standards.

### Essential Performance Features

- The programmed angles are accurately maintained with a tolerance of +/- 2°.
- The programmed speeds are accurately maintained with a tolerance of +/- 5%.
- The selected mode and mechanical settings do not change during operation.

### Frequent functions:

- Unpacking (device and accessories)
- Inserting the motion element and locking it at the system rack (chair)
- Mounting the armrest
- Device connections (mains power connection, connecting the control pendant and the motion element)
- Mechanical adjustments of the motion element on the rotational axes
- Programming the control pendant (range of motion, speed, functions)
- Storage

## 2.1 Description of the device components

**Note: Please fold out page III !**

1. Programming unit
2. Patient chip card
3. Compartment for storage of programming unit
4. Armrest for healthy arm
5. Wheels
6. Connection for power cord
7. Power switch (ON/OFF)
8. Instrument fuse
9. Connection for motion element
10. Thumbscrew for adjustment of the backrest angle
11. Folds down the backrest (transport position)
12. Adjustment for horizontal adduction/abduction
13. Holding tube for height adjustment
14. Locking screw for height adjustment
15. Insertion tube for height adjustment
16. Motor A
17. Holding tube for length adjustment (upper arm)
18. Clamping lever for length adjustment (upper arm)
19. Insertion tube for length adjustment (upper arm)
20. Locking screw for elbow range
21. Motor B
22. Elbow support
23. Forearm support
24. Strap for forearm restraint
25. Holding tube for length adjustment (forearm)
26. Clamping lever for length adjustment (forearm)
27. Insertion tube for length adjustment (forearm)
28. Locking screw for swivel mechanism
29. Handle
30. Headrest adjustment<sup>1</sup>
31. Height adjustment for shoulder restraint<sup>1</sup>
32. Height adjustment for shoulder restraint<sup>1</sup>
33. Ball joint for headrest<sup>1</sup>
34. Locking pin


Subject to technical modifications

---

<sup>1</sup> Only with **ARTROMOT®-S3 comfort**

## 2.2 Description of the programming unit

### 2.2.1 Programming unit in normal mode

## 2.2.2 Programming unit in MENU selection mode

## 2.2.3 Programming unit set to ROM programming

## 2.2.4 Programming unit in general programming mode

## 2.3 Explanation of symbols

Also refer to symbol overview on page V.

### Standard protocols:

	Adduction
	Abduction
	Internal rotation
	External rotation
	Pause
	Timer
	Speed
	New patient
	Reverse on load - motor A
	Reverse on load - motor B
	Motor A On/Off
	Motor B On/Off
	Transport setting
A+B	Synchronized/ non-synchronized mode
	Total therapy time
	Service menu

### Comfort protocols:

	Warm-up protocol
AIII B	Isolation program
	Therapy documentation Abduction/adduction
	Therapy documentation Internal/external rotation
	Oscillation
	Stretching in abduction
	Stretching in internal rotation
	Stretching in external rotation

## 2.4 Explanation of symbols (connections and nameplate)

Alternating current

For serial numbers < 20,000:  
protection class I equipment. The medical device must be connected to a system with protective earth conductor!

For serial numbers > 20,000:  
protection class II equipment. The medical device has a double or reinforced insulation.

Type B applied part

Power switch OFF

Power switch ON

The year and month of manufacture are indicated next to this factory symbol.

The name next to this factory symbol is the manufacturer.

The number next to this symbol is the article reference number.

Device complies with Council Directive 93/42/EEC about medical devices, tested and approved by DQS Medizinprodukte GmbH.

The number next to this symbol is the serial number.

Caution! Observe warnings set forth in operation manual!

Refer to instruction manual! ISO 7010-M002

Refer to instruction manual! ISO 7000-1641

Do not dispose with unsorted municipal waste.

Protect from moisture.

## IP21

The IP rating indicates the level of protection and thus the suitability of the device for use under different ambient conditions.

The rating IP21 means:

2 is the level of protection against contact and solid objects.

The digit 2 means:

- Protection from contact: protected from contact with a finger.
- Protection against foreign object: protected against solid foreign object (diameter of 12.5 mm and greater).

1 indicates the degree of protection against water.

The digit 1 means: protection against vertically falling water drops.

Observe temperature limits (storage)! ISO 7000-0632

The overall weight of the device is indicated next to this symbol.

Do not rest or lean on or push the assembled device because it might tip over.

## 3. Safety information

### Definitions

Read the safety statements before use of the CPM device. The safety statements are classified as follows:

#### **Danger!**

This term indicates an imminent hazard. If not avoided, this hazard will result in death or serious injury.

#### **Warning!**

This term indicates a hazard. If not avoided, this hazard can result in death or serious injury.

#### **Caution!**

Indicates a potential hazard. If not avoided, this hazard can result in minor personal injury and/or product/property damage.

### Safety information

#### **Danger!**

Explosion hazard –

**ARTROMOT®-S3** is not designed for use in areas of medical locations where an explosion hazard may occur. An explosion hazard may result from the use of flammable anesthetics, skin cleansing agents and disinfectants.

#### **Warning!**

Patient hazard –

- Only **authorized individuals** are allowed to operate the **ARTROMOT®-S3** device. Individuals are authorized after receiving training in the operation of the device and reading this operation manual.
- Before using the device, the operator must ascertain that it is in correct working order and operating condition. In particular, the cables and connectors must be checked for signs of damage. Damaged parts must be replaced immediately, before use.
- **Before therapy**, a **test run** consisting of several exercise cycles must be completed, first without and then with the patient. Check that all setting screws are tightened.
- Stop therapy immediately, when you have doubts about the device settings and/or the therapy protocol.

## Warning!

Patient hazard –

- Ensure an **anatomically correct setup of the CPM device** suitable for the patient to be treated. For this purpose, check the following settings/positions (see numbers on device):
  1. Horizontal adduction/abduction
  2. Height adjustment
  3. Adjustment of upper arm length
  4. Adjustment of elbow angle
  5. Adjustment of forearm length
  6. Adjustment of headrest.  
When the optional patient restraint is used<sup>1</sup>
  7. Adjustment of headrest and shoulder restraint<sup>1</sup>
- It is not permitted to change the adjustment of 1 through 7 while a patient is sitting in the CPM device.
- Movements must **not cause any pain or irritation**.
- Patients must be **fully conscious** while being instructed in the use of the CPM device and during therapy.
- The choice of the **therapy parameters** to program and of the **therapy protocols** to use is restricted to the responsible physician or therapist. It is the physician's or therapist's decision whether or not to use the CPM device on a specific patient.
- The patient must be familiar with the functions of the **ARTROMOT®-S3 programming unit** and the unit must be **within easy reach** of the patient, allowing him or her to stop therapy, if needed. **Patients unable to operate the programming unit**, e.g. paralytic patients, must never be left unattended during therapy.
- Write the patient's name on the **patient chip card**. The card should only be used for this patient. If the patient chip card is used for another patient, be sure to **delete the previous patient's data** from the card first (see sections 4.1 and 5.3, paragraph "New Patient"). Use **original chip cards only**.

- All accessories used with the **ARTROMOT®-S3** device must first be approved by DJO.
- Modifications to the medical device described in this document without the manufacturer's written consent is prohibited.
- Ensure that the medical device is securely set up on solid ground during use.
- Do not allow **parts of the body or any objects** (such as blankets, cushions or cables) to get caught **in the moving parts** of the CPM device.

## Warning!

- Extreme caution should be taken when in use around small children and babies! Sufficient distance to the device is mandatory for their safety!
- Never leave the device unattended when it is switched on! Switch the device off and disconnect the power cord from the wall outlet!
- After use, store the device in a safe place! Ensure device stability also during storage!  
Do not rest or lean on or push the assembled device because it might tip over.

## Warning!

Shock hazard –

Strictly observe the following warnings. Failure to do so endangers the lives of the patient, the user and other persons involved.

- Before use allow the **ARTROMOT®-S3** to reach room temperature. If the device has been transported at **temperatures below 0 °C (32 °F)**, leave it to dry at room temperature for about 2 hours, until any condensation has disappeared.
- The **ARTROMOT®-S3** must only be operated in dry rooms.

<sup>1</sup> Only with **ARTROMOT®-S3 comfort**

- When disconnecting the device from the power line, remove the plug from the wall outlet first, before disconnecting the cable from the device.
- When connecting the device to other equipment or when creating a medical system, check that the sum of leakage currents will not cause any hazard. Please contact DJO, if you have questions in this matter.
- Do not use multiple portable socket outlets (MPSO) to connect the device to the power line. **ARTROMOT®-S3** must be connected to a **properly installed wall outlet with a non-fused earthed wire**. Before connecting the power cord, it must be completely unrolled and placed such that it will not get caught in the moving parts of the device.
- Before cleaning and service interventions, **disconnect the device from the power line** by removing the power cord from the wall outlet.
- **Liquids must not be allowed to enter the CPM device** or the programming unit. If liquids have entered into the devices, **ARTROMOT®-S3** must be immediately checked by a service technician, before it can be reused.

### **Warning!**

Equipment malfunction –

- Magnetic and electrical fields are capable of interfering with the proper performance of the device. For this reason make sure that all external devices operated in the vicinity of the CPM device comply with the relevant EMC requirements. X-ray equipment, MRI devices and radio systems are possible sources of interference as they may emit higher levels of electromagnetic radiation. Keep the CPM device away from these devices and verify its performance before use.
- Refer **repair and maintenance to authorized** persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.

- **Route all cables such** that they will not get caught in the moving parts during operation and that they do not present a stumbling hazard.
- **Inspect** the ARTROMOT®-S3 for damage and loose connections at least **once a year**.

### **Caution!**

Preventing chafing and pressure sores –  
When your patient is **adipose**, very **tall** or very **short**, be sure to prevent chafing and pressure sores.  
Patient hazard, damage to the carriage –  
Do not use the carriage to transport patients.

### **Caution!**

Equipment damage –

- **Check that the voltage and frequency ratings of your local** power line are those indicated on the nameplate.
- The seat withstands a **maximum continuous load** of 150 kg, the motion element of 12 kg.
- Do not allow **any objects** (such as blankets, cushions, or cables) to get caught in the **moving parts** of the CPM device.
- Do not expose the **ARTROMOT®-S3** device to direct sunlight, because some of the components may reach inadmissibly high temperatures.
- Be aware that the connectors can only be inserted in the correct orientation and secure all connections with the locks. Furthermore, the device must be set up at a safe distance from radiators to prevent excessive temperature rises.
- The presence of children, pets and rodents does not normally impair the functioning of the device. However, avoid contamination of the device by children or animals, from dust and lint, and keep them at a safe distance from the device. The safety statements set forth apply.

## 4. Device setup

**Note:** For a better understanding of each step, please fold out pages III and VIII.

### 4.1 Connecting the device, performance check

The equipment supplied includes these items:

base unit (chair), motion element, programming unit (1), patient chip card (2), arm rest (4), shoulder and head rest (ARTROMOT-S3 comfort only), power cord (not shown), operating instructions

1. Set up the device and connect the required cables by reversing the transport instructions given in section 6.3.
2. Follow these steps to set the carriage to the home position:

### Initial adjustment for new patients

Write the patient's name down on the back of the chip card. Insert the original patient chip card (2) into the programming unit (1).

Press the MENU key on the programming unit two times to access programming level 2 (with each key press, you advance one level).

Press the "new patient" parameter key **→0←** and select this function (a check mark appears in the circle next to the function).

Press the **START** key. The CPM device automatically enters the **home position**.

### Setup with programmed chip card

Insert the original patient chip card (2) into the programming unit (1).

Press the **START** key.

The device automatically enters the **starting position** (maximum adduction, half-way between internal / external rotation).

### Performance check:

If the programming unit can be operated as described above and **ARTROMOT®-S3** enters the home position (for home position values, refer to section 5.3), the device has passed the performance check.

The device also runs performance checks regularly during operation. This is what happens, if a problem is identified:

- An audio signal sounds.
- The device switches off immediately.
- The message "ERROR" and an error code (e.g. ERROR 5) appear on the display.

In this situation, you may attempt to restart the unit by turning it briefly off and on again with the power switch. With the unit switched off, check that all plugs are correctly connected. Switch the unit on again: if the error message persists, have the unit inspected by a Service technician, before using it again.

If it has been determined that the **ARTROMOT®-S3** device operates perfectly, ask the patient to sit down on the **ARTROMOT®-S3** chair.

## 4.2 Adjusting the device to the patient

### Note!

For the following adjustments, the patient's arm should not yet be placed on the armrest. After you have adjusted the device to the approximate patient measurements, the patient can place his/her arm on the armrest for a check of the setup and for the fine adjustment.

The settings are numbered 1 through 5. For easy orientation, you will find the same numbers on the CPM device. Always set up the CPM device in this order.

Write the settings down on the back of the patient's chip card.

Before you adjust the der ARTROMOT®-S3 to the patient, you may have to convert the device for use on the left or right shoulder. (see section 6.4 Conversion)

Before any treatment is possible, adjust the device as follows:

### Backrest, headrest, shoulder restraint, armrest

Before you start setting up the device for the exercise protocol, you must adjust it to a position that is anatomically correct for the respective patient.

For this purpose, the **ARTROMOT®-S3** comes with various padded and exchangeable supports with the corresponding restraints where applicable.

- Open thumbscrew (10) and adjust the **backrest** to a position that is comfortable for the patient. (Fig. F)
- In the next step adjust the **headrest**<sup>1</sup>: Height adjustment with locking screw (32), headrest adjustment with locking screw (30), fine adjustment with ball joint (33).
- Use locking screw (31) to position the **shoulder restraint**<sup>1</sup> approx. 1 cm above the shoulders.

- Set the **armrest** for the healthy arm to a height that allows the patient to sit up straight (locking screw 14).

Anatomically correct adjustment

<sup>1</sup> Only with ARTROMOT®-S3 comfort

### ① Horizontal adduction/abduction (Fig. A)

It is the purpose of the adjustment procedures to accommodate the patient in the most comfortable position possible.

- Press on the indexing knob (12) and make sure that it locks into place after the adjustment.

### ② Height adjustment (Fig. B)

Before opening locking screw (14), hold the moving part of the CPM device by grasping the holding tube for the length adjustment (17), to prevent the part from falling down.

- Adjust the height so that the axis of motor A is level with the pivot of the shoulder joint (see illustration on page 56). The pivot of motor A and the pivot of the shoulder joint must be on the same level.
- Tighten the locking screw.

### ③ Upper arm length (Fig. C)

The device for adjustment of the upper arm length is self-locking.

During adjustment, slightly lift motor B and ensure that the insertion tube does not wedge in the holding tube.

- Loosen clamping lever (18) and slightly lift the motor during adjustment.
- Tighten the clamping lever.

### ④ Elbow angle (Fig. D)

In most cases an angle of 90° to 60° is adjusted for elbow flexion.

- Loosen locking screw (20). To facilitate the adjustment, slightly lift motor B.
- Complete the adjustment and tighten the locking screw.

#### Note!

If the elbow angle is changed to a flexion value of greater than or less than 90°, the setting for the upper arm length will have to be changed, too.

### ⑤ Forearm length (Fig. E)

- Loosen clamping lever (26) and pull out the handle until the forearm can easily be accommodated between elbow support and handle.
- Tighten the clamping lever.

#### Note!

The armrest inclination can be adjusted until optimally suited to the arm.

To do so, loosen the locking screw (28) below the armrest. Adjust the inclination and tighten the locking screw.

### Checking the setup, fine adjustment

- Place the patient's arm on the armrest.
- Check adjustments 1 through 5 to ensure that
  - the pivot of motor A and the pivot of the shoulder joint are on the same level
  - the axis of motor B, the center of the elbow joint and the center of the shoulder joint form a straight line.
- To adapt the armrest even better to the patient, you can also adjust its height and inclination. To do so, loosen the locking screw (28) below the armrest.
- Check that all locking screws are tightened and all clamping levers properly closed.

## 5. Setting the treatment values

### Warning!

Patient hazard –

**Before therapy, a test run** consisting of several exercise cycles must be completed without the patient. Then repeat the test run with the patient and check that the movement does not cause any pain.

**Note: See also 2.2 and 2.3 as well as page V!**

### Note!

Programming is only possible when the patient's chip card has been inserted.

For information about therapy parameters and programming of special functions, please refer to sections 5.1 to 5.5.

See section 5.6 for programming examples.

### Important!

The programming unit **ARTROMOT®-S3** "Graphics" can be connected to all products of the **ARTROMOT®** shoulder series with serial numbers 3000 and higher.

It is possible to exchange chipcards between the "Text" and "Graphics" versions.

#### Please note:

1. If you use a "Graphics" programming unit with a formatted chipcard of the "Text" version, please note the following
  - the selected language of the "Text" version remains activated"; it is irrelevant for the "Graphics" version
  - the **reverse on load** setting saved in the "Text" version is automatically used for both motors
  - the **key lock** is disabled.
2. If you use a "Text" programming unit with a formatted chipcard of the "Graphics" version, please note the following
  - the selected language automatically changes to the default setting: German
  - the reverse on load setting saved for motor A in the "Graphics" version is automatically used for both motors
  - a key lock set in the "Graphics" version will remain active, because it is irrelevant for the "Text" version.

### 5.1 General information on programming ARTROMOT®-S3

1. You activate the programming mode by briefly pressing the **MENU** key on the programming unit.
2. The treatment parameters and functions are allocated to four (Standard

model) or six (Comfort model) programming levels (four per level). To be able to program a parameter you will have to access the corresponding programming level. This is also done with the **MENU** key. With each key press you advance one level. The code M1, M2, etc. that appears in the middle of the display indicates the programming level.

- You activate the treatment parameters and functions with the **four parameter keys** below the display. The symbols above the four parameter keys indicate the assigned parameters and functions.

This is what happens when you press one of the parameter keys to select a parameter:

- The corresponding symbol appears on the display in a larger format.
- The set value is displayed.
- The symbol above the parameter key appears in reverse video.

- With the + / - keys (plus/minus) you change the displayed value. When you press and hold the key, the value will change at a faster rate.

Some of the (special) functions can only be enabled and disabled. This is done by pressing the corresponding parameter key or with the + / - keys. Active parameters are identified with a check mark in the circle next to the symbol.

- Having programmed all parameters, press the **STOP** key to save the values.
- Then press the **START** key: ARTROMOT®-S3 checks the set values, moves to the position halfway between the set internal and external rotation values as well as to the set maximum adduction value and stops.
- Press the **START** key again to start therapy.

Next the carriage will move the maximum internal rotation value in the **synchronized mode**. The carriage will then simultaneously move to the maximum external rotation position and the maximum abduction position and subsequently to the maximum adduction position and the maximum internal rotation position. After reaching this position, the cycle restarts: maximum abduction and maximum external rotation.

After activation of the **START** key in the **non-synchronized mode** the motors will perform random movements, each motor reversing after reaching the maximum values.

## Note!

- Refer to sections 5.3 and 5.5 for a description of the parameters.
- To **view the set parameter values**, press the corresponding parameter key. Before, however, you have to press the **STOP** key and access the correct menu level.
- To prevent accidental changes of the parameter settings, **lock the keys** by simultaneously pressing the **+** (**plus**) and **-** (**minus**) keys for approx. 4 seconds.

Press both keys again for approx. 4 seconds to unlock.

- Selecting the "New Patient" function will automatically delete the data on the patient chip card. When you have finished programming the unit and press the **STOP** key, the settings will automatically also be saved to the patient chip card.
- Emergency stop function:** ARTROMOT®-S3 will stop immediately, when any of the keys is pressed during therapy.

Patient treatment can be resumed by pressing the **START** key. The device will automatically change the direction.

- The manufacturer recommends an application period of up to 1 hour per therapy session.
- For a uniform, smooth and gentle transition between the two directions of motion, the speed is automatically reduced before reaching the reversal point and, on reaching the reversal point, it is continuously increased until it reaches the set speed.

## Patients with a programmed chip card

- First complete the mechanical adjustments.
- Then insert the chip card (the patient is not yet positioned on the CPM device).
- Press the **START** key: the device will move to the starting position as specified by the parameters stored on the chip card and stops.
- Position the patient on the CPM device and press the **START** key to initiate therapy.

## 5.2 Programming ARTROMOT®-S3 Standard

Different programming levels are provided to program the ARTROMOT®-S3 Standard models.

You change between levels by pressing the **MENU** key repeatedly.

The display always indicates the currently selected level.

The following **treatment values, settings** and **information** can be entered/viewed on the programming unit (1):

### LEVEL 1:

- Abduction 
- Adduction 
- Internal rotation 
- External rotation 

### MENU

### LEVEL 2:


- Pause 
- Therapy timer 
- Speed 
- New patient 

### LEVEL 3:

- Reverse on load - motor A 
- Reverse on load - motor B 
- Motor A ON/OFF 
- Motor B ON/OFF 

### MENU

### LEVEL 4:

- Transport setting 
- Synchronized/ non-synchronized mode 
- Total therapy time 
- Service menu 

## Note!

- **While you adjust the values**, the carriage will move to the set range. This allows you to easily and quickly determine the ROM where the patient does not experience pain.
- **Special functions** can be programmed and retrieved with **ARTROMOT®-S3 Comfort** devices (see sections 5.4 and 5.5).
- The angle entered last for the respective direction of motion will be saved in each case.

## 5.3 Treatment value details – Standard models

- To access the programming level you need, press the **MENU** key repeatedly.
- The treatment parameters are selected with the different **parameter keys**.
- To change the treatment values, press the **+ / – keys**.
- To enable/disable a function, press the respective **parameter key** again.
- To save the entered values, press the **STOP** key.

**LEVEL 1:**

- **Abduction** 
  - Maximum value: 175 degrees
- **Adduction** 
  - Maximum value: 30 degrees


**⚠ Caution!**

Patient hazard –

When using the shoulder restraint, be sure not to program abduction values greater than 80°.

- **Internal rotation** 
  - Maximum value: - 90 degrees
- **External rotation** 
  - Maximum value: 90 degrees

**Note!**

- The programmed values and the values (degrees) measured on the patient may deviate slightly.
- To ensure a physiological motion, the carriage simultaneously moves to the following positions in the synchronized mode:
  - Maximum abduction value simultaneous with maximum external rotation value
  - Maximum adduction value simultaneous with maximum internal rotation value
- To ensure that the movements you program are physiological in nature, we recommend to program the carriage in the following order:
  1. Maximum adduction value 
  2. Maximum internal rotation value 
  3. External rotation (to relieve the shoulder)  = 0°
  4. Maximum abduction value 
  5. Maximum external rotation value 

- In the synchronized mode, after activation of **START** the carriage will first move to the maximum adduction position and to the position halfway between internal and external rotation. Then the carriage will stop, allowing the patient to be positioned. When the **START key** is pressed again (start of therapy), the carriage will first move to the maximum internal rotation position. Motor A (abduction/adduction) stops during this procedure. Once the maximum internal rotation position has been reached, both motors (A and B) will simultaneously move to the maximum abduction and external rotation position. After reaching this position, the cycle restarts: maximum adduction and maximum internal rotation.
- The elevation movement (flexion) is programmed by means of the abduction/adduction values (for mechanical adjustment instructions, refer to section 5.6 "Examples of use").
- Horizontal adduction/abduction is adjusted manually. These are possible settings:
 

Horizontal adduction:	<b>120°</b>
Horizontal abduction:	<b>0°</b>

**LEVEL 2:**■ **Pauses** 

Pauses occur at each programmed maximum value.

These are the two pause points:

- Maximum adduction value simultaneous with maximum internal rotation value
- Maximum abduction value simultaneous with maximum external rotation value

Pauses are adjustable in steps of 1 second between 0 and 30 seconds.

**Default: no pause**

Follow these steps to program pauses:

- Access menu level 2, then press the parameter key to select the special function .


- The display indicates the symbol for the special function and the selected pause durations.

top: abduction/external rotation pause

bottom: adduction/internal rotation pause (see display image)

The pause duration for abduction/external rotation appears in a dotted frame.


- To change the abduction/external rotation pause, press the "+" or "-" key. Select a duration of 10 seconds, for instance.

- Press the **parameter key** again when done. The dotted frame automatically moves to the line below, indicating that the adduction/internal rotation pause can be adjusted.

The **parameter key** is only used to toggle between these two functions.

- To modify the value, press the "+" (plus) or "-" (minus) key. Select a duration of 5 seconds, for instance.

- Save the programmed values by pressing the **STOP** key, then press the **START** key to begin the treatment.

## ■ Therapy timer

**Default setting is continuous operation.**

A clock symbol in the upper right-hand corner of the display identifies the continuous mode of operation. The clock indicates the elapsed therapy time.

In the continuous mode, the device must be stopped with the **STOP** key.

However, the therapy timer can be set in steps of **1 minute to any value between 1 and 300 minutes.**

At the end of the therapy time, the device automatically stops in the starting position (maximum adduction position / halfway between internal and external rotation).

In this case, a circle replaces the clock symbol. The circle fills as the therapy time progresses.

## ■ Speed

The speed can be adjusted between 1 % and 100 % in steps of 1 %. 100 % is equivalent to 230°/minute

**Default: 100 %**

## ■ New patient

With this function, the CPM device will move to the home position.

Select the function and press the **START** key:

- The carriage moves to the home position.
- The programmed treatment parameters are deleted.
- All values stored on the chip card are deleted.
- The device stops halfway between the abduction/adduction and internal/external rotation position.

The "New patient" function (home position) selects the following settings:

- Abduction: 41°
- Adduction: 39°
- Internal rotation: 1°
- External rotation: -1°
- Pauses: 0
- Timer: continuous operation
- Speed: 100 %
- Reverse on load - motor A: 25
- Reverse on load - motor B: 25
- Motor A: ON
- Motor B: ON
- Synchronized mode: ON
- Total therapy time: 0
- Special functions: disabled

**LEVEL 3:**

■ **Reverse on load - motor A**  **(safety circuit)**

The motors automatically switch to opposite directions of movement when the patient's resistance (load) exceeds the set value.

Adjustable levels for reverse on load feature: 1 - 25. At level 1, very low resistance will cause the device to reverse; at level 25, a high resistance is required to initiate the reversal.

**Default: level 25**

■ **Reverse on load - motor B**  **(safety circuit)**

The motors automatically switch to opposite directions of movement when the patient's resistance (load) exceeds the set value.

Adjustable levels for reverse on load feature: 1 - 25. At level 1, very low resistance will cause the device to reverse; at level 25, a high resistance is required to initiate the reversal.

**Default: level 25**

 **Caution!**

Patient hazard –

The reverse on load feature is a safety measure to protect the patient in the event of cramps, spasms, locked joints and similar situations. The manufacturer cannot be held liable for misuse of this feature.

■ **Motor A ON/OFF**  **A**

To permit fully isolated movements, the motors can individually be turned on and off. Motor A effects adduction and abduction, motor B effects internal rotation and external rotation.

For an isolated internal / external rotation movement, set motor A to the programmed position (abduction/adduction), then turn the motor off.

In normal mode the display indicates the symbol "OFF" during therapy, instead of the programmed adduction/abduction angles.

**Default: motor A ON**

■ **Motor B ON/OFF**  **B**

To permit fully isolated movements, the motors can individually be turned on and off. Motor A effects adduction and abduction, motor B effects internal rotation and external rotation.

For an isolated adduction/abduction movement, set motor B to the programmed position (internal / external rotation), then turn the motor off.

In normal mode the display indicates the symbol "OFF" during therapy, instead of the programmed internal/external rotation angles.


**Default: motor B ON**

**Note!**

Please note that one motor (A or B) must always be turned on. Otherwise the following symbol will be displayed when you press **START** :  OFF 

**LEVEL 4:**

■ **Transport setting** 

With this function, the carriage will move to a position optimally suited for packing the CPM device. Select the function and press the START key. The carriage moves to the transport position. In the top left-hand corner of the display you will see  (also refer to section 6.3 Transport).

■ **Synchronized/ A+B non-synchronized mode**

Motors A and B can run in synchronized or non-synchronized mode.

**Synchronized:**

Motors A and B perform a synchronized movement according to the shoulder joint's physiological motion pattern: Starting from the initial position (maximum adduction / halfway between internal and external rotation position), the carriage will first move to the maximum internal rotation position. The carriage will then simultaneously move to the maximum abduction position and the maximum external rotation position and subsequently to the maximum adduction position and the maximum

internal rotation position. After reaching this position, the cycle restarts: maximum abduction and maximum external rotation.

During operation the synchronized mode is indicated in the top left-hand corner by the symbol .

#### Non-synchronized:

Both motors run independently of each other within the selected ranges of motion.

To enable the non-synchronized mode, disable the synchronized mode.

During operation the non-synchronized mode is indicated in the top left-hand corner by the symbol .

**Default: synchronized mode enabled**

### Caution!

Patient hazard –

The synchronized operation is recommended as a general rule. However, medical and/or therapeutic indications may call for the non-synchronized mode. In the non-synchronized mode, particular caution and attention on the part of the physician/therapist is required to rule out any risk for the patient.

#### ■ Total therapy time

Under menu item "total therapy time" you can view each patient's total therapy time (duration of all the patient's therapy sessions).

#### To delete the stored therapy time:

Press and hold the **parameter key** for 5 seconds or select the "New patient" function.

#### ■ Service menu

For service purposes only, refer to Service Manual.

#### Reminder:

You save the selected parameter values by pressing the **STOP** key.

## 5.4 Programming ARTROMOT®-S3 Comfort

ARTROMOT®-S3 devices of the Comfort series offer two more programming levels for additional functions.

The programming levels are selected in the same way as on Standard models.

Programming levels 1, 2, 3 and 6 are identical with programming levels 1 to 4 of the Standard models.

All special functions are disabled in the standard models.

The following treatment values, settings and information can be entered/viewed on the programming unit (1):

#### LEVEL 1:

- Abduction 
- Adduction 
- Internal rotation 
- External rotation 

#### LEVEL 2:

- Pause 
- Therapy timer 
- Speed 
- New patient 

#### LEVEL 3:

- Reverse on load - motor A 
- Reverse on load - motor B 
- Motor A ON/OFF 
- Motor B ON/OFF 

#### LEVEL 4:

- Warm-up protocol 
- Isolation program 
- Therapy documentation Abduction/adduction 
- Therapy documentation Internal/external rotation 

**LEVEL 5:**

- Oscillation
- Stretching in abduction
- Stretching in internal rotation
- Stretching in external rotation

**MENU****MENU****LEVEL 6:**

- Transport setting
- Synchronized/non-synchronized mode
- Total therapy time
- Service menu

A+B

**5.5 Protocol details – Comfort model**

- To access the programming level you need, press the MENU key repeatedly.
- The treatment parameters are selected with the different parameter keys.
- To change the treatment values, press the + / – keys.
- To enable/disable a function, press the respective parameter key again.
- To save the entered values, press the STOP key.

All special functions are disabled in the standard models.

**Note!**

**LEVEL 1:** equivalent to level 1 of the Standard model (see: 5.3)

**LEVEL 2:** equivalent to level 2 of the Standard model (see: 5.3)

**LEVEL 3:** equivalent to level 3 of the Standard model (see: 5.3)

**LEVEL 6:** equivalent to level 4 of the Standard model (see: 5.3)

**LEVEL 4:****■ Warm-up protocol**

During warm-up, the patient will slowly become used to the set maximum values for adduction/abduction and internal/external rotation.

After activation of the warm-up protocol the carriage moves between the maximum adduction and abduction positions and the center position of the programmed range of motion (adduction/internal rotation and abduction/external rotation).

With each cycle, 3° is added to the abduction/external rotation range of motion until the maximum abduction/external rotation positions are reached.

The adduction and internal rotation ranges are completed in each cycle. Once the maximum abduction and external rotation values have been attained, the carriage switches to the normal mode.

If a therapy time is activated, the system restarts with the warm-up protocol when the therapy time has elapsed and the devices is started again.

During warm-up, the display will show the symbol in the upper, left-hand corner.

**Default: disabled**

**■ Isolation protocol AIIIB**

For this special function, both motors are turned on, but they will never move simultaneously.

These are the steps of the special function:

- At first motor A moves to the programmed maximum values for adduction and abduction for ten cycles, then it stops. In this phase, motor B is disabled (symbol for motor B: **OFF**)
- Then motor B will move to the programmed maximum values for internal and external rotation for 10 cycles and stops, after this the full cycle will start over with motor A. Motor A is disabled while motor B is operating (symbol for motor A: **OFF**)

- The stop position at the end of each tenth cycle (for both adduction/abduction and internal/external rotation) can be selected in steps of 25% between 0% and 100% of the programmed maximum range of motion.
- Steps 1 and 2 can be repeated as often as needed. You can stop treatment with the **STOP** key; after the programmed therapy duration, it will stop automatically.

During the isolation protocol, the display will show the symbol  in the upper, left-hand corner.


**Default: disabled**

Follow these steps to program the special function:

- Access menu level 4, then press the parameter key to select the special function. The display will indicate:
  - the symbol representing the special function
  - the status of the special function (circle with or without check mark)
  - two boxes with the current percentage of the stop positions for motor A and B
  - a check mark in a circle
- You activate the function with the "+" key and deactivate it with the "-" key.

In this case it is not possible to activate/deactivate the function by pressing the **parameter key** again. A check mark in the circle identifies the activated function.

- Next press the parameter key of the isolation protocol. The dotted frame moves to the upper percentage value, allowing you to enter the stop position for motor A.
- To modify the value, press the "+" (plus) or "-" (minus) key. Select 75%, for instance.
- Press the parameter key again when done. The dotted frame automatically moves to the lower value, allowing you to adjust the stop position of motor B.
- To modify the value, press the "+" (plus) or "-" (minus) key. Select 25%, for instance.
- Save the programmed values by pressing the **STOP** key, then press the **START** key to begin the treatment.

**This is what happens:**

Abduction/adduction is treated ten times first (motor A). Motor B is positioned at 25% of the programmed range of motion for internal and external rotation.

Subsequently, internal / external rotation is treated ten times (motor B). Motor A is positioned at 75% of the programmed range of motion for adduction and abduction.

**Note!**

The percentage values can be modified only when the special function is enabled (check mark).

■ **Therapy documentation adduction/abduction**

This is a special function of ARTROMOT®-S3 Comfort which allows the entire therapy documentation to be represented.

The carriage run times as well as the range of motion of the sessions are recorded.

The collected data are presented graphically in the form of a coordinate system (X-axis = range of motion / Y-axis = time) where the upper curve illustrates the trend of the abduction movement and the lower curve the trend of the adduction movement.

■ **Therapy documentation internal / external rotation**

This is a special function of ARTROMOT®-S3 Comfort which allows the entire therapy documentation to be represented.

The carriage run times as well as the range of motion of the sessions are recorded.

The collected data are presented graphically in the form of a coordinate system

(X-axis = range of motion / Y-axis = time) where the upper curve illustrates the trend of the external rotation movement and the lower curve the trend of the internal rotation movement.

## LEVEL 5:

### ■ Oscillation

The „oscillation“ function permits a more efficient therapy in the last 10° before the maximum abduction and adduction values are reached.

For this exercise, the carriage begins at the starting position (maximum adduction, halfway between internal and external rotation). First of all the carriage moves to the maximum internal rotation position. Then it will move to the maximum abduction position, at the same time as the maximum external rotation value.

When the abduction/external rotation value has been reached, the carriage reverses 10° toward adduction/internal rotation and then moves back again to the abduction/external rotation position.

The movement through the final 10° is repeated five times at a slow speed. At the end of the cycle, the carriage will again move to the maximum adduction position, at the same time as the maximum internal rotation value, and then starts another cycle with five repetitions through the last 10° of abduction/external rotation.

This sequence can be repeated as often as needed. You can stop treatment with the **STOP** key; after the programmed therapy duration, it will stop automatically.

**Default: disabled**

### ■ Stretching in abduction

With the special „Stretching in abduction“ function the carriage gently stretches the joint by lifting the arm. Only the abduction movement is treated here. Motor B is automatically disabled. The position of motor B cannot be changed once this special function has been activated.

Starting at the middle position the carriage will first move to the programmed

adduction position and then to the programmed abduction position.

Subsequently the carriage reverses 5° towards adduction and then moves very slowly back again to the programmed abduction position (display ).

After that it attempts to attain another 5°, moving even slower than before (display ).

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the carriage moves in the opposite direction.

This stretch cycle is repeated ten times. After that the carriage moves to the programmed maximum adduction position and restarts the abduction stretch cycle.

This sequence can be repeated as often as needed. You can stop treatment with the **STOP** key; after the programmed therapy duration, it will stop automatically.

**Default: disabled**

## Note!

- If a pause is programmed when you enable the function, the carriage will stop for the duration of the pause each time the maximum stretching value is attained.
- Enabling the stretching in abduction function automatically disables motor B. The only movement performed is adduction / abduction.
- It is not possible to combine a second stretching function with the special function „stretching in abduction“.
- During operation, the display will show the symbol  in the upper, left-hand corner.

### ■ Stretching in internal rotation

With the special „Stretching in internal rotation“ function the carriage gently stretches the joint towards the inside.

Starting at the middle position the carriage will first move to the programmed external rotation position and then to the programmed internal rotation position.

Subsequently the carriage reverses 5° towards external rotation and then moves very slowly back again to the programmed internal rotation position (display ). After that it attempts to attain another 5°, moving even slower than before (display ).

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the carriage moves in the opposite direction.

This stretch cycle is repeated ten times.

After that the carriage moves to the programmed maximum external rotation position and restarts the internal rotation stretch cycle.

This sequence can be repeated as often as needed. You can stop treatment with the **STOP** key; after the programmed therapy duration, it will stop automatically.

**Default: disabled**

## Note!

- If a pause is programmed when you enable the function, the carriage will stop for the duration of the pause each time the maximum stretching value is attained.
- Enabling the stretching in internal rotation function automatically disables motor A.

The only movement performed is internal / external rotation.

- It is not possible to combine a second stretching function with the special function „stretching in internal rotation“.
- During operation, the display will show the symbol  in the upper, left-hand corner.

## ■ Stretching in external rotation

With the special „Stretching in external rotation“ function the carriage gently stretches the joint towards the outside.

Starting at the middle position the carriage will first move to the programmed internal rotation position and then to the programmed external rotation position.

Subsequently the carriage reverses 5° towards internal rotation and then moves very slowly back again to the programmed external rotation position (display ). After that it attempts to attain another 5°, moving even slower than before (display ).

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the carriage moves in the opposite direction.

This stretch cycle is repeated ten times.

After that the carriage moves to the programmed maximum internal rotation position and restarts the external stretch cycle.

This sequence can be repeated as often as needed. You can stop treatment with the **STOP** key; after the programmed therapy duration, it will stop automatically.

**Default: disabled**

## Note!

- If a pause is programmed when you enable the function, the carriage will stop for the duration of the pause each time the maximum stretching value is attained.
- Enabling the stretching in external rotation function automatically disables motor A. The only movement performed is internal / external rotation.
- It is not possible to combine a second stretching function with the special function „stretching in external rotation“.
- During operation, the display will show the symbol  in the upper, left-hand corner.


## Reminder:

You save the selected parameter values by pressing the **STOP** key.

## 5.6 Examples of use/ sample protocols

### 5.6.1 Isolated adduction/abduction


1. Complete the mechanical settings and adjust the carriage to the patient as described in section 4.2.
2. Press the **MENU key** on the programming unit to access level 1 (M1).
3. Press the **parameter key for internal rotation**  or **external rotation**  and press the +/- keys to set the carriage to the desired rotation position to be maintained during the isolated abduction/adduction movements.
4. Next set the desired range of motion for adduction  and abduction by first pressing the **parameter key for adduction** and set the values with the +/- keys.
5. Set the **abduction** value in the same way .

#### Note!

- For an exclusive adduction/abduction protocol, motor B for the rotation movement must be turned off as described in steps 6 to 8 below.
- Having programmed the ranges of motion, you can program additional options, such as pause or speed.

6. Press the **MENU key** on the programming unit to access level 3 (M3).
7. Press the Motor B ON/OFF parameter key  B, to activate the parameter.
8. Deactivate motor B by pressing the parameter key **Motor B ON/OFF** again or the „-“ (minus) key. The check mark in the circle next to the symbol must have disappeared.
9. Save the settings with the **STOP key**, set the carriage to the starting.

### 5.6.2 Isolated internal/external rotation

1. Complete the mechanical settings and adjust the carriage to the patient as described in section 4.2.
2. Press the **MENU key** on the programming unit to access level 1 (M1).
3. Press the **parameter key for adduction**  or **abduction**  and press the +/- keys to set the carriage to the desired position to be maintained during the isolated rotation movements.
4. Next set the desired range of motion for internal/external rotation by first pressing the **parameter key for internal rotation**  and set the values with the +/- keys.
5. Set the **external rotation** value in the same way .

#### Note!

- For an exclusive rotation protocol, motor A for adduction/abduction must be turned off as described in steps 6 to 8 below.
- Having programmed the ranges of motion, you can program additional options, such as pause or speed.

6. Press the **MENU key** on the programming unit to access level 3 (M3).
7. Press the **Motor A ON/OFF**  A, parameter key to activate the parameter.
8. Deactivate motor A by pressing the parameter key **Motor A ON/OFF** again or the „-“ (minus) key. The check mark in the circle next to the symbol must have disappeared.
9. Save the settings with the **STOP key**, set the carriage to the starting position with the **START key** and initiate therapy by pressing the **START key** again.

### 5.6.3 Isolated elevation (flexion)

1. First adjust an angle of 90° for external rotation and deactivate motor B as follows:
  - a. Press the MENU key on the programming unit to access level 1 (M1).
  - b. Press the **parameter key** for **external rotation**  and move the carriage to the rotation position of 90°, using the + / - keys.

#### Note!

- For an exclusive elevation protocol, motor B for the rotation movement must be turned off as described in steps c to e below.
- Having programmed the ranges of motion, you can program additional options, such as pause or speed.

- c. Press the **MENU key** on the programming unit to access level 3 (M3).
  - d. Press the **Motor B ON/OFF**  parameter key to activate the parameter.
  - e. Deactivate motor B by pressing the parameter key Motor B ON/OFF again or the „-“ (minus) key. The check mark in the circle next to the symbol must have disappeared.
  - f. Save the entered values by press the **STOP key**.
2. Complete the mechanical settings and adjust the carriage to the patient as follows:
  - a. Adjust an angle of 90° for horizontal adduction/abduction ①, using indexing knob (12).
  - b. Secure and hold the motion element at the holding tube for length adjustment and open the locking screw (14) for height adjustment ②. Select a height at which the pivot points of motor A and the shoulder are at the same level.
  - c. Now open clamping lever (18) for adjustment of the upper arm length ③, locking screw (20) for adjustment of the elbow angle ④ and locking screw (28) for adjustment of the forearm rest inclination and set the motion

element to a position that is comfortable to the patient and suitable for the therapy session. When done, tighten all screws and levers.

- d. Open clamping lever (26), adjust the forearm length ⑤ to the patient and close the clamping lever.
3. In the next step, program the elevation range of motion as follows:
  - a. Press the **MENU key** on the programming unit to access level 1 (M1).
  - b. Next set the desired range of motion for adduction and abduction by first pressing the **parameter key** for **adduction**  and set the values with the +/- keys.
  - c. Set the **abduction**  value in the same way.

#### Note!

Having programmed the ranges of motion, you can program additional options, such as pause or speed.

4. Save the settings with the **STOP key**, set the carriage to the starting position with the **START key** and initiate therapy by pressing the **START key** again.

## 6. Care, Maintenance, Transport, Conversion

The ARTROMOT®-S3 is suitable for re-use if the following points are observed.

### 6.1 Care/re-use

#### **Warning!**

Shock hazard –

Unplug the device from the power line before cleaning.

Shock hazard, equipment damage –

Liquids must not enter the device or the programming unit.

- The ARTROMOT®-S3 must be disinfected by wiping down with a disinfectant. Thus, it complies with the special hygiene standards for medical technical equipment.
- **The enclosure** and removable **armrest assemblies** can be cleaned with **commonly used disinfectants** and **mild household detergents**.
- Only use a **damp cloth** to wipe the CPM device down.

#### **Warning!**

Patient hazard –  
patient contamination

- Before using the device on another patient, be sure to clean and disinfect it according to the instructions given here.

#### **Note!**

The manufacturer recommends using only special disinfectants approved for medical devices and with the characteristics specified below under "Caution".

Disinfection with the following disinfectants and in compliance with the instructions for use has been evaluated and approved by the manufacturer in the risk assessment:

- DESCOSEPT AF lemon  
(Part no.: 00-311L-xxx)  
BAuA Reg.No.: N-55153, CE-0482  
Manufacturer: Dr. Schumacher GmbH  
Contact time: 2 minutes minimum

#### **Caution!**

Equipment damage –

- The plastic material used is not resistant to mineral acids, formic acid, phenols, cresols, oxidants and strong organic or inorganic acids with a pH value below 4.
- Use only clear disinfectants to prevent discoloration of the device.
- Do not expose the CPM device to strong ultraviolet radiation (sunlight) and fire.

### 6.2 Maintenance (fuse replacement)

#### **Check before each use**

Visually inspect the device for signs of mechanical damage before each use.

If you detect damage or malfunctions that may impair the safety of the patient or of the operator, have the device repaired before using it.

## Technical inspections

Due to their design, the bearings and joints of the ARTROMOT®-S3 do not require maintenance and all materials are protected from corrosion, BUT: For safety, the devices require regular maintenance. To maintain the functional and operational safety, check all components for damage and loose connections at least once a year.

These checks should be performed by persons with adequate training and experience. Damaged and worn parts must immediately be replaced with original spare parts by authorized, specially trained staff.

### ⚠ Warning!

Patient hazard, equipment malfunction and damage

- Refer repair and maintenance to authorized persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.

ORMED GmbH will make all documents required for servicing, such as circuit diagrams, parts lists, descriptions or calibration instructions, available to authorized experts.

These checks should be performed by persons with adequate training and experience. Damaged and worn parts must immediately be replaced with original spare parts by authorized staff.

**Other than that, the manufacturer does not require any other regular maintenance.**

### Note!

Regarding further technical or other inspections or inspection intervals, observe applicable local requirements, such as IEC 62353, DGUV3 or similar regulations to be observed by users of medical equipment or electric devices.

## Fuse replacement

### ⚠ Warning!

Patient hazard, equipment malfunction and damage

The replacement of fuses must be referred to specialists as defined in IEC 60364 or other applicable standards (e.g. biomedical technicians, electricians, electronics installers).

**Fuses used must be**

**S# < 20.000: T1A L250Vac**

**S# > 20.000: T2A H250Vac fuses.**

- Before replacing fuses, turn off the ARTROMOT®-S3 and disconnect the device from the power line.
- Use an appropriate tool to remove the fuse holder (8) situated between the power switch (7) and the power connector (6) (Fig. 1).
- Replace the fuses and reinsert the fuse holder (Fig. 2). Check that the fuse holder locks properly into place.


Fig. 1


Fig. 2


## 6.3 Transport

The following operating steps must be completed before transporting the ARTROMOT®-S3:


1. Activate the „transport setting“ →  and start the ARTROMOT®-S3.
2. Turn off the ARTROMOT®-S3 power switch (7).
3. Remove the power cord (6) and disconnect the motion element and the programming unit (9).
4. Open locking screw (14) and pull out the armrest (4) for the healthy arm.
5. Grasping the holding tube (17) of the motion element, open locking screw (14) and pull out the motion element.
6. Set horizontal adduction/abduction to 0° (indexing knob 12).
7. Open the locking screw for the headrest adjustment (30) and remove the headrest<sup>1</sup>.
8. Open the locking screw for the shoulder restraint (31) and remove<sup>1</sup>.
9. Open the locking screw for the height adjustment of the patient restraint (32) and remove<sup>1</sup>.
10. Open the locking screw for adjustment of the backrest (10), fold the backrest to the front all the way until flat on the seat and tighten the screw.
11. Remove the two locking pins (34). Pull out the legs towards the sides, reverse and reinsert. Then put the locking pins back in place.

12. Only use the original shipping box for transporting the device. DJO can not be held liable for transport damage, if the device was not shipped in its original shipping box.
13. Place the ARTROMOT®-S3 – legs down – on the bottom of the box. The markings on the bottom of the box indicate the correct position.

14. Then put the enclosed polystyrene block on the chair, observing the cut-outs in the block.
15. Put the programming unit (1) in the supplied box. Place the motion element, the armrest, the components of the neck and shoulder restraint and the power cord in the respective cut-outs in the polystyrene block.

16. Then close the box.

Reverse the above steps to reassemble the device after transport.

### **Warning!**

Shock hazard –

**Before use** allow the ARTROMOT®-S3 to reach room temperature. If the device has been transported at **temperatures below 0 °C (32 °F)**, leave it to dry at room temperature for about 2 hours, until any condensation has disappeared.


<sup>1</sup>Only with **ARTROMOT®-S3 comfort**

## 6.4 Conversion

The ARTROMOT®-S3 can be used on the left and on the right shoulder joint. However, the device must be converted first. This is done very quickly.

1. Activate the „transport setting“ →  function in the menu (see also 5.3) and start the ARTROMOT®-S3.
2. On both sides adjust an angle of 90° for horizontal adduction/abduction (indexing knob 12) (see conversion drawing ①).
3. Open locking screw (14), remove the armrest for the healthy arm and place it on the seat (see conversion diagram ① a).
4. Grasping the holding tube (17) of the motion element, open locking screw (14).
5. Remove the motion element and insert it on the opposite side. Tighten locking screw (14) (see conversion diagram ① b, c).
6. Open locking screw (20) and rotate the forearm rest 180°.  
Tighten locking screw (20) (see conversion diagram ② d).
7. Hold the forearm rest and open locking screw (28).
8. Swivel the forearm rest to the other side around motor B and tighten locking screw (28) (see conversion diagram ③ e).
9. Insert the armrest for the healthy arm into the holding tube and tighten locking screw (14).
10. Set horizontal adduction/abduction to the desired value on the side of the motion element (normally a value between 0° and 40°, as needed).

### Conversion diagram:

## 7. Environmental Protection Statement

The product described in this operation manual must not be disposed of with unsorted household or municipal waste. It requires separate disposal. Please contact DJO for information about the possible recycling of the product.

The expected service life of the device and all supplied parts and accessories is 6 years. The application beyond this period is the responsibility of the user.

## 8. Specifications

<b>Model:</b>	ARTROMOT®-S3
<b>Device name:</b>	ARTROMOT®-S3 Standard, Part no.: 80.00.023 ARTROMOT®-S3 Comfort, Part no.: 80.00.024
<b>Input ratings:</b>	100 – 240 Vac / 50 – 60 Hz tolerance -15% to +10%
<b>Current consumption:</b>	100V      240V
standby (ON):	5VA      5VA
operation (maximum):	40VA      40VA
<b>Fuses:</b>	serial numbers < 20.000: 2x T1A L250 Vac serial numbers > 20.000: 2x T2A H250 Vac according to IEC 60127-2/3, UL 248-14, CSA C22.2 no. 248.14
<b>Protection class:</b>	serial numbers <20,000: I serial numbers >20,000: II
<b>IP class:</b>	serial numbers <20,000: IPX0 serial numbers >20,000: IP21
<b>Applied part:</b>	type B
<b>Max. continuous load</b>	
seat:	150 kg
motion element:	12 kg
<b>Dimensions (transport):</b>	
length:	87.5 cm
width:	57.5 cm
height:	58 cm

<b>Adjustment ranges (min./max.):</b>	
height adjustment: (measured from seat)	35 – 71 cm
upper arm length:	20 – 32 cm
forearm length:	29 – 46 cm
<b>Seat height:</b>	48 cm
<b>Accuracy of measured values:</b>	
goniometer in the measuring range:	from 90° to 90° in internal/external rotation
accuracy:	+/- 2°
<b>weight:</b>	25 kg
<b>materials used:</b>	ABS, POM (Delrin 100), PUR, PA, FR4, aluminum, stainless steel, brass
<b>MDD:</b>	class IIa
<b>Standards compliance:</b>	93/42/EEC IEC 60601-1 IEC 60601-1-6 IEC 60601-1-9 IEC 60601-1-11 (S# >20,000) IEC 62366 IEC 62304 EN ISO 14971 ANSI/UL 60601-1 (S# <20,000) CAN/CSA C22.2 No. 601.1-M90 (S# <20,000) ANSI AAMI ES 60601-1 1st Edition (S# >20,000) CAN CSA 22.2 No. 60601-1-08 (S# >20,000)
<b>EMC</b> (electromagnetic compatibility)	IEC 60601-1-2
<b>Manufactured in compliance with:</b>	EN ISO 13485
<b>Ambient conditions (storage, transport)</b>	
temperature:	-25 °C to +70 °C
relative humidity:	at 70 °C up to 93 % no condensation
atmospheric pressure:	500 hPa to 1060 hPa
<b>Ambient conditions (operation)</b>	
temperature:	+5 °C to +40 °C
relative humidity:	15 % to 93 %
atmospheric pressure:	700 hPa to 1060 hPa

## 9. Contact

We would be happy to answer any questions you may have about our products and services.

### **ARTROMOT® International:**

Please contact your local dealer, the DJO Headquarters in the USA, the DJO International Headquarters in England or DJO in Germany.

### **DJO International**

Headquarters: DJO UK Ltd.  
1a Guildford Business Park  
Guildford  
Surrey, GU2 8XG  
United Kingdom  
Telephone: +44 (0)1483 459659  
Fax: +44 (0)1483 459470  
e-mail: [info@DJOglobal.eu](mailto:info@DJOglobal.eu)  
Web: [www.DJOglobal.eu](http://www.DJOglobal.eu)

### **DJO Headquarters**

DJO, LLC  
1430 Decision Street  
92081 Vista – California / USA  
Telephone: +1 760 727 1280  
Fax: +1 800 936 6569  
e-mail: [webmaster@DJOglobal.com](mailto:webmaster@DJOglobal.com)  
Web: [www.DJOglobal.com](http://www.DJOglobal.com)

### **Manufacturer / Headquarters Germany:**

DJO  
ORMED GmbH  
Merzhauser Strasse 112  
79100 Freiburg, Germany  
Telephone: +49 (0) 761 4566 01  
Fax: +49 (0) 761 456655-01  
e-mail: [infoservice@DJOglobal.com](mailto:infoservice@DJOglobal.com)  
Web: [www.DJOglobal.de](http://www.DJOglobal.de)

### **Warranty:**

2 years (mechanical parts)  
2 years (electronics)

## 10. Technical Service

### 10.1 Technical hotline

Do you have any technical questions?  
Do you need technical service?

Telephone: +49-180-5-1 ormed de  
              +49-180-5-1-67 63 33  
Fax: +49-180-5-3 ormed de  
      +49-180-5-3-67 63 33

### 10.2 Shipment

To prevent damage during transport, only use the original shipping box. These boxes can be obtained from DJO.

Before packing the CPM device, set it to the transport position (see chapter 5 and 6).

### 10.3 Spare parts

Refer to the Service Manual for the most recent list of spare parts.

When ordering spare parts, always specify:

- item
- description
- part number
- quantity
- serial number of the CPM device

#### Note!

Refer repairs to authorized, specially trained staff.

DJO offers service training for your personnel.

#### Note!

Country-specific power cords are available from DJO. Please contact DJO or your DJO distributor.

Surcharges may apply in certain cases to spare parts ordered in low quantities.

Item	Description	Part No.	Qty
1.	patient chip card	o.0034.048	<input type="checkbox"/>
2.	Patient chip card (protocol)	o.0037.035	<input type="checkbox"/>
3.	marker pen for patient chip card	o.0031.006	<input type="checkbox"/>
4.	Power cord EU version H05VV-F3G 3 × 1 mm <sup>2</sup> length: 4 m	o.0034.118	<input type="checkbox"/>

## 11. CE-Declaration of Conformity

DECLARATION OF CONFORMITY

In compliance with the provisions of the Council Directive 93/42/EEC, Annex II, of 14 June 1993 / Amendment of 5 September 2007, concerning medical devices, the company

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declares under its sole responsibility that the products of the product line

**ARTROMOT®** (see Annex)

fulfill the essential requirements of Annex I  
of the Council Directive 93/42/EEC.

With reference to Rule 9 of the Directive 93/42/EEC, Annex 9,  
the product is a device of risk class IIa.

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Freiburg, 17 August 2018

-----  
- QA Management Representative -  
-Bernhard Krohne-

This certificate is valid until expiry of the certificate referred to.  
(The certificate for the year of manufacture can be downloaded from:  
<https://www.djoglobal.de/arzt/qualitaet.html>)

**Annex:**  
ARTROMOT®-S3  
ARTROMOT®-S4  
ARTROMOT® ACTIVE-K  
ARTROMOT®-K1  
ARTROMOT®-SP3  
ARTROMOT®-E2

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**Bijlage:**  
ARTROMOT®-S3  
ARTROMOT®-S4  
ARTROMOT® ACTIVE-K  
ARTROMOT®-K1  
ARTROMOT®-SP3  
ARTROMOT®-E2

## 12. IEC 60601-1-2 and IEC 60601-1-11

The **ARTROMOT®-S3** device is subject to particular precautions regarding electromagnetic compatibility (EMC). The device must be installed and put into service strictly in compliance with the EMC directives set forth in the accompanying documents.

Portable and mobile RF communication systems may affect the **ARTROMOT®-S3**.

The **ARTROMOT®-S3** should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, **ARTROMOT®-S3** should be observed to verify normal operation in the configuration in which it will be used.

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If it is necessary to replace assemblies or cables, only the manufacturer's original parts must be used to ensure continued compliance with EMC requirements after repair. This requirement applies to the power supply unit, cables and cable lengths, drive unit consisting of the motor and the control system, the programming unit incl. the coiled cable and the connector.

Regarding the electromagnetic compatibility, there is no maintenance work required during the expected service life of the product.

The product designation **ARTROMOT®-S3** used in the text below includes all product variants.

### **WARNING!**

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

## 12.1 Electromagnetic emissions

### Guidance and Manufacturer's Declaration – Electromagnetic Emissions

**ARTROMOT®-S3** is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-S3** unit is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions to CISPR 11	Group 1	<b>ARTROMOT®-S3</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions to CISPR 11	Class B	<b>ARTROMOT®-S3</b> is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions to IEC 61000-3-2	complies with class A requirements	
Voltage fluctuations/flicker emissions to IEC 61 000-3-3	complies	

## 12.2 Electromagnetic immunity

### Guidance and Manufacturer's Declaration – Electromagnetic Immunity

**ARTROMOT®-S3** is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-S3** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment – guidance
Electrostatic discharge (ESD) to IEC 61 000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/bursts to IEC 61 000-4-4	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power should be that of a typical commercial or hospital environment.
Surges to IEC 61 000-4-5	± 1 kV line-to-line voltage	± 1 kV line-to-line voltage	Mains power should be that of a typical commercial or hospital environment.
Voltage dips To IEC 61000-4-11	0% $U_T$ ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°  0% $U_T$ ; 1 cycle and 70% $U_T$ ; 25/30 cycles  Single phase at 0°	0% $U_T$ ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°  0% $U_T$ ; 1 cycle and 70% $U_T$ ; 25/30 cycles  Single phase at 0°	Mains power should be that of a typical commercial or hospital environment. If the user of the <b>ARTROMOT®-S3</b> unit requires continued operation during power mains interruptions, it is recommended that the <b>ARTROMOT®-S3</b> unit be powered from an uninterruptible power supply or a battery.
Voltage interruptions to IEC 61000-4-11	0% $U_T$ ; 250/300 cycles	0% $U_T$ ; 250/300 cycles	
Power frequency (50/60 Hz) magnetic field to IEC 61 000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.

**Note:**  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and Manufacturer's Declaration – Electromagnetic Immunity

**ARTROMOT®-S3** is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-S3** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment – guidance
Conducted RF to IEC 61000-4-6	6 V 0,15 MHz-80 MHz 80% AM at 1 kHz	6 V 0,15 MHz-80 MHz 80% AM at 1 kHz	
Radiated RF to IEC 61000-4-3	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communications equipment to IEC 61000-4-3	27 V/m 380-390 MHz PM 18 Hz  28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine  9 V/m 704-787 MHz PM 217 Hz  28 V/m 800-960 MHz PM 18 Hz  28 V/m 1,7-1,99 GHz PM 217 Hz  28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz  9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	27 V/m 380-390 MHz PM 18 Hz  28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine  9 V/m 704-787 MHz PM 217 Hz  28 V/m 800-960 MHz PM 18 Hz  28 V/m 1,7-1,99 GHz PM 217 Hz  28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz  9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ARTROMOT®-S3** exceeds the applicable **ARTROMOT®-S3** unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **ARTROMOT®-S3** unit.

## 12.3 RF (Radio Frequency) Interference

EMC disturbance can affect the performance of the **ARTROMOT®-S3**. To prevent the user from unacceptable risks the device runs performance checks regularly during operation. This is what happens, if a problem is identified:

- an audio signal sounds
- the unit stops immediately
- the message ERROR and an error code appear on the display

In this case try to quit the error with the STOP key. If the error persists you may attempt to restart the unit by turning it briefly off and on with the power switches.

### **WARNING!**

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the **ARTROMOT®-S3**, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

### **Notice!**

(for home use)

Portable and mobile RF communication systems may affect the ARTROMOT®-S3 device.

For this reason, make sure that wireless communication equipment, such as wireless home network devices, cell phones, cordless telephones and their base stations and walkie-talkies operate at a minimum distance of 3,3 m from the device. (Calculated on the basis of the maximum power output of a typical cell phone of 2 W).

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### Symbol Overview standard models ARTROMOT®-S3

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## Symbol Overview comfort models ARTROMOT®-S3

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## AMENDMENT TO THE INSTRUCTION FOR USE (IFU)

This document outlines the amendment to the Instruction for Use (IFU) for necessary updates and changes.

### List of the affected devices:

Reference Number	Article Name	GTIN
80.00.040	<b>ARTROMOT®</b> -K1 Standard	00888912403191
80.00.041	<b>ARTROMOT®</b> -K1 Standard Chip	00888912403207
80.00.042	<b>ARTROMOT®</b> -K1 Comfort	00888912403214
80.00.043	<b>ARTROMOT®</b> -K1 Comfort Chip	00888912403221
80.00.044	<b>ARTROMOT®</b> -K1 Standard Basic	00190446751823
80.00.045	<b>ARTROMOT®</b> -K1 Classic	00888912403238
80.00.047	<b>ARTROMOT®</b> -K1 Classic Basic	00190446751830
80.00.023	<b>ARTROMOT®</b> -S3 Standard	00888912403108
80.00.024	<b>ARTROMOT®</b> -S3 Comfort	00888912403115
80.00.033	<b>ARTROMOT®</b> -E2 Compact	00888912403146
80.00.031	<b>ARTROMOT®</b> -E2	00888912403139
80.00.035	<b>ARTROMOT®</b> -SP3 Standard	00888912403153
80.00.036	<b>ARTROMOT®</b> -SP3 Standard Chip	00888912403160
80.00.037	<b>ARTROMOT®</b> -SP3 Comfort	00888912403177
80.00.038	<b>ARTROMOT®</b> -SP3 Comfort Chip	00888912403184
80.00.070	<b>ARTROMOT®</b> ACTIV E-K	00888912403283
80.00.084	<b>ARTROMOT®</b> -S4 Comfort	00190446278764
80.00.022	OptiFlex®-S	00190446794783
80.00.046	OptiFlex®-K1 Classic	00190446796015
80.00.048	OptiFlex®-K1 Standard	00190446796039
80.00.049	OptiFlex®-K1 Comfort Chip	00190446796046

### First change: Warning

New warnings have been implemented due to labeling changes that add or strengthen a warning, precaution, or include additional information.

Safety warnings:

#### **Warning!**

"The treatment of minors or persons with limited physical, sensory or mental abilities (e.g. movement disorders, epilepsy, seizures, mental disorders, behavioural disorders or neurodevelopmental disorders) may only be carried out under the constant observation of attendants (e.g. specialist staff, legal guardians or custodians)

– After data storage, write the patient's name on the patient chip card. The card should only be used for this particular patient. If the patient chip card is used for another patient, be sure to delete the previous patient data from the card first. Use original chip cards only."

## Second change: Label

### Symbols

Table 1: Symbols changes

Symbol change	Old version	New version	Description of symbol acc. ISO 15223 + 2014/
UDI	New requirement		Indicates a carrier that contains unique device identifier information. This symbol identifies the UDI carrier, including the AIDC and human readable information. (01) GTIN code of device (11) Date of production (21) serial number of device
Manufacturer + Date of production + country of manufacture	 Name Adress  2015	 YYYY-MM	Indicates the medical device manufacturer Indicates the date when the medical device was manufactured DE is the two or three letter country code as defined in ISO 3166-1. Indicates the country of manufacture
Manufacturer's name	<b>ORMED GmbH</b> A DJO Company	<b>ORMED GmbH</b> a company of Enovis	Name of the manufacturer
Manufacturer's address spelling	Bötzinger Straße 90 79111 Freiburg / Germany	Bötzinger Strasse 90 79111 Freiburg / Germany	Address of the manufacturer
Norm/Standard ETL classification	Conforms to ANSI/ AAMI ES STD 60601-1 Certified CAN/CSA C22.2 No. 60601-1	 InterTek 5029964 MEE	Medical electrical equipment relevant standards according ETL classification – control number above

### Other changes

Manufacturers contact was changed to: Tel. +49 761 4566-01 · Fax +49 761 4566-5501, medizintechnik@enovis.com, enovis-medtech.de

### Third change: Electromagnetic emissions, Electromagnetic immunity and RF (radio frequency) interference

The EMC table provides a detailed comparison of electromagnetic compatibility as per the latest standards.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
ARTROMOT® CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the ARTROMOT® CPM device unit is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions to CISPR 11: 2015+A1:2016+A2:2019	Group 1	<b>ARTROMOT®</b> CPM device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions to CISPR 11: 2015+A1:2016+A2:2019	Class B	<b>ARTROMOT®</b> CPM device is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions to IEC 61000-3-2: 2005+A1:2008 +A2:2009	complies with class A requirements	
Voltage fluctuations/flicker emissions to IEC 61000-3-3: 2013	complies	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

**ARTROMOT®** CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®** CPM device unit is used in such an environment.

Immunity test	IEC 60601-1-2:2014 +A1:2020		Electromagnetic environment – guidance
	test level Home healthcare	Compliance level Home healthcare	
Electrostatic discharge (ESD) to IEC 61000-4-2:2008	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/bursts to IEC 61000-4-4:2012	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power should be that of a typical commercial or hospital environment.
Surges to IEC 61000-4-5:2014	± 1 kV line-to-line voltage	± 1 kV line-to-line voltage	Mains power should be that of a typical commercial or hospital environment.
Voltage dips To IEC 61000-4-11:2017	0% Ut ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% Ut ; 1 cycle and 70% Ut ; 25/30 cycles Single phase at 0°	0% Ut ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% Ut ; 1 cycle and 70% Ut ; 25/30 cycles Single phase at 0°	Mains power should be that of a typical commercial or hospital environment. If the user of the <b>ARTROMOT®</b> CPM device unit requires continued operation during power mains interruptions, it is recommended that the <b>ARTROMOT®</b> CPM device unit be powered from an uninterruptible power supply or a battery.
Voltage interruptions to IEC 61000-4-11:2017	0% Ut; 250/300 cycles	0% Ut; 250/300 cycles	
Power frequency (50/60 Hz) magnetic field to IEC 61000-4-8:2009	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical commercial or hospital environment.

Note: Ut is the a.c. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

**ARTROMOT®** CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®** CPM device unit is used in such an environment.

Immunity test	IEC 60601-1-2:2014 +A1:2020		Electromagnetic environment – guidance
	test level Home healthcare	Compliance level Home healthcare	

Conducted RF to IEC 61000-4-6:2013	3 V 0.15 MHz-80 MHz 80% AM at 1kHz 6V 0.15 MHz-80 MHz 80% AM at 1kHz at ISM frequency bands and Amateur radio bands	3 V 0.15 MHz-80 MHz 80% AM at 1kHz 6V 0.15 MHz-80 MHz 80% AM at 1kHz at ISM frequency bands and Amateur radio bands	
Radiated RF to IEC 61000-4-3: 2006+A1:2007 +A2:2010	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communica- tions equipment to IEC 61000-4-3: 2006+A1:2007 +A2:2010	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	
Radiated fields in close proximity IEC 61000-4-39:2017	8 A/m 30 kHz No modulation (CW) 65 A/m 134.2 kHz PM 2.1 kHz 7,5 A/m 13.56 MHz PM 50 kHz	8 A/m 30 kHz No modulation (CW) 65 A/m 134.2 kHz PM 2.1 kHz 7,5 A/m 13.56 MHz PM 50 kHz	

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ARTROMOT®** CPM device exceeds the applicable **ARTROMOT®** CPM device unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **ARTROMOT®** CPM device unit.

The **ARTROMOT®** CPM movementsplint is intended for operation in the specified electromagnetic environment. The customer or user of the **ARTROMOT®** CPM device should ensure that it is used in such an environment. The required interference levels of the specified standards are met.

Hersteller/Manufacturer:

ORMED GmbH  
A company of Enovis  
Bötzingen Straße 90 · 79111 Freiburg · Germany  
Tel. +49 761 4566 -01 · Fax +49 761 4566 -5501  
medizintechnik@enovis.com  
www.enovis-medtech.de