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**Anhang:**  
ARTROMOT®-S3  
ARTROMOT®-S4  
ARTROMOT® ACTIVE-K  
ARTROMOT®-K1  
ARTROMOT®-SP3  
ARTROMOT®-E2

# Contents

<b>Figures ARTROMOT®-E2</b>	<b>III</b>
<b>Device description ARTROMOT®-E2</b>	<b>IV</b>
<b>Symbol overview ARTROMOT®-E2/-E2 compact</b>	<b>VI</b>
<b>Device description ARTROMOT®-E2 compact</b>	<b>IX</b>
<b>Figures ARTROMOT®-E2 compact</b>	<b>XI</b>
<b>1. How to use the CPM device</b>	<b>46</b>
1.1 Fields of application	46
1.2 Therapy objectives	46
1.3 Indications	46
1.4 Contraindications	46
<b>2. Description of the ARTROMOT®-E2/-E2 compact</b>	<b>47</b>
2.1 Description of the ARTROMOT®-E2 device components	48
2.2 Description of the ARTROMOT®-E2 compact device components	48
2.3 Description of the programming unit	49
2.4 Explanation of symbols	52
2.5 Explanation of symbols (connections and nameplate)	53
<b>3. Safety information</b>	<b>55</b>
<b>4. Device setup</b>	<b>59</b>
4.1 Connecting the ARTROMOT®-E2/-E2 compact, performance check	59
4.2 Adjusting the device to the patient	60
<b>5. Setting the treatment values</b>	<b>64</b>
5.1 General information on programming the ARTROMOT®-E2/-E2 compact	65
5.2 Programming the ARTROMOT®-E2/-E2 compact	66
5.3 Therapy parameter details	67
5.4 Application and programming examples	74
<b>6. Care, Maintenance, Transport, Conversion</b>	<b>76</b>
6.1 Care/Re-use	76
6.2 Maintenance (fuse replacement)	77
6.3 Transport	78
6.4 Conversion	81
<b>7. Environmental Protection Statement</b>	<b>83</b>
<b>8. Specifications</b>	<b>83</b>
<b>9. Contact</b>	<b>85</b>
<b>10. Technical Service</b>	<b>86</b>
10.1 Technical hotline	86
10.2 Shipment	86
10.3 Spare parts	86
<b>11. CE-Declaration of Conformity</b>	<b>87</b>
<b>12. IEC 60601-1-2 and IEC 60601-1-11</b>	<b>267</b>
12.1 Electromagnetic emissions	267
12.2 Electromagnetic immunity	268
12.3 Recommended separation distances	270

# 1. How to use the CPM device

## 1.1 Fields of application

The **ARTROMOT®-E2/-E2 compact** is a motor-operated **Continuous Passive Motion (CPM)** device providing motion to the elbow joint.

Suitable for use in hospitals, clinics, general practices and rental services as well as in the patient's home, it is an important supplement to medical and therapeutic treatment.

## 1.2 Therapy objectives

CPM therapy with the **ARTROMOT®-E2/-E2 compact** is mainly used to prevent the negative effects of immobilization, to allow patients to regain painless mobility of joints at an early stage and to promote healing and achieve a positive functional result.

Other objectives of therapy include:

- improvement of joint metabolism
- prevention of joint stiffness (arthrofibrosis)
- promotion of the regeneration and healing of cartilage and damaged ligaments
- faster hematoma/fluid resorption
- improved lymph and blood circulation
- thrombosis and embolism prophylaxis

## 1.3 Indications

The CPM device is indicated in the treatment of most injuries and postoperative conditions and diseases of the elbow joint. Examples:

- joint distortion and contusion
- arthrotomy and arthroscopy procedures in combination with synovectomy, arthrolysis or other intra-articular interventions
- all types of arthroplasty
- mobilization of joints in anesthetized patients
- exercise-stable fractures after surgery and pseudoarthrosis
- endoprosthetic implants
- myoplasty
- corrective osteotomy

## 1.4 Contraindications

**Do NOT use** the **ARTROMOT®-E2/-E2 compact** on patients with:

- acute inflammatory processes in the joints, unless on the order of a physician
- spastic paralysis
- unstable osteosynthesis

## 2. Description of the ARTROMOT®-E2/-E2 compact

The motorized CPM device supports the following movements of the elbow joint:

**Extension/flexion** 5° - 0° - 140°

**Pronation/supination** 90° - 0° - 90°

It can be reconfigured for use on either side.

### Note!

To unambiguously represent the current position of the CPM device, pronation and extension values below 0° are marked with the symbol "-" both on the display and in this document.

### These are some of the outstanding ARTROMOT®-E2/-E2 compact features:

- anatomically correct adjustment
- physiological movements
- maximum possible ranges of motion
- programming unit for precise adjustment of patient-specific therapy values
- chip card for storage of the programmed therapy parameters
- easy transport

### Biocompatibility

Those parts of the ARTROMOT®-E2/-E2 compact that come into contact with the patient when the device is used as intended, are designed to fulfil the biocompatibility requirements of the applicable standards.

### Essential Performance Features

- The programmed angles are accurately maintained with a tolerance of +/- 2°.
- The programmed speeds are accurately maintained with a tolerance of +/- 5%.
- The selected mode and mechanical settings do not change during operation.

### Frequent functions:

- Unpacking (device and accessories)
- Inserting the motion element and locking it at the system rack
- Mounting the armrest
- Device connections (mains power connection, connecting the control pendant and the motion element)
- Mechanical adjustments of the motion element on the rotational axes
- Programming the control pendant (range of motion, speed, functions)
- Storage

## 2.1 Description of the ARTROMOT®-E2 device components

**Note:** See device description on pages IV and V.

1. Knob for horizontal adduction/abduction (horizontal extension/flexion)
2. Screw for double joint
3. Double joint
4. Screw for height adjustment
5. Upper arm support
6. Locking pin for height adjustment of upper arm support
7. Screw for elbow angle adjustment
8. Cam lever for forearm length adjustment
9. Slider
10. Screw for adjustment of the backrest angle
11. Pin for swivel motion of upper arm support
12. Screw for forearm module swivel motion
13. Locking screw for right/left positioning
14. Armrest for healthy arm
15. Programming unit
16. Compartment for storage of programming unit
17. Patient chip card
18. Connection for power cord
19. Power switch (ON/OFF)
20. Fuse
21. Connectors for motion element
22. Locking pins
23. Straps for forearm support
24. Motor A
25. Motor B
26. Palmar hand support
27. Dorsal hand support
28. Folds down the backrest (transport position)
29. Wheels
30. Bracket for adjustment of the positioning angle
31. Strap for upper arm support
32. Connector for programming unit

Subject to technical modifications

## 2.2 Description of the ARTROMOT®-E2 compact device components


**Note:** See device description on pages IX and X.

1. Screw for height adjustment
2. Screw for elbow angle adjustment
3. Cam lever for forearm length adjustment
4. Slider
5. Pin for swivel motion of upper arm support
6. Screw for forearm module swivel motion
7. Locking screw for right/left positioning
8. Fuse
9. Connectors for motion element
10. Straps for forearm support
11. Motor A
12. Motor B
13. Palmar hand support
14. Dorsal hand support
15. Programming unit
16. Connector for programming unit
17. Patient chip card
18. Connection for power cord
19. Power switch (ON/OFF)
20. Upper arm support
21. Strap for upper arm support
22. Wheel locks
23. Compartment for storage of programming unit
24. Locking screw
25. Wheels
26. Pin for height adjustment of upper arm support
27. Bracket for adjustment of the positioning angle

Subject to change without notice

## 2.3 Description of the programming unit


### 2.3.1 Programming unit in normal mode

### 2.3.2 Programming unit in MENU selection mode

### 2.3.3 Programming unit in ROM programming mode (range of motion)


### 2.3.4 Programming unit in general programming mode

## 2.4 Explanation of symbols

Also refer to symbol overview on page VI and VII.

	extension
	flexion
	pronation
	supination
	pause
	timer
	speed
	warm-up protocol
	load reversal motor A
	load reversal motor B
	motor A ON/OFF
	motor B ON/OFF
	isolation program
	total therapy time
	extension/flexion therapy documentation
	pronation/supination therapy documentation

	stretching in extension
	stretching in flexion
	stretching in pronation
	stretching in supination
	EROM repeat extension/ pronation
	EROM repeat flexion/ supination
	transport setting
	new patient
	synchronized/ non-synchronized mode
	service menu

## 2.5 Explanation of symbols (connections and nameplate)

Alternating current

for serial numbers < 20,000:  
protection class I equipment. The medical device must be connected to a system with protective earth conductor!

for serial numbers > 20,000:  
protection class II equipment. The medical device has a double or reinforced insulation.

Type B applied part

Power switch OFF

Power switch ON

The year and month of manufacture are indicated next to this factory symbol.

The name next to this factory symbol is the manufacturer.

The number next to this symbol is the article reference number.

Device complies with Council Directive 93/42/EEC about medical devices, tested and approved by DQS Medizinprodukte GmbH.

The number next to this symbol is the serial number.

Caution! Observe warnings set forth in operation manual!

Observe Operation Manual! ISO 7010-M002

Observe Operation Manual! ISO 7000-1641

Do not dispose with unsorted municipal waste.

Protect from moisture.

## IP21

The IP rating indicates the level of protection and thus the suitability of the device for use under different ambient conditions.

The rating IP21 means:

2 is the level of protection against contact and solid objects.

The digit 2 means:

- Protection from contact: protected from contact with a finger.
- Protection against foreign object: protected against solid foreign object (diameter of 12.5 mm and greater).

1 indicates the degree of protection against water.

The digit 1 means: protection against vertically falling water drops.

Observe temperature limits (storage)! ISO 7000-0632

The overall weight of the device is indicated next to this symbol.

Do not rest or lean on or push the assembled device because it might tip over.

### 3. Safety information

#### Definitions

Read the safety statements before use of the CPM device. The safety statements are classified as follows:

#### **Danger!**

Indicates an imminent hazard. If not avoided, this hazard will result in death or serious injury.

#### **Warning!**

Indicates a hazard. If not avoided, the hazard can result in death or serious injury.

#### **Caution!**

Indicates a potential hazard. If not avoided, the hazard may result in minor injury and/or product/property damage.

#### Safety information

#### **Danger!**

Explosion hazard –

The **ARTROMOT®-E2/-E2 compact** is not designed for use in areas of medical locations where an explosion hazard may occur. An explosion hazard may result from the use of flammable anesthetics, skin cleansing agents and disinfectants.

#### **Warning!**

Patient hazard –

- **Only authorized individuals** are allowed to operate the ARTROMOT®-E2/-E2 compact. Individuals are authorized after receiving training in the operation of the device and reading this operation manual.
- Before using the device, the operator must ascertain that it is in correct working order and operating condition. In particular, the cables and connectors must be checked for signs of damage. Damaged parts must be replaced immediately, before use.
- **Before therapy, a test run** consisting of several exercise cycles must be completed, first without and then with the patient. Check that all setting screws are tightened.
- Stop therapy immediately, when you have doubts about the device settings and/or the therapy protocol.

## **Warning!**

Patient hazard –

- Ensure an **anatomically correct setup of the CPM device** suitable for the patient to be treated. Therefore, carefully verify the following settings/positions:  
(also refer to the numbers on the device):
  1. Horizontal adduction/abduction (horizontal extension/flexion)
  2. Height adjustment
  3. Adjustment of the positioning angle
  4. Adjustment of the forearm length
  5. Axis adjustment motor A and motor B
  6. Backrest adjustment
  7. ROM adjustment
- It is not permitted to change the adjustment of 1 through 6 while the CPM device is applied on a patient.
- Movements must **not cause any pain or irritation**.
- Patients must be **fully conscious** while being instructed in the use of the CPM device and during therapy.
- The **choice of the therapy parameters** to program and of the **therapy protocols** to use is restricted to the responsible **physician or therapist**. It is the physician's or therapist's decision whether or not to use the CPM device on a specific patient.
- The patient must be familiar with the functions of the ARTROMOT®-E2/-E2 compact programming unit and the unit must be **within easy reach** of the patient, allowing him or her to stop therapy, if needed. **Patients unable to operate the programming unit**, e.g. paralytic patients, must never be left unattended during therapy.
- Write the patient's name on the **patient chip card**. The card should only be used for this patient. If the **patient chip card** is used for another patient, be sure to delete the **previous patient's data** from the card first (see sections 4.1 and 5.3, paragraph "New Patient"). Use **original chip cards only**.

- All **accessories** used with the ARTROMOT®-E2/-E2 compact must be approved by DJO.
- Modifications to the medical device described in this document without the manufacturer's written consent is prohibited.
- Ensure that the medical device is securely set up on solid ground during use.
- Do not allow **parts of the body or any objects** (such as blankets, cushions or cables) to get caught **in the moving parts** of the CPM device.
- Always seat the patient on a chair with four legs and without armrests when using the **ARTROMOT®-E2 compact**.

## **Warning!**

- Extreme caution should be taken when in use around small children and babies! Sufficient distance to the device is mandatory for their safety!
- Never leave the device unattended when it is switched on! Switch the device off and disconnect the power cord from the wall outlet!
- After use, store the device in a safe place! Ensure device stability also during storage!

## **Warning!**

Shock hazard –

Strictly observe the following warnings. Failure to do so endangers the lives of the patient, the user and other persons involved.

- **Before use** allow the ARTROMOT®-E2/-E2 compact to reach room temperature. If the device has been transported at **temperatures below 0 °C (32°F)**, leave it to dry at room temperature for about 2 hours, until any condensation has disappeared.
- The ARTROMOT®-E2/-E2 compact must only be operated in **dry rooms**.
- When disconnecting the device from the power line, remove the plug from the wall outlet first, before disconnecting the cable from the device.
- When connecting the device to other equipment or when creating a medical system, check that the sum of leakage currents will not cause any hazard. Please contact DJO, if you have questions in this matter.
- Do not use multiple portable socket outlets (MPSO) to connect the device to the power line. The ARTROMOT®-E2/-E2 compact must be connected to a **properly installed wall outlet with a non-fused earthed wire**. Before connecting the power cord, it must be completely unrolled and placed such that it will not get caught in the moving parts of the device.
- Before cleaning and service interventions, **disconnect the device from the power line by removing the power cord from the wall outlet**.
- **Liquids must not be allowed to enter the CPM device** or the **programming unit**. If liquids have entered into the device, the ARTROMOT®-E2/-E2 compact must be checked by a service technician, before it can be reused.

## **Warning!**

Equipment malfunction –

- Magnetic and electrical fields are capable of interfering with the proper performance of the device. For this reason make sure that all external devices operated in the vicinity of the CPM device comply with the relevant EMC requirements. X-ray equipment, MRI devices, radio systems, cell phones, etc. are possible sources of interference as they may emit higher levels of electromagnetic radiation. Keep the CPM device away from these devices and verify its performance before use.
- Refer **repair** and **maintenance** to authorized persons.
- **Route all cables** such that they will not get caught by the moving parts during operation and that they do not present a stumbling hazard.
- **Inspect** the ARTROMOT®-E2/-E2 compact for damage and loose connections at least once a year.

## **Caution!**

Preventing chafing and pressure sores –

If your patient is **adipose**, very **tall** or very **short**, be sure to prevent chafing and pressure sores.

Patient hazard, damage to the CPM device –

The CPM device **ARTROMOT®-E2** must not be used to transport people.

## **Caution!**

Equipment damage –

- Check that the voltage and frequency ratings of your local **power line** are those indicated on the nameplate.
- The **seat** withstands a **maximum continuous load** of **150 kg / 330 lb**.
- The **arm support element** withstands a **maximum continuous load** of **9 kg / 20 lb**.
- Do not allow **any objects** (such as blankets, cushions, or cables) to get caught in the **moving parts** of the CPM device.
- Do not expose the ARTROMOT®-E2/-E2 compact to direct sunlight, because some of the components may reach inadmissibly high temperatures. Furthermore, the device must be set up at a safe distance from radiators to prevent excessive temperature rises.
- The presence of children, pets and rodents does not normally impair the functioning of the device. However, avoid contamination of the device by children or animals, from dust and lint, and keep them at a safe distance from the device. The safety statements set forth apply.
- Be aware that the connectors can only be inserted in the correct orientation and secure all connections with the locks.

## 4. Device setup

**Note:** See device description on pages IV/V and IX/X.

### 4.1 Connecting the ARTROMOT®-E2/-E2 compact, performance check

The equipment supplied includes these items:

Base unit (chair), motion element, programming unit (15), patient chip card (17), arm rest (4), power cord (not shown), operating instructions

1. Set up the device and connect the required cables by reversing the transport instructions given in section 6.3.
2. Follow these steps to set the carriage to the **home position**:

#### Initial adjustment for new patients

Write the patient's name down on the back of the chip card. Insert the original patient chip card (17) into the programming unit (15).

Press the **MENU** key on the programming unit seven times in rapid succession or give it one long press to access programming level seven (with each key press, you advance one level).

Press the "new patient" parameter key **→0←** and select this function (a check mark appears in the circle next to the function).

Press the **START** key. The CPM device automatically enters the **home position**.

#### Adjustment with programmed chip card

Insert the original patient chip card (17) into the programming unit (15).

Press the **START** key.

The CPM device automatically moves to the **starting position** (middle position of the set extension/flexion and pronation/supination values).

#### Performance check

If the programming unit can be operated as described above and the ARTROMOT®-E2/-E2 compact moves to the home position (for home position values, please refer to sections 5.3), the CPM device has passed the performance check.

The device also runs performance checks regularly during operation. This is what happens, if a problem is identified:

- An audio signal sounds.
- The device switches off immediately.
- The message "ERROR" and an error code (e.g. ERROR 5) appear on the display.

In this situation, you may attempt to restart the device by turning it briefly off and on again with the power switch. If the error message persists, have the device inspected by a Service technician, before using it again.

If it has been determined that the ARTROMOT®-E2/-E2 compact operates perfectly, ask the patient to sit down on the ARTROMOT®-E2/-E2 compact chair.

## 4.2 Adjusting the device to the patient

### Note!

For the following adjustments, the patient's arm should not yet be placed on the armrest. After you have adjusted the device to the approximate patient measurements, the patient can place his/her arm on the armrest for a check of the setup and for the fine adjustment.

The settings are numbered 1 through 5. For easy orientation, you will find the same numbers on the CPM device. Always set up the CPM device in this order.

Write the settings down on the back of the patient's chip card.

Before you adjust the ARTROMOT®-E2/-E2 compact to the patient, you may have to convert the device for use on the left or right elbow joint. (see section 6.4 Conversion)

Before any treatment is possible, adjust the device as follows:

### ARTROMOT®-E2: backrest, armrest

Before you start setting up the device for the exercise protocol, you must adjust it to a position that is anatomically correct for the respective patient.

- Open thumbscrew (10) and adjust the backrest to a position that is comfortable for the patient (Fig. F).
- Set the armrest for the healthy arm to a height that allows the patient to sit up straight (locking screw 14).

### ARTROMOT®-E2 compact: positioning the chair and the CPM device

Before you start setting up the device for the exercise protocol, you must adjust it to a position that is anatomically correct for the respective patient.

- Only use chairs with four legs and without armrests. The patient must be seated safely and straight on the chair.
- Once you have set up the ARTROMOT®-E2 compact, lock the brakes (22) of the castor wheels (25).

### 4.2.1 Adjusting the ARTROMOT®-E2 to the patient

**Note:** The individual steps are illustrated on page III.

#### 1 a,b and c Horizontal Extension/Flexion (Fig. A)

Horizontal adduction/abduction is set manually. It is the purpose of the adjustment to achieve congruent axes of the motor axis and the elbow joint.

- To begin with, press knob (1). Leave it engaged in the open position.
- Now open screw (2). The double joint now moves freely.
- Set the motion element with motor A to the appropriate position. Make sure that a congruence of the axis of motor A and the elbow joint is achieved.
- Push knob (1) in the opposite direction.
- Tighten screw (2).

#### 2 Height adjustment (Fig. B)

Before opening locking screw (4), hold the moving part of the CPM device by grasping the double joint (3), to prevent the part from falling down.

- Open screw (4) and set the motion element to the appropriate height. Ensure that the entire upper arm lies on the support (5). The patient should sit up straight and be relaxed.
- Tighten screw (4).
- Now check the agreement of the axis of motor A and the elbow joint. If required, make a fine adjustment of the vertical orientation, using locking pin (6).

#### 3 Adjustment of the positioning angle (Fig. C)

##### **Caution!**

Patient hazard / equipment damage –

**Hold motor B for this adjustment to secure the device.**

- Hold the motion element at bracket (30).
- Now open screw (7) and set the appropriate angle.
- Tighten screw (7).

#### 4 Adjustment of the forearm length (Fig. D)

##### **Caution!**

Patient hazard / equipment damage –

**Hold motor B for this adjustment to secure the device.**

- Open the cam lever (8) and set the appropriate length. Make sure that the slider (9) is free to move at least 2 cm in either direction. There must be enough space for the fingers between the hand support and motor B.
- Tighten the cam lever (8).

#### 5 Axis adjustment motor B (Fig. E)

In most cases the setting is 0 on the scale. In cases of a deformed forearm, the setting may have to be changed.

- Open the locking screw (13) and perform the adjustment.
- Tighten the locking screw (13).

#### 6 Backrest adjustment (Fig. F)

Set the backrest to the foremost position before setting the motion element to an angle of 0° for horizontal adduction/abduction. To do so, open adjusting screw

(10), set the backrest to the foremost position and tighten the screw.

To achieve optimal congruence between the axis of motor A and the patient's rotational axis of the elbow joint, the backrest angle can be individually adjusted in all other horizontal adduction/abduction positions.

### Checking the setup, fine adjustment

Please check the following before using the CPM device on a patient:

- Check settings 1 to 6 and make sure that the pivot of motor A agrees with the pivot of the elbow joint, i.e., that the axis of motor B vertically passes through the elbow joint.
- Check that all locking screws and clamping levers are tight.

### 4.2.2 Adjusting the ARTROMOT®-E2 compact to the patient

**Note:** The individual steps are illustrated on page XI.

**Note:** Step 1 does not exist for the ARTROMOT®-E2 compact.

#### 2 Height adjustment (Fig. A)

Before opening locking screw (1), hold the moving part of the CPM device by grasping motor A (11), to prevent the part from falling down.

- Open screw (1) and set the motion element to the appropriate height. Ensure that the entire upper arm lies on the support (20). The patient should sit up straight and be relaxed.
- Tighten screw (1).
- Now check the agreement of the axis of motor A and the elbow joint. If required, make a fine adjustment of the vertical orientation, using locking pin (26).

#### 3 Adjustment of the positioning angle (Fig. B)

##### ⚠ Caution!

Patient hazard / equipment damage –

**Hold motor B for this adjustment to secure the device.**

- Hold the motion element at bracket (27).
- Now open screw (2) and set the appropriate angle.
- Tighten screw (2).

#### 4 Adjustment of the forearm length (Fig. C)

##### ⚠ Caution!

Patient hazard / equipment damage –

**Hold motor B for this adjustment to secure the device.**

- Open the cam lever (3) and set the appropriate length. Make sure that the slider (4) is free to move at least 2 cm in either direction. There must be enough space for the fingers between the hand support and motor B.
- Tighten the cam lever (3).

#### 5 Axis adjustment motor B (Fig. D)

In most cases the setting is 0 on the scale. In cases of a deformed forearm, the setting may have to be changed.

- Open the locking screw (7) and perform the adjustment.
- Tighten the locking screw (7).

### Checking the setup, fine adjustment

Please check the following before using the CPM device on a patient:

- Check settings 2 to 5 and make sure that

the pivot of motor A agrees with the pivot of the elbow joint, i.e., that the axis of motor B vertically passes through the elbow joint.

- Check that all locking screws and clamping levers are tight.

## 5. Setting the treatment values

### **Warning!**

Patient hazard –

**Before therapy, a test run** consisting of several exercise cycles must be completed without the patient. Then repeat the test run with the patient and check that the movement does not cause any pain.

**Note:** See also 2.2 and 2.3, as well as page VI/VII.

### **Note!**

Programming is only possible when the patient's chip card has been inserted.

For information about therapy parameters and programming of special functions, please refer to sections 5.1 to 5.3.

See section 5.4 for programming examples.

### **Important!**

The programming unit **ARTROMOT®-E2/-E2 compact "Graphics"** can be connected to all products of the **ARTROMOT® family of elbow CPM machines**.

It is possible to exchange chip cards between the "Text" and "Graphics" versions.

#### **Please note:**

1. If a "Graphics" programming unit is operating with a formatted "Text" chip card
  - the **selected language** of the "Text" version remains activated; it is irrelevant for the "Graphics" version.
  - the **reverse on load** setting saved in the "Text" version is automatically used for both motors.
  - the **key lock** is disabled.
2. If a "**Text**" programming unit is operating with a formatted "Graphics" chip card
  - the **selected language** automatically changes to the default setting: German.
  - the **reverse on load** setting saved for motor A in the "Graphics" version is automatically used for both motors.
  - a **key lock** set in the "Graphics" version will remain active, only when no changes were made with the "Text" version; the key lock function is **not an option with the "Text" version**.

## 5.1 General information on programming the ARTROMOT®-E2/-E2 compact

1. You activate the programming mode by briefly pressing the **MENU** key on the programming unit.
2. The treatment parameters and functions are allocated to seven programming levels (four functions per level).

To be able to program a parameter you will have to access the corresponding programming level. This is also done with the **MENU** key. With each brief key press you advance one level. The code M1, M2, etc. that appears in the middle of the display indicates the programming level. To back up one level, e.g. from level 2 to 1 or from level 1 to 7, hold the MENU key down for some seconds.

3. You activate the treatment parameters and functions with the **four parameter keys** below the display. The symbols above the four parameter keys indicate the assigned parameters and functions.

This is what happens when you press one of the parameter keys to select a parameter:

- The corresponding symbol appears on the display in a larger format.
- The set value is displayed.
- The symbol above the parameter key appears in reverse video.

4. With the + / - keys (plus/minus) you change the displayed value. When you press and hold the key, the value will change at a higher rate.

Some of the (special) functions can only be enabled and disabled. This is done by pressing the corresponding parameter key or with the + / - keys. Active parameters are identified with a check mark in the circle next to the symbol.

5. Having programmed all parameters, press the **STOP** key to save the values.
6. Then press the **START** key: the ARTROMOT®-E2/-E2 compact checks the set values, moves to the position halfway between the set extension and flexion values as well as to the set supination and pronation values and stops.

7. Press the **START** key again to start therapy.

The CPM device will now begin with the set physiological motion pattern in the **synchronized mode**. The motor with the larger range of motion operates at the programmed speed, while the speed of the second motor is adapted. This ensures that both motors will reach their target positions at almost the same time.

After activation of the **START** key in the **non-synchronized mode** the motors will perform random movements, each motor reversing after reaching the maximum values. The programmed speed is valid for both motors.

### Note!

- Refer to section 5.3 for a description of the parameters.
- To **view the set parameter values**, press the corresponding parameter key. Before, however, you have to press the **STOP** key and access the correct menu level.
- To prevent accidental changes of the parameter settings, you can **lock the keys**. To do so, simultaneously press **keys + and - for approx. 3 seconds**.

Press both keys again for approx. 3 seconds to unlock.

- Selecting the "new patient" function will automatically delete the data on the patient chip card. When you have finished programming the unit and press the **STOP** key, the settings will automatically also be saved to the patient chip card.
- **Emergency stop function:** When any of the keys is pressed during therapy, the ARTROMOT®-E2/-E2 compact switches immediately off.
- When the **START** key is pressed again in the synchronized mode and in the non-synchronized mode, the carriage will reverse its direction.
- The manufacturer recommends an application period of up to 1 hour per therapy session.

- For a uniform, smooth and gentle transition between the two directions of motion, the speed is automatically reduced before reaching the reversal point and, on reaching the reversal point, it is continuously increased until it reaches the set speed.

### Patients with a programmed chip card

- First complete the mechanical adjustments.
- Then insert the chip card (the patient is not yet positioned on the CPM device).
- Press the **START** key: the device will move to the starting position as specified by the parameters stored on the chip card and stops.
- Position the patient on the CPM device and press the **START** key to initiate therapy.

### 5.2 Programming the ARTROMOT®-E2/-E2 compact

Different programming levels are provided to program the ARTROMOT®-E2/-E2 compact. You change between levels by pressing the **MENU key** repeatedly.

The currently selected level is indicated on the display.

The following **treatment values, settings and information** can be entered/viewed on the programming unit (15):

#### LEVEL 1:

- extension
- flexion
- pronation
- supination

#### LEVEL 2:

- pause
- therapy timer

- speed
- warm-up protocol

#### LEVEL 3:

- reverse on load - motor A
- reverse on load - motor B
- motor A ON/OFF
- motor B ON/OFF

#### LEVEL 4:

- isolation protocol
- total therapy time
- therapy documentation: extension/flexion
- Therapy documentation: pronation/supination

#### LEVEL 5:

- stretching in extension
- stretching in flexion
- stretching in pronation
- stretching in supination

#### LEVEL 6:

- EROM repeat extension/ pronation
- EROM repeat flexion/supination

#### LEVEL 7:

- transport setting
- new patient
- operating mode synchronized/ non-synchronized
- Service menu

## Note!

- **While you adjust the values**, the carriage will move to the set range. This allows you to easily and quickly determine the ROM where the patient does not experience any pain.
- The angle entered last for the respective directions of motion will be saved in each case.

## 5.3 Therapy parameter details

- To access the programming level you need, repeatedly press the **MENU key**.
- The treatment parameters are selected with the different **parameter keys**.
- To change the treatment values, press the **+ / - keys**.
- To enable/disable a function, press the respective **parameter key** again.
- To save the entered values, press the **STOP key**.

### LEVEL 1:

- **extension**  maximum value: - 5 degrees
- **flexion**  maximum value: 140 degrees
- **pronation**  maximum value: - 90 degrees
- **supination**  maximum value: 90 degrees

## Note!

- The programmed values and the values (degrees) measured on the patient may deviate slightly.
- To ensure a physiological motion, the carriage simultaneously moves to the following positions in the synchronized mode:
  - maximum extension value and maximum pronation value
  - maximum flexion value and maximum supination value.

- To ensure that the movements you program are physiological in nature, we recommend to program the carriage in the following order:

1. maximum extension value

2. maximum pronation value

3. supination  = 0°  
(to relieve the elbow)

4. maximum flexion value

5. maximum supination value

- When the **START** key is pressed in the synchronized mode, the carriage will first move to the position halfway between the set maximum extension/flexion values as well as the pronation/supination values. The carriage will stop in this position. Treatment begins when you press the **START key** again.
- Horizontal adduction/abduction is adjusted manually. These are possible settings:

**horizontal adduction :** 120°

**horizontal abduction:** 0°

### LEVEL 2:

#### ■ Pauses

Pauses occur at each programmed maximum value.

These are the two pause points:

- maximum extension value and maximum pronation value

- maximum flexion value and maximum supination value.

Pauses are separately adjustable for motor A and motor B between 1 second and 60 minutes. Pauses can be set in steps of 1 second to any value between 0 and 59 seconds, and in steps of 1 minute to any value between 1 and 60 minutes.

**Default: no pause**

Follow these steps to program pauses:

- Access menu level 2, then press the parameter key to select the special function .
- The display shows the symbol for the special function and the selected pause durations.

The display shows the symbol for the special function and the selected pause durations.


top: flexion/supination pause .

bottom: extension/pronation pause .

The flexion/supination pause duration appears in a box.

- To change the flexion/supination pause, press the "+" or "-" key.

Select a duration of 10 seconds, for instance.


- Press the **parameter key** again when done. The frame automatically moves to the line below, indicating that the extension/pronation pause can be adjusted .

The **parameter key** is only used to toggle between these two functions.

- To modify the value, press the "+" (plus) or "-" (minus) key.

Select a duration of 5 seconds, for instance.

- Save the programmed values by pressing the **STOP** key, then press the **START** key to initiate the treatment.

## ■ Therapy timer

**Default setting is continuous operation.**

A clock symbol in the upper right-hand corner of the display identifies the **continuous mode of operation**. The clock indicates the elapsed therapy time.

In the **continuous mode**, the device must be stopped with the **STOP** key.

However, the therapy timer can be set in **steps of 1 minute to any value between 1 and 59 minutes and in steps of 30 minutes to any value between 1 and 24 hours**. When the time has elapsed, the device switches automatically off and stops in the middle between the set limit values.

In this case, a circle replaces the clock symbol. The circle fills as the therapy time progresses.

## ■ Speed

The speed can be adjusted between 1 % and 100 % in steps of 1 %.

1 % is equivalent to 14°/minute  
100 % is equivalent to 230°/minute

**Default: 100 %**


## ■ Warm-up protocol

During warm-up, the patient will slowly become used to the set maximum values for extension/flexion and pronation/supination.

The warm-up protocol starts in the middle between the two maximum values set for extension/flexion and pronation/supination. The range of motion increases with each cycle, until the programmed maximum values are reached after a total of 15 cycles.

After reaching the maximum values, the device automatically enters the normal mode.

If a therapy time is activated, the system restarts with the warm-up protocol when the therapy time has elapsed and the devices is started again.

During warm-up, the display will show the symbol  in the upper, left-hand corner.

**Default: disabled**

### LEVEL 3:

#### ■ Reverse on load - motor A A (safety feature)

The device automatically reverses both motors when the patient's resistance (load) exceeds the set value.

Adjustable levels for reverse on load feature: 1 - 25. At level 1, very low resistance will cause the device to reverse; at level 25, a high resistance is required to initiate the reversal.

**Default: level 25**

#### ■ Reverse on load - motor B B (safety setting)

The device automatically reverses both motors when the patient's resistance (load) exceeds the set value.

Adjustable levels for reverse on load feature: 1 - 25.  
At level 1, very low resistance will cause the device to reverse; at level 25, a high resistance is required to initiate the reversal.

**Default: level 25**

### **Caution!**

Patient hazard –

The reverse on load feature is a safety measure to protect the patient in the event of cramps, spasms, locked joints and similar situations. The manufacturer cannot be held liable for misuse of this feature.

#### ■ Motor A ON/OFF A

To permit fully isolated movements, the motors can individually be turned on and off. Motor A effects the extension and flexion movements, while motor B effects the pronation and supination movements.

For an isolated pronation/supination movement, set motor A to the appropriate position (extension/flexion), then turn the motor off.

In normal mode the display indicates the symbol "OFF" during therapy, instead of the programmed extension/flexion angles.

**Default: motor A ON**

#### ■ Motor B ON/OFF B

To permit fully isolated movements, the motors can individually be turned on and off. Motor A effects the extension and flexion movements, while motor B effects the pronation and supination movements.

For an isolated extension/flexion movement, set motor B to the appropriate position (pronation/supination), then turn the motor off.

In normal mode the display indicates the symbol "OFF" during therapy, instead of the programmed pronation/supination angles.

**Default: motor B ON**

### Note!

Please note that one motor (A or B) must always be turned on. Otherwise the following symbol will be displayed when you press **START**:  OFF 

### LEVEL 4:

#### ■ Isolation protocol AIII B

For this special function, both motors are turned on, but they will never move simultaneously.

These are the steps of the function:

- To begin with, motor A moves to the programmed maximum values for extension and flexion for three cycles and then stops. In this phase, motor B is disabled (symbol for motor B: **OFF**)
- Then motor B will move to the programmed maximum values for pronation and supination for one cycle and stops, then the full cycle will start over with motor A. Motor A is disabled while motor B is operating (symbol for motor A: **OFF**)
- The stop position at the end of the last cycle (for both extension/flexion and pronation/supination) can be selected in steps of 25% between 0% and 100% of the programmed maximum range of motion.
- Steps 1 and 2 can be repeated as often as needed. You can stop treatment with the STOP key; after the pro-

grammed therapy duration, it will stop automatically.

During warm-up, the display will show the symbol  in the upper, left-hand corner.

**Default: disabled**

Follow these steps to program the special function:

- Access menu level 4, then press the parameter key to select the special function **A|||B**.


The display will indicate:

- the symbol representing the special function
  - the status of the special function (circle with or without check mark)
  - two boxes with the current percentage of the stop positions for motor A and B
  - a check mark in a circle.
- You activate the function with the "+" key and deactivate it with the "-" key.

In this case it is not possible to activate/deactivate the function by pressing the parameter key again.

A check mark in the circle identifies the activated function.

- Next press the parameter key of the isolation protocol. The marking moves to the upper percentage value, allowing you to enter the stop position for motor A.
- To modify the value, press the "+" (plus) or "-" (minus) key.  
Select 75%, for instance.
- Press the **parameter key** again when done. The marking automatically moves to the lower box, allowing you to adjust the stop position of motor B.
- To modify the value, press the "+" (plus) or "-" (minus) key.  
Select 25%, for instance.
- Save the programmed values by pressing the **STOP** key, then press the **START** key to begin the treatment.

**This is what happens:**

Extension/flexion is treated three times first (motor A). Motor B is positioned at 25% of the programmed range of motion for pronation and supination.

Subsequently, pronation/supination is treated once (motor B). Motor A is positioned at 75% of the programmed range of motion for extension and flexion.

### Note!

The percentage values can be modified only when the special function is enabled (check mark).

#### ■ Total therapy time

Under menu item "total therapy time" you can view each patient's total therapy time (duration of all the patient's therapy sessions).

#### ■ Deleting the stored therapy time:

Press and hold the **parameter key** for 5 seconds or select the "new patient" function.

#### ■ Extension/flexion therapy documentation

This is a special function of the ARTROMOT®-E2/-E2 compact which allows the entire therapy documentation to be represented.

The carriage run times as well as the range of motion of the sessions are recorded.

The collected data are presented graphically in the form of a coordinate system (X-axis = range of motion / Y-axis = time) where the upper curve illustrates the trend of the flexion movement and the lower curve the trend of the extension movement.

### ■ Pronation/supination therapy documentation

This is a special function of the ARTROMOT®-E2/-E2 compact which allows the entire therapy documentation to be represented.

The carriage run times as well as the range of motion of the sessions are recorded.

The collected data are presented graphically in the form of a coordinate system (X-axis = range of motion / Y-axis = time) where the upper curve illustrates the trend of the supination movement and the lower curve the trend of the pronation movement.

## LEVEL 5:

### ■ Stretching in extension

With the special "stretching in extension" function the joint will be gently stretched (arm stretching movement). Only the extension/flexion movement is treated here. Motor B is automatically disabled. The position of motor B cannot be changed once this special function has been activated.

Starting at the middle position the carriage will first move to the programmed flexion limit and then to the programmed extension limit.

Subsequently the carriage reverses 5° towards flexion and then moves very slowly back to the programmed extension value (symbol: ). Subsequently the carriage attempts to attain another 5° at an even slower speed (symbol: ).

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the carriage moves in the opposite direction.

This stretch cycle is repeated ten times.

After that the carriage moves to the programmed flexion limit and restarts the stretching in extension cycle.

This sequence can be repeated as often as needed. You can stop treatment with the STOP key; after the programmed therapy duration, it will stop automatically.

**Default: disabled**

## Note!

- If a pause has been programmed, the carriage will stop for the pause each time the maximum stretching value is attained.
- Enabling the stretching in extension function automatically disables motor B. The only movement performed is extension/flexion.
- It is not possible to combine a second stretching function with the special function "stretching in extension".
- During operation, the display will show the symbol  in the upper, left-hand corner.

### ■ Stretching in flexion

With the special "stretching in flexion" function the joint will be gently stretched (arm bending movement). Only the extension/flexion movement is treated here. Motor B is automatically disabled. The position of motor B cannot be changed once this special function has been activated.

Starting at the middle position the carriage will first move to the programmed extension limit and then to the programmed flexion limit.

Subsequently the carriage reverses 5° towards extension and then moves very slowly back to the programmed flexion value (symbol: ). Then the carriage attempts to attain another 5° at an even slower speed (symbol: ).

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the carriage moves in the opposite direction.

This stretch cycle is repeated ten times.

After that the carriage moves to the programmed extension limit and restarts the stretching in flexion cycle.

This sequence can be repeated as often as needed. You can stop treatment with the STOP key; after the programmed therapy duration, it will stop automatically.

**Default: disabled**

## Note!

The percentage values can be modified only when the special function is enabled (check mark).


**Note!**

- If a pause has been programmed, the carriage will stop for the pause each time the maximum stretching value is attained.
- Enabling the stretching in flexion function automatically disables motor B. The only movement performed is extension/flexion.
- It is not possible to combine a second stretching function with the special function "stretching in extension".
- During operation, the display will show the symbol  in the upper, left-hand corner.

### ■ Stretching in pronation

With the special "stretching in pronation" function the joint will be gently stretched in the pronation direction. Only the pronation/supination movement is treated here. Motor A is automatically disabled. The position of motor A cannot be changed once this special function has been activated.

Starting at the middle position the carriage will first move to the programmed supination position and then to the programmed pronation position.

Subsequently the carriage reverses 5° towards supination and then moves very slowly back to the programmed pronation value (symbol: ). Subsequently the carriage attempts to attain another 5° at an even slower speed (symbol: .

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the carriage moves in the opposite direction.

This stretch cycle is repeated ten times.

After that the carriage moves to the programmed maximum supination position and restarts the stretching in pronation cycle.

This sequence can be repeated as often as needed. You can stop treatment with the **STOP** key; after the programmed therapy duration, it will stop automatically.

**Default: disabled**


**Note!**

- If a pause has been programmed, the carriage will stop for the pause each time the maximum stretching value is attained.
- Enabling the stretching in pronation function automatically disables motor A. The only movement performed is pronation/supination.
- It is not possible to combine a second stretching function with the special function "stretching in pronation".
- During operation, the display will show the symbol  in the upper, left-hand corner.

### ■ Stretching in supination

With the special "stretching in supination" function the joint will be gently stretched in the supination direction. Only the pronation/supination movement is treated here. Motor A is automatically disabled. The position of motor A cannot be changed once this special function has been activated.

Starting at the middle position the carriage will first move to the programmed pronation position and then to the programmed supination position.

Subsequently the carriage reverses 5° towards pronation and then moves very slowly back to the programmed supination value (symbol: ). Subsequently the carriage attempts to attain another 5° at an even slower speed (symbol: .

If a high resistance toward the additional 5° is sensed, the reverse on load function is automatically activated and the carriage moves in the opposite direction.


This stretch cycle is repeated ten times.

After that the carriage moves to the programmed maximum pronation position and restarts the stretching in supination cycle.

This sequence can be repeated as often as needed. You can stop treatment with the **STOP** key; after the programmed therapy duration, it will stop automatically.

**Default: disabled**

## Note!

- If a pause has been programmed, the carriage will stop for the pause each time the maximum stretching value is attained.
- Enabling the stretching in supination function automatically disables motor A. The only movement performed is pronation/supination.
- It is not possible to combine a second stretching function with the special function "stretching in supination".
- During operation, the display will show the symbol  in the upper, left-hand corner.

## LEVEL 6:

### ■ EROM repeat extension/pronation

The special "EROM repeat extension/pronation" function allows a more efficient exercise in the last 10° before the set maximum extension and pronation values.

For this protocol, the CPM device starts in the middle between the set maximum extension/flexion and pronation/supination values. First of all the carriage moves simultaneously to the maximum flexion angle and to the maximum supination angle. Then it moves simultaneously to the maximum extension angle and to the maximum pronation angle.

When the set extension/pronation value has been reached, the carriage reverses 10° toward flexion/supination and then moves back again to the maximum extension/pronation angle. The movement through the final 10° is repeated five times at a slow speed.

At the end of the cycle, the carriage will again move to the maximum flexion value, at the same time as the maximum supination value, and then starts another cycle with five repetitions through the last 10° of extension/pronation.

This sequence can be repeated as often as needed. You can stop treatment with the **STOP** key; after the programmed therapy duration, it will stop automatically.

**Default: disabled**

### ■ EROM repeat flexion/supination

The special "EROM repeat flexion/supination" function allows a more efficient exercise in the last 10° before the set maximum flexion and supination values.

For this protocol, the CPM device starts in the middle between the set maximum extension/flexion and pronation/supination values. First of all the carriage moves simultaneously to the maximum extension angle and to the maximum pronation angle. Then it moves simultaneously to the maximum flexion angle and to the maximum supination angle.

When the flexion/supination value has been reached, the carriage reverses 10° toward extension/pronation and then moves back again to the maximum flexion/supination angle. The movement through the final 10° is repeated five times at a slow speed.


At the end of the cycle, the carriage will again move to the maximum extension value, at the same time as the maximum pronation value, and then starts another cycle with five repetitions through the last 10° of flexion/supination.

This sequence can be repeated as often as needed. You can stop treatment with the **STOP** key; after the programmed therapy duration, it will stop automatically.

**Default: disabled**

## LEVEL 7:

### ■ Transport setting

With this function, the carriage will move to a position optimally suited for packing the CPM device. Select the function and press the **START** key. The carriage moves to the transport position. In the top left-hand corner of the display you will see  (also refer to chapter 6 Transport).

### ■ New patient

With this function, the CPM device will move to the home position.

- Select the function and press the **START** key:
  - The carriage moves to the home position.

- The programmed treatment parameters are deleted.
- All values stored on the chip card are deleted.
- The carriage stops in the middle position of the extension/flexion and pronation/supination angles.

The "new patient" function (home position) selects the following settings:

- extension: 25°
- flexion: 35°
- pronation: -5°
- supination: 5°
- pauses: 0
- timer: continuous operation
- speed: 100 %
- reverse on load motor A: 25
- reverse on load motor B: 25
- motor A: ON
- motor B: ON
- synchronized mode: ON
- total therapy time: 0
- special functions: disabled

#### ■ Operating mode A+B synchronized/non-synchronized

Motors A and B can run in synchronized or non-synchronized mode.

##### synchronized:

Motors A and B perform a synchronized movement according to the elbow joint's physiological motion pattern:

Starting at the middle position of the set extension/flexion and pronation/supination angles, the carriage first moves simultaneously to the maximum extension value and the maximum pronation value. Afterwards it moves simultaneously to the maximum flexion and maximum supination values. After reaching this position, the cycle restarts: maximum extension and maximum pronation.

During operation the synchronized mode is indicated in the top left-hand corner by the symbol .

##### non-synchronized:

Both motors run independently of each other within the selected ranges of motion.

To enable the non-synchronized mode, disable the synchronized mode.

During operation the non-synchronized mode is indicated in the top left-hand corner by the symbol .

**Default: synchronized mode enabled**

### Caution!

Patient hazard –

The synchronized operation is recommended as a general rule. However, medical and/or therapeutic indications may call for the non-synchronized mode. In the non-synchronized mode, particular caution and attention on the part of the physician/therapist is required to rule out any risk for the patient.

#### ■ Service MENU

For service purposes only, refer to Service Manual.

#### Reminder:

You save the selected parameter values by pressing the **STOP** key.

## 5.4 Application and programming examples

### 5.4.1 Isolated extension/flexion

1. Complete the mechanical settings and adjust the carriage to the patient as described in section 4.2.
2. Press the **MENU key** on the programming unit to access level 1 (M1).
3. Press the **pronation** () or **supination parameter key** () , the press the **+ / - keys** to set the carriage to the appropriate rotation position to be maintained during the isolated extension/flexion movements.
4. Next set the desired range of motion for extension/flexion by first pressing the **extension parameter key** () , then set the values with the **+ / - keys**.
5. Set the flexion value () in the same way.

## Note!

- For an exclusive extension/flexion protocol, motor B for the rotation movement must be turned off as described in steps 6 through 8 below.
- Having programmed the ranges of motion, you can program additional options, such as pause or speed.

6. Press the **MENU key** on the programming unit to access level 3 (M3).
7. Press the **Motor B ON/OFF** parameter key (⊙ B) to activate the parameter.
8. Deactivate motor B by pressing the parameter key **Motor B ON/OFF** again or the - (minus) key. The check mark in the circle next to the symbol must have disappeared.
9. Save the settings with the **STOP key**, set the carriage to the starting position with the **START key** and initiate therapy by pressing the **START key** again.

### 5.4.2 Isolated pronation/supination

1. Complete the mechanical settings and adjust the carriage to the patient as described in section 4.2.
2. Press the **MENU key** on the programming unit to access level 1 (M1).
3. Activate the **parameter key extension** (↔) or **flexion** (↕), then press the **+ / - keys** to set the carriage to the appropriate position to be maintained during the isolated pronation/supination movements.
4. Next set the desired range of motion for pronation/supination by first pressing the **pronation parameter key** (↶), then set the values with the **+ / - keys**.
5. Set the supination value (↷) in the same way.

## Note!

- For an exclusive rotation protocol, motor A for extension/flexion must be turned off as described in steps 6 through 8 below.
- Having programmed the ranges of motion, you can program additional options, such as pause or speed.

6. Press the **MENU key** on the programming unit to access level 3 (M3).
7. Press the **Motor A ON/OFF** parameter key (⊙ A) to activate the parameter.
8. Deactivate motor A by pressing the parameter key **Motor A ON/OFF** again or the - (minus) key. The check mark in the circle next to the symbol must have disappeared.
9. Save the settings with the **STOP key**, set the carriage to the starting position with the **START key** and initiate therapy by pressing the **START key** again.

## 6. Care, Maintenance, Transport, Conversion

The ARTROMOT®-E2/-E2 compact is suitable for re-use if the following points are observed.

### 6.1 Care/Re-use

#### **Warning!**

Shock hazard – Unplug the device from the power line before cleaning.

Shock hazard, equipment damage – Liquids must not enter the device or the programming unit.

- The ARTROMOT®-E2/-E2 compact must be **disinfected by wiping down** with a disinfectant. Thus, it complies with the special hygiene standards for medical technical equipment.
- **The enclosure** and removable **armrest assemblies** can be cleaned with **commonly used disinfectants** and **mild household detergents**.
- Only use a **damp cloth** to wipe the CPM device down.

#### **Warning!**

Patient hazard – patient contamination

- Before using the device on another patient, be sure to clean and disinfect it according to the instructions given here.

#### **Note!**

The manufacturer recommends using only special disinfectants approved for medical devices and with the characteristics specified below under "Caution".

Disinfection with the following disinfectants and in compliance with the instructions for use has been evaluated and approved by the manufacturer in the risk assessment:

- DESCOSEPT AF Lemon  
(Art.-Nr: 00-311L-xxx)  
BAuA Reg.No.: N-55153, CE-0482  
Manufacturer: Dr. Schumacher GmbH  
Contact time: 2 minutes minimum

#### **Caution!**

Equipment damage –

- The plastic material used is not resistant to mineral acids, formic acid, phenols, cresols, oxidants and strong organic or inorganic acids with a pH value below 4.
- Use only clear disinfectants to prevent discoloration of the device.
- Do not expose the CPM device to strong ultraviolet radiation (sunlight) and fire.

## 6.2 Maintenance (fuse replacement)

### Check before each use

Visually inspect the device for signs of mechanical damage before each use.

If you detect damage or malfunctions that may impair the safety of the patient or of the operator, have the device repaired before using it.

### Technical inspections

Due to their design, the bearings and joints of the ARTROMOT®-E2/-E2 compact do not require maintenance and all materials are protected from corrosion, BUT: For safety, the devices require regular maintenance. To maintain the functional and operational safety, check all components for damage and loose connections at least once a year.

These checks should be performed by persons with adequate training and experience. Damaged and worn parts must immediately be replaced with original spare parts by authorized staff.

## Warning!

Patient hazard, equipment malfunction and damage

- Refer repair and maintenance to authorized persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.
- ORMED GmbH will make all documents required for servicing, such as circuit diagrams, parts lists, descriptions or calibration instructions, available to authorized experts.

The inspections can be carried out by DJO Technical Service within the framework of a service agreement. Contact DJO for details.

**Other than that, the manufacturer does not require any other regular maintenance.**

## Note!

Regarding further technical or other inspections or inspection intervals, observe applicable local requirements, such as IEC 62353, DGUV3 or similar regulations to be observed by users of medical equipment or electric devices.

## Fuse replacement

### **⚠ Warning!**

Patient hazard, equipment malfunction and damage

The replacement of fuses must be referred to specialists as defined in IEC 60364 or other applicable standards (e.g. bio-medical technicians, electricians, electronics installers).

#### **Fuses used must be**

**S# < 20.000: T1A L250Vac**

**S# > 20.000: T2A H250Vac fuses.**

- Before replacing fuses, turn off the ARTROMOT®-E2/-E2 compact and disconnect the device from the power line.
- Use an appropriate tool to remove the fuse holder (8) situated between the power switch (7) and the power connector (6) (Fig. 1).
- Replace the fuses and reinsert the fuse holder (Fig. 2). Check that the fuser holder locks properly into place.


Fig. 1

Fig. 2

## 6.3 Transport


### 6.3.1 ARTROMOT®-E2

The following operating steps must be completed before transporting the ARTROMOT®-E2:


1. Activate the "transport setting" function →  in the menu (see also 5.3) and start the ARTROMOT®-E2.
2. Turn off the ARTROMOT®-E2 power switch (19).
3. Remove the power cord (18) and disconnect the motion element (21) and the programming unit (32).
4. Open locking screw (4) and pull out the armrest (14) for the healthy arm.
5. Grasping the motion element at the bracket (30) for adjustment of the positioning angle, open the locking screw (4) and pull out the motion element.
6. Set horizontal adduction/abduction to 0° (indexing knob 1).
7. Open the locking screw for adjustment of the backrest (10), fold the backrest to the front all the way until flat on the seat and tighten the screw.
8. Remove the two locking pins (22). Pull out the legs towards the sides, reverse and reinsert. Then put the locking pins back in place.

9. Only use the original shipping box for transporting the device. Ormed GmbH cannot be held liable for transport damage, if the device was not shipped in its original shipping box.
10. Place the ARTROMOT®-E2 – legs first – on the bottom of the box. The markings on the bottom of the box indicate the correct position.

11. Then put the enclosed polystyrene block on the chair, observing the cut-outs in the block.
12. Put the programming unit (15) in the supplied box. Place the motion element, the armrest and the power cord in the respective cut-outs in the polystyrene block of the package.

13. Then close the box.

**Reverse the above steps to reassemble the device after transport.**

### **⚠ Warning!**


Shock hazard –

**Allow the ARTROMOT®-E2 to reach room temperature** before use. If the device has been transported at **temperatures below 0 °C (32°F)**, leave it to dry at room temperature for about 2 hours, until any condensation has disappeared.

## **6.3.2 ARTROMOT®-E2 compact**


The following operating steps must be completed before transporting the ARTROMOT®-E2 compact:

1. Activate the “transport setting” function →  in the menu (see also 5.3) and start the ARTROMOT®-E2 compact.
2. Turn off the ARTROMOT®-E2 compact power switch (19).
3. Remove the power cord (18) and disconnect the motion element (9) and the programming unit (16).
4. Remove pin (26) on the underside and take off the upper arm support.
5. Remove pin (5) and take off the swiveling mechanism of the upper arm support. Set the guide of the upper arm support to the transport setting (parallel with motor A).
6. Grasping the motion element at the bracket (27) for adjustment of the positioning angle, open the locking screw (1) and pull out the motion element.
7. Open the locking screw for adjustment of the elbow angle adjustment (2) and separate the oval tube from the motion element.
8. Open locking screw (24) on the underside of the stand, then remove the column with the electronics box from the star-shaped stand.
9. Disassemble the stand by loosening the five legs with the special tool supplied and remove them.
10. Unscrew the locking screw (7) from motor B (12), remove the entire forearm support including the slider (4) and the hand support (13 + 14). Afterwards screw the locking screw (7) back into motor B.
11. Only use the original shipping box for transporting the device. Ormed GmbH & Co. KG cannot be held liable for transport damage, if the device was not shipped in its original shipping box.
12. Put the programming unit (15) in the supplied box, using sufficient padding.


13. Pack the following items in the longer box provided, again using sufficient padding:

- Upper arm support
- Oval tube of the motion element
- Column with electronics box
- Entire forearm support with slider and hand support


14. Pack the following items in the shorter box provided, again using sufficient padding:

- Star-shaped stand
- 5 legs with castors
- 5 screws to attach the legs
- 5 washers to attach the legs
- 1 Allen key

15. Please put the motion element, the power cord and the three smaller boxes in the respective cut-outs of the polystyrene block.

16. Then close the box.

### **Warning!**

Shock hazard –

**Allow the ARTROMOT®-E2 compact to reach room temperature** before use. If the device has been transported at **temperatures below 0 °C (32°F)**, leave it to dry at room temperature for about 2 hours, until any condensation has disappeared.

## 6.4 Conversion

### Note!

On the display the respective sides are indicated as follows:

-  The CPM device (upper arm support) is prepared for the left elbow.
-  The CPM device (upper arm support) is prepared for the right elbow.
-  This symbol appears during conversion of the CPM device.

### 6.4.1 Reconfiguration of the ARTROMOT®-E2

The ARTROMOT®-E2 can be used on the left and on the right elbow joint. However, the device must be converted first. This is done very quickly.

1. Activate the "transport setting" function  in the menu (see also 5.3) and start the ARTROMOT®-E2.
2. Adjust an angle of 90° for horizontal adduction/abduction (indexing knob 1) (see conversion drawing ①).
3. Open locking screw (4) for the height adjustment, remove the armrest for the healthy arm and place it on the seat (see conversion diagram ①).


4. Hold the double joint (3) of the motion element and open the locking screw for height adjustment (4) on this side as well.
5. Remove the motion element and insert it on the opposite side. Tighten the locking screw (4). (see conversion diagram ① b, c).

6. Pull the pin for the swiveling movement of the upper arm support (11) downward and swing the support 180° towards the inside.

Attention: Release the pin (11) while the support is swinging; you will hear it lock into place on the other side (see conversion diagram ② d)

7. Open the locking screw for the swiveling movement of the forearm support (12) and swing the support with motor B 180°. Tighten locking screw (12) (see conversion diagram ② e).
8. Open the locking screw for right/left positioning (13) two revolutions. Move the forearm support to the lower recess in the slot and turn it 180° in this position. Then move the forearm support back to the middle position (setting 0) and tighten locking screw (13) again (see conversion diagram ③ f).

Conversion diagram ARTROMOT®-E2:

## 6.4.2 Reconfiguration of the ARTROMOT®-E2 compact


The ARTROMOT®-E2 compact can be used on the left and on the right elbow joint. However, the device must be converted first. This is done very quickly.

1. Activate the "transport setting" function →  in the menu (see also 5.3) and start the ARTROMOT®-E2 compact.
2. Release the wheel brakes (22) and roll the ARTROMOT®-E2 compact to the side of the elbow to be treated.
3. Pull the pin for the swiveling movement of the upper arm support (5) downward and swing the support 180° towards the inside.

Attention: Release the pin (5) while the support is swinging; you will hear it lock into place on the other side (see conversion diagram ① a).

4. Open the locking screw for the swiveling movement of the forearm support (6) and swing the support with motor B 180°. Tighten locking screw (6) (see conversion diagram ① b).
5. Open the locking screw for right/left positioning (7) two revolutions. Move the forearm support to the lower recess in the slot and turn it 180° in this position. Then move the forearm support back to the middle position (setting o) and tighten locking screw (7) again (see conversion diagram ② c).

Conversion diagram ARTROMOT®-E2 compact:

## 7. Environmental Protection Statement

The product described in this operation manual must not be disposed of with unsorted household or municipal waste. It requires separate disposal. Please contact DJO for information about the possible recycling of the product.

The expected service life of the device and all supplied parts and accessories is 6 years. The application beyond this period is the responsibility of the user.

## 8. Specifications

<b>Model:</b>	ARTROMOT®-E2
<b>Device name:</b>	ARTROMOT®-E2, Part no.: 80.00.031 ARTROMOT®-E2 compact, Part no.: 80.00.033
<b>Input ratings:</b>	100 – 240 Vac / 50 – 60 Hz tolerance -15% to +10%
<b>Current consumption:</b>	100V      240V
standby (ON):	5VA      5VA
operation (maximum):	33VA      33VA
	330mA    140mA
<b>Fuses:</b>	serial numbers < 20.000: 2x T1A L250 Vac serial numbers > 20.000: 2x T2A H250 Vac according to IEC 60127-2/3, UL 248-14, CSA C22.2 no. 248.14
<b>Protection class:</b>	serial numbers < 20,000: I serial numbers > 20,000: II
<b>IP class:</b>	serial numbers < 20,000: IPX0 serial numbers > 20,000: IP21
<b>Applied part:</b>	Type B
<b>Max. continuous load</b>	
seat:	150 kg
motion element:	9 kg
<b>Dimensions (transport):</b>	
ARTROMOT®-E2	
length:	87.5 cm
width:	57.5 cm
height:	58.0 cm
ARTROMOT®-E2 compact	
length:	87.5 cm
width:	57.5 cm
height:	29.0 cm

<b>Adjustment ranges (min./max.):</b>	
height adjustment: (ARTROMOT®-E2 measured from seat)	35 – 71 cm
upper arm length:	29 – 46 cm
<b>Seat height (ARTROMOT®-E2):</b>	48 cm
<b>Accuracy of measured values:</b>	
goniometer in the measuring range:	from -90° to +90° in pronation/supination
accuracy:	+/- 2°
<b>weight:</b>	
ARTROMOT®-E2	25 kg
ARTROMOT®-E2 compact	17 kg
<b>materials used:</b>	ABS, POM (Delrin 100), PUR, PA, FR4, aluminum, stainless steel, brass
<b>MDD:</b>	class IIa
<b>Standards compliance:</b>	93/42/EEC IEC 60601-1 IEC 60601-1-6 IEC 60601-1-9 IEC 60601-1-11 (S# >20,000) IEC 62366 IEC 62304 EN ISO 14971 ANSI/UL 60601-1 (S# <20,000) CAN/CSA C22.2 No. 601.1-M90 (S# <20,000) ANSI AAMI ES 60601-1 1st Edition (S# >20,000) CAN CSA 22.2 No. 60601-1-08 (S# >20,000)
<b>EMC (electromagnetic compatibility)</b>	IEC 60601-1-2
<b>Manufactured in compliance with:</b>	EN ISO 13485
<b>Ambient conditions (storage, transport)</b>	
temperature:	-25 °C to +70 °C
relative humidity:	at 70 °C up to 93 % no condensation
atmospheric pressure:	500 hPa to 1060 hPa
<b>Ambient conditions (operation)</b>	
temperature:	+5 °C to +40 °C
relative humidity:	15 % to 93 %
atmospheric pressure:	700 hPa to 1060 hPa

## 9. Contact

We would be happy to answer any questions you may have about our products and services.

### **ARTROMOT® International**

Please contact your local dealer, the DJO Headquarters in the USA, the DJO International Headquarters in England or DJO in Germany.

### **DJO International**

Headquarters: DJO UK Ltd.  
1a Guildford Business Park  
Guildford  
Surrey, GU2 8XG, United Kingdom  
Telephone: +44 (0) 1483 459659  
Fax: +44 (0) 1483 459470  
Email: [info@DJOglobal.eu](mailto:info@DJOglobal.eu)  
Web: [www.DJOglobal.eu](http://www.DJOglobal.eu)

### **DJO Headquarters**

DJO, LLC  
2900 Lake Vista Drive  
Lewisville, TX 75067, USA  
Telephone: +1 800 321 9549  
Fax: +1 800 936 6569  
Email: [webmaster@DJOglobal.com](mailto:webmaster@DJOglobal.com)  
Web: [www.DJOglobal.com](http://www.DJOglobal.com)

### **Manufacturer / Headquarters Germany**

ORMED GmbH  
A DJO Company

Serial numbers < 30,000:  
Merzhauser Straße 112  
79100 Freiburg, Germany

Serial numbers > 30,000:  
Bötzingen Straße 90  
79111 Freiburg, Germany

Telephone: +49 (0) 761 4566 01  
Fax: +49 (0) 761 456655-01  
Email: [medizintechnik@DJOglobal.com](mailto:medizintechnik@DJOglobal.com)  
Web: [www.DJOglobal.de](http://www.DJOglobal.de)

### **Quarters Austria**

DJO  
ORTHOMED Medizintechnik GmbH  
Annagasse 5/1/Top 4  
1010 Vienna, Austria  
Telephone: +43 1 5320834  
Fax: +43 1 532083431  
Email: [orthomed@orthomed.cc](mailto:orthomed@orthomed.cc)  
Website: [www.orthomed.cc](http://www.orthomed.cc)

### **Warranty:**

2 years (mechanical parts)  
2 years (electronics)

## 10. Technical Service

### 10.1 Technical hotline

Do you have any technical questions?  
Do you need technical service?

Telephone: +49-180-5-1 ormed de  
                  +49-180-5-1-67 63 33  
Fax:           +49-180-5-3 ormed de  
                  +49-180-5-3-67 63 33

### 10.2 Shipment

To prevent damage during transport, only use the original shipping box. These boxes can be obtained from DJO.

Before packing the CPM device, set it to the transport position (see chapter 5 and 6).

### 10.3 Spare parts

Refer to the Service Manual for the most recent list of spare parts.

When ordering spare parts, always specify:

- item
- description
- part number
- quantity
- serial number of the CPM device

#### Note!

Refer repairs to authorized, specially trained staff.

DJO offers service training for your personnel.

#### Note!

Country-specific power cords are available from DJO. Please contact DJO or your DJO distributor.

Surcharges may apply in certain cases to spare parts ordered in low quantities.

Item	Description	Part No.	Qty
1.	patient chip card	0.0034.048	<input type="checkbox"/>
2.	Patient chip card (protocol)	0.0037.035	<input type="checkbox"/>
3.	marker pen for patient chip card	0.0031.006	<input type="checkbox"/>
4.	Power cord EU version HO5VV-F3G 3 × 1 mm <sup>2</sup> length: 4 m	0.0034.118	<input type="checkbox"/>

## 11. CE-Declaration of Conformity

DECLARATION OF CONFORMITY

In compliance with the provisions of the Council Directive 93/42/EEC, Annex II, of 14 June 1993 / Amendment of 5 September 2007, concerning medical devices, the company

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**declares under its sole responsibility that the products of the product line**

**ARTROMOT® (see Annex)**

**fulfill the essential requirements of Annex I  
of the Council Directive 93/42/EEC.**

**With reference to Rule 9 of the Directive 93/42/EEC, Annex 9,  
the product is a device of risk class IIa.**

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**Freiburg, 23 October 2020**

-----  
**- QA Management Representative -  
-Bernhard Krohne-**

This certificate is valid until expiry of the certificate referred to.  
(The certificate for the year of manufacture can be downloaded from:  
<https://www.djoglobal.de/unternehmen/qualitaet-und-verantwortung>)

**Annex:**  
ARTROMOT®-S3  
ARTROMOT®-S4  
ARTROMOT® ACTIVE-K  
ARTROMOT®-K1  
ARTROMOT®-SP3  
ARTROMOT®-E2

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**Bijlage:**  
ARTROMOT®-S3  
ARTROMOT®-S4  
ARTROMOT® ACTIVE-K  
ARTROMOT®-K1  
ARTROMOT®-SP3  
ARTROMOT®-E2

## 12. IEC 60601-1-2 and IEC 60601-1-11

The **ARTROMOT®-E2/-E2 compact** device is subject to particular precautions regarding electromagnetic compatibility (EMC). The device must be installed and put into service strictly in compliance with the EMC directives set forth in the accompanying documents.

Portable and mobile RF communication systems may affect the **ARTROMOT®-E2/-E2 compact**.

The **ARTROMOT®-E2/-E2 compact** should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, **ARTROMOT®-E2/-E2 compact** should be observed to verify normal operation in the configuration in which it will be used.

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If it is necessary to replace assemblies or cables, only the manufacturer's original parts must be used to ensure continued compliance with EMC requirements after repair. This requirement applies to the power supply unit, cables and cable lengths, drive unit consisting of the motor and the control system, the programming unit incl. the coiled cable and the connector.

Regarding the electromagnetic compatibility, there is no maintenance work required during the expected service life of the product.

The product designation **ARTROMOT®-E2/-E2 compact** used in the text below includes all product variants.

### **WARNING!**

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

## 12.1 Electromagnetic emissions

### Guidance and Manufacturer's Declaration – Electromagnetic Emissions

**ARTROMOT®-E2/-E2 compact** is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-E2/-E2 compact** unit is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions to CISPR 11	Group 1	<b>ARTROMOT®-E2/-E2 compact</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions to CISPR 11	Class B	<b>ARTROMOT®-E2/-E2 compact</b> is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions to IEC 61000-3-2	not applicable	
Voltage fluctuations/flicker emissions to IEC 61000-3-3	not applicable	

## 12.2 Electromagnetic immunity

### Guidance and Manufacturer's Declaration – Electromagnetic Immunity

**ARTROMOT®-E2/-E2 compact** is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-E2/-E2 compact** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment – guidance
Electrostatic discharge (ESD) to IEC 61 000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/bursts to IEC 61 000-4-4	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power should be that of a typical commercial or hospital environment.
Surges to IEC 61 000-4-5	± 1 kV line-to-line voltage	± 1 kV line-to-line voltage	Mains power should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines to IEC 61000-4-11	0% $U_T$ ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% $U_T$ ; 1 cycle and 70% $U_T$ ; 25/30 cycles Single phase at 0°	0% $U_T$ ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% $U_T$ ; 1 cycle and 70% $U_T$ ; 25/30 cycles Single phase at 0°	Mains power should be that of a typical commercial or hospital environment. If the user of the <b>ARTROMOT®-E2/-E2 compact</b> unit requires continued operation during power mains interruptions, it is recommended that the <b>ARTROMOT®-E2/-E2 compact</b> unit be powered from an uninterruptible power supply or a battery.
Voltage interruptions to IEC 61000-4-11	0% $U_T$ ; 250/300 cycles	0% $U_T$ ; 250/300 cycles	
Power frequency (50/60 Hz) magnetic field to IEC 61 000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

**Note:**  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and Manufacturer's Declaration – Electromagnetic Immunity

**ARTROMOT®-E2/-E2 compact** is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-E2/-E2 compact** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment – guidance
Conducted RF to IEC 61000-4-6	6 V 0,15 MHz-80 MHz 80% AM at 1 kHz	6 V 0,15 MHz-80 MHz 80% AM at 1 kHz	
Radiated RF to IEC 61000-4-3	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communications equipment to IEC 61000-4-3	27 V/m 380-390 MHz PM 18 Hz  28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine  9 V/m 704-787 MHz PM 217 Hz  28 V/m 800-960 MHz PM 18 Hz  28 V/m 1,7-1,99 GHz PM 217 Hz  28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz  9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	27 V/m 380-390 MHz PM 18 Hz  28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine  9 V/m 704-787 MHz PM 217 Hz  28 V/m 800-960 MHz PM 18 Hz  28 V/m 1,7-1,99 GHz PM 217 Hz  28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz  9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ARTROMOT®-E2/-E2 compact** exceeds the applicable **ARTROMOT®-E2/-E2 compact** unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **ARTROMOT®-E2/-E2 compact** unit.

## 12.3 RF (Radio Frequency) Interference

EMC disturbance can affect the performance of the **ARTROMOT®-E2/-E2 compact**. To prevent the user from unacceptable risks the device runs performance checks regularly during operation. This is what happens, if a problem is identified:

- an audio signal sounds
- the unit stops immediately
- the message **ERROR** and an error code appear on the display

In this case try to quit the error with the **STOP** key. If the error persists you may attempt to restart the unit by turning it briefly off and on with the power switches.

### **WARNING!**

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the **ARTROMOT®-E2/-E2 compact**, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

### **Notice!**

(for home use)

Portable and mobile RF communication systems may affect the **ARTROMOT®-E2/-E2 compact** device.

For this reason, make sure that wireless communication equipment, such as wireless home network devices, cell phones, cordless telephones and their base stations and walkie-talkies operate at a minimum distance of 3,3 m from the device. (Calculated on the basis of the maximum power output of a typical cell phone of 2 W).

## Piktogrammübersicht ARTROMOT®-E2/-E2 compact

### Symbol Overview ARTROMOT®-E2/-E2 compact

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