

ENGLISH

BREAK-IN PERIOD FOR SHOES WITH HEAT MOLDED O-B2R CUSTOM INSERTS:

To ensure that your shoes become a functioning extension of your pedorthic medical care, please follow these instructions:

1. Place your new shoes (with the inserts in them) on your feet (with socks) and wear them for 30 to 60 minutes – only on carpeting at first. 2. Look for any areas of redness on your foot. 3. Once you have verified that the shoes do not rub your skin (absence of redness), wear your shoes around your home for a day or two; check again for areas of redness. 4. Once you have verified that there are no issues with rubbing, you are ready to wear them outside the home. 5. Continue to check your shoes and feet each day looking for anything out of the ordinary. Before wearing the shoe, please read the following instructions completely to ensure keeping your shoes in usable condition.

OVERALL CARE OF SHOES (EXCEPT WOOL): 1. Never put in a washing machine. 2. Clean shoes regularly. 3. Use a fabric protector or leather protector before wearing to make the shoe last.

LEATHER Care: 1. If your shoes are dirty, first clean them with a damp cloth to get them ready for an application of shoe cream. 2. Use a cleaning cream for leather shoes. This cream will keep the leather clean and supple. Apply the cream with a clean dry cloth and work it into the leather. Buff or brush out to provide a fresh finish.

NUBUCK Care: 1. Use a small suede brush to work away the dirt. 2. There are sprays designed to refurbish suede or nubuck materials. Try the spray in a small area first to test for any discoloration. 3. Do not immerse the nubuck shoes in water. Shoe polish should not be used.

LYCRA® Care: 1. Use a fabric protector on this material to help with the dirt. 2. If the shoes do get soiled, use a small amount of mild soap and water to remove the dirt.

CANVAS Care: 1. Make sure shoes are dry before you begin. 2. Remove laces and soak in solution of mild detergent and water for 15 minutes. 3. Do not soak the shoes in water, instead rinse the shoes under running water. 4. Using a cleaning solution made up of 3-4 cups of water and 1 cup of liquid soap, use a scrub brush with light pressure. 5. Air dry until completely dry.

MESH® Care: 1. Using a soft brush, brush away the dirt. 2. Combine warm water and small amount of laundry detergent mixture in a small bowl. 3. Using a soft cloth apply the mixture in circular motions. 4. Using a second cloth, wet it then wipe the exterior of the shoes to remove the mixture you put on. 5. Place them in a cool, dry spot out of sunlight to allow for shoes to air dry.

WOOL Care: 1. Brush of excess dirt with towel, remove laces and insoles. 2. Use a neutral detergent, avoid bleach. 3. Put shoes in linen bag and wash using cold water and a wool or delicate cycle. 4. Washing insoles and laces, use a mild wool wash detergent and cold water you can scrub them by hand. If you want to wash your shoes via hand, use this solution as well. 5. Do not tumble dry the shoes, air dry them for best result. 6. Once dry reinsert laces and insoles.

INTENDED USER: The intended user should be the patient, a licensed medical professional or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use. Dr. Comfort Diabetic Shoes: are intended to provide support and protection and help in the prevention of diabetic foot ulcers. Dr. Comfort Flex-OA Shoes are intended to provide support and stabilisation and help in the management of knee osteoarthritis.

CONTRAINDICATIONS: None

WARNINGS AND PRECAUTIONS: Correct sizing is very important. The shoe should not be too tight or too loose. If you experience any pain, swelling, callus, ulceration, sensation changes or any unusual reactions while using this product, or if the fit feels too tight stop using this product and consult your medical professional immediately. Check your feet each day for redness. Check your footwear every three months to ensure it still fits adequately, protects, and supports the foot. Always wear insoles in your shoes. Always wear socks with your shoes. Do not wear over open wounds. Do not use this device if it is damaged and/or the packaging has been opened.

Note: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

INTENDED FOR SINGLE PATIENT USE.

CAUTION: THIS PRODUCT IS MADE WITH NATURAL RUBBER LATEX WHICH MAY CAUSE ALLERGIC REACTION.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.



GENUINE
Dr. Comfort®
EST. USA

INSTRUCTIONS FOR USE GEBRAUCHSANWEISUNG ISTRUZIONI PER L'USO BRUKSANVISNINGAR INSTRUCTIONS D'UTILISATION



CE



MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany

©2021 DJO, LLC



DJO, LLC
1430 Decision Street
Vista CA 92081-8553
U.S.A.

IFU-010-EU REV_1 2021-MAY-18

