Figure 1


Figure 2


Figure 3

Figure 4


Figure 5

Figure 6

Figure 7

ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE: This brace is designed to alleviate pain by supporting the body in guiding and stabilizing the knee joint.

INDICATIONS: For patients with symptoms, physical findings or radiographic evidence of moderate to severe unicompartmental knee osteoarthritis, post-operative meniscal repair only where unicompartmental unloading is indicated, or varus-valgus conditions that may benefit from unicompartmental unloading.

CONTRAINDICATIONS: Do not use if you are allergic to any of the materials contained within this product.

WARNING AND PRECAUTIONS: The brace was designed to complement the variety of medical treatments common to the above afflictions and should be used under the supervision of a medical professional who will determine the frequency, duration and directions of use. This device is to be prescribed and fitted by a qualified healthcare professional. ROM settings should not be changed without supervision of a medical professional. If you experience any pain, discomfort, swelling, sensation changes, circulation problems, or any unusual reactions while using the brace, consult a medical professional immediately. If you develop an allergic reaction and/or experience itching, red skin after coming into contact with any part of this device, please stop using the brace and contact your healthcare professional immediately. The brace should be worn directly against dry skin. The brace is intended for use on a single patient and is not intended for use in contact sports. Braces should not be worn on both right and left legs during activities that may cause the braces to contact each other resulting in potential loss of control and injury.

The brace should be snug but not impair circulation. Consult a medical professional for your activity limitations. Do not use over open wounds. Do not use the device if it is damaged and/or the packaging has been opened.

Note: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

PRE-APPLICATION INFORMATION: (ONLY DONE AT INITIAL FITTING)

- To begin application of the brace open all buckles and disconnect magnetic latches on the Thigh and Calf **(Figure 1)**
- While sitting on the edge of a chair, flex the leg at 45 degrees, position hinge at top of the knee cap and slightly in front of midline of the leg on the affected compartment **(Figure 1)**
- Connect the magnetic latches on Thigh and Calf and ensure buckles are closed in position **(Figure 2)**
- Tighten Thigh and Calf circumferential straps to snug fit, trim excess material **(Figure 3)**
- Straighten leg, position central harness on unaffected side of knee joint with leading edge centered at midline of leg **(Figure 4)**
- While holding the position of the central harness, separate posterior Thigh load strap from allogator tab of central harness, pull tight, re-attach to allogator tab, trim excess material. Repeat for Calf **(Figure 4)**
- Ensure BOA™ dials are in fully extended position by rotating dial counter-clockwise and pulling straps towards back of knee **(Figure 5)**
- While holding the position of the central harness, separate anterior Thigh load strap from allogator tab of central harness, pull tight, re-attach to allogator tab, trim excess material. Repeat for Calf **(Figure 6)**
- For high tensioning, move leg back into flexed 45 degree position, rotate BOA™ dials clockwise to tighten as necessary, note location on printed graduated scale. Repeat for Calf **(Figure 7)**
- To remove brace, open both buckles and disconnect the magnetic latches by sliding them forward

APPLICATION INFORMATION:

- Disconnect magnetic latches on Thigh and Calf
- While sitting on the edge of a chair, flex the leg at 45 degrees, position hinge at top of the knee cap and slightly in front of midline of the leg on the affected compartment **(Figure 1)**
- With buckles in the open position, connect the magnetic latches on Thigh and Calf
- Bring Thigh and Calf buckles to closed position
- Adjust Thigh and Calf circumferential straps as necessary
- To remove brace, open both buckles and disconnect the magnetic latches by sliding them forward

ADJUSTING TELESCOPING CUFFS:

- Remove either Thigh or Calf Cuff lining pad by pulling pad away from cuff
- Using a standard flat head screwdriver, rotate screw head counterclockwise to remove. Hold on to screw
- Repeat for second screw
- Lift bar up away from Cuff and re-position to desired location
- Insert screw through cutout on upright bar and tighten in to threaded insert until 2 clicks are reached indicating full engagement
- Repeat for second screw
- Secure lining pad by placing over inside surface of cuff, ensure full coverage, and lightly press into place over velcro engagement

HINGE STOP INSTALLATION INSTRUCTIONS: Detach condyle pad. For extension stops, loosen, and remove 4 screw from existing stop. Insert desired extension stop into front of hinge, aligning single notch, re-insert and tighten screw. Similarly for flexion stops, loosen, and remove 4 screw, insert desired flexion stop into back of hinge aligning both notches, re-insert and tight screw.

Extension Stops: 0°, 10°, 20°, 30°, 45°
Flexion Stops: 20°, 30°, 60°, 75°, 90°

CLEANING INSTRUCTIONS: Hand wash in warm water (30°C) using mild soap. Rinse thoroughly. Air dry. Do not machine dry. Do not dry clean. Do not iron. **Note:** If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months for soft goods and two years for the Custom ROAM frame and the Schamier ab dem Kaufdatum.

Intended for use on a single patient. Rx Only

LATEX NOT MADE WITH NATURAL RUBBER LATEX

COMPOSITION

Zinc, Aluminum, Stainless Steel, Steel Brass, Polyamide, Polycarbonate, Polyoxymethylene, Polypropylene, Polyurethane, Acrylon.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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