ENGLISH

Before using the device, please read the following instructions completely and carefully. Correct application is vital to the proper functioning of the device.

INTENDED USER PROFILE

The intended user should be a licensed medical professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

INTENDED USE/INDICATIONS

DonJoy Oostestrap is designed to provide support, stabilization and partial immobilization (restriction of movement) of the spine. Postural dysfunction associated with Parkinson's disease. • Traumatic fractures • Osteoporosis • Stable vertebral fractures • Osteoarthritis • Dorsal hyper-kypnosis • Muscular weakness

PERFORMANCE CHARACTERISTICS

Compression: Providing mild compression of limb or body segment through elastic construction. Support: Soft-good/Semi-rigid designed to restrict motion through elastic or semi-rigid construction. Immobilization: Providing immobilization or controlled movement of the body segment. Protection: Providing mild protection of the body segment for acute and prophylactic care.

CONTRAINDICATIONS

• Pregnancy and circulation, pulmonary, cardiovascular or skeletal conditions which pose risk to be made worse as a result of compression of the abdomen and/or torso, unstable, displaced fractures. • Do not use if you are allergic to any of the materials contained within this product.

WARNINGS AND PRECAUTIONS

• This device is prescribed, adjusted and explained by a health care professional. • The product is intended for single patient use only. • We recommend that the first application be made under supervision of a health care professional. • Contact your health care professional if you are in doubt about the application. • Do not wear the device near to open flames. • The frequency and duration of use should be determined by your healthcare professional. • Direct contact of the device to very sensitive skin may cause redness or irritation. Immediate medical attention is needed in case of intolerance. • Do not fit the Oostestrap too tightly, as this could lead to local pressure symptoms. • Do not apply to open wounds. • If patient experiences any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your healthcare professional immediately. • NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to use of this device.

ASSEMBLY, FITTING AND APPLICATION INSTRUCTIONS

This device is designed to be worn during normal everyday activities. It can be worn over a T-shirt. The medical advisor/orthopaedic technician will explain the application of this back brace and make all adjustments necessary for a proper fit. • Both shoulder straps need length adjustment after identifying the ideal length correspondence with your body dimensions. • Cut both shoulder straps to length when you are confident you have identified the right length for ease of use. Wear the back brace a couple of days before you fit the straps to length. Hook fasteners are marked with symbols to help fasten the straps in the right order.

ASSEMBLING THE MODULAR BACK BRACE

This back brace is formed of 2 parts: the lumbar belt with tension strap and the dorsal frame with shoulder straps. Both parts are fixed together to accommodate different body dimensions, spine length and waist circumference. Your medical advisor/orthopaedic technician will put both part together before application and cutting the shoulder straps to length.

HOW TO ASSEMBLE THE DORSAL FRAME WITH SHOULDER STRAPS TO THE LUMBAR BELT

Open the lumbar belt and lay the belt flat on the table with the inside upwards. • The vertical hook strap (1) on the lumbar belt fixes to the vertical loop strip (2) of the dorsal frame. The dorsal frame is positioned in the middle of the lumbar belt with the distal border aligned with the distal border of the lumbar belt (3). • Open the distal border of the pocket in which the dorsal frame is placed and fix the horizontal loop strap (4) of the lumbar belt to the hook strip (5) in the pocket of the dorsal frame. In this way the dorsal frame is kept in the pocket (6) and fixed firmly to the lumbar belt. • Fix the YYY-tabs (7) to the crossing shoulder strap ends (8).

FOR AN INDIVIDUALIZED FIT

For an individualized fit, support and functioning, the medical advisor/orthopaedic technician can bend the aluminum dorsal frame to suit spine curvature. The aluminum stay can be taken out of the pocket, and after bending to suit the patient, placed back in the pocket.

APPLICATION OF THE BACK BRACE

1. Open the lumbar belt and loosen all straps. 2. Put both hook fixation ends of the tension straps (**) and both YYY-tabs of the shoulder straps left and right on the lumbar belt. 3. Put your arms as a back support through the loops made by the crossing shoulder straps. 4. Close the lumbar belt around your waist (*). 5. Secure the lumbar belt firmly around your waist by pulling on both lumbar belt tension straps (**) and fixing them anteriorly. 6. Grip both support strap fixations (YYY-tab) and pull at the same time on both straps until you feel enough backward traction on your shoulders. 7. If necessary, reposition both shoulder straps in the middle of your shoulder for a better fit and improved pulling action. 8. Repeat Step 5 until the desired backward traction on both shoulders is achieved. 9. Fix both shoulder strap fixations (YYY-tab) in the front area of the lumbar belt, within reach for easy adjustments. The pulling tension of the shoulder straps can be easily adjusted to patient needs during daily activities.

CLEANING INSTRUCTIONS

• The aluminum dorsal frame should be taken out of the brace during wash. Please check you place the stay back after washing in the pocket, in a way so that the shape follows the curvature of your spine. • Hook fasteners of the straps and belt should be closed for washing to avoid any damage. • Hand wash the product with lukewarm water (30°) and neutral soap, then rinse carefully with clean water. • Air dry and keep it away from heat sources. Do not machine dry. Do not dry clean. Do not iron. • Contact your medical advisor/orthopaedic technician in case of problems or worn components which need replacement.

COMPOSITION

Polyurethane fibre, Polyester, Polyamide, Cotton, Elastane, Carbon, Aluminium

NOT MADE WITH NATURAL RUBBER LATEX

WARRANTY D.J.O. LLC will repair or replace all or parts of the unit and its accessories for material or workmanship defects for a period of 12 months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

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