## ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

### INTENDED USER PROFILE:

The intended user should be a licensed medical professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable to perform all the instructions, warnings and cautions provided in the information for use.

### INTENDED USE/INDICATIONS:

The Donjoy ACL Everyday is designed to support the knee and provide protected range of motion for those with mild ACL/ MCL/ LCL instabilities. Soft-good/ semi-rigid designed to restrict motion through elastic or semi-rigid construction. Providing mild protection to a limb or body segment for acute and prophylactic care.

### CONTRAINDICATIONS:

None

### WARNINGS AND PRECAUTIONS:

- Loss of circulation, patient discomfort and patient re-injury are potential effects caused by device failure.
- ROM settings should not be changed without supervision of a medical professional.
- This product must be prescribed and fitted by a healthcare professional.
- The frequency, duration of use and directions for use should be determined by your healthcare professional.
- If you develop an allergic reaction and/or experience itching, red skin after coming into contact with any part of this device, please stop using it and contact your healthcare professional immediately.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
- The support should be snug but not impair circulation.
- Do not use over open wounds.
- Do not use the device if it is damaged and/or the packaging has been opened.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

### APPLICATION INFORMATION:

- While sitting on the edge of a chair, step through the brace. Flexing the knee to 45 degrees, position the brace as shown in figure 1, centering the hinge about 1" (2.54 cm) above the top of the kneecap and slightly past midline of the leg. Make sure the brace is not twisted on the leg. (Fig. 1)
- Secure the bottom calf strap first (strap #1), then the top thigh strap (strap #2).
- Keeping the knee brace bent at 45 degrees, push both condyles (hinges) toward the back of the leg. (Fig. 2)
- Maintaining the position of the hinges, secure lower thigh strap (strap #3), then upper calf strap (strap #4) on the back of the leg. (Fig. 3)
- Draw the hinges forward slightly before fastening the front tibial (shin) strap (strap #4). Do not pull the hinges beyond midline of the leg. (Fig. 3)
- Adjust the front tibial strap (strap #5). Be careful not to pull this strap too tight and risk circulation problems.
- When all strapping is complete, the hinge should center just above the top of the kneecap and slightly behind midline. To ensure that no slippage occurs during activity, it is recommended to do a few knee bends or walk for a few minutes, then re-tighten the straps following steps 2-6. (Fig. 4)

### INSTRUCTIONS FOR HINGE ADJUSTMENT:

Your brace comes with 10 degree extension stops installed. To change the installed stops, first completely remove the screw. (Fig. 5) Then remove the stops. Insert the new stops, along with the holes in the hinge plate, then reinsert the screw. Using a Phillips head screwdriver, ensure that the screws are seated in the recess and completely tightened. Both medial and lateral stops must be identical in degree. (Fig. 6) Extension Stops: 0°, 10°, 20°, 30°, and 40°. Flexion Stops: 0°, 45°, 60°, 75°, and 90° (no flexion stops are installed upon shipment)

### USE AND CARE:

- During activity, re-tightening is necessary, re-tighten all straps following steps in the application instructions.
- All straps and pads may be trimmed to fit individual variations in limb circumference.
- The brace may be used in fresh or salt water. Following use in water activities, drain thoroughly and rinse with clean, fresh water and then air dry.
- Hand wash your liners and condyle pads in cold water with mild detergent. AIR DRY ONLY. DO NOT HEAT DRY LINERS OR PADS.
- The hinge mechanism requires little maintenance in normal use. However, you may apply a dry lubricant such as Teflon™ spray to the internal mechanism of the hinge. Dry lubrication is recommended following use in water.
- Periodically inspect hinge screws for tightness and re-tighten as needed.

**WARRANTY:** DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

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NOTE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PROFESSIONAL.

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- Go up to the edge of a chair and pass the leg through the brace. Bend the knee at 45 degrees, position the brace as shown in figure 1, centering the hinge about 1" (2.54 cm) above the top of the kneecap and slightly past midline of the leg. Make sure the brace is not twisted on the leg. (Fig. 1)
- Secure the bottom calf strap first (strap #1), then the top thigh strap (strap #2).
- Keeping the knee brace bent at 45 degrees, push both condyles (hinges) toward the back of the leg. (Fig. 2)
- Maintaining the position of the hinges, secure lower thigh strap (strap #3), then upper calf strap (strap #4) on the back of the leg. (Fig. 3)
- Draw the hinges forward slightly before fastening the front tibial (shin) strap (strap #4). Do not pull the hinges beyond midline of the leg. (Fig. 3)
- Adjust the front tibial strap (strap #5). Be careful not to pull this strap too tight and risk circulation problems.
- When all strapping is complete, the hinge should center just above the top of the kneecap and slightly behind midline. To ensure that no slippage occurs during activity, it is recommended to do a few knee bends or walk for a few minutes, then re-tighten the straps following steps 2-6. (Fig. 4)

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