

# AirSelect™ Elite Walker

# AIRCAST®

Put the walker on only in a seated position and use the included socks for optimal hygiene and the best possible comfort.



## 1. Preparing the walker

- The air chambers must not be filled before the walker is put on
- Loosen straps, remove front plate, open toe covering and foot flaps




## 2. Place the foot in the walker

- Position the heel at the rear
- Close the foot flaps and toe covering



## 3. Put on the front plate

- The AIRCAST® logo is over the toes
- The sides of the lower front plate (foot part) are on the inside of the walker.
- The sides of the upper front plate are on the outside of the walker



## 4. Close straps

- Close straps from bottom to top (A, B & C)
- Tighten straps so that the walker fits snugly but comfortably



## 5. Adjusting the air chamber pressure

- First fill air chambers on the injured side
- Rotate the knob on the right side of the walker to the corresponding number:  
Number 1: fills the right air chamber  
Number 2: fills the left air chamber  
Number 3: fills the rear air chamber
- After all air chambers have been inflated, turn the knob to the lock symbol



## 6. Inflating the air chambers

- Press the center of the large pump (top right) and slowly release it again
- Pump until the air chambers fit snugly but comfortably
- Repeat the procedure for the other air chambers
- If the air chambers are too full, release some air

## Taking off the walker

**For a proper fit when the walker is put back on, the air should be released before the walker is taken off**

- Rotate the large selector wheel to the corresponding number of the air chamber (Fig. 5)
- Press the small release button on the upper right of the walker
- Repeat the procedure for the other air chambers
- Loosen straps, remove front plate, and open toe covering and foot flaps

## Size chart

SIZE	SIDE	SHOE SIZE
XS	UNIVERSAL	32-34
S	UNIVERSAL	35-38
M	UNIVERSAL	39-42
L	UNIVERSAL	43-45
XL	UNIVERSAL	46-48



## Precautionary measures

- Do not overexpand the air chambers
- Check the skin regularly for pressure marks and skin irritation. In the event of skin changes or discomfort, reduce the air chamber pressure
- Do not use the walker if there is no feeling in the foot or lower leg due to postoperative anesthesia or nerve, skin or other disorders

**enovis™**

*Creating Better Together™*

**enovis-medtech.de**

ORMED GmbH, Bötzingen Straße 90, 79111 Freiburg, GERMANY

## Important note!

Please be aware that neither DJO, LLC nor any of its subsidiaries can provide medical advice. The contents of these application instructions do not represent any medical, legal or other type of professional advice. Information with regard to different medical conditions, diseases and physical conditions and their treatment is not intended as a replacement for a consultation by a physician or other healthcare professional. The results may vary individually. These application instructions contain general information for guidance. However, no liability can be assumed for the accuracy of all information and no legal claims can be asserted against the author of the application instructions or against the party who transferred these application instructions.