## ENGLISH

**BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.**

**INTENDED USER PROFILE:** The intended user should be a licensed medical professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

**INTENDED USE/INDICATIONS:** The Procure Nextep Contour Shortie walker is designed to provide support and immobilization to the ankle and foot. It may be suitable for use following stress fractures of the lower limb, stable ankle and/or foot fractures and acute lateral ankle sprains. Providing immobilization or controlled movement of a limb or body segment.

**CONTRAINDICATIONS:** This device is contraindicated for tibial fractures and unstable fractures of the lower leg, ankle and foot.

### WARNINGS AND PRECAUTIONS:

- This device is to be used under the supervision of a healthcare professional. The determination of when to apply the brace and the frequency and duration of use should be strictly at the discretion of the treating physician. Always consult with your physician before making changes to the brace.
- Application of this device is recommended only when the fracture is demonstrably stable and there are acceptable limits of angular and rotational deformity.
- Do not use this device on patient's incapable of communicating physical discomfort.
- Check your skin frequently for skin damage and skin irritation.
- Do not use this device if you cannot feel sensations whether due to post-op anaesthesia, nerve, skin or other conditions.
- Use caution when walking on slippery or wet surfaces to avoid injury.
- Do not use over open wounds.
- Accessory pads may be used where necessary to aid in patient comfort.
- Do not over tighten straps. This may result in reduction of blood flow or sensation.
- Do not modify the device or use the device other than intended.
- If you develop an allergic reaction and/or experience itchy, red skin after coming into contact with any part of this device, please stop using it and contact your healthcare professional immediately.
- Do not use this device if it is damaged and/or the packaging has been opened.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, consult your medical professional immediately.

**NOTE:** Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

### APPLICATION INFORMATION:

- 1) Select correct size of walker for the patient. The end of the patient's toes should be within 1 in (2.5 cm) from , without extending past, the front end of the walker.
- 2) Remove the wrap from the walker and place the patient's foot in the wrap. (Figure 1)
- 3) Close the wrap gently but firmly around the foot, ankle, and lower leg, using the hook and loop strips to secure the wrap in the desired position. (Figure 2)
- 4) Place the foot, within the wrap, firmly into the rear shell with the heel as far back as possible. (Figure 3)
- 5) Secure the patient's foot, ankle, and lower leg in place by looping the straps through the buckles and fastening the hook and loop closures. Begin with the instep strap first, then proceed in any order until all three straps are secured. Adjust the tightness of straps until the walker is secure and comfortable. (Figure 4)

### USE AND CARE:

Hand wash liners in cold water using mild soap, rinse thoroughly. AIR DRY.  
NOTE: If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

**WARRANTY:** DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

### INTENDED FOR SINGLE PATIENT USE.

NOT MADE WITH NATURAL RUBBER LATEX.

Rx only.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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