
AIRCAST®

DEUTSCH

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ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE

The intended user should be a licensed medical professional, the patient, the caregiver, or a family member providing assistance. The user should be able to read, understand and be physically capable of performing all the instructions, warnings and cautions provided in the information for use.

INTENDED USE/INDICATIONS

- Pain relief in patients with symptomatic Hallux Valgus (bunion)
- Mild or moderate mal-positioning of the big toe
- Maintaining position after Hallux Valgus (bunion) surgery

Performance characteristics:

Soft/good/semi-rigid designed to restrict motion through elastic or semi-rigid construction. Providing immobilization or controlled movement of the limb or body segment. Providing mild protection of a limb or body segment for acute and prophylactic (preventive) care.

KONTRAINDICATIONS

- If it is not possible to correct the angle (< 20°) of the base joint of the big toe passively, we advise against wearing the splint and recommend consulting your medical professional.
- If you have the following conditions please consult a physician before using the device: gout attack, peripheral arterial disease or impaired blood circulation, rheumatic diseases, hypersensitivity to pressure, banded joint arthritis, polyneuropathy or soft tissue swelling.
- Do not use this device if you have skin diseases or injuries in the affected area of the body.
- Any known allergies to the medical devices materials.
- No to utilize si tiene cualquier alergia conocida a los materiales del dispositivo médico.

WARNING & PRECAUTIONS

Post-operative treatment: Use only as directed by your health care professional. To achieve maximum benefit, the splint should be worn for several hours. It is recommended that the brace is worn through the night, or in the evening whilst at rest. You may need to gradually build up wear time over a few days. The product is intended for single patient use only. The product must be worn directly on the skin. The splint is not suitable for use while walking. It should be used during periods of rest and/or at night. Do not apply any lotions or ointments on the relevant part of the body while wearing the splint. If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact a medical professional immediately. Note: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INSTRUCTIONS

The splint set consists of (Fig. A):
1 - Big toe strap
2 - Mid-foot strap
3 - Removable side cushion (fixed to the mid-foot part)

The splint can be worn on the right or left foot. For application to the left foot, the straps must be on the left side of the splint. For application to the right foot, the straps must be on the right side of the splint.

- Loosen both big toe strap and mid-foot strap.
- Slide foot through the mid-foot strap.
- Guide the big toe through the big toe strap.
- Adjust the fit of the splint as necessary.
- Tighten both straps firmly to ensure the splint stays in place. The splint should fit snugly against the side of the foot without any discomfort.

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CLEANING INSTRUCTIONS

Straps and cushion: machine washable at gentle cycle (max 30°C) with a mild detergent. Hook and loop of the straps should be closed for washing to avoid any damage. Do not bleach. Rinse thoroughly. Air dry. Do not machine dry. Do not dry clean. Do not iron.

Rigid splint: hand wash with mild detergent. Do not bleach. Rinse thoroughly. Air dry. Do not machine dry. Do not dry clean. Do not iron.

COMPOSITION

Polyamide, polypropylene, Polyvinyl Chloride PVC, Polyurethane

NOT MADE WITH NATURAL RUBBER LATEX

WARRANTY DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

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Treatment post-operative: Use only as directed by your health care professional. To achieve maximum benefit, the brace should be worn for several hours. It is recommended that the brace is worn through the night, or in the evening whilst at rest. You may need to gradually build up wear time over a few days. The product is intended for single patient use only. The product must be worn directly on the skin. The splint is not suitable for use while walking. It should be used during periods of rest and/or at night. Do not apply any lotions or ointments on the relevant part of the body while wearing the splint. If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact a medical professional immediately. Note: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

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El kit de inmovilización se compone de (Fig. A):
1 - Cinturino per l'alluce
2 - Cinturino per il mesopiede
3 - Imbottitura interna rimovibile (fissata alla parte metatarsale)

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