**DONJOY®**

# Immo AT4 Everest

## IMMO AT4 EVEREST

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
# IMMO AT4 EVEREST

# IMMO AT4 EVEREST

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

#### INTENDED USER PROFILE

The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance. The user should be able to:

- Read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

#### INTENDED USE/INDICATIONS

- Immobilisation of the knee in defined position
- Impairment of Bone / Joint capsule / Ligament-structure (following mostly after trauma and/or surgery)
- Temporary stabilization of the knee in a defined position
- All indications which need an immobilization in a defined position, e.g.:
  - Pre-operative
  - Post-operative
  - Post-traumatic
- Providing immobilization or controlled movement of the limb or body segment.

#### OPTIONAL COLD THERAPY INDICATIONS:

- Localized external application of cold therapy and compression after surgery or injury, to help minimize inflammation, swelling, and pain.

#### CONTRAINDICATIONS

Do not use if you are allergic to any of the materials contained in this product. We recommend that the first application be made by a professional.

#### WARNINGS AND PRECAUTIONS

- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
- Do not use over open wounds.
- Do not use this device if it was damaged and/or package has been opened.
- **NOTE:** Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.
- **CAUTION:** This brace is intended to be used for immobilization of a patient's lower limb following injury or surgery. It is NOT intended to support the full body weight of a patient while standing or walking.

#### APPLICATION INSTRUCTIONS

Attach the 3 Panels in a way that the two lateral stabilizers are positioned right and left in the middle of the leg (Fig. 1).

- 1) With the patient lying down, open-up the brace and place it under the patient's extended leg. Adjust the internal anti-slip strap, under the patella (Fig. 2).
- 2) Wrap the brace over the leg, taking care to place the front opening right in the middle of the knee (Fig. 3).
- 3) Close the immobilizer and tighten the closure straps according to the circumference of the thigh and calf. Once fitted, the immobilizer can be easily put into place and removed by using the buckles (Fig. 4). The Immobilizer is delivered with a 0° degree immobilization angle, but it can be adapted to 20° by exchanging the lateral stays and by bending the posterior stay to the desired angle. Therefore open the contact closure tab on top of the casings and remove the stays, exchange and/or bend them and reinsert them into their casings (Fig. 5).

#### CARE

Wash in lukewarm water (30°C) with a mild detergent. Do not bleach. Rinse thoroughly. Air dry. Do not machine dry. Do not dry clean. Do not iron.

#### CARE - Cold Pack

Get inserts may be cleaned with water, moist towel or sponge.

NOT MADE WITH NATURAL RUBBER LATEX.

#### COLD THERAPY & COMPRESSION (OPTIONAL)

For optional application combining Cold therapy and Compression, the Knee Cuff contains 1 Inflatable Pad, 1 Hand Pump and 2 Cold Packs (Fig. 6a & 6b).

#### INTENDED USE/INDICATIONS

- Localized external application of cold therapy and compression after surgery or injury, to help minimize inflammation, swelling, and pain.

#### CONTRAINDICATIONS - COLD COMPRESSION

DO NOT use this device on patients with:  
Raynaud's or other vasospastic diseases; diabetics; sickle-cell anemia, arteriosclerosis, cryoglobulinemia, hypercoagulable clotting disorders or other peripheral vascular diseases causing ischemia or poor local circulation; compromised local circulation; local tissue infection; cold allergy, cold hypersensitivity, previous cold injury.

#### APPLICATION INSTRUCTIONS - COLD COMPRESSION

- 1) Prior to use place the cold pack in freezer for 4 hours minimum. (Domestic deep freeze temperature between -12°C & -20°C). An hygienic bag is supplied for this purpose.
- 2) With the brace settled on the knee, open the 2 middle straps (Fig. 7).
- 3) Cover the Cold-Pack with the cover provided and attach the cold-Pack onto the knee Cuff (Fig. 8a & 8b). Feel Cold Pack and ensure temperature is adequate before use.
- 4) Place Knee Cu with Cold-Pack on the middle of the knee and attach it with the lateral hook straps to the brace (Fig. 9).
- 5) Close the 2 middle straps of the brace (Fig. 10).
- 6) Overlaying straps can be cut.
  - **CAUTION:** It needs to be applied with Inflatable Pad and Cold Pack before cutting the straps (Fig. 11a & 11b)
- 7) For compression: Screw the pump to the Knee Cu, turn the valve to 'ON' position. In ate the Air Cell to the desired level of compression (Fig. 12a). Once the desired inflation level has been reached turn the valve to 'OFF' and remove the pump. (Fig. 12b)
- 8) To de ate, screw the pump to the Knee Cu, position the valve on ON and press the button at the base of the pump to release the air (Fig. 13).
  - With the two removable straps provided the Inflatable Pad and Cold Pack can still be used once the immobilizing brace is no longer needed (Fig. 14).
  - **WARNING!** Be sure to place a barrier/cloth between your skin and the cold pack (See WARNING & PRECAUTIONS).
  - It is critical that your Skin and soft-tissue temperature is not lower than 5°C, so make sure to wear the cover between the Cold Pad and skin!
  - **WARNING!** Do not inflate and tighten the air liner excessively. Ensure level of compression is not too tight so as not too tight so as not to restrict blood or fluid flow (see WARNINGS & PRECAUTIONS).
  - **WARNING!** Upon instruction of a Health Care Professional: Apply the Cold Pack into the knee for time recommended by a Health Care Professional. Do not exceed 20-30 minutes period of application. Alternate with adequate rest periods. If needed, renew application. Replace the Cold Pack in freezer until needed. Do not sleep with the Cold Pack applied.

#### WARNINGS AND PRECAUTIONS - COLD COMPRESSION

- **WARNING!** THIS DEVICE CAN BE COLD ENOUGH TO CAUSE SERIOUS INJURY. SERIOUS ADVERSE REACTIONS AND SAFETY HAZARDS MAY OCCUR WHEN USING THIS DEVICE.

**WARRANTY** DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

**NOTICE:** WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

#### USE OF KNEE CUFF WITH INFLATABLE PAD & COLD PACK :

- A health care practitioner must state how long and how often the device should be used and the length of breaks between uses.
- Extreme care must be taken when using any cold therapy system. Any cold products may cause injuries, e.g. burns, blisters, frostbite, swelling, discoloration and skin or tissue necrosis if improperly used or monitored.
- **WARNING!** DO NOT inflate and tighten the Knee Cuff excessively. Ensure level of compression is not too tight so as not to restrict blood or fluid flow.
- **WARNING!** This device can be cold enough to cause serious injury, including tissue necrosis. You must be able to check your skin condition under the Cold Pack frequently (at least every 15 minutes). DO NOT use this device if you cannot check your skin condition frequently (at least every hour). Check for increased pain, burning, numbness, tingling, increased redness, discoloration, itching, increased swelling, blisters, irritation or other changes in skin condition under the Cold Pack or around the treatment area. If you experience any of these conditions, immediately discontinue use of this device and contact your physician.
- **WARNING!** Application of the Cold Packs directly on the skin may result in serious injury, including tissue necrosis. DO NOT let any part of the Cold Pack touch your skin. ALWAYS use with a barrier/cloth between your skin and the Cold Pack.
- Keep the product away from children.
- All Cold Pack ingredients are non-toxic. They are, however, not intended for ingestion.
- Do not heat the Cold Pack (microwave, oven, hot water, etc..)
- Do not use on open wound(s) and/or sensitive skin.
- Do not use this product if it is damaged or defective.
- Do not use this product in any way that may damage its integrity.
- Do not sit or lay over the Cold Pack or fall asleep with the Cold Pack on the skin.
- If contents get in contact with eye or skin, rinse with water immediately.

#### LIMITATIONS ON SPECIAL PATIENT POPULATIONS

Limit the use of this device with patients who are unresponsive, incapacitated, have altered mental status or altered pain perception. Post-surgical patients under sedation or on analgesics or anesthetics, as well as patients taking hypnotics, anxiolytics, or antidepressants, must be monitored frequently during use of this device. These patients may not be able to perceive pain, burning, numbness, tingling or decreased sensation and may be susceptible to injury. Discontinue cold therapy immediately at the first sign of cold injury.

#### COMPOSITION

##### Knee Immobilizer -

Fabrics: Polyamide; Polyurethane; Polyester; elasthane, silicone, polyethersulfone, silicone

Buckles : acetal

Stays : Aluminium

##### Knee Cuff -

Inflatable Pad: Polyester, nylon, polyethylene

Pump: Silicone, Polyethylene, PVC, ABS Resin (thermoplastic polymer), synthetic rubber and steel wire, polypropylene, polycarbonate

##### Cold Pack -

Gel: propylene glycol, sodium carboxymethyl cellulose, water.

Fabric: cotton, nylon, polyester

Protective cloth : Polyester, cotton, polyurethane, polyethersulfone

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MDSS GmbH
DJO FRANCE S.A.S.