### DEUTSCH

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### ENGLISH

**BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.**

#### INTENDED USER PROFILE

The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance. The user should be able to:

- Read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

#### INTENDED USE/INDICATIONS

The ankle brace is designed to provide support and stabilization for a wide range of conditions including, ankle ligament tears, sprains, chronic instabilities and post-operative use. For maximum comfort and support, wear an absorbent sock and a laced athletic shoe when wearing the brace. Take care not to damage clothing when attaching/detaching hook and loop straps. Soft/good/semi-rigid designed to restrict motion through elastic or semi-rigid construction.

#### CONTRAINDICATIONS

Do not use if you are allergic to any of the materials contained in this product.

#### WARNINGS AND PRECAUTIONS

1. The brace should only be worn under supervision of a health care professional.
2. If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
3. Do not use over open wounds.
4. Do not use this device if it was damaged and/or package has been opened.

**NOTE:** Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

#### BOA® LACING SYSTEM (FIG. 3)

##### To Loosen BOA® (Fig. 3a):

- Lift BOA Knob.
- Pull the lacing

##### To Close and Tighten BOA® (Fig. 3b):

- Push BOA Knob down and then turn BOA Knob clockwise.
- BOA Knob must be in the down locked position when wearing brace.

#### APPLICATION INSTRUCTIONS

##### PHASE I – IMMOBILISATION (Strong Stability)

- Loosen the laces, unfasten straps and open brace.
  - Slide foot into the back of brace. Pull the lateral shells upwards and align with the sides of the ankle (Fig. 0).
  - Close the straps according to numbering:
1. Pull the top strap around the back of the calf and fix to the external side of the brace without tightening too much (Fig. 1). This helps position the lateral shells correctly.
  2. Pull the bottom Achilles strap around the back of ankle and fix to the external side of the brace without tightening too much (Fig. 2).
  3. Tighten the brace by turning the BOA Lacing System clockwise until the desired level of support (Fig. 3).

Note: If the ankle brace is too tight, open the BOA system, loosen the laces, push the knob down and re-tighten.

4. Pull the front cross-strap tightly through the buckle and fix firmly onto front of brace (Fig. 4).

##### 5. FIGURE-8 STRAP (Fig. 5a-e)

- a- To increase stability, you can extend the front cross-strap by using the figure-8 strap. Attach it to the cross-strap. The arrow on the strap should be visible.
  - b- Wrap the extended strap, centered around your Achilles tendon.
  - c- Bring the strap forward, crossing over your forefoot.
  - d- Run the strap under your foot.
  - e- Fasten the strap across the front.
- If the figure-8 strap is too long, shorten it.

##### PHASE II – STABILISATION (Early mobilization and Moderate Stability)

- Warning:
- Flexible shells should be fitted by a qualified healthcare professional.
- Do not heat flexible shells directly on the body.
- Ensure flexible shells are not too hot prior to application.
- Do not heat flexible shells to temperatures higher than 550°C.

##### Heating and Application of the flexible shells

6. Remove the brace and remove the lateral rigid shells by opening the hook & loop closure of the internal shell pockets at the top end of the brace (Fig. 6).

7. Set a heat gun to 350°C. Heat up the flexible shells for 20 seconds on each side.

8. Remove the Figure-8 strap. Pull the anterior cross strap tightly through the buckle and close using the hook and loop (Fig. 8).

Apply the brace in the usual manner (Fig. 0-5e). After approximately 5 minutes the flexible shells will have cooled down and are adjusted to the ankle.

##### PHASE III – MOBILISATION (Moderate Stability)

8. Remove the Figure-8 strap. Pull the anterior cross strap tightly through the buckle and close using the hook and loop (Fig. 8).

#### CARE

- Remove the shells.
- Hand wash in lukewarm water (30° C) with a mild detergent
- Air dry (do not place on radiators/heaters or other heat sources)
- Do not use bleach or chlorine, do not iron.












#### COMPOSITION

Plastic, Silicone, PES (polyether sulfone), Polyester, Polyamide, Elastane, natural rubber (latex), Polyurethane, Nylon

**FOR SINGLE PATIENT USE ONLY.**  
**CAUTION: THIS PRODUCT CONTAINS NATURAL RUBBER LATEX WHICH MAY CAUSE ALLERGIC REACTION.**  
(< 0.1% - NOT IN DIRECT CONTACT WITH SKIN)

**WARRANTY** DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

**NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.**

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