ARTROMOT®-S4

SW-Version: 1.X

- DE · GEBRAUCHSANLEITUNG
- GB · OPERATING INSTRUCTIONS
- FR · MODE D'EMPLOI
- ES · MODO DE EMPLEO
- IT · ISTRUZIONI PER L'USO
- NL · GEBRUIKSAANWIJZING

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Anhang:
ARTROMOT®-S3
ARTROMOT®-S4
ARTROMOT® ACTIVE-K
ARTROMOT®-K1
ARTROMOT®-SP3
ARTROMOT®-E2

Contents

Description of the ARTROMOT®-S4	II
Photos illustrating the setup of the ARTROMOT®-S4	III
1. How to use the CPM device (Intended Use)	59
1.1 Fields of application	59
1.2 Therapy objectives	59
1.3 Indications	59
1.4 Contraindications	59
1.5 Secondary effects	59
1.6 Biocompatibility	60
1.7 Essential Performance	60
1.8 Frequently used functions	60
1.9 User profiles	60
1.10 Intended patient population	61
1.11 Intended contact points of the device with the patient's tissue	61
1.12 Intended application setting	61
1.13 Expected service life	61
1.14 Ambient conditions during operation	61
2. Description of the ARTROMOT®-S4	62
2.1 Description of the device components	63
2.2 Description of the programming unit	64
2.3 Explanation of symbols	68
2.4 Explanation of motor controls	71
2.5 Explanation of symbols (connections and rating plate)	72
3. Safety information	74
4. Device setup	78
4.1 Scope of delivery, device connection, performance check	78
4.2 Adjusting the device to the patient	80
4.3 Switching the device off	82
5. Setting the treatment values	83
5.1 General information on programming the ARTROMOT®-S4	83
5.2 Overview of the parameters available for the different motion patterns	86
5.3 Other functions/settings	95
5.4 Application and programming examples	100
6. Care, Maintenance, Transport, Conversion	103
6.1 Care / re-use	103
6.2 Maintenance (fuse replacement)	104
6.3 Transport	105
6.4 Conversion	108
7. Environmental Protection Statement	109
8. Specifications	109
9. Contact	111
10. Technical service	111
10.1 Technical hotline	111
10.2 Shipment	111
10.3 Spare parts	112
11. Declaration of conformity	113
12. IEC 60601-1-2 and IEC 60601-1-11	338

1. How to use the CPM device (Intended Use)

1.1 Fields of application

The ARTROMOT®-S4 is a motor-operated **Continuous Passive Motion (CPM)** device used to mobilize of the shoulder joint.

Here, the patient sits on the device chair during therapy. The arm to be treated is rested on the arm support assemblies (19, 20, 25) of the motion element (32), which is mechanically adjusted so that the joint axes of the motion element coincide with the joint axes of the shoulder joint. After programming the range of motion for the individual motors (A, B and C), the motion element moves the arm and thus the patient's shoulder joint slowly and continuously in the set, pain-free range of motion, without the patient having to actively contract the muscles. Movement of the shoulder joint takes place exclusively passively through the movement element.

Suitable for use in hospitals, clinics, general practices and rental services as well as in the patient's home, it is an important supplement to medical and therapeutic treatment.

1.2 Therapy objectives

CPM therapy with the ARTROMOT®-S4 is mainly used to prevent the negative effects of immobilization, to allow patients to regain painless mobility of joints at an early stage and to promote healing and achieve a positive functional result.

Other objectives of therapy include:

- improvement of joint metabolism
- prevention of joint stiffness (arthrofibrosis)
- promotion of the regeneration and healing of cartilage and damaged ligaments
- faster hematoma/fluid resorption
- improved lymph and blood circulation
- thrombosis and embolism prophylaxis

1.3 Indications

The CPM device is indicated in the treatment of most injuries and postoperative conditions and diseases of the shoulder joint. Examples:

- joint distortion and contusion
- arthrotomy and arthroscopy procedures in combination with synovectomy, arthrolysis, or other intra-articular interventions
- all types of arthroplasty
- mobilization of joints in anesthetized patients
- surgically treated fractures, pseudoarthrosis, if exercise-stable
- decompression surgery (acromioplasty)
- endoprosthetic implants
- soft tissue surgery in the armpit and the shoulder girdle
- tumor surgery in the shoulder region

1.4 Contraindications

Do NOT use the ARTROMOT®-S4 on patients with:

- acute inflammatory processes in the joints, unless on the order of a physician
- spastic paralysis
- unstable osteosynthesis

1.5 Secondary effects

Currently, there is no evidence of desired or undesired secondary effects resulting from CPM therapy.

1.6 Biocompatibility

Those parts of the ARTROMOT®-S4 that come into contact with the patient when the device is used as intended, are designed to fulfill the biocompatibility requirements of the applicable standards.

1.7 Essential Performance

- The programmed angles are accurately maintained with a tolerance of $\pm 2^\circ$.
- The programmed speeds are accurately maintained with a tolerance of $\pm 5\%$.
- The selected mode and mechanical settings do not change during operation.
- The device can be stopped with the STOP button and all other buttons.

1.8 Frequently used functions

- Unpacking (device and accessories)
- Inserting the system rack and locking it at the chair and backrest
- Mounting the armrest
- Inserting the motion element and locking it at the system rack
- Device connections (mains power connection, connecting the control unit and the motion element)
- Mechanical adjustments of the motion element on the rotational axes
- Programming the control unit (range of motion, speed, functions)
- Storage

1.9 User profiles

The ARTROMOT®-S4 is intended for use by medical specialists ("professional users") properly trained in the settings and use of the device, such as: specialist physicians, physiotherapists, occupational therapists and private care nurses.

Users must be able to do the following:

- Read and understand the instructions for use, warnings and safety notes.
This requires at least 8 years of reading experience
- Reading and understanding the precautions and contraindications of the device
- Observing and understanding the safety labelling
- Observing the information on the display

Furthermore, the medical device is also intended for use by patients ("lay users"). To this end, the patient must be trained in the operation of the device, and must at least understand how to start and stop the device.

A patient advisor or an equivalent person from the rental service organisation is required for maintaining the device and acting as a contact for the patient. This person is in charge of assembling the device, adjusting the mechanical settings according to the patient's size, initial programming of the device according to medical specifications and instruction of the patient in the programming.

The patient caregiver must explain the operation of the ARTROMOT®-S4 to the patient and inform the patient about the warnings, precautions and hazard information contained in the instructions for use.

1.10 Intended patient population

The ARTROMOT®-S4 is intended for passive mobilisation of the shoulder joint. To this end, the device settings should be optimally adapted to the patient. The following requirements apply to the patient:

- At least 8 years of reading experience is required to read and understand the operating instructions.
- The patient weight should not exceed 175 kg
- To enable an anatomically correct adjustment, the patient should not be taller than 210 cm
- To enable safe operation, the patient must be awake, alert and mentally competent

1.11 Intended contact points of the device with the patient's tissue

The patient sits on the device seat, touching the surface of the seat (31), the backrest (30) and the headrest (28). The to be treated arm of the patient rests on the arm support assemblies (19, 20, 25) of the motion element (32). The healthy arm of the patient rests on the armrest (4).

The device may only come into contact with healthy skin. A closed wound treatment for the affected arm or shoulder must therefore be ensured.

1.12 Intended application setting

Professional setting:

The ARTROMOT®-S4 is intended for use in a clinical setting (hospital/rehabilitation facility) or medical practice and can be moved from room to room.

Home environment:

On the basis of a doctor's prescription, the ARTROMOT®-S4 can also be used in the patient's home environment (e.g. as part of a rental service model).

1.13 Expected service life

The "expected service life" of the medical device, including all parts and accessories supplied, is 10 years under normal use and correct maintenance. Any use beyond this period is the responsibility of the operator.

1.14 Ambient conditions during operation

See Chapter 8. "Technical data" – "Ambient conditions (operation)"

2. Description of the ARTROMOT®-S4

Operating principle

The motorized CPM device supports the following movements of the shoulder joint:

Adduction/abduction 0° - 30° - 160°

Internal/external rotation 90° - 0° - 90°

Extension/Flexion (Elevation) 0° - 30° - 160°

Anteversion/retroversion 110° - 0° - 10°

It can be reconfigured for use on either side.

Note!

To unambiguously represent the current position of the CPM device, the values for internal rotation and retroversion are marked with the symbol "–" both on the display and in this document.

Note!

Elevation denotes the movement of raising the arm above 90° in abduction. Extension and flexion are adjusted via the values entered for adduction and abduction (see chapter 5.4)

These are some of the features of the ARTROMOT®-S4:

- anatomically correct adjustment
- physiological movements
- extensive ranges of movement
- programming unit for precise adjustment of patient-specific therapy values
- chip card for storage of the programmed therapy parameters
- easy transport

2.1 Description of the device components

Note: Please fold out page III!

1. Programming unit (type B applied part)
2. Patient chip card
3. Compartment for storage of programming unit
4. Armrest for healthy arm
5. Height adjustment of armrest for healthy arm
6. Width adjustment of armrest for healthy arm
7. Wheels
8. Crank handle for height adjustment (rotational axis motor A and B)
9. Crank handle for adjustment of shoulder width (rotational axis motor A and B)
10. Crank handle for adjusting the depth of the backrest (rotational axis motor A)
11. Motor A (with control keys)
12. Motor B (with control keys)
13. Motor C (with control keys)
14. Swivel arm – motor A
15. Receiving tube for length adjustment (upper arm)
16. Clamping lever for length adjustment (upper arm)
17. Insertion tube for length adjustment (upper arm)
18. Locking screw for elbow angle range
19. Elbow support
20. Forearm support (type B applied part)
21. Strap for forearm restraint
22. Receiving tube for length adjustment (forearm)
23. Clamping lever for length adjustment (forearm)
24. Insertion tube for length adjustment (forearm)
25. Hand rest (type B applied part)
26. Clamping lever for length adjustment (forearm support)
27. Length adjustment of forearm support
28. Headrest with ball joint
29. Adjustment for depth of headrest
30. Backrest
31. Seat
32. Motion element
33. Thumbscrew for mounting the motion element on the system rack
34. Width adjustment and lateral mounting of motion element
35. System rack
36. Connection backrest/system rack
37. Control cable motor C
38. Control cable motor B
39. Control cable motor A
40. Connection for motor control cable or control pendant
41. Connection for motor control cable or control pendant
42. Connection for power cord
43. Power switch (ON/OFF)
44. Instrument fuse
45. Thumbscrew for mounting the system rack


Subject to change without notice.

Note!

- **The motion element** (32), with the integrated motors A-C (11-13), the length settings (15-17, 27) and the arm support assemblies (19, 20, 25), is **the moving application part** within the meaning of IEC 80601-2-78.
- As of serial number: 2,000
The legs of the device chair are removable to enable space-saving shipping. For assembly and disassembly of the legs of the device chair, see Chapter 6.3 „Transport“ – 6.3.2. „Packed – dismantled state“.

2.2 Description of the programming unit

2.2.1 Programming unit in normal mode

Note!

When tapping on the display during operation, the question "STOP?" will be asked and the session can be terminated. By tapping "STOP?", you stop the CPM device.


⚠ Warning!

Patient hazard –

To operate the ARTROMOT®-S4, patients must at least know how to start and stop the device.

The MENU button on the remote control can be locked to prevent unintentional change of the settings. For this purpose, please consult the instructions in chapter 5.1 „General programming instructions for the ARTROMOT®-S4“.


2.2.2 Display in normal mode, standard view


Note!

When tapping the display during operation, the question "STOP?" will be asked and the session can be terminated. By tapping "STOP?", you stop the CPM device.


2.2.3 Display – session stopped

2.2.4 Display in MENU selection mode


2.2.5 Display in ROM programming mode


2.2.6 Display in general programming mode

selected function
(here: speed)

set value of selected
function (here: 50%)

2.2.7 Display in Expert programming mode

2.3 Explanation of symbols

2.3.1 Parameters available for configuration in each menu

Main Menu	
Parameter	Selectable functions / parameters
 Anteversion/Retroversion	 Anteversion  Retroversion
 Adduction/Abduction	 Adduction  Abduction
 Rotation	 Internal Rotation  External Rotation
 Expert	 Therapy  Warm-up  New Patient
	 Documentation  Ri/Le Conversion  Settings
 Speed	
 Timer	

Documentation	
Parameter	Selectable functions / parameters
 Documentation	 Therapy Time  Therapy Trend


Therapy				
Parameter	Selectable functions / parameters			
 PNF	 Pause	 Speed	 Timer	
 Isolation	 Cycles	 Standby		
	 Pause	 Speed	 Timer	
 Oscillation	 Cycles	 Pause		
	 Speed	 Timer		
 Stretching	 Retroversion	 Abduction	 Internal Rotation	
	 External Rotation	 Cycles	 Pause	
	 Speed	 Timer		
 Wave	 Left	 Right		
	 Speed	 Timer		

Special treatment parameters available for configuration				
Parameter	Selectable functions / parameters			
 Pause	 Pause up	 Pause down		
 Standby	 Anteversion/Retroversion	 Adduction/Abduction	 Rotation	

Settings	
Parameter	Selectable functions / parameters
 Settings	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">  Load Reversal </div> <div style="text-align: center;">  Transport </div> <div style="text-align: center;">  Language </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div style="text-align: center;">  Time/Date </div> <div style="text-align: center;">  Unit Run Time </div> <div style="text-align: center;">  More </div> </div>

More	
Parameter	Selectable functions / parameters
 More	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">  Brightness </div> <div style="text-align: center;">  Volume </div> <div style="text-align: center;">  Factory Defaults </div> </div>

2.3.2 Additional symbols on the control pendant

Settings	
Parameter	Selectable functions / parameters
 Lock	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">  Menu </div> <div style="text-align: center;">  Angle </div> <div style="text-align: center;">  Displayed View </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div style="text-align: center;">  Expert </div> <div style="text-align: center;">  Motor On/Off </div> <div style="text-align: center;">  Motor Controls </div> </div>
 Factory Defaults	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">  Software Information </div> <div style="text-align: center;">  Operating Hours </div> <div style="text-align: center;">  Error Log </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div style="text-align: center;">  Calibration </div> <div style="text-align: center;">  Load Display </div> <div style="text-align: center;">  SW Update </div> </div>
 Displayed View	

2.4 Explanation of motor controls

Each motor of the ARTROMOT®-S4 has controls for direct ROM programming.

To program the ranges of motion with the motor controls, use

- the arrow buttons to adjust values
- the OK button to confirm and save the last input, to switch to the next parameter of the motor concerned

Note!

- When you program the ranges of motion with the motor controls, the CPM device moves to the respective settings at the same time.
- When the motor controls are pressed, the following parameters are always adjusted first:
 - Motor A: Anteverision – OK – Retroversion – OK
 - Motor B: Adduction – OK – Abduction – OK
 - Motor C: Internal Rotation – OK – External Rotation – OK

If you want to adjust only one of the maximum values, it is sufficient to set this value at the motor concerned and to press the OK button twice afterwards. The CPM device automatically identifies the maximum value to adjust and saves the value. The second value of the respective plane of motion remains unchanged.

The direct adjustment via the motor controls can be enabled and disabled in the lockout menu of the control pendant (see section: 5.3).

Default setting: enabled

- OK button to
- confirm and save the last input,
 - advance to the next parameter of the selected motor

arrow buttons to adjust values

2.5 Explanation of symbols (connections and rating plate)

Symbol	Description
	Double or reinforced insulation
	Alternating current
	Power switch OFF
	Power switch ON
	The name next to this factory symbol is the manufacturer
	The year and month of manufacture are indicated next to this factory symbol
	The number next to this symbol is the article reference number
	The number next to this symbol is the serial number
	This is a medical device EN ISO 15223-1:2016
	Device complies with Council Directive 93/42/EEC about medical devices, tested and approved by DQS Medizinprodukte GmbH
	Device complies with ETL requirements for the USA and Canada
	Caution! Observe warnings set forth in instruction manual!
	Refer to instruction manual! ISO 7010-M002
	Refer to instruction manual! ISO 7000 -1641
	Do not dispose with unsorted municipal waste

Symbol	Description
	Protect from moisture
	Type B applied part IEC 60417-5840
	No pushing ISO 7010-P017
	Warning "Crushing of hands" ISO 7010-W024 Warning symbol pinch point hazard! Warning! Depending on the device settings, the moving parts of the device present pinch points. Pay particular attention to small children and babies!
	Motor A swivel arm With decreasing retroversion angles, the swivel arm of motor A moves below the frame. There is a higher risk of crushing injuries at small retroversion angles.
	3V type CR2032 (in control pendant: not user-replaceable! Contact Technical Service)
IP21	The IP rating indicates the level of protection and thus the suitability of the device for use under different ambient conditions. IP21 means: 2 is the level of protection against contact and solid objects. The digit 2 means: – Protection against contact: protected against contact with a finger – Protection against solid object: protected against solid objects (diameter of 12.5 mm and greater) 1 is the level of protection against moisture – The digit 1 means: protection against vertically falling water drops
	The overall weight of the device is indicated next to this symbol.
	Limit values for temperatures (labelling on packaging for storage and transport) – values indicate the temperature limit values to which the medical device can be safely exposed. ISO 7000-0632
	Limit values for air pressure (labelling on packaging for storage and transport) – values indicate the air pressure range to which the medical device can be safely exposed. ISO 7000-2621
	Limit values for humidity (labelling on packaging for storage and transport) – values indicate the humidity range to which the medical device can be safely exposed. ISO 7000-2620

3. Safety information

Definitions

Be sure to read the safety information before using the CPM device. The safety information is classified as follows:

Danger!

indicates an imminent hazard. If not avoided, this hazard will result in death or serious injury.

Warning!

indicates a hazard. If not avoided, the hazard can result in death or serious injury.

Caution!

indicates a potential hazard. If not avoided, this hazard may result in minor injury and/or product/property damage.

Safety information

Danger!

Explosion hazard –

The **ARTROMOT®-S4** is not designed for use in areas of rooms used for medical purposes where an explosion hazard may occur. An explosion hazard may result from the use of flammable anesthetics, skin cleansing agents, and disinfectants.

Warning!

Patient hazard –

- The **ARTROMOT®-S4** is **not** designed for use in an oxygen-enriched atmosphere.
- The **ARTROMOT®-S4** is **not** designed for use in combination with flammable substances.
- Only **authorized individuals** are allowed to operate the **ARTROMOT®-S4**. Individuals are authorized after receiving training in the operation of the device and reading these operating instructions. All persons who have undergone the appropriate training can instruct users in the correct operation; alternatively, DJO offers training worldwide; please contact your DJO office or your local dealer for details (see chapter 9 "Contact"). The training takes about 2 hours. Equipment training should be repeated every 2 years (this is a recommendation, not a requirement).

- Before using the CPM device, the operator must ascertain that it is in correct working order and operating condition. The cables and connectors, in particular, must be checked for signs of damage. Damaged parts must be replaced immediately, before use.
- **Before therapy, a test run** consisting of several exercise cycles must be completed, first without and then with the patient. Check that all setting screws are tightened.
- Stop therapy immediately when you have doubts about the device settings and/or the therapy protocol.
- Ensure an **anatomically correct setup of the CPM device** suitable for the patient to be treated. For this purpose, check the following settings/positions (see numbers on device):
 1. Height adjustment of motion element
 2. Shoulder width adjustment of motion element
 3. Backrest
 4. Length adjustment of upper arm
 5. Adjustment of elbow angle
 6. Forearm guide
 7. Forearm length adjustment
- The CPM device must be stopped before any of the adjustments 1 through 7 is modified.
- Movement must always be **free of pain and irritation**.
- The patient must be **fully conscious** during instruction and when using the device.
- The **choice** of the **therapy parameters** to program and of the therapy protocols to use is **restricted to** the responsible **physician or therapist**. The doctor or therapist must decide on a case-by-case basis whether the device can be used with each particular patient.
- The patient must be familiar with the functions of the ARTROMOT®-S4 **programming unit** and the unit must be **within easy reach** of the patient, allowing him or her to stop therapy if needed.
- **Minors or persons with reduced physical, sensory or mental capabilities** must **never be left unattended** during therapy sessions; specialist staff, a parent or legal guardian must always be present.
- Write the patient's name on the **patient chip card**. The card should only be used for this patient. If the **patient chip card** is used for another patient, be sure to delete the **previous patient data** from the card first (see sections 4.1 and 5.3, paragraph "New Patient"). Use **original chip cards only**.
- The ARTROMOT®-S4 must only be used with **accessories** approved by ORMED.
- Modifications to the medical device described in this document without the manufacturer's written consent is prohibited.
- Do not allow **parts of the body or objects** (such as blankets, cushions, or cables) to get caught in the **moving parts** of the CPM device.

Warning!

Shock hazard –

Strictly observe the following warnings. Failure to do so endangers the lives of the patient, the user, and other persons involved.

- Allow the ARTROMOT®-S4 to reach room temperature **before use**. If the device was transported at **temperatures below 0°C (32°F)**, allow it to stand at room temperature for about 2 hours, until all condensation has disappeared. If the devices was stored at temperatures in the maximum storage temperature range, also allow it to stand for 2 hours to reach room temperature.
- The ARTROMOT®-S4 can only be operated in **dry rooms**.
- When disconnecting the device from the power line, remove the plug from the wall outlet first, before disconnecting the cable from the device.
- When connecting the device to other equipment or when creating a medical system, check that the sum of leakage currents will not cause any hazard. Please contact DJO if you have questions in this matter.
- Do not use multiple portable socket outlets (MPSO) to connect the device to the power line. The ARTROMOT®-S4 must be connected to a **properly installed wall outlet with a non-fused earthed wire**. Before connecting the power cord, it must be completely unrolled and placed such that it will not get caught in the moving parts of the device.
- Before cleaning and service interventions, **disconnect the device from the power line** by removing the power cord from the wall outlet.
- **Liquids must not be allowed to enter the CPM device** or the **programming unit**. If liquids have entered into the devices, the ARTROMOT®-S4 must be immediately checked by a service technician, before it can be reused.

Warning!

Equipment malfunction –

- **Magnetic and electrical fields** are capable of interfering with the proper performance of the device. For this reason make sure that all external devices operated in the vicinity of the device comply with the relevant EMC requirements. X-ray equipment, MRI devices, radio systems, cell phones, etc. are possible sources of interference as they may emit higher levels of electromagnetic radiation. Keep the device away from such equipment and verify its performance before use.
- **Refer repair and maintenance** to authorized persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.
- **Route all cables** such that they will not get caught by the moving parts during operation and that they do not present a stumbling hazard.
- **Inspect** the ARTROMOT®-S4 for damage and loose connections at least **once a year**.

Caution!

Patient hazard –

Preventing chafing and pressure sores. When your patient is **adipose**, very **tall** or **very short**, be sure to prevent chafing and pressure sores

Caution!

Patient hazard, equipment damage –

- Do not use the CPM device to transport persons.
- Extreme caution should be taken when in use around small children and babies or pets! Sufficient distance to the device is mandatory for their safety!
- Never leave the device unattended when it is switched on! Switch the device off and disconnect the power line from the wall outlet!
- After use, store the device in a safe place! Ensure device stability also during storage!
- **This device is not a toy!**

Caution!

Equipment damage –

- Check that the voltage and frequency ratings of your local power line are those indicated on the nameplate.
- The set is designed for a maximum continuous load of 175 kg.
- Do not allow any objects (such as blankets, cushions, or cables) to get caught in the moving parts of the CPM device.
- Do not expose the ARTROMOT®-S4 to direct sunlight, because some of the components may reach inadmissibly high temperatures.
- Be aware that the connectors can only be inserted in the correct orientation and secure all connections with the locks.

Caution!

Patient hazard –

Under severe ambient conditions, the following applied parts may reach temperatures of up to 46 °C:

- TOUCH display control unit (45.8 °C)
- Membrane keypad motors (42.0 °C)
- Drive housing motors (42.6 °C)
- Power supply housing (41.2 °C)

Avoid touching these areas for more than 1 minute.

Children and people with skin and tissue defects or impaired perception of temperatures should avoid touching these areas.

4. Device setup

Note: For a better understanding of each step, please fold out pages II and III.

4.1 Scope of delivery, device connection, performance check

The equipment supplied includes these items:

Seat, programming unit (1), patient chip card (2), power cord, system rack (35), motion element (32), forearm support (20), armrest for healthy arm (4), element for height adjustment of armrest (5), operating instructions

1. Mount the system rack onto the chair by inserting it in the holder on the back until only the OK label is visible. Secure the rack by tightening the thumbscrew (45).
2. Fasten the motion element (32) on the width adjustment/lateral adjustment unit (34) of the system rack. Fix it by fully tightening the thumbscrew (33). Check that the red mark points to the recess.
3. Connect the motor control cable (39) to one of the corresponding sockets (40 or 41) and screw it tight.
4. Connect the programming unit (1) to one of the corresponding sockets (40 or 41) and by plugging the cable in and screwing it tight.
5. Connect the backrest (30) to the crank handle for adjusting the depth of the backrest (10) by introducing the end of the crank handle into the guide rail in the backrest (36).
6. Connect the power cord to socket (42) of the device and connect the power plug to a wall outlet with a non-fused earthed wire (100 to 240 Volt, 50/60 Hz).
7. Turn the power switch (43) on.
8. The display of the programming unit prompts you to indicate whether the motion element is mounted on the right or on the left. Indicate the appropriate side: the motion element will automatically move to the position in which the forearm support can be mounted.

9. Secure the forearm support by inserting the receiving tube for length adjustment (22) into the guide and fastening it with the clamping lever for length adjustment of the forearm support (26).
10. Press the START button: the device automatically moves to the starting position (similar to "New Patient" function).
11. Finally, mount the armrest for the healthy arm on the motion element on the opposite side. For this, open the thumbscrew for width adjustment (6) and insert the fastening element of the height adjustment unit, then tighten the thumbscrew again. Open the clamping lever for height adjustment of the armrest (5) and insert the armrest (4), then close the clamping lever again.

Caution!

Personal injury –

Only the supplied control pendant and the motion element may be connected to sockets (40) and (41). It is not permitted to connect any other devices to these sockets.

Initial adjustment for new patients

- Write the patient's name on the back of the chip card. Insert the original patient chip card (2) into the control pendant (1).
- Press the **MENU** button on the programming unit.
- Select the "Expert" menu.
- On the display, tap the "New Patient" field and confirm by tapping "OK".
- Press the **START** button. The CPM device automatically enters the **home position**. Any existing parameter settings will be deleted.

Adjustment with programmed chip card

- Insert the original patient chip card (2) into the programming unit (1).
- Press the **START** button.
- The device automatically enters the **starting position** (maximum adduction value, midway of internal / external rotation settings, midway of anteversion/retroversion settings).
- Perform the mechanical adjustments according to the notes on the back of the chip card.

Performance check:

If the programming unit can be operated as described above and the **ARTROMOT®-S4** enters the home position (for home position values, refer to section 5.3), the device has passed the performance check and functions correctly.

The device also runs performance checks regularly during operation. This is what happens, if a problem is identified:

- an audio signal sounds
- the device switches off immediately
- a symbol indicating the location of the problem, the message ERROR and an error code (e.g. ERROR 5) appear on the display.

In this situation, you may attempt to restart the device by turning it briefly off and on again with the power switch. If the error message persists, have the device inspected by a Service technician, before using it again.

If it has been determined that the ARTROMOT®-S4 operates perfectly, ask the patient to sit down on the ARTROMOT®-S4 chair.

Note!

Refer to the ARTROMOT®-S4 Service Manual for a list of error messages and troubleshooting hints. The Service Manual can be ordered from DJO, or contact your local dealer or Technical Service.

You can find all contact addressed in chapter 9 "Contact".

4.2 Adjusting the device to the patient

Note!

For the following adjustments, the patient's arm should not yet be placed on the forearm support of the motion element. After you have adjusted the device to the approximate patient measurements, the patient can place his/her arm on the armrest for a check of the setup and for the final adjustment.

The settings are numbered 1 through 7. For easy orientation, you will find the same numbers on the CPM device. Always set up the CPM device in this order.

1. Height adjustment of the motion element (rotational axis motor B – height)
2. Adjustment of the shoulder width of the motion element (rotational axis motors A and B – width)
3. Adjustment of the depth of the backrest (rotational axis motor A – depth)
4. Upper arm length
5. Elbow positioning angle (rotational axis motor C)
6. Forearm support (rotational axis motor C)
7. Forearm length

Then, position the patient's arm and perform the final adjustment.

On completion of the final adjustment, note down the mechanical adjustment values on the back of the patient chip card.

Before you adjust the ARTROMOT®-S4 to the patient, you may have to convert the device for use on the left or right shoulder. (see section 6.4 "Conversion")

Before any treatment is possible, adjust the device as follows and check the settings:

4.2.1 Patient positioning

Before you start setting up the motion element, you must adjust the CPM device so it is comfortable and anatomically correct for the respective patient.

For patient positioning, the ARTROMOT®-S4 provides a number of padded, exchangeable supports which can be attached and secured where required.

1. Ask the patient to sit down on the seat and to lean back.
2. Next, adjust the **headrest (28)**. This ensures that the patient is sitting in the middle of the chair: The height is adjusted by means of the height adjustment unit of the motion element (8); the depth is adjusted by changing the position of the guide rods (29), the ball joint of the headrest can be used for fine adjustment.
3. Adjust the width (6) and height (5) of the **armrest** for the healthy arm to a position where the patient sits up straight (shoulders at the same height) and comfortable.

4.2.2 Adjusting the motion element

Standard settings such as for a normal lateral adduction/abduction movement combined with anteversion/retroversion and internal/external rotation.

1 Height adjustment (Fig. A)

The purpose of this adjustment is to set the rotational axis of motor B to the shoulder level. To achieve this, the pivot of motor B must be at the same level as the pivot of the shoulder joint (see illustration in 4.2.3).

This adjustment is performed with the crank handle for height adjustment (8).

2 Shoulder width (Fig. B)

The purpose of this adjustment is to laterally line up the rotational axes of motor A and motor B with the shoulder joint. To achieve this, the pivot of motor A and the pivot of motor B must be set to the same width as the pivot of the shoulder joint (see illustration in 4.2.3).

This adjustment is performed with the crank handle for width adjustment (9).

3 Depth adjustment of the backrest (Fig. C)

The purpose of this adjustment is to set the rotational axis of motor A to the same depth as the shoulder joint. To achieve this, the pivot of motor A must be at the same depth as the pivot of the shoulder joint (see illustration in 4.2.3). Furthermore, this adjustment guarantees a more comfortable position for the patient on the chair during the treatment.

This adjustment is performed with the crank handle for depth adjustment of the backrest (10).

4 Upper arm length (Fig. D)

The device for adjustment of the upper arm length is self-locking. During adjustment, slightly lift motor C and ensure that the insertion tube does not jam.

- Loosen clamping lever (16) and slightly lift motor C during adjustment.
- Perform the desired adjustment.
- Tighten the clamping lever.

5 Elbow angle (Fig. E)

The rotational axis of motor C must be aligned with the shoulder joint such that the rotational axis of motor C goes through the upper arm to the pivot of the shoulder joint in rotation (see illustration in 4.2.3).

In most cases an angle of 90° to 60° is adjusted for elbow flexion.

- Loosen locking screw (18). To facilitate the adjustment, slightly lift motor C.
- Complete the necessary adjustment and tighten the locking screw.

Note!

If the elbow angle is changed to a flexion of less than 90°, the setting for the upper arm length needs to be changed, too.

The maximum adjustment range is from 0° to 90°.

6 Forearm support (Fig. F)

The rotational axis of motor C must be aligned with the shoulder joint in rotation (see illustration in 4.2.3).

In most cases, the forearm guide (red mark) is set to the adjustment point 15.

- To change the adjustment, loosen clamping lever (26) and adjust a length where the rotational axis of motor C goes through the upper arm to the pivot of the shoulder joint in rotation.
- Tighten the clamping lever.

7 Forearm length (Fig. G)

- Loosen clamping lever (23) and pull out the hand rest (25) until the forearm can easily be accommodated between elbow support and hand rest.

- Tighten the clamping lever.

Note!

The hand rest inclination can be adjusted until optimally suited to the arm. For this purpose, it can be set to three different positions.

4.2.3 Checking the setup, final adjustment

1. Check adjustments 1 through 7 to ensure that:
 - a) the pivot of motor A and the pivot of the patient's shoulder joint correspond for anteversion/retroversion
 - b) the pivot of motor B and the pivot of the shoulder correspond for abduction/adduction
 - c) the rotational axis of motor C, the center of the elbow joint, and the center of the shoulder joint form a straight line.
2. Check that all locking screws are tightened and all clamping levers properly closed.

Anatomically correct adjustment

4.3 Switching the device off

1. To completely switch off the device, turn off the **power switch** (43).
2. To completely disconnect the device from the power supply (all poles), remove the power cord from the wall outlet.

Note!

- Set up the device in a location where the power switch and the power plug are **readily accessible at all times**.
- **Stop system (ending the therapy session):**
The ARTROMOT®-S4 will stop immediately, when any of the three buttons is pressed during the therapy session. As an alternative, you can tap the display during operation. The question "STOP?" will appear in the header of the display. Tap the field to end the session. When the START button is pressed afterwards, the CPM device automatically moves in the opposite direction of the last movement.

5. Setting the treatment values

Warning!

Patient hazard –

Before therapy, a test run consisting of several exercise cycles must be completed without the patient. Then, repeat the test run with the patient and check that the movements do not cause any pain.

Note: See also 2.2 and 2.3.

Note!

Programming is only possible when the patient's chip card has been inserted.

For information about therapy parameters and programming of standard and special functions, please refer to sections 5.1 through 5.3.

See section 5.4 for programming examples.

The ranges of motion can also be directly programmed at the respective motors, see section 2.4 for details.

5.1 General information on programming the ARTROMOT®-S4

The ARTROMOT®-S4 comes with a manual control pendant with touch screen for control and input of the treatment values.

Three buttons on the control pendant provide additional functions:

- START button to start the session
- STOP button to stop the session
- MENU button to call up the programming menu

Different motion patterns such as PNF, Wave, or Isolated motion patterns can be selected during operation.

All treatment parameters can be preconfigured for each patient.

1. You activate the programming mode by briefly pressing the **MENU** button on the programming unit.

As an alternative, you can directly select the motion patterns on the display by tapping the symbol.

2. The treatment parameters and functions are allocated to different programming levels. The status bar of the associated symbol indicates the status (activated/deactivated/current setting) of each parameter.

Simply **tap** the appropriate **parameter symbol** to select a parameter and view the corresponding functions.

The currently displayed menu is indicated in the header of the display on the right side.

3. Simply **tap** the appropriate **function symbol** to select a function. The selected function will be enlarged on the display.

4. Some of the functions can only be activated and deactivated, while values can be entered for other functions.

Functions that can be activated: Some of the functions can only be activated and deactivated; no other adjustments are required (e.g., warm-up). They are activated and deactivated by tapping. When activated, the status field is green; when deactivated, it is gray.

Functions that can be activated and have additional adjustment options:

Toggle the on/off button. When the field is green, the corresponding function is activated; when it is black, the function is deactivated. If the function is activated, the adjustment options are displayed in color and can be selected. If the function is deactivated, the adjustment options are gray and cannot be selected.

Adjustable functions: Select the function and, using the arrow left/arrow right buttons, change the value displayed between the arrow buttons.

5. To return to the last menu item/parameter, tap its name in the header of the display on the left.
6. Having programmed all parameters, press the **STOP** button to save the values.
7. Next, press the **START** button: The ARTROMOT®-S4 checks the set values, positions the motion element midway between the set internal and external rotation values, midway between the set anteversion and retroversion values as well as to the set maximum adduction value and stops.
8. Press the **START** button again to start therapy.

PNF motion pattern:

Then, the device will first move to the maximum internal rotation position and simultaneously to the maximum anteversion position. The device will then simultaneously move to the maximum external rotation position, the maximum abduction position and the maximum retroversion position, and subsequently to the maximum adduction position, the maximum internal rotation position and the maximum anteversion position. After reaching this position, the cycle restarts: maximum abduction, maximum external rotation and maximum retroversion at the same time.

Wave motion pattern:

Next, the carriage will first move the maximum anteversion position. All 3 motors are synchronized in this position and the resulting pattern resembles a sideways figure eight.

Patients with programmed chip card

- First, complete the mechanical adjustments.
- Then, insert the chip card (the patient is not yet positioned on the CPM device).
- Press the **START** button: the device will move to the starting position as specified by the parameters stored on the chip card and stops.
- Position the patient on the CPM device and press the **START** button to initiate therapy.

Note!

- Refer to sections 5.3 and 5.5 for a description of the parameters.
- To prevent accidental changes of the parameter settings, you can **lock the buttons**. To do so, stop the CPM device (no menu is open) and simultaneously press the **MENU and STOP buttons for approx. 3 seconds**: the lockout menu displays in which individual functions can be locked (for details, see: section 5.3.2). Press both buttons again for approx. 3 seconds to unlock and disable the lockout function in the lockout menu.
- Selecting the "New Patient" function will automatically delete the data on the patient chip card.

Note!

- When you have finished programming the device and press the **STOP** button, the settings will automatically also be saved to the patient chip card.
- **Stop system (ending the therapy session):**
The ARTROMOT®-S4 will stop immediately, when any of the three buttons is pressed during the therapy session. As an alternative, you can tap the display during operation. The question "STOP?" will appear in the header of the display. Tap the field to end the session. When the START button is pressed afterwards, the CPM device automatically moves in the opposite direction of the last movement.

Note!

- **While you adjust the values**, the CPM device moves to the set range. This allows you to easily and quickly determine the ROM where the patient does not experience pain.
- The angle entered last for the respective directions of motion will be saved in each case.
- The individual settings of the ARTROMOT®-S4 are always programmed in context-related function levels.
- You can navigate between levels and functions by tapping the appropriate symbol on the control pendant.

5.2 Overview of the parameters available for the different motion patterns

Motion patterns / adjustable parameters	PNF	Wave	Isolation	Oscillation	Stretching
Anteversion (forward)					
Retroversion (backward)					
Adduction (downward)					
Abduction (upward)					
Int. Rotation					
Ext. Rotation					
Timer (therapy timer)					
Speed					
Warm-up					
Pause up					
Pause down					
Cycles					
Standby					
Wave counter-clock-wise (to the left)					
Wave clockwise (to the right)					

5.2.1 Description of the selectable motion patterns

The motion patterns can be selected in the "Expert" menu  under "Therapy" .

PNF

The PNF motion pattern is based on the normal pattern of motion of the shoulder joint according to PNF (proprioceptive neuromuscular facilitation). This means that movements toward the maximum abduction, external rotation and retroversion positions are performed at the same time. When these positions have been reached, the simultaneous movements to the maximum adduction, internal rotation and anteversion positions are performed. After reaching this position, the cycle restarts: simultaneous maximum abduction, maximum external rotation and maximum retroversion.

The ranges of motion are programmed as described in section 5.1.

Wave

Up to three motors can be used for the Wave motion pattern. Motor A, performing the anteversion / retroversion movement, is required in any case. All motors are synchronized too in this mode, and the resulting motion pattern resembles a sideways figure eight.

This means that not the maximum positions are entered at the same time but that each motor, after reaching the set maximum position, individually changes its direction.

The ranges of motion are programmed as described in section 5.1. Additionally, the direction of the wave

can be selected: counter-clockwise  or clockwise .

Isolation program

For this special function, all motors are switched on but will never move simultaneously.

These are the steps of the special function:

- At first, motor A moves to the programmed maximum values for anteverision and retroversion. This cycle is performed 5 to 20 times as configured. Then, the motor stops. Meanwhile, motors B and C are inactive.
- Then, motor B moves to the programmed maximum values for adduction and abduction for the set number of cycles, after which it stops. Meanwhile, motors A and C are inactive.
- Then, motor C moves to the programmed maximum values for internal and external rotation for the set number of cycles and stops also. Finally, the full cycle will start over with motor A. Motors A and B are inactive while motor C is operating.
- The stop position at the end of the last cycles (for adduction/abduction as well as for anteverision and retroversion and internal/external rotation) can be selected in steps of 25% between 0% and 100 % of the programmed maximum range of motion.
- Steps 1 and 2 can be repeated as often as needed. You can stop treatment with the **STOP** button; after the programmed therapy duration, it will stop automatically.

Default setting: disabled

To program the isolation pattern, proceed as follows:

1. Program the range of motion as described in section 5.1.
2. Select the motion pattern in the "Expert" menu under "Therapy" by tapping the appropriate symbol. The sub-menu for activation / deactivation of the function will be displayed. The possible settings are displayed in gray.
3. Activate the function. The possible subsequent functions are displayed in color.
4. Set the desired number of cycles (5 to 20).
5. Set the respective standby positions for the motors (parked position while another motor is operating)
Default setting: 25 %
6. If appropriate, set pauses, speed, or timer.
7. Press **START** to save the settings and start the treatment.

Oscillation

The Oscillation pattern permits a more efficient therapy in the last 10° before the maximum abduction, retroversion, and external rotation positions are reached.

For this exercise, the device begins at the starting position (maximum adduction, midway between internal and external rotation, and midway between anteversion and retroversion).

First, the device moves simultaneously to the programmed maximum internal rotation position and the programmed maximum anteversion position. Then, it will move to the programmed maximum abduction position, at the same time as the maximum external rotation and maximum retroversion position.

When the programmed abduction/external rotation position has been reached, the device reverses 10° toward adduction as well as internal rotation and anteversion, and then moves back again to the maximum abduction/external rotation/retroversion position. The movement through the final 10° is performed for the set number of cycles (5 to 20 times) at a slow speed.

At the end of the cycle, the device will again move to the programmed maximum adduction position, at the same time as the programmed maximum internal rotation and anteversion positions, and then starts another cycle for the set number of cycles through the last 10° of abduction.

This sequence can be repeated as often as needed. You can stop treatment with the **STOP** button; after the programmed therapy duration, it will stop automatically.

Default setting: disabled

To program the oscillation pattern, proceed as follows:

1. Program the range of motion as described in section 5.1.
2. Select the motion pattern in the "Expert" menu under "Therapy" by tapping the appropriate symbol. The sub-menu for activation / deactivation of the function will be displayed. The possible settings are displayed in gray.
3. Activate the function. The possible subsequent functions are displayed in color.
4. Set the desired number of cycles (5 to 20).
5. Set the respective standby positions for the motors (parked position while another motor is operating).
Default setting: 25%
6. If appropriate, set pauses, speed, or timer.
7. Press **START** to save the settings and start the treatment.

Stretching

With the stretching pattern, the joints are gently stretched in the set direction.

Options are:

1. Stretching in retroversion

2. Stretching in abduction

3. Stretching in internal rotation

4. Stretching in external rotation

The treatment is only performed in the selected direction of movement. The other motors are automatically deactivated.

Starting in the middle position, the device first moves to the programmed maximum position in the direction opposite to the stretching direction and then to the programmed maximum position in the stretching direction.

Subsequently, the device reverses 5° in the opposite direction and then moves very slowly back again to the set maximum position in the stretching direction. After that, it attempts to attain another -5°, moving even slower than before.

If a high resistance toward the additional 5° is sensed, the reverse-on-load function is automatically activated and the device moves in the opposite direction.

This stretching cycle is repeated for the set number of cycles (5 to 20 times).

Then, the device moves to the programmed maximum position in the direction opposite to the stretching direction before it starts over performing a stretching cycle in the direction of the selected stretching movement.

This sequence can be repeated as often as needed.

You can stop treatment with the **STOP** button; after the programmed therapy duration, it will stop automatically.

Default setting: disabled

To program the stretching pattern, proceed as follows:

1. Program the range of motion as described in section 5.1. Please make sure that the motors for the planes of motion not currently needed are in the position in which the stretching movement in the selected plane of motion (abduction, retroversion, internal rotation or external rotation) will occur.

Important!

Please remember that, when performing the stretching movements, the motors not currently needed are automatically deactivated and remain in the current position.

2. Select the motion pattern in the "Expert" menu under "Therapy" by tapping the appropriate symbol. The sub-menu for selection of the stretching direction will be displayed.
3. Select the desired stretching direction by tapping.
4. Activate the function. The possible subsequent functions are displayed in color.
5. Set the desired number of cycles (5 to 20).
6. If appropriate, set pauses, speed, or timer.
7. Press START to save the settings and start the treatment.

Note!

- If a pause has been programmed, the device will stop for the pause each time the maximum stretching position is reached.
- When the "Stretching" function is activated, the motors not needed will automatically be deactivated. The only movement performed is stretching in the set direction of movement.
- Only one stretching pattern can be activated at a time.

5.2.2 Description of the available parameters of the motion patterns

	Anteversion (moving the arm forward) maximum value: 110°
	Retroversion (moving the arm backward) maximum value: -10°
	Adduction (lowering the arm) maximum value: 0°
	Abduction (raising the arm) maximum value: 160°
	Internal rotation (rotating the arm towards the body) maximum value: -90°
	External rotation (rotating the arm away from the body) maximum value: 90°

Note!

- The programmed values and the values (degrees) measured on the patient may deviate slightly.
- To ensure physiological motion in the PNF motion pattern, the device simultaneously moves to the following positions:
 - Maximum abduction simultaneous with maximum external rotation and retroversion
 - Maximum adduction simultaneous with maximum internal rotation and anteversion
- To ensure that the movements you program are physiological in nature, we recommend to program the device in the following order:

1. Maximum adduction 
2. Maximum internal rotation 
3. Maximum anteversion 
4. External rotation (to relieve the shoulder)  = 0°
5. Maximum retroversion 
6. Maximum abduction 
7. Maximum external rotation 

Note!

- After activation of the **START** button, the device will first enter a neutral position based on the programmed values. In each case, this is the maximum adduction position, the position midway between anteversion and retroversion and midway between internal rotation and external rotation. Then, the device will stop, allowing the patient to be positioned. When the **START button** is pressed again (start of therapy), the device first moves to the maximum position of internal rotation at the same time as to the maximum anteversion position in the PNF motion pattern. Motor B (abduction/adduction) is inactive during this procedure. Once the maximum internal rotation and anteversion positions have been reached, all three motors (A, B and C) will simultaneously move to the maximum abduction, retroversion and external rotation positions. After reaching this position, the cycle restarts: maximum adduction at the same time as maximum anteversion and maximum internal rotation.
- The elevation movement (flexion) is programmed by means of the abduction/adduction values (for mechanical adjustment instructions, refer to section 5.6 "Examples of use").
- During programming of the ranges of motion (maximum values), the device moves to the set positions and the display indicates the movement in steps of 1 degree. This allows you to determine easily, safely and quickly the ROM where the patient does not experience pain.
- The difference between the set values of one direction of motion must always be 10° at minimum. Therefore, it is not possible to select values where the difference would be less.
- For a uniform, smooth and gentle transition between the two directions of motion, the speed is automatically reduced before reaching the reversal point and, on reaching the reversal point, it is continuously increased until it reaches the set speed.

Speed

The speed can be adjusted in steps of 5 % between 5 % and 100 %.

Default setting: 50 %

Timer (therapy timer)

The therapy time is always indicated at the top right of the display (in the header). When the arrow points upward, it is the elapsed time (continuous operation); when the arrow points downward, it is the remaining therapy time (timer activated).

Default setting: Continuous operation (∞).

In **continuous operation**, you must press the **STOP** button to end the therapy session.

However, you can also select therapy durations of **1 to 59 minutes in steps of 1 minute** and of **1 to 12 hours in steps of 30 minutes**.

When the time has elapsed, the device switches **automatically** off and stops in the neutral/starting position.

Note!

The manufacturer recommends an application period of up to 1 hour per therapy session.

Warm-up protocol

The warm-up protocol allows the patient to slowly become used to the set maximum positions.

The device starts the warm-up protocol at the set maximum adduction position and midway between the maximum anteversion/retroversion positions and internal/external rotation positions. The range of motion increases with each cycle, until all set maximum values are reached after a total of 15 cycles. The warm-up protocol can be selected to precede any motion pattern.

Default setting: disabled

**Pauses:**

Pauses occur at each programmed maximum value. These are the two pause points:

**Pause up:**

Maximum adduction value simultaneous with maximum internal rotation value and maximum anteversion value

**Pause down:**

Maximum abduction value simultaneous with maximum external rotation value and maximum retroversion value
They are adjustable in steps of 1 second between 0 and 30 seconds.

Default setting: no pause

Note!

The pause functions are not available for the Wave motion pattern.

**Standby**

The standby function is only available in the Isolation motion pattern. With the standby function, you can determine the parking position of a motor by specifying the percentage of the respective range of motion in a given plane (motor A, B or C) while movements are performed in the currently active plane.

The stop position at the end of the last cycle (for all three planes of motion) can be selected in steps of 25% between 0% and 100 % of the programmed maximum range of motion.

Options are:

**Standby Ante / Retro:**

While abduction/adduction or internal/external rotation exercises are performed, the stop position for motor A (in the anteversion/retroversion plane of motion) can be specified as a percentage of the set range of motion for motor A.

**Standby Ab-/Adduction:**

While ante-/retroversion or internal/external rotation exercises are performed, the stop position for motor B (in the adduction/abduction plane of motion) can be specified as a percentage of the set range of motion for motor B.

**Standby Internal/External Rotation:**

While ante-/retroversion or adduction/abduction exercises are performed, the stop position for motor C (in the internal/external rotation plane of motion) can be specified as a percentage of the set range of motion for motor C.

Default: 25% each

**Cycles:**

The number of cycles defines the number of repetitions of a function. This parameter is available for motion patterns: isolation, oscillation, stretching. The number of cycles is adjustable from 5 to 20 repetitions.

Default setting: 5

5.3 Other functions/settings

New Patient (Reset)

With this function, the CPM device will move to the home position.

Select the function and activate it by tapping "OK"; then press **START**:

- The device moves to the home position.
- The programmed treatment parameters are deleted.
- All values stored on the chip card are deleted.
- The device stops midway between the angles set for ante-/retroversion, abduction/adduction and internal/external rotation.

The "New Patient" function (home position) selects the following settings:

- Anteversion: 30°
- Retroversion: 20°
- Abduction: 50°
- Adduction: 40°
- Internal rotation: 5°
- External rotation: -5°
- Motion pattern: PNF
- Pauses: 0
- Timer: continuous operation
- Speed: 50 %
- Load reversal: 100%
- Motor A: ON
- Motor B: ON
- Motor C: ON
- Total therapy time: 0
- Special functions: OFF

Conversion (right/left conversion):

With this function, the CPM device can easily be converted for treatment of the right arm or of the left arm.

Select the function and tap "OK".

Then, you will be guided interactively through the conversion procedure by on-screen instructions (for details, also refer to section: 6.4). Each time you confirm completion of a step with "OK", the next step follows automatically.

Displayed view:

With the View setting, you toggle between the displayed views during operation of the device. The symbol is shown in the upper left corner of the status line when the device has stopped; it can only be selected at this point.

Tap the symbol and the view changes automatically.

Available views are: "Expert view" and "Standard view".

Documentation

The Documentation function allows you to reproduce the therapy data saved by the carriage.

The following documented data can be reviewed:

Total patient therapy time

The total patient therapy time is the added sum of operating hours stored on the chip card.

(added sum of all therapy sessions stored on the chip card).

Deleting the stored therapy time on the chip card:

Select the function with the "OK" button: the stored period of time will be displayed. Tap the "Reset" field for 3 seconds: a count-down timer is displayed to reset the time; or select the "New Patient" function.

Documentation of the therapy trend

This is a special function of ARTROMOT®-S4 that allows the entire therapy documentation to be reviewed on the display.

The device run times as well as the range of motion of the sessions are recorded.

Three displays are available, representing the three planes of motion. Each of them shows two curves in a coordinate system (X-axis = ROM, Y-axis = time) and a green area between them.

The first display shows the anteversion (upper curve) and the retroversion trend (lower curve).

The second display shows the abduction (upper curve) and the adduction trend (lower curve).

The third display shows the external rotation (upper curve) and the internal rotation trend (lower curve).

The different displays are selected by tapping the field below.

Deleting the stored documentation of the therapy trend:

The documentation is deleted with the "total patient therapy time" (see above) or by selection of the "New Patient" function.

Settings

With "Settings" you access the setup menu. From this menu, a number of basic settings for the carriage can be performed.

Adjustable parameters:

Adaptive load reversal (reverse-on-load feature for patient safety/spasm protection)

The reverse-on-load feature is a safety function.

The device automatically starts moving in the opposite direction when the patient's resistance (load) in the current direction of motion is too high (e.g., due to a spasm).

The reverse-on-load function is automatically active when the device operates. Therefore, it can be selected in all operating modes.

The motors continuously measure the average energy necessary to carry out a movement (depending on the load of the motion element) and provide this energy.

If a significant rise in energy demand which exceeds the set tolerance range is identified, the reverse-on-load function is activated, the device stops and reverses its direction of movement. Then, it starts over with a regular cycle.

The tolerance range between the average energy necessary to carry out a movement and the threshold for activation of the reverse-on-load function can be adjusted from 5 to 100% in small steps of 5%.

- 5 % = narrow tolerance range
This means: Very sensitive response of the reverse-on-load function; a minimal load increase is enough to activate the function.
- 100 % = wide tolerance range
This means: Less sensitive response of the reverse-on-load function; a substantial load increase is necessary to activate the function.

Default setting: 100 %

Caution!

Personal injury –

The reverse-on-load feature is only a safety measure in the event of cramps, spasms, locked joints and similar conditions. The manufacturer cannot be held liable for misuse of this feature.

Transport setting

With this function, the device will move to a position optimally suited for packing. Select the function and tap "OK". Then, you will be guided interactively through the procedure by on-screen instructions.

(for more details, see also section: 6.3)

Language selection

You can choose the following languages:

- Danish
- German
- English
- French
- Italian
- Dutch
- Portuguese
- Polish
- Russian
- Spanish
- Swedish
- Czech
- Turkish
- Hungarian
- Chinese
- Japanese
- Korean

Time/date

Adjustment of the internal device clock.

Adjustable parameters:

- Time  in: hh/mm
- Date  in dd/mm/yy

Total unit run time

The total unit run time is the added sum of operating hours stored **on the control pendant**.

(added sum of all stored therapy sessions)

Brightness

Adjustment of the display brightness.

Default setting: 100 %

Volume

Adjustment of the volume of audio signals.

0 = no audio signal

Default setting: 100 %

Service menu

For service purposes only, refer to Service Manual.

Lockout function

With the lockout function, access to individual operating modes or functions can be blocked. This reduces the scope of programming options and makes operation easier.

To activate the lockout function, simultaneously press the MENU and STOP buttons for 3 seconds while the device is switched on. The device must be in the stop mode. The display automatically shows the lockout menu.

The individual functions/menus can be locked and unlocked by tapping the slider below the appropriate symbol.

When the field is green, the function/menu is enabled and NOT locked. When the slider is touched, the color changes to gray and the corresponding function/menu is locked.

It is possible to lock more than one function/menu at a time.

A padlock in the status bar indicates the locked status.

Default setting: all functions/menus active (not locked)

The following functions/menus can be locked:

Menu

If deactivated, the main menu level is locked.

The ranges of motion can still be selected by tapping the appropriate field on the display while the device is in stop mode; the displayed view can also be changed. No other modifications are possible.

Angle

If deactivated, angles cannot be set; the same applies if a motor is deactivated.

All other functions/menus are not affected.

Displayed View

If deactivated, the view currently selected for the display is locked.

Expert

If deactivated, the "Expert" menu is locked.

The ranges of motion, timer, speed and displayed view can still be changed.

Motor ON/OFF

If deactivated, activation and deactivation of the motors is locked.

All other functions/menus are not affected.

Motor Controls

If deactivated, adjustment of the ranges of motion by means of the motor controls is locked.

All other functions/menus are not affected.

5.4 Application and programming examples

5.4.1 Isolated adduction/abduction

1. Complete the mechanical settings and adjust the device to the patient as described in section 4.2.
2. Then, program the fixed positions for motors A (anteversion/retroversion) and C (internal/external rotation) and the range of motion for motor B (abduction/adduction) as described in 5.1 in the following order:
 - First, program the fixed position for execution of abduction/adduction in the anteversion/retroversion plane of motion (motor A), and then deactivate the motor.
 - Subsequently, program the fixed position for execution of abduction/adduction in the internal/external rotation plane of motion (motor C), and then deactivate the motor.
 - Finally, program the range of motion for abduction/adduction (motor B).
3. Save the settings with the **STOP button**, set the device to the starting position with the **START button** and initiate therapy by pressing the **START button** again.

Note!

- For an exclusive adduction/abduction protocol, motors C (rotation) and A (anteversion/retroversion) must be switched off.
- Having programmed the ranges of motion, you can program additional options, such as pause or speed.
- The stop positions of the different motors can also be adjusted when they are deactivated.

5.4.2 Isolated internal/external rotation

1. Complete the mechanical settings and adjust the device to the patient as described in section 4.2.
2. Then, program the fixed positions for motors A (anteversion/retroversion) and B (abduction/adduction) and the range of motion for motor C (internal/external rotation) as described in 5.1 in the following order:
 - First, program the fixed position for execution of internal/external rotation in the anteversion/retroversion plane of motion (motor A), and then deactivate the motor.
 - Then, program the fixed position for execution of internal/external rotation in the abduction/adduction plane of motion (motor B), and then deactivate the motor.
 - Finally, program the range of motion for internal/external rotation (motor C).
3. Save the settings with the **STOP button**, set the device to the starting position with the **START button**, and initiate therapy by pressing the **START button** again.

Note!

- For an exclusive internal/external rotation protocol, motors B (abduction/adduction) and A (anteversion/retroversion) must be switched off.
- Having programmed the ranges of motion, you can program additional options, such as pause or speed.
- The fixed stop positions of the different motors can also be adjusted when they are deactivated.

5.4.3 Isolated anteversion/ retroversion

1. Complete the mechanical settings and adjust the device to the patient as described in section 4.2.
2. Then, program the fixed positions for motors B (abduction/adduction) and C (internal/external rotation) and the range of motion for motor A (anteversion/retroversion) as described in 5.1 in the following order:
 - First, program the fixed position for execution of anteversion/retroversion in the abduction/adduction plane of motion (motor B), and then deactivate the motor.
 - Subsequently, program the fixed position for execution of anteversion/retroversion in the internal/external rotation plane of motion (motor C), and then deactivate the motor.
 - Finally, program the range of motion for anteversion/retroversion (motor A).
3. Save the settings with the **STOP button**, set the device to the starting position with the **START button**, and initiate therapy by pressing the **START button** again.

Note!

- For an exclusive anteversion/retroversion protocol, motors B (abduction/adduction) and C (internal/external rotation) must be switched off.
- Having programmed the ranges of motion, you can program additional options, such as pause or speed.
- The fixed stop positions of the different motors can also be adjusted when they are deactivated.

5.4.4 Isolated elevation (flexion)

5.4.4.1 For extension/flexion (elevation) with flexion of elbow

1. Complete the mechanical settings and adjust the device to the patient as described in section 4.2.
2. Program the desired position for motor A (anteversion/retroversion) where the isolated elevation (flexion) movement is to be carried out (usually 90°), and then deactivate motor A.
3. Program the desired position for motor C (internal/external rotation) to be maintained during the isolated elevation (flexion) movement (usually 90°), and then deactivate motor C.
4. Program the desired range of motion for the isolated elevation (flexion) movement at motor B (abduction/adduction).
5. Save the settings with the **STOP button**, set the device to the starting position with the **START button**, and initiate therapy by pressing the **START button** again.

Note!

- For an exclusive flexion (elevation) protocol, motors A (anteversion/retroversion) and C (internal/external rotation) must be switched off.
- Having programmed the ranges of motion, you can program additional options, such as pause or speed.
- The fixed stop positions of the different motors can also be adjusted when they are deactivated.

5.4.4.2 For elevation with extended elbow

1. Complete the mechanical settings and adjust the device to the patient as described in section 4.2 – settings **1** to **3**.
2. Perform adjustment **4** and set the upper arm length to the minimum.
3. Perform adjustment **5** and adjust an elbow angle of 0°.
4. Turn the elbow support (19) by 90°, normally outward, toward motor C.
5. Perform adjustment **6** and adjust the forearm guide as required (arm length), but usually to the minimum position 0, i.e., the elbow should rest comfortably in the forearm support (19).
6. Perform adjustment **7** and adjust the forearm length as required (arm length), i.e., the hand should rest comfortably in the hand rest (25).
7. Program the desired position for motor A (anteversion/retroversion) where the isolated elevation (flexion) movement is to be carried out (usually 90°), and then deactivate motor A.
8. Program the desired position for motor C (internal/external rotation) to be maintained during the isolated elevation (flexion) movement (usually 0°), and then deactivate motor C.
9. Program the desired range of motion for the isolated elevation (flexion) movement at motor B (abduction/adduction).
10. Save the settings with the **STOP button**, set the device to the starting position with the **START button**, and initiate therapy by pressing the **START button** again.

Note!

- For an exclusive flexion (elevation) protocol, motors A (anteversion/retroversion) and C (internal/external rotation) must be switched off.
- Having programmed the ranges of motion, you can program additional options, such as pause or speed.
- The fixed stop positions of the different motors can also be adjusted when they are deactivated.

6. Care, Maintenance, Transport, Conversion

6.1 Care / re-use

The ARTROMOT®-S4 is suitable for re-use if the following points are observed.

Warning!

Shock hazard –
Always remove the power cord from the wall outlet before cleaning.

Shock hazard, equipment damage –
Liquids must not enter the device or the programming unit.

- The ARTROMOT®-S4 must be disinfected by **wiping down with a disinfectant**. Thus, it complies with the special hygiene standards for medical technical equipment.
- **The enclosure** and removable **arm support assemblies** can be cleaned with commonly used **disinfectants** and **mild household detergents**.
- Use a **damp cloth** to wipe the CPM device down.

Warning!


Patient hazard –
Patient contamination

Before using the device on another patient, be sure to clean and disinfect it according to the instructions given here.

Note!

The manufacturer recommends using only special disinfectants approved for medical devices and with the characteristics specified below under "Caution".

Disinfection with the following disinfectants and in compliance with the instructions for use has been evaluated and approved by the manufacturer in the risk assessment:

- DESCOSEPT AF lemon
(Part no: 00-311L-xxx)
BAuA Reg.Nr. N-55153, -0482
Manufacturer: Dr. Schumacher GmbH.
Contact time: 2 minutes minimum

Caution!

Equipment damage –

- The plastic material used is not resistant to mineral acids, formic acid, phenols, cresols, oxidants and strong organic or inorganic acids with a pH value below 4.
- Use only clear disinfectants to prevent discoloration of the material.
- Do not expose the CPM device to strong ultraviolet radiation (sunlight) or fire.

6.2 Maintenance (fuse replacement)

Check before each use

Visually inspect the device for signs of mechanical damage before each use.

If you detect damage or malfunctions that may impair the safety of the patient or of the operator, have the device repaired before using it.

Technical inspections

Due to their design, the bearings and joints of the ARTROMOT®-S4 do not require maintenance and all materials are protected from corrosion, BUT: For safety, the devices require regular maintenance. To maintain the functional and operational safety, check all components for damage and loose connections at least once a year.

These checks should be performed by persons with adequate training and experience. Damaged and worn parts must immediately be replaced with original spare parts by authorized staff.

Warning!

Patient hazard –

Refer repair and maintenance to authorized persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.

DJO will make all documents required for servicing, such as circuit diagrams, parts lists, descriptions or calibration instructions, available to authorized experts.

In Germany, the inspections can be carried out by DJO Technical Service within the framework of a service agreement.

Other than that, the manufacturer does not require any other regular maintenance.

Note!

Regarding further technical or other inspections or inspection intervals, observe applicable local requirements, such as IEC 62353, DGUV3 or similar regulations to be observed by users of medical equipment or electric devices.

Fuse replacement

⚠ Warning!

Patient hazard, equipment malfunction and damage –

The replacement of fuses must be referred to specialists as defined in IEC 60364 or other applicable standards (e.g. biomedical technicians, electricians, electronics installers).

Use only fuses of the type T2A H250Vac.

- Before replacing fuses, turn off the power switch (43) of the ARTROMOT®-S4 and disconnect the device from the power line.
- Use an appropriate tool to release and remove the fuse holder (44) situated between the power switch (43) and the power connector (42) (Fig. 1).
- Replace the fuses and reinsert the fuse holder (Fig. 2). Check that the fuse holder locks properly into place.

Fig. 1

Fig. 2

6.3 Transport

The following adjustments need to be carried out to prepare the ARTROMOT®-S4 for transport:

6.3.1 Within buildings – assembled

1. Set motor A to 30° in anteversion (similar to "New Patient" setting).
2. Turn off the ARTROMOT®-S4 power switch (43).
3. Disconnect the power cord (42) from the wall outlet and from the device, roll it up and place it on the seat.
4. Using the height adjustment unit (8), set the height of the motion element to 12 at minimum.
5. Open the clamping levers for adjustment of the upper arm length (16) and forearm length (23) and adjust the minimum length. Close the clamping levers again.
6. Turn the headrest (28) by 90° so that it points downward. Open the thumb-screw that secures the motion element on the system rack (33), turn the motion element by about 45° toward the inside, then tighten the thumbscrew to fix the motion element in this position.
7. Set the width of the armrest for the healthy arm (5) to the minimum position.
8. Rotate the armrest toward the inside over the seat of the chair as far as possible.
9. Holding the system rack in the upper section, tip the chair towards yourself until the wheels touch the floor and get a grip.
10. Push the CPM device forward and transfer it to its new location.
11. Put it into service again by performing steps 2 to 8 in reverse order (also refer to: 4.2).

6.3.2 Packed – disassembled

6.3.2.1 Disassembly

1. Activate the "transport setting" function →  in the menu (also refer to 5.3). The display will show the preparations for transport in an interactive manner. Follow the instructions and set the CPM device to the transport position. Each time you confirm completion of a step with "OK", the next step is automatically displayed.

These are the steps:

- a. Open the clamping lever for adjusting the length of the **forearm support** (26) and **remove** it. Leave the clamping lever open.
 - b. Perform adjustment **5** and adjust an elbow angle of **90°**.
 - c. Perform adjustment **4** and set the upper arm length to **4**.
 - d. You will be prompted to start the function. Press **START**: the motors of the ARTROMOT®-S4 move automatically to the transport position.
2. Turn off the ARTROMOT®-S4 power switch (43).
 3. Remove the power cord (42), disconnect the motion element (40/41) and the programming unit (40/41).
 4. Open the clamping lever for adjusting the height of the armrest for the healthy arm (5) and remove the part.
 5. Open the locking screw for adjusting the width of the armrest for the healthy arm (6) and remove.
 6. Open the thumbscrew (33) that secures the motion element on the system rack and remove.
 7. Unwind the crank handle for adjusting the depth of the backrest (19) and remove it from the system rack: first, it separates from the backrest (36 – clip closure), then you can remove it from the rack. Fold the backrest forward.
 8. Open the thumbscrew that secures the system rack (45) on the chair and remove the system rack.

As of Serial number 2,000: Disassembling the legs of the device chair

9. Remove the two locking pins.


10. Pull the chair legs out to the side.

11. Turn the seat around once.

12. Fasten the red strap around the seat and the folded backrest so that the backrest is secured to the seat, and put the chair legs back on.

13. Reinsert the two locking pins.

Note!

- Reverse the above steps to reassemble the device after transport.
- When the CPM device is reassembled after the transport (only step 1 is left) and switched on, the display prompts you to indicate whether the motion element is mounted on the right or on the left.
 1. Indicate the appropriate side: the motion element will automatically move to the position in which the forearm support can be mounted.
 2. Install the forearm support.
 3. Press the START button: the device automatically moves to the starting position (similar to "New Patient" function).

6.3.2.2 Packaging

The ARTROMOT®-S4 is transported in 3 shipping boxes.

Box 1

Content:

- Chair
- Armrest for healthy arm
- Power cord

Packaging:

1. Put the power cord and the components of the armrest for the healthy arm in the cut-outs provided in the box.
2. Place the ARTROMOT®-S4 – legs first – into the box. The markings on the bottom of the box indicate the correct position.
3. Then close the box.

Box 2

Content:

- System rack
- Headrest
- Crank handle for adjusting the depth of the backrest

Packaging:

1. Push the system rack into the box.
2. Place the headrest and the crank handle in the cut-outs provided for these parts.
3. Then close the box.

Box 3

Content:

- Motion element
- Programming unit
- Forearm support
- Operating instructions/chip cards

Packaging:

1. Place the motion element, the programming unit, and the operating instructions in the cut-outs in the foam pad provided for them.
2. Then close the box.

Only use the original shipping box to transport the device. Ormed GmbH (DJO) cannot be held liable for transport damage if the device was not shipped in its original shipping box.

Warning!

Shock hazard –

Allow the ARTROMOT®-S4 to reach room temperature **before use**. If the device was transported at **temperatures below 0°C (32°F)**, allow it to stand at room temperature for about 2 hours, until all condensation has disappeared. The same applies if the devices was stored at temperatures in the maximum storage temperature range.

6.4 Conversion

The ARTROMOT®-S4 can be used for the left and for the right shoulder joint. However, the device must be converted first. This is done very quickly.

Select the "Ri/LE Conversion" function  in the menu (also refer to 5.3). The display will show the conversion steps in an interactive manner. Follow the instructions to convert the CPM device. Each time you confirm completion of a step with "OK", the next step is automatically displayed.

These are the steps:

1. Open the clamping lever for adjusting the height of the armrest for the healthy arm (5) and remove.

Open the locking screw for adjusting the width of the armrest for the healthy arm (6) and remove.

Mount the armrest for the healthy arm on the opposite side by carrying out the two preceding steps in reverse order.

2. Open the clamping lever for adjusting the length of the **forearm support** (26) and **remove**.
3. Perform adjustment  and adjust an elbow angle of **0°**.
4. Open the thumbscrew (33) that secures the motion element on the system rack and remove. Then mount and secure it on the opposite side.
5. You will be prompted to press START. Press the START button: the motors will automatically move to the appropriate position on the opposite side.
6. Install the forearm support and close clamping lever (26).
7. Perform adjustment  and adjust an elbow angle of **90°**.
8. Press START: the device will move to the starting position stored on the chip card.

7. Environmental Protection Statement

The product described in this operation manual must not be disposed of with unsorted household or municipal waste. It requires separate disposal. Please contact DJO for information about the possible recycling of the product.

8. Specifications

Model:	ARTROMOT®-S4
Version:	COMFORT
Part No.	80.00.084
Input ratings:	100–240 V _{AC} / 50–60 Hz tolerance -15 % to +10 %
Current consumption:	100 V 240 V
Standby (ON):	5.5 VA 5.5 VA
Operation (maximum):	80 VA 80 VA
Fuses:	2 × T2A H250Vac breaking capacity 1500 A according to IEC 60127-2/3, UL 248-14, CSA C22.2 no. 248.14
Battery (control pendant):	type CR2032, 3 V, 230 mAh, cannot be replaced by the user!
Battery life (control pendant):	approx. 5 years
Protection class:	II
IP class:	IP21
Applied part:	type B
Max. load on seat:	175 kg
Max. load on motion element:	13 kg
Physical	
Adjustment ranges (min./max.):	The adjustment is done continuously, the display indicates the movement in steps of 1 degree.
Adduction/Abduction	0° - 30° - 160°
Internal rotation/External rotation	90° - 0° - 90°
Extension/Flexion (Elevation)	0° - 30° - 160°
Retroversion/Anteversion	110° - 0° - 10°

Accuracy of measured values	
Goniometer in the measuring range:	Retroversion – Anteversion –10° – 110° Adduction – Abduction 30° – 160° Rotation –90° – 90°
Accuracy:	±2°
Weight:	SN < 2.000: 28 kg SN ≥ 2.000: 33 kg
Materials used:	PA, POM, PUR, FR4, aluminum, stainless steel
MDD:	class IIa
Standards compliance:	93/42/EEC IEC 60601-1 IEC 60601-1-6 IEC 60601-1-11 IEC 62304 IEC 62366 (S4 < 2.000) IEC 62366-1 (S4 ≥ 2.000) IEC 80601-2-78 (SN ≥ 2.000) ANSI/AAMI ES60601-1 ANSI/AAMI HA60601-1-11 CAN CSA 22.2 No. 60601-1:14
EMC (electromagnetic compatibility)	IEC 60601-1-2
Manufactured in compliance with:	EN ISO 13485
Ambient conditions (storage, transport)	
Ambient temperature range and respective relative humidity:	-25 °C to +5 °C and +5 °C to 35 °C at a relative humidity of up to 90 %, no condensation, >35 °C to 70 °C at a water vapor pressure of up to 50 hPa
Atmospheric pressure:	500 hPa to 1060 hPa
Ambient conditions (operation)	
Temperature:	+5 °C to +40 °C
Relative humidity:	15 % to 90 %, no condensation, but without required partial pressure of water vapor greater than 50 hPa
Atmospheric pressure:	700 hPa to 1060 hPa
Warranty:	2 years (mechanical parts) 2 years (electronics)

Subject to change without notice.

9. Contact

Do not hesitate to contact us with any questions about the product or if you need help with the setup, adjustments or use of the product, about service and maintenance, or if you would like to report unexpected behavior or incidents. We will be happy to help.

ARTROMOT® International:

Please contact your local dealer, the DJO Headquarters in the USA, the DJO International Headquarters in the United Kingdom, or DJO in Germany.

DJO International Headquarters:

DJO UK Ltd.
1a Guildford Business Park
Guildford
Surrey, GU2 8XG
United Kingdom
Telephone: +44 (0)1483 459659
Facsimile: +44 (0)1483 459470
E-Mail: info@DJOglobal.eu
Web: www.DJOglobal.eu

DJO Headquarters

DJO, LLC
2900 Lake Vista Drive
Lewisville, TX 75067 - USA
Telephone: +1 800 321 9549
Facsimile: +1 800 936 6569
E-Mail: chattproductsupport@djoglobal.com
Web: www.DJOglobal.com

Manufacturer / Headquarters Germany:

ORMED GmbH
A DJO Company

Serial number < 2000:
Merzhauser Strasse 112
79100 Freiburg – Germany

Serial number ≥ 2000:
Bötzingen Straße 90
79111 Freiburg – Germany

Telephone: +49 (0) 761 4566 01
Facsimile: +49 (0) 761 456655-01
E-Mail: medizintechnik@DJOglobal.com
Web: www.DJOglobal.de

Branch office Austria:

ORTHOMED Medizintechnik GmbH
A DJO Company
Annagasse 5/1/Top 4
1010 Vienna – Austria
Telephone: 0043 1 5320834
Facsimile: 0043 1 532083431
E-Mail: orthomed@orthomed.cc
Web: www.orthomed.cc

10. Technical service

10.1 Technical hotline

Do you have any technical questions?
Do you need technical service?
Telephone: +49-180-5-1 ormed de
+49-180-5-1 67 63 33
Facsimile: +49-180-5-3 ormed de
+49-180-5-3 67 63 33

10.2 Shipment

To prevent damage during transport, only use the original shipping box. These boxes can be obtained from DJO.

Before packing the CPM device, set it to the transport position (see section 6.3.2).

10.3 Spare parts

Refer to the Service Manual for the most recent list of spare parts.

When ordering spare parts, always specify:

- Item
- Description
- Part number
- Quantity
- Serial number of the device

Note!

Refer repairs to authorized, specially trained staff.

DJO offers service training for your staff.

Note!

Country-specific power cords are available from DJO. Please contact DJO or your DJO distributor.

Item	Description	Part No.	Qty
1.	Patient chip card (protocol)	0.0050.330	<input type="checkbox"/>
2.	Marker pen for patient chip card	0.0031.006	<input type="checkbox"/>
3.	Power cord EU version H05VV-F3G 3 × 1 mm ² , length 4 m	0.0034.118	<input type="checkbox"/>

11. Declaration of conformity

In compliance with the provisions of the Council Directive 93/42/EEC, Annex II, of 14 June 1993 / Amendment of 5 September 2007, concerning medical devices, the company

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declares under its sole responsibility that the products of the product line

ARTROMOT® (see Annex)

**fulfill the essential requirements of Annex I
of the Council Directive 93/42/EEC.**

**With reference to Rule 9 of the Directive 93/42/EEC,
the product is a device of risk class IIa.**

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Freiburg, 15 April 2021

**- QA Management Representative ARTROMOT -
-Bernhard Krohne-**

This certificate is valid until expiry of the certificate referred to.
(The certificate for the year of manufacture can be downloaded from:
<https://www.djoglobal.de/unternehmen/qualitaet-und-verantwortung/>)

Annex:

ARTROMOT®-S3
ARTROMOT®-S4
ARTROMOT® ACTIVE-K
ARTROMOT®-K1
ARTROMOT®-SP3
ARTROMOT®-E2

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Menu de configuration Expert

The screenshot shows a menu interface with a blue header containing '< Traitement' and 'Menu'. Below the header are six settings cards arranged in a 2x3 grid. The top row includes 'Anté/Rétro' (50..6°), 'Ad/Abduction' (35..124°), and 'Rotation' (-20..25°). The bottom row includes 'Vitesse' (100%), 'Minuterie' (∞), and 'Expert' (with a doctor icon). A dashed box highlights the 'Ad/Abduction' card, and a toggle switch labeled 'B' is shown below it.

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The screenshot shows a programming screen for 'Ad/Abduction'. The header has '< Menu' and 'Ad/Abduction'. The main display shows a large '71°' for the current motor angle. Below it is a toggle switch labeled 'B' which is turned on. To the left and right are two 'Adduction' cards with icons of a person moving their arm, showing '35°' and '124°' respectively. At the bottom is a horizontal slider with a blue bar and a white knob, with tick marks at 10°, 50°, 90°, 130°, and 170°. Green arrows point left and right from the slider.

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12. IEC 60601-1-2 and IEC 60601-1-11

The product designation **ARTROMOT®-S4** used in the text below includes all their variants.

The **ARTROMOT®-S4** device is subject to particular precautions regarding electromagnetic compatibility (EMC). The device must be installed and put into service strictly in compliance with the EMC directives set forth in the accompanying documents.

Portable and mobile RF communication systems may affect the **ARTROMOT®-S4**.

The **ARTROMOT®-S4** should not be used adjacent to or stacked with other equipment.

If adjacent or stacked use is necessary, **ARTROMOT®-S4** should be observed to verify normal operation in the configuration in which it will be used.

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If it is necessary to replace assemblies or cables, only the manufacturer's original parts must be used to ensure continued compliance with EMC requirements after repair. This requirement applies to the power supply unit, cables and cable lengths, drive unit consisting of the motor and the control system, the programming unit incl. the coiled cable and the connector.

With regards to the electromagnetic compatibility, there is no maintenance operation required during the expected product life time!

Warning!

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

12.1 Electromagnetic emissions

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

ARTROMOT®-S4 is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-S4** unit is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions to CISPR 11	Group 1	ARTROMOT®-S4 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions to CISPR 11	Class B	ARTROMOT®-S4 is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions to IEC 61000-3-2	complies with class A requirements	
Voltage fluctuations/flicker emissions to IEC 61 000-3-3	complies	

12.2 Electromagnetic immunity

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

ARTROMOT®-S4 is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-S4** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment – guidance
Electrostatic discharge (ESD) to IEC 61 000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/bursts to IEC 61 000-4-4 Surges to IEC 61 000-4-5	± 2 kV 100 kHz repetition frequency ± 1 kV line-to-line voltage	± 2 kV 100 kHz repetition frequency ± 1 kV line-to-line voltage	Mains power should be that of a typical commercial or hospital environment. Mains power should be that of a typical commercial or hospital environment.
Voltage dips To IEC 61000-4-11	0% U_T ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T ; 1 cycle and 70% U_T ; 25/30 cycles Single phase at 0°	0% U_T ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T ; 1 cycle and 70% U_T ; 25/30 cycles Single phase at 0°	Mains power should be that of a typical commercial or hospital environment. If the user of the ARTROMOT®-S4 unit requires continued operation during power mains interruptions, it is recommended that the ARTROMOT®-S4 unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field to IEC 61 000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.
Voltage interruptions to IEC 61000-4-11	0% U_T ; 250/300 cycles	0% U_T ; 250/300 cycles	

Note: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

ARTROMOT®-S4 is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-S4** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment – guidance
Conducted RF to IEC 61000-4-6	6V 0,15–80 MHz 80% AM at 1kHz	6V 0,15–80 MHz 80% AM at 1 kHz	
Radiated RF to IEC 61000-4-3	10 V/m 80 MHz–2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz–2,7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communications equipment to IEC 61000-4-3	27 V/m 380–390 MHz PM 18 Hz 28 V/m 430–470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704–787 MHz PM 217 Hz 28 V/m 800–960 MHz PM 18 Hz 28 V/m 1,7–1,99 GHz PM 217 Hz 28 V/m 2,4–2,57 GHz at PM 217 Hz 9 V/m 5,1–5,8 GHz at PM 217 Hz	27 V/m 380–390 MHz PM 18 Hz 28 V/m 430–470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704–787 MHz PM 217 Hz 28 V/m 800–960 MHz PM 18 Hz 28 V/m 1,7–1,99 GHz PM 217 Hz 28 V/m 2,4–2,57 GHz at PM 217 Hz 9 V/m 5,1–5,8 GHz at PM 217 Hz	

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ARTROMOT®-S4** exceeds the applicable **ARTROMOT®-S4** unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **ARTROMOT®-S4** unit.

12.3 RF (Radio Frequency) Interference

EMC disturbance can affect the performance of the **ARTROMOT®-S4**. To prevent the user from unacceptable risks the device runs performance checks regularly during operation.

This is what happens, if a problem is identified:

- an audio signal sounds
- the unit stops immediately
- the message ERROR and an error code appear on the display

In this case try to quit the error with the STOP key. If the error persists you may attempt to restart the unit by turning it briefly off and on with the power switch.

Warning!

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the **ARTROMOT®-S4**, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Notice!

(for home use)

Portable and mobile RF communication systems may affect the **ARTROMOT®-S4** device. For this reason, make sure that wireless communication equipment, such as wireless home network devices, cell phones, cordless telephones and their base stations and walkie-talkies operate at a minimum distance of 3.3 m from the device. (Calculated on the basis of the maximum power output of a typical cell phone of 2 W).

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