ARTROMOT®-K1 classic

SW-Version: 1.X

- DE · Gebrauchsanleitung
- GB · Operating Instructions
- FR · Mode d'emploi
- ES · Instrucciones de uso
- IT · Istruzioni per l'uso
- NL · Gebruiksaanwijzing

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<https://www.djoglobal.de/unternehmen/qualitaet-und-verantwortung>

Contents

Description of the ARTROMOT®-K1 classic	III
ARTROMOT®-K1 classic setup illustrations	VI
1. How to use the CPM device	27
1.1 Fields of application	27
1.2 Therapy objectives	27
1.3 Indications	27
1.4 Contraindications	27
2. Description of the ARTROMOT®-K1 classic	28
2.1 Description of the device components	29
2.2 Description of the control pendant	30
2.3 Explanation of symbols	32
2.4 Explanation of symbols (connections and nameplate)	33
3. Safety information	35
4. Device setup	38
4.1 Connecting the device, performance check	38
4.2 Adjusting the device to the femur length	38
4.3 Adapting the leg support assemblies/footplate	39
4.4 Connecting the potential equalization conductor (SN > 35,000)	39
5. Setting the treatment values	40
5.1 General information on programming ARTROMOT®-K1 classic	40
5.2 Therapy parameter details for the ARTROMOT®-K1 classic	41
6. Care, maintenance, transport	42
6.1 Care/re-use	42
6.2 Maintenance (fuse replacement)	43
6.3. Transport	44
7. Environmental Protection Statement	45
8. Specifications	45
9. Contact	47
10. Technical service	48
10.1 Technical hotline	48
10.2 Shipment	48
10.3 Spare parts	48
11. CE-Declaration of Conformity	49
12. IEC 60601-1-2 and IEC 60601-1-11	154
12.1 Electromagnetic emissions	154
12.2 Electromagnetic immunity	155
12.3 RF (Radio Frequency) Interference	157

1. How to use the CPM device

1.1 Fields of application

ARTROMOT®-K1 classic is a motor-operated Continuous Passive Motion (CPM) device used to mobilize knee and hip joints.

Suitable for use in hospitals, clinics, general practices and rental services as well as in the patient's home, it is an important supplement to medical and therapeutic treatment.

1.2 Therapy objectives

CPM therapy with the **ARTROMOT®-K1 classic** is mainly used to prevent the negative effects of immobilization, to allow patients to regain painless mobility of joints at an early stage and to promote healing and achieve a positive functional result.

Other objectives of therapy include:

- improvement of joint metabolism
- prevention of joint stiffness
- promotion of the regeneration and healing of cartilage and damaged ligaments
- faster hematoma/fluid resorption
- improved lymph and blood circulation
- thrombosis and embolism prophylaxis

1.3 Indications

The CPM device is indicated in the treatment of most injuries and diseases of the knee and hip joints as well as in postoperative treatment after knee and hip joint surgery. Examples:

- joint distortion and contusion
- arthrotomy and arthroscopy procedures in combination with synovectomy, arthrolysis, or other intra-articular interventions
- mobilization of joints in anesthetized patients
- operative treatment of fractures, pseudoarthrosis, and osteotomy
- cruciate ligament replacement or reconstruction
- endoprosthetic implants

1.4 Contraindications

Do NOT use the **ARTROMOT®-K1 classic** on patients with:

- acute inflammatory processes in the joints, unless on the order of a physician
- spastic paralysis
- unstable osteosynthesis

2. Description of the ARTROMOT®-K1 classic

The motorized CPM device permits extension and flexion of the knee joint in the range of $-10^{\circ}/0^{\circ}/120^{\circ}$, and of the hip joint in the range of $0^{\circ}/7^{\circ}/115^{\circ}$.

It can be used on either side and requires no configuration change.

These are some of the outstanding **ARTROMOT®-K1 classic** features:

- anatomically correct adjustment
- physiological movements
- control pendant for precise adjustment of patient-specific therapy values
- symbols for easy operation of the control pendant

Biocompatibility

Those parts of the **ARTROMOT®-K1 classic** that come into contact with the patient when the device is used as intended, are designed to fulfill the biocompatibility requirements of the applicable standards.

Essential Performance Features

- The programmed angles are accurately maintained with a tolerance of $\pm 2^{\circ}$.
- The programmed speeds are accurately maintained with a tolerance of $\pm 2\%$.
- The selected mode and mechanical settings do not change during operation.

Frequent functions:

- a) Unpacking (device and accessories)
- b) Device connections (mains power connection, connecting the control pendant and the motion element)
- c) Mechanical adjustments of the motion element on the rotational axes
- d) Programming the control pendant (range of motion, speed, functions)
- e) Storage

Warning!

Patient hazard –

- If the Patient operates the **ARTROMOT®-K1 classic**, the Patient needs to know how to **START** or **STOP** the device.
- To prevent the parameters from being changed unintentionally, lock the adjustment controls on the remote control. Contact your local DJO office for assistance.

2.1 Description of the device components


Note: Please fold out page III !

1. Thigh support assembly (type B applied part)
2. Clamping lever to adjust the height of the thigh support assembly
3. Thigh length scale (femur length scale)
4. Thigh length fixation screws (femur length)
5. Knee hinge
6. Calf support assembly (type B applied part)
7. Clamping levers to adjust the height of the calf support assembly
8. Calf length fixation screws (tibia length)
9. Calf length scale (tibia length scale)
10. Strap to secure the foot to the footplate
11. Footplate (type B applied part)
12. Footplate angle fixation screw, also allows swiveling the footplate
13. Clamping lever to adjust footplate rotation and level and to remove the footplate
14. Connection for control pendant
15. Connection for power cord
16. Fuse cap
17. Power switch (ON/OFF)
18. Rating plate
19. Control pendant (type B applied part)
20. Compartment for storage of control pendant
21. Potential equalization terminal¹

¹ Applies only to **ARTROMOT®-K1** devices with serial numbers > 35,000

2.2 Description of the control pendant

2.2.1 Control pendant in normal mode:

2.2.2 Control pendant in speed or pause programming mode

2.3 Explanation of symbols

Symbols on the control pendant


 Extension
(stretching the knee)

 Flexion
(bending the knee)

 Speed

 Pause (extension
and flexion)

Symbols that may appear in the display

 Go to start position
(see Notes in 4.1)

 Controls on
control pendant
locked (see Notes in 5.1)

 Controls on control
pendant unlocked
(see Notes in 5.1)

 Service menu activated,
for service purposes only
(also refer to Service
Manual)

2.4 Explanation of symbols (connections and nameplate)

-  Alternating current
-  For serial numbers < 20,000:
Protection class I equipment. The medical device must be connected to a system with protective earth conductor!
-  For serial numbers > 20,000:
Protection class II equipment. The medical device has a double or reinforced insulation.
-  Type B applied part
-  Power switch OFF
-  Power switch ON
-  The year and month of manufacture are indicated next to this factory symbol.
-  The name next to this factory symbol is the manufacturer.
-  The number next to this symbol is the article reference number.
-  Device complies with Council Directive 93/42/EEC about medical devices, tested and approved by DQS Medizinprodukte GmbH.
-  The number next to this symbol is the serial number.
-  Caution! Observe warnings set forth in instructions for use!
-  "Refer to instruction manual!" ISO 7010-M002
-  "Refer to instruction manual!" ISO 7000 -1641
-  Do not dispose with unsorted municipal waste.

Protect from moisture.

IP21

The IP rating indicates the level of protection and thus the suitability of the device for use under different ambient conditions.

The rating IP21 means:

2 is the level of protection against contact and solid objects.

The digit 2 means:

- Protection against contact: protected against contact with a finger.
- Protection against solid objects: protected against solid objects (diameter of 12.5 mm and greater).

1 is the level of protection against moisture.

The digit 1 means: protection against vertically falling water drops.

Warning! Depending on the device settings, the moving parts of the device present pinch points!

Pay particular attention to small children and babies!

"Observe temperature limits "(storage)! ISO 7000-0632

The overall weight of the device is indicated next to this symbol.

For serial numbers > 35,000: Potential equalization terminal.

3. Safety information

Definitions

Be sure to read the safety information before using the CPM device. The safety information is classified as follows:

Danger!

Indicates an imminent hazard. If not avoided, this hazard will result in death or serious injury.

Warning!

Indicates a hazard. If not avoided, the hazard can result in death or serious injury.

Caution!

Indicates a potential hazard. If not avoided, the hazard may result in minor injury and/or product/property damage.

Safety information

Danger!

Explosion hazard –

The **ARTROMOT®-K1 classic** is not designed for use in areas of rooms used for medical purposes where an explosion hazard may occur. An explosion hazard may result from the use of flammable anesthetics, skin cleansing agents, and disinfectants.

Warning!

Patient hazard –

- Only authorized individuals are allowed to operate the **ARTROMOT®-K1 classic**. Individuals are authorized after receiving training in the operation of the device and reading these instructions for use.
- Before using the CPM device, the operator must ascertain that it is in correct working order and operating condition.
The cables and connectors, in particular, must be checked for signs of damage.
Damaged parts must immediately be replaced with original spare parts.
- **Before therapy, a test run** consisting of several exercise cycles must be completed, first without and then with the patient. Check that all setting screws are tightened.
- Stop therapy immediately when you have doubts about the device settings and/or the therapy protocol.

Warning!

Patient hazard –

- It is important that the patient's **position** is anatomically **correct**. Therefore, carefully verify the following settings/positions:

1. femur length
2. knee joint axis
3. tibia length and leg rotation
4. leg support assemblies

- Movements must **not** cause **any pain or irritation**.
- The patient must be **fully conscious** while being instructed in the use of the CPM device and when using the device.
- The **choice** of the **therapy parameters** to program is restricted to the responsible physician or therapist. The **physician or therapist** must decide on a case-by-case basis whether the device can be used with a particular patient.
- The patient must be familiar with the functions of the **ARTROMOT®-K1 classic control pendant** and the control pendant must be **within easy reach** of the patient, allowing him or her to stop therapy if needed. **Patients unable to operate the control pendant**, e.g. paralytic patients, must never be left unattended during therapy.

The **ARTROMOT®-K1 classic** must only be used with **accessories** approved by DJO.

- Do not allow **parts of the body or any objects** (such as blankets, cushions, or cables) to **get caught in the moving parts of the CPM device**.
- Modifications to the medical device described in this document without the manufacturer's written consent are prohibited.
- The simultaneous treatment of both legs by simultaneous use of two CPM devices is not permitted because the motion elements might interfere with each other.
- Stability of the CPM device must always be ensured while it is in use. The **ARTROMOT®-K1 classic** must only be set up on surfaces that guarantee its stability. Very soft or instable surfaces (such as waterbeds) are NOT suitable.

Warning!

- Special care must be taken when small children and babies are present while the device is in use! Sufficient distance to the device is mandatory for their safety!
- Never leave the device unattended when it is switched on! Switch the device off and disconnect the power line from the wall outlet!
- After use, store the device in a safe place! Ensure device stability also during storage!

Warning!

Shock hazard –

Strictly observe the following warnings. Failure to do so endangers the lives of the patient, the user, and other persons involved.

- Allow the **ARTROMOT®-K1 classic** to reach room temperature **before use**.
- The time required for from the minimum storage temperature between uses until **ARTROMOT®-K1 classic** is ready for its intended use when the ambient temperature is 20°C is 5 hours.
- The time required for **ARTROMOT®-K1 classic** to cool from the maximum storage temperature between uses until Intellect RPW 2 is ready for its Intended use when the ambient temperature is 20 °C is 5 hours
- The **ARTROMOT®-K1 classic** device must only be operated in **dry rooms**.
- When disconnecting the device from the power line, remove the plug from the wall outlet first, before disconnecting the cable from the device.
- When connecting the device to other equipment or when creating a medical system, check that the sum of leakage currents will not cause any hazard. Please contact DJO, if you have questions in this matter.
- Do not use multiple portable socket outlets (MPSO) to connect the device to the power line. The **ARTROMOT®-K1 classic** must be connected to a **properly installed wall outlet with a non-fused earthed wire**.

Before connecting the power cord, it must be completely unrolled and placed such that it will not get caught in the moving parts of the device.

- Before cleaning and service interventions, **disconnect the unit from the power line** by removing the power cord from the wall outlet.
- **Liquids must not be allowed to enter the CPM device** or the control pendant.
If liquids have entered into the devices, the **ARTROMOT®-K1 classic** must be immediately checked by a service technician, before it can be reused.

Warning!

Equipment malfunction –

- Magnetic and electrical fields are capable of interfering with the proper performance of the unit. For this reason make sure that all external devices operated in the vicinity of the device comply with the relevant EMC requirements. X-ray equipment, MRI devices, radio systems, cell phones, etc. are possible sources of interference as they may emit higher levels of electromagnetic radiation.
Keep the device away from such equipment and verify its performance before use.
- Refer **repair and maintenance** to authorized persons.
Persons are authorized after training by a specialist trained and commissioned by the manufacturer.
- **Route all cables below the device frame to either side**, ensuring that they cannot get caught by the moving parts during operation.
- **Check the ARTROMOT®-K1 classic** for damage and loose connections at least **once a year**. Damaged and worn parts must immediately be replaced with original spare parts by authorized staff.

Caution!

Preventing chafing and pressure sores –
When your patient is **adipose**, very **tall** or **very short**, be sure to prevent chafing and pressure sores. Place the leg concerned in a moderate abduction position if deemed appropriate.

Caution!

Equipment damage –

- Check that the voltage and frequency ratings of your local **power line** are those indicated on the nameplate.
- The **leg support element** withstands a maximum continuous load of:
20 kg serial number < 20,000 /
25 kg serial number > 20,000
- Do not allow **any objects** (such as blankets, cushions, or cables) to get caught in the **moving parts** of the unit.
- Do not expose the **ARTROMOT®-K1 classic** to direct sunlight, because some of the components may warm up excessively.

Furthermore, the device must be set up at a safe distance from radiators to prevent excessive temperature rises.

- The presence of children, pets and rodents does not normally impair the functioning of the device. However, avoid contamination of the device by children or animals, from dust and lint, and keep them at a safe distance from the device. The safety statements set forth apply.

Warning!

Patient hazard –

- While the **ARTROMOT®-K1 classic** is in use, do not repair and maintenance the device! There are no parts that can be repaired during use.
- There are no user-serviceable parts inside the device. If a malfunction occurs, discontinue use immediately and contact your local DJO office for assistance and service information.

4. Device setup

Note: For a better understanding of each step, please fold out pages III and VI!

4.1 Connecting the device, performance check

The equipment supplied includes these items:

base unit, control pendant (19), power cord (not shown), instructions for use

1. Connect the control pendant (19) to the provided socket (14) by plugging it in and closing the bayonet lock.
2. Connect the **power cord** to socket (15) of the device and connect the **power plug** to a wall outlet with a non-fused earthed wire (100–240 Volt, 50/60 Hz).
3. Turn the **power switch** (17) on.
4. Follow these steps to set the carriage to the **home position**:

Press the **Extension** control,

and, holding it depressed, rotate until **30°** appears in the display above the control.

In the same manner, select a **Flexion value of 35°**.

Then press the **START/STOP** button. When the carriage has reached this range and does not stop automatically, press the **START/STOP** button again to stop any movement.

Note!

The **ARTROMOT®-K1 classic** classic will stop automatically in the home position range, only if it was positioned outside this range (30° to 35°) at the time of programming (also refer to 5.2)

Performance check

If the control pendant can be operated as described above and the **ARTROMOT®-K1 classic** enters the home position (for home position values, refer to section 4.1), the device has passed the performance check and functions correctly.

The device also runs performance checks regularly during operation. This is what happens, if a problem is identified:

- An audio signal sounds.
- The device switches off immediately.
- The message **ERROR** and an error code (e.g. **ERROR 5**) appear on the display.

In this situation, you may attempt to restart the device by turning it briefly off and on again with the power switch. With the CPM device switched off, check that all plugs are correctly connected. Switch the CPM device on again: if the error message persists, have the device inspected by a Service technician, before using it again.

4.2 Adjusting the device to the femur length

1. Measure the **length of the patient's thigh (femur)** from the greater trochanter to the lateral knee joint cavity (Fig. A).
2. Set the carriage to the home position (see 4.1).
3. Adjust the measured value at the **femur scale** (3) of the carriage (Fig. B):
 - Loosen the two fixation screws (4).
 - Extend the scale (3) to the required length.
 - Tighten the fixation screws (4) to set the scale to the new length.

⚠ Caution!

Equipment damage –

Do not attempt to extend the femur scale beyond the stop.

4.3 Adapting the leg support assemblies/footplate

1. Set the **leg support assemblies and the footplate** (1, 6 and 11) to the expected positions before accommodating the patient.

- Loosen fixation screws (8) to adjust the **footplate** (11) to the length of the patient's lower leg (Fig. C).

Loosen clamping lever (13) and adapt the footplate's rotation and height to the patient (Fig. D).

Loosen fixation screw (12) and adapt the angle to the patient's foot (turn the screw a few revolutions until the footplate can be easily adjusted).

For **short patients** you can reverse the footplate bracket 180° (Fig. H) to adapt the footplate to shorter calves:

- Loosen clamping lever (13) and remove the footplate (11).
- Loosen the fixation screws (12).
- Reverse the bracket 180°.
- Screw the footplate to the bracket and tighten the clamping lever.

Note!

When reversing the footplate, ensure that the pins below the clamping lever engage with the recesses in the bracket.

- To adjust the height of the **support assemblies for calf** (6) and **thigh** (1), loosen clamping levers (2) and (7) (Figs. E/F).
2. Place the patient's leg on the carriage and repeat the steps outlined under 1 above to adjust the CPM device to the patient.

⚠ Caution!

Equipment damage –

Cover the **leg support assemblies** with **disposable tissues**, when using the **ARTROMOT®-K1 classic** immediately after surgery. This helps prevent discoloration.


⚠ Caution!

Patient hazard –

Ensure that the rotational axes of the CPM device and of the knee joint coincide both in the vertical and in the horizontal plane (Fig. G).


Symbol 1:

Measurement of the patient's femur length from the greater trochanter to the knee joint cavity.


Symbol 2:

Set the carriage to the home position (see 4.1) and adjust it to the measured femur length.

Symbol 3:

Adjust height of calf and thigh support assemblies. Adjust the footplate to the height and length of the lower leg.

4.4 Connecting the potential equalization conductor

To prevent interference with monitoring devices when used in treatment rooms, the CPM device must be connected to a potential equalization system according to DIN 42801.

For this purpose, when using devices with serial numbers > 35,000, connect the equalization conductor to the potential equalization terminal (21) of **ARTROMOT®-K1 classic** classic and to the potential equalization connection of the treatment room.

Observe the requirements of IEC 60601-1 for potential equalization connections.

5. Setting the treatment values

Warning!

Patient hazard –

Before **therapy**, a **test run** consisting of several exercise cycles must be completed without the patient. Then, repeat the test run with the patient and check that the movements do not cause any pain.

Note!

See also 2.2 and 2.3 as well as page VI!

5.1 General information on programming ARTROMOT®-K1 classic

1. You activate a function by briefly pressing a control on the control pendant.
2. You select a treatment value by pressing the respective control and turning it in either direction.

You increase a value by turning the control clockwise towards the + (plus) symbol and you decrease a value by turning the control counterclockwise towards the – (minus) symbol. For the first 5° the values in the display change in steps of 1°, then the interval changes to 5° so that you reach the target value faster.

During adjustment of the extension and flexion angles, the values in the display change as you turn the control after pressing it.

When adjusting speed or pause, the information shown in the display changes automatically as you press the control.

The selected parameter (speed or pause) is immediately represented by its symbol in a large format plus the current value. (Also refer to 2.2.2).

The current value can be changed by pressing and turning the control. When you have set the new value, simply release the control and the standard display reappears automatically after approx. 5 seconds (see also 2.2.1).

3. Subsequently press the **START/STOP** button to start therapy.

Note!

Refer to section 5.2 for a description of the parameters.

- To prevent accidental changes of the parameter settings, you can **lock the buttons**. To do so, **simultaneously** press controls “**Extension**” and “**Speed**” for approx. **4 seconds**.

Press both controls again for approx. 4 seconds to unlock.

- **Emergency stop function:** The **ARTROMOT®-K1 classic** will stop immediately, when any of the controls is pressed during therapy. Patient treatment can be resumed by pressing the **START/ STOP** button. The device will automatically change the direction of movement.
- If the carriage is positioned within the programmed range of motion at the time therapy starts, the therapy session will start immediately. If the carriage is positioned outside the programmed range of motion at the time therapy starts, it will first enter the position “extension +10°”. The carriage will stop in this position and you can initiate the therapy session by pressing the **START/STOP** button again.
- The manufacturer recommends an application period of up to 1 hour.
- For a uniform, smooth and gentle transition between the two directions of motion, the speed is automatically reduced before reaching the reversal point and, on reaching the reversal point, it is continuously increased until it reaches the set speed.

5.2 Therapy parameter details for the ARTROMOT®-K1 classic

- You select a function by pressing the corresponding control.
- You change the treatment values by turning the depressed control.
- You initiate the treatment session by pressing the START/STOP button.

■ Extension (stretching)

- Maximum **knee extension**: -10°
- Maximum **hip extension**: 7°

■ Flexion (bending)

- Maximum **knee flexion**: 120°
- Maximum **hip flexion**: 115°

Note!

The programmed value and the value measured at the patient's knee may differ slightly.

■ Speed

The speed can be adjusted in steps of 5% between 5% and 100%.

Default setting: 50 %

■ Pauses

Pauses occur at the selected limits where stretching turns into bending and bending into stretching (selected extension and flexion values). The configured value applies to both extension and flexion pause.

Pauses can be set to any value between 0 and 59 seconds in steps of 1 second, and to values between 1 and 59 minutes in steps of 1 minute.

When selecting the seconds, the value changes in 1-second steps for the first 5 seconds. Subsequently the interval changes to 5-second steps.

The minutes are always adjusted in 1-minute steps.

Default setting: no pause

6. Care, maintenance, transport

The **ARTROMOT®-K1 classic** is suitable for re-use if the following points are observed.

6.1 Care/re-use

Warning!

Shock hazard –

Unplug the device from the power line before cleaning.

Shock hazard, equipment damage –
Liquids must not enter the device or the control pendant.

- The **ARTROMOT®-K1 classic** must be **disinfected by wiping down** with a disinfectant. Thus, it complies with the special hygiene standards for medical technical equipment.
- The **enclosure** and removable **leg support assemblies** can be cleaned with commonly used **disinfectants** and mild **household detergents**.
- Use only a **damp cloth** to wipe the CPM device down.
- Clean the armour of the control pendant on a regular basis. Remove it from the control pendant before cleaning it. Make sure that it is completely dry before re-attaching it to the control pendant.

Warning!

Patient hazard – patient contamination

- Before using the device on another patient, be sure to clean and disinfect it according to the instructions given here.

Note!

The manufacturer recommends using only special disinfectants approved for medical devices and with the characteristics specified below under "Caution":

Disinfection with the following disinfectants and in compliance with the instructions for use has been evaluated and approved by the manufacturer in the risk assessment:

- DESCOSEPT AF lemon
(Part no.: 00-311L-xxx)
BAuA Reg.No. N-55153, C€-0482
Manufacturer: Dr. Schumacher GmbH
Contact time: 2 minutes minimum

Caution!

Equipment damage –

- The plastic material used is not resistant to mineral acids, formic acid, phenols, cresols, oxidants and strong organic or inorganic acids with a pH value below 4.
- Use only clear disinfectants to prevent discoloration of the material.
- Do not expose the CPM device to strong ultraviolet radiation (sunlight) or fire.

6.2 Maintenance (fuse replacement)

Check before each use

Visually inspect the device for signs of mechanical damage before each use.

If you detect damage or malfunctions that may impair the safety of the patient or of the operator, have the device repaired before using it.

Technical inspections

Due to their design, the bearings and joints of the **ARTROMOT®-K1 classic** do not require maintenance and all materials are protected from corrosion, BUT: For safety, the devices require regular maintenance. To maintain the functional and operational safety, check all components for damage and loose connections at least **once a year**.

These checks should be performed by persons with adequate training and experience. Damaged and worn parts must immediately be replaced with original spare parts by authorized staff.

Warning!

Patient hazard –
equipment malfunction and damage

- Refer repair and maintenance to authorized persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.

ORMED GmbH will make all documents required for servicing, such as circuit diagrams, parts lists, descriptions or calibration instructions, available to authorized experts.

The inspections can be carried out by DJO Technical Service within the framework of a service agreement. Contact DJO for details.

Other than that, the manufacturer does not require any other regular maintenance.

Note!

Regarding further technical or other inspections or inspection intervals, observe applicable local requirements, such as IEC 62353, DGUV3 or similar regulations to be observed by users of medical equipment or electric devices.

Fuse replacement

Warning!

Patient hazard –
equipment malfunction and damage




The replacement of fuses must be referred to specialists as defined in DIN VDE 0105 or IEC 60364 or other applicable standards (e.g. biomedical technicians, electricians, electronics installers).

Fuses used must be
S# < 20,000: T1A L250Vac
S# > 20,000: T2A H250Vac.

6.3. Transport

The following adjustments need to be carried out to prepare the **ARTROMOT®-K1 classic** for transport:

1. Adjust the femur length to 49 cm and the tibia length to 42 cm.
2. Set the carriage to 0° by pressing the START/STOP button when the carriage is in the 0° position. This will stop the carriage in that position.
3. Turn off the **ARTROMOT®-K1 classic** power switch.
4. Disconnect the power cord and the control pendant.
5. Only use the original shipping box to transport the device. Ormed GmbH cannot be held liable for transport damage, if the device is not shipped in its original shipping box.
6. Set the ankle joint to the horizontal position.
7. Put the foam parts ① and the power cord ② in the box.
8. Then lower the **ARTROMOT®-K1 classic** into the box.
9. Fold down the block ③ that secures the frame and insert the foam part ④ for the angle joint.
10. Then adjust the carriage to a tibia length of 50 cm so that the ankle joint will be fixed by the foam part ④. Insert the brace ⑤.
11. Put the control pendant (19) in the supplied small box.
12. Fold down the fixation element ⑥ for the control pendant and store the control pendant in the box.

7. Environmental Protection Statement

The product described in these instructions for use must not be disposed of with unsorted household or municipal waste. It requires separate disposal. Please contact DJO for information about the possible recycling of the product.

The expected service life of the device and all supplied parts and accessories is 6 years. The application beyond this period is the responsibility of the user.

8. Specifications

Model:	ARTROMOT®-K1
Device name:	ARTROMOT®-K1 classic, Part no.: 80.00.045
Input ratings:	100–240 Vac / 50–60 Hz tolerance –15 % bis +10 %
Current consumption:	100 V 240 V Standby (ON): 5 VA 5 VA Operation (maximum): 40 VA 40 VA
Fuses:	serial numbers < 20,000: 2 × T1A L250 Vac serial numbers > 20,000: 2 × T2A H250 Vac according to IEC 60127-2/3, UL 248-14, CSA C22.2 no. 248.14
Protection class:	serial numbers < 20,000: I serial numbers > 20,000: II
IP class:	serial numbers < 20,000: IPX0 serial numbers > 20,000: IP21
Applied part:	Type B
Max. load on carriage:	serial numbers < 20,000: 20 kg serial numbers > 20,000: 25 kg
Dimensions	
length:	96 cm
width:	35 cm
height:	min. 23 cm to max. 56 cm

Adjustment ranges (min./max.)	
Femur range:	approx. 31 to 49 cm
Tibia range:	approx. 25 to 57 cm
Accuracy of measured values	
Goniometer in the measuring range:	from -10° to +120°
Accuracy:	+/- 2°
Weight:	11 kg
Materials used:	ABS, POM, PUR, PA, FR 4, aluminum, stainless steel, brass
MDD:	class IIa
Standards compliance:	93/42/EEC IEC 60601-1 IEC 60601-1-6 IEC 60601-1-11 (S# > 20,000) IEC 62366 IEC 62304 EN ISO 14971 ANSI/UL 60601-1 (S# < 20,000) CAN/CSA C22.2 No. 601.1-M90 (S# < 20,000) ANSI AAMI ES 60601-1 (S# > 20,000) CAN CSA 22.2 No. 60601-1 (S# > 20,000)
EMC (electromagnetic compatibility)	IEC 60601-1-2
Manufactured in compliance with:	EN ISO 13485
Ambient conditions (storage, transport)	
Temperature:	-25 °C to +70 °C
Relative humidity:	at 70 °C to 93 % no condensation
Atmospheric pressure:	500 hPa to 1060 hPa
Ambient conditions (operation)	
Temperature:	+5 °C to +40 °C
Relative humidity:	15 % to 93 %
Atmospheric pressure:	700 hPa to 1060 hPa

9. Contact

We would be happy to answer any questions you may have about our products and services.

ARTROMOT® International

Please contact your local dealer, the global DJO Headquarters in the USA, the DJO International Headquarters in England, or DJO Germany.

DJO International

Headquarters: DJO UK Ltd.
1a Guildford Business Park
Guildford
Surrey, GU2 8XG, United Kingdom
Telephone: +44 (0) 1483 459659
Fax: +44 (0) 1483 459470
Email: info@DJOglobal.eu
Website: www.DJOglobal.eu

Global DJO Headquarters

DJO, LLC
2900 Lake Vista Drive
Lewisville, TX 75067, USA
Telephone: +1 800 321 9549
Fax: +1 800 936 6569
Email: webmaster@DJOglobal.com
Website: www.DJOglobal.com

Manufacturer / Headquarters Germany

ORMED GmbH
A company of Enovis

Serial numbers < 45,000:
Merzhauser Straße 112
79100 Freiburg, Germany

Serial numbers > 45,000:
Bötzingen Straße 90
79111 Freiburg, Germany

Telephone: +49 (0) 761 4566 01
Fax: +49 (0) 761 456655-01
Email: medizintechnik@enovis.com
Website: www.enovis-medtech.de

Quarters Austria

ORTHOMED Medizintechnik GmbH
Annagasse 5/1/Top 4
1010 Vienna, Austria
Telephone: +43 1 5320834
Fax: +43 1 532083431
Email: orthomed@orthomed.cc
Website: www.orthomed.cc

Warranty:

2 years (mechanical parts)
2 years (electronics)

10. Technical service

10.1 Technical hotline

Do you have any technical questions?
Do you need technical service?

Telephone: +49-180-5-1 ormed de
 +49-180-5-1 67 63 33
Fax : +49-180-5-3 ormed de
 +49-180-5-3 67 63 33

10.2 Shipment

To prevent damage during transport, only use the original shipping box. These boxes can be obtained from DJO.

Before packing the CPM device, set it to the transport position (see section 6.3).

10.3 Spare parts

Refer to the Service Manual for the most recent list of spare parts.

When ordering spare parts, always specify:

- Item
- Description
- Part number
- Quantity
- Serial number of the device

Note!

Refer repairs to authorized, specially trained staff.

DJO offers service training for your staff.

Note!

Country-specific power cords are available from DJO. Please contact DJO or your DJO distributor.

Surcharges may apply in certain cases to spare parts ordered in low quantities.

Item	Description	Part no.	Qty
1.	Patient chip card	0.0034.035	<input type="checkbox"/>
2.	Patient chip card (protocol)	0.0031.006	<input type="checkbox"/>
3.	Marker pen for patient chip card	2.0037.024	<input type="checkbox"/>
4.	Power cord EU-Version H05VV-F3G 3 × 1 mm ² length 4 m	0.0034.118	<input type="checkbox"/>
5.	Fuse T2A H250Vac 1500A Breaking Capacity	0.0000.107	<input type="checkbox"/>

11. CE-Declaration of Conformity

The certificate can be download from:

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<https://www.djoglobal.de/unternehmen/qualitaet-und-verantwortung>

12. IEC 60601-1-2 and IEC 60601-1-11

For **ARTROMOT®-K1** or **ARTROMOT®-K1 Classic** devices with serial numbers **>#35.000**.

The product designation **ARTROMOT®-K1** used in the text below includes the **ARTROMOT®-K1** and **ARTROMOT®-K1 Classic** products with all their variants.

The **ARTROMOT®-K1** device is subject to particular precautions regarding electromagnetic compatibility (EMC). The device must be installed and put into service strictly in compliance with the EMC directives set forth in the accompanying documents.

Portable and mobile RF communication systems may affect the **ARTROMOT®-K1**. The **ARTROMOT®-K1** should not be used adjacent to or stacked with other equipment.

If adjacent or stacked use is necessary, **ARTROMOT®-K1** should be observed to verify normal operation in the configuration in which it will be used.

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If it is necessary to replace assemblies or cables, only the manufacturer's original parts must be used to ensure continued compliance with EMC requirements after repair. This requirement applies to the power supply unit, cables and cable lengths, drive unit consisting of the motor and the control system, the programming unit incl. the coiled cable and the connector.

With regards to the electromagnetic compatibility, there is no maintenance operation required during the expected product life time!

WARNING!

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

12.1 Electromagnetic emissions

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

ARTROMOT®-K1 is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-K1** unit is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions to CISPR 11	Group 1	ARTROMOT®-K1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions to CISPR 11	Class B	ARTROMOT®-K1 is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions to IEC 61000-3-2	complies with class A requirements	
Voltage fluctuations/flicker emissions to IEC 61 000-3-3	complies	

12.2 Electromagnetic immunity

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

ARTROMOT®-K1 is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-K1** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment – guidance
Electrostatic discharge (ESD) to IEC 61 000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/bursts to IEC 61 000-4-4	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power should be that of a typical commercial or hospital environment.
Surges to IEC 61 000-4-5	± 1 kV line-to-line voltage	± 1 kV line-to-line voltage	Mains power should be that of a typical commercial or hospital environment.
Voltage dips To IEC 61000-4-11	0% U_T ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T ; 1 cycle and 70% U_T ; 25/30 cycles Single phase at 0°	0% U_T ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T ; 1 cycle and 70% U_T ; 25/30 cycles Single phase at 0°	Mains power should be that of a typical commercial or hospital environment. If the user of the ARTROMOT®-K1 unit requires continued operation during power mains interruptions, it is recommended that the ARTROMOT®-K1 unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field to IEC 61 000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.
Voltage interruptions to IEC 61000-4-11	0% U_T ; 250/300 cycles	0% U_T ; 250/300 cycles	

Note: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

ARTROMOT®-K1 is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-K1** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment - guidance
Conducted RF to IEC 61000-4-6	6 V 0,15 MHz-80 MHz 80% AM at 1 kHz	6 V 0,15 MHz-80 MHz 80% AM at 1 kHz	
Radiated RF to IEC 61000-4-3	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communications equipment to IEC 61000-4-3	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ARTROMOT®-K1** exceeds the applicable **ARTROMOT®-K1** unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **ARTROMOT®-K1** unit.

12.3 RF (Radio Frequency) Interference

EMC disturbance can affect the performance of the **ARTROMOT®-K1**. To prevent the user from unacceptable risks the device runs performance checks regularly during operation. This is what happens, if a problem is identified:

- an audio signal sounds
- the unit stops immediately
- the message ERROR and an error code appear on the display

In this case try to quit the error with the STOP key. If the error persists you may attempt to restart the unit by turning it briefly off and on with the power switch.

WARNING!

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the **ARTROMOT®-K1**, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Notice!

(for home use)

Portable and mobile RF communication systems may affect the **ARTROMOT®-K1 classic** device.

For this reason, make sure that wireless communication equipment, such as wireless home network devices, cell phones, cordless telephones and their base stations and walkie-talkies operate at a minimum distance of 3.3 m from the device. (Calculated on the basis of the maximum power output of a typical cell phone of 2 W).

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


ARTROMOT®-K1 classic setup
illustrations

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AMENDMENT TO THE INSTRUCTION FOR USE (IFU)

This document outlines the amendment to the Instruction for Use (IFU) for necessary updates and changes.

List of the affected devices:

Reference Number	Article Name	GTIN
80.00.040	ARTROMOT® -K1 Standard	00888912403191
80.00.041	ARTROMOT® -K1 Standard Chip	00888912403207
80.00.042	ARTROMOT® -K1 Comfort	00888912403214
80.00.043	ARTROMOT® -K1 Comfort Chip	00888912403221
80.00.044	ARTROMOT® -K1 Standard Basic	00190446751823
80.00.045	ARTROMOT® -K1 Classic	00888912403238
80.00.047	ARTROMOT® -K1 Classic Basic	00190446751830
80.00.023	ARTROMOT® -S3 Standard	00888912403108
80.00.024	ARTROMOT® -S3 Comfort	00888912403115
80.00.033	ARTROMOT® -E2 Compact	00888912403146
80.00.031	ARTROMOT® -E2	00888912403139
80.00.035	ARTROMOT® -SP3 Standard	00888912403153
80.00.036	ARTROMOT® -SP3 Standard Chip	00888912403160
80.00.037	ARTROMOT® -SP3 Comfort	00888912403177
80.00.038	ARTROMOT® -SP3 Comfort Chip	00888912403184
80.00.070	ARTROMOT® ACTIV E-K	00888912403283
80.00.084	ARTROMOT® -S4 Comfort	00190446278764
80.00.022	OptiFlex®-S	00190446794783
80.00.046	OptiFlex®-K1 Classic	00190446796015
80.00.048	OptiFlex®-K1 Standard	00190446796039
80.00.049	OptiFlex®-K1 Comfort Chip	00190446796046

First change: Warning

New warnings have been implemented due to labeling changes that add or strengthen a warning, precaution, or include additional information.

Safety warnings:

Warning!

"The treatment of minors or persons with limited physical, sensory or mental abilities (e.g. movement disorders, epilepsy, seizures, mental disorders, behavioural disorders or neurodevelopmental disorders) may only be carried out under the constant observation of attendants (e.g. specialist staff, legal guardians or custodians)

– After data storage, write the patient's name on the patient chip card. The card should only be used for this particular patient. If the patient chip card is used for another patient, be sure to delete the previous patient data from the card first. Use original chip cards only."

Second change: Label

Symbols

Table 1: Symbols changes

Symbol change	Old version	New version	Description of symbol acc. ISO 15223 + 2014/
UDI	New requirement		Indicates a carrier that contains unique device identifier information. This symbol identifies the UDI carrier, including the AIDC and human readable information. (01) GTIN code of device (11) Date of production (21) serial number of device
Manufacturer + Date of production + country of manufacture	 Name Adress 2015	 YYYY-MM	Indicates the medical device manufacturer Indicates the date when the medical device was manufactured DE is the two or three letter country code as defined in ISO 3166-1. Indicates the country of manufacture
Manufacturer's name	ORMED GmbH A DJO Company	ORMED GmbH a company of Enovis	Name of the manufacturer
Manufacturer's address spelling	Böszinger Straße 90 79111 Freiburg / Germany	Böszinger Strasse 90 79111 Freiburg / Germany	Address of the manufacturer
Norm/Standard ETL classification	Conforms to ANSI/ AAMI ES STD 60601-1 Certified CAN/CSA C22.2 No. 60601-1	 InterTek 5029964 MEE	Medical electrical equipment relevant standards according ETL classification – control number above

Other changes

Manufacturers contact was changed to: Tel. +49 761 4566-01 · Fax +49 761 4566-5501, medizintechnik@enovis.com, enovis-medtech.de

Third change: Electromagnetic emissions, Electromagnetic immunity and RF (radio frequency) interference

The EMC table provides a detailed comparison of electromagnetic compatibility as per the latest standards.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
ARTROMOT® CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the ARTROMOT® CPM device unit is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions to CISPR 11: 2015+A1:2016+A2:2019	Group 1	ARTROMOT® CPM device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions to CISPR 11: 2015+A1:2016+A2:2019	Class B	ARTROMOT® CPM device is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions to IEC 61000-3-2: 2005+A1:2008 +A2:2009	complies with class A requirements	
Voltage fluctuations/flicker emissions to IEC 61000-3-3: 2013	complies	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

ARTROMOT® CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®** CPM device unit is used in such an environment.

Immunity test	IEC 60601-1-2:2014 +A1:2020		Electromagnetic environment – guidance
	test level Home healthcare	Compliance level Home healthcare	
Electrostatic discharge (ESD) to IEC 61000-4-2:2008	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/bursts to IEC 61000-4-4:2012	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power should be that of a typical commercial or hospital environment.
Surges to IEC 61000-4-5:2014	± 1 kV line-to-line voltage	± 1 kV line-to-line voltage	Mains power should be that of a typical commercial or hospital environment.
Voltage dips To IEC 61000-4-11:2017	0% Ut ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% Ut ; 1 cycle and 70% Ut ; 25/30 cycles Single phase at 0°	0% Ut ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% Ut ; 1 cycle and 70% Ut ; 25/30 cycles Single phase at 0°	Mains power should be that of a typical commercial or hospital environment. If the user of the ARTROMOT® CPM device unit requires continued operation during power mains interruptions, it is recommended that the ARTROMOT® CPM device unit be powered from an uninterruptible power supply or a battery.
Voltage interruptions to IEC 61000-4-11:2017	0% Ut; 250/300 cycles	0% Ut; 250/300 cycles	
Power frequency (50/60 Hz) magnetic field to IEC 61000-4-8:2009	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.

Note: Ut is the a.c. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

ARTROMOT® CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®** CPM device unit is used in such an environment.

Immunity test	IEC 60601-1-2:2014 +A1:2020		Electromagnetic environment – guidance
	test level Home healthcare	Compliance level Home healthcare	

Conducted RF to IEC 61000-4-6:2013	3 V 0.15 MHz-80 MHz 80% AM at 1kHz 6V 0.15 MHz-80 MHz 80% AM at 1kHz at ISM frequency bands and Amateur radio bands	3 V 0.15 MHz-80 MHz 80% AM at 1kHz 6V 0.15 MHz-80 MHz 80% AM at 1kHz at ISM frequency bands and Amateur radio bands	
Radiated RF to IEC 61000-4-3: 2006+A1:2007 +A2:2010	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communica- tions equipment to IEC 61000-4-3: 2006+A1:2007 +A2:2010	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	
Radiated fields in close proximity IEC 61000-4-39:2017	8 A/m 30 kHz No modulation (CW) 65 A/m 134.2 kHz PM 2.1 kHz 7,5 A/m 13.56 MHz PM 50 kHz	8 A/m 30 kHz No modulation (CW) 65 A/m 134.2 kHz PM 2.1 kHz 7,5 A/m 13.56 MHz PM 50 kHz	

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ARTROMOT®** CPM device exceeds the applicable **ARTROMOT®** CPM device unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **ARTROMOT®** CPM device unit.

The **ARTROMOT®** CPM movement splint is intended for operation in the specified electromagnetic environment. The customer or user of the **ARTROMOT®** CPM device should ensure that it is used in such an environment. The required interference levels of the specified standards are met.

Hersteller/Manufacturer:

ORMED GmbH
A company of Enovis
Bötzingen Straße 90 · 79111 Freiburg · Germany
Tel. +49 761 4566-01 · Fax +49 761 4566-5501
medizintechnik@enovis.com
www.enovis-medtech.de