DONJOY®

ISOFORM™ TLSO








ENGLISH.....	2
ESPAÑOL.....	4
DEUTSCH.....	6
ITALIANO.....	8
FRANÇAIS.....	10
NEDERLANDS.....	12
PORTUGUÊS.....	14
ČEŠTINA.....	16
SLOVENČINA.....	18
DANSK.....	20
SVENSKA.....	22
SUOMI.....	24
NORSK.....	26
简体中文.....	28
日本語.....	30

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DJO, LLC
5919 Sea Otter Place,
Suite 200
Carlsbad, CA 92010 · USA

ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE: The intended user should be a suitably trained professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

INTENDED USE/INDICATIONS: The TLSO brace is designed to support and stabilize the spine for acute pain relief through postural extension alignment and may be suitable for people with compression fractures of the thoracic spine, osteoporosis, Thoracic mechanical back pain and Kyphosis.

PERFORMANCE CHARACTERISTICS: Soft-good/semi-rigid designed to restrict motion through elastic or semi-rigid construction. Providing immobilization or controlled movement of the limb or body segment.

CONTRAINDICATIONS: Pregnancy and circulation, pulmonary, cardiovascular or skeletal conditions which have risk to be made worse as a result of compression and/or pressure, unstable, displaced fractures. Do not use if you are allergic to any of the materials contained within this product.

WARNINGS AND CAUTIONS: If you experience any pain, swelling, sensation changes, excessive discoloration, irritation to the skin or any unusual reaction while using this product, consult your medical professional. If any distortion of the support panels, fraying or deterioration of lace cables is noticed, discontinue use. Brace is intended to be worn over clothing. Brace should be appropriately sized to the patient to ensure proper support and function. If the brace does not fit properly, it may be adjusted by a healthcare professional. Care should be taken with patients, especially the elderly, not to apply the brace too tightly, which MAY cause excessive pressure. Patients should be instructed on appropriate tensioning during ongoing use. Patient should see a healthcare professional for a replacement if a product quality issue occurs within the warranty timeframe. Brace should not be exposed to moisturizers, chemicals, or solvents, as they may affect durability.

APPLICATION INFORMATION:

SIZING:

- Begin by taking a circumferential measurement of the patient approximately 1" (2.5 cm) above the widest part of the hips. For patients whose circumferential measurement is between sizes, select the larger brace.
- Match the patient's size to the Sizing Chart below.

SIZE	CIRCUMFERENTIAL MEASUREMENT
S/M	Less than 48" (122 cm)
L/XL	49" - 60" (125 cm - 152 cm)

NOTE: It is not uncommon for the patient to experience post-surgical swelling, followed by weight loss (10-20+ lbs/4.5-9.1 kg). The TLSO brace is designed with two (2) adjustable belts which can be easily sized to achieve optimal fit.

PREPARE BRACE:

- Open both of the anterior/side Velcro® panels. Keep the shoulder straps attached and extend to maximum length. (Figure A)

APPLICATION INSTRUCTIONS:

- With patient standing, lift the TLSO over the patient's head. Lower the TLSO into position, keeping the patient's head between the two shoulder straps. Let the front/ anterior panel and back/posterior panel hang over the torso.
- Close both of the compression side panels so that they are even and comfortably snug on the Anterior panel (Figure B).
- Adjust the sternal Y-Bar while the patient is seated in prescribed posture. Pull up on the black release button (on top portion of Anterior panel) to release the sternal bar lock. Slide the sternal Y-bar into position, so that the top of the sternal pad fits 1/2" below the sternal notch. Be sure to keep the bottom edge of the Anterior Panel aligned with the pubic symphysis (Figure C).
- To adjust the shoulder straps, release and adjust the posterior Velcro® straps first. Adjust to desired length. Be sure to maintain the bottom edge of front panel with the pubic bone and the back panel between shoulder blades and lower back (Figure D).

Underarm shoulder strap Option: Release the posterior velcro from the back panel and insert into lateral strap slot. Adjust to desired length.

- To activate compression, grasp pull-handle and pull straight out from the torso. After the desired level of compression has been achieved, fasten the pull handle to anterior

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portion of the compression panels. Adjust compression after five minutes, then re-adjusting as needed (*Figure E*).

REMOVAL OF BRACE:

1. Release the compression handles and remove one of the shoulder straps.
2. Release one side of the side panel from front panel and slide the orthosis off the other shoulder like removing a jacket or backpack.

CUSTOMIZING TLSO: (To be performed by a licensed medical professional only)

- **Waist Belt:** If needed, to modify the waist belt, adjust the side panels by un-hooking the Velcro on the backpanel and moving panels more center for a smaller waist, or closer to the edge of the back panel for a larger waist.
- **Shoulder Straps:** If needed, the extra length in the straps can be trimmed with scissors.
- **Sternal Bar:** May be bent as needed by using a stable edged of a table. Using both hands, bend the bar in a slow and controlled fashion to adjust to desired degree of offset. The stock sternal Y-bar is set at 15°.

HOT/COLD PACK (SOLD SEPARATELY)

1. Follow instructions on hot/cold pack. Insert hot/cold pack into pocket and adhere to back panel.

ADJUSTING PULL TAB LENGTH (OPTIONAL)

1. Remove lace spool from the pull tab pocket and wind lace around the spool to shorten the length of the pull tab. Unwind lace from the spool to lengthen the pull tab. (*Figure F*)
2. Insert lace spool back into pull tab pocket.

INSERT/REMOVE FRONT PANEL

1. Open inside panel pocket. Remove semi-rigid panel. (*Figure G*)
2. Open inside panel pocket and insert semi-rigid panel.

CLEANING INSTRUCTIONS: Separate the brace's removable REAR PANEL before cleaning.

Hand wash in water (86°F/30°C), using mild soap. Air dry.

NOTE: If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

Do NOT iron.

Do NOT tumble dry.

Do NOT bleach.

MATERIAL CONTENTS: Nylon, Polyester, Polyurethane, Polyethylene, Acrylic, Polypropylene, PVC, Synthetic Rubber, EVA, Aluminum, Steel, Glass, Acetal, Calcium Carbonate

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations; the provisions of such local regulations will apply.

**FOR SINGLE PATIENT USE ONLY
RX ONLY**

NOT MADE WITH NATURAL RUBBER LATEX.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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