ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER

The intended user should be a licensed healthcare professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

INTENDED USE/INDICATIONS

- The device is designed to relieve pain, improve proprioception and knee function by supporting the body in guiding and stabilizing the knee joint.
 - Mild to moderate single compartment osteoarthritis (medial or lateral)
 - General instability of the knee joint by pain
 - Various valgus knee conditions that may benefit from unicompartamental unloading
- Performance characteristics:**
Provide mild compression of limb or body segment to help improve proprioception through elastic construction. Provide support and stabilization through elastic or semi rigid construction.

CONTRAINDICATIONS

- Any known allergies to the medical device materials.
- Peripheral arterial disease.
- Impaired blood flow in the leg (e.g. poor circulation, diabetes mellitus).
- Impaired lymphatic drainage or soft tissue swelling of unknown origin.

WARNING AND PRECAUTIONS

- The product is intended for single patient use only.
- The product is to be prescribed and fitted by a qualified healthcare professional.
- The frequency, duration of use and directions for use should be determined by your prescribing healthcare professional.
- If patients experience pain, swelling, redness or any other unusual reactions while using this product, they should contact their doctor immediately.
- If you have skin conditions/trauma in the relevant part of the body, the product should be fitted and worn only after consultation with a medical professional.
- Do not apply on open wounds.
- The brace should not be worn on both left and right legs during activities that may cause the braces to contact each other resulting in potential loss of control and injury.
- If patients experience pain, swelling, redness or any other unusual reactions while using this product, they should contact their doctor immediately.
- Do not use this device if it is damaged and/or packaging has been opened.
- Do not use this product in any way that may damage its integrity.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

FITTING, APPLICATION, AND REMOVAL INSTRUCTIONS

* PUTTING ON THE DONJOY OA GO

- NOTE:** The hinge bar should be located on the unaffected side of the knee:
 - for medial relief, the hinge is positioned laterally/on the outside of the knee
 - for lateral relief, the hinge is positioned medially/on the inside of the knee
- Open the braces at the top using the following steps:
 - 1. Pull the BOA* closure system to the left to open the brace (Fig. 1)
 - 2. Undo the upper hook by pulling the tab, then remove it from the baseplate (Fig. 2, 3 and 4)
 - 3. Undo the hook/loop fastening on the brace (Fig. 5 and 6)
- While seated, slip on and pull the brace up until you ensure the kneecap is centered in the middle of the opening (Fig. 7 and 8). Stand up to complete the fitting.
- While standing, wrap the flaps around your thigh so the brace fits properly and fix the hook/loop closure (Fig. 9). Replace the upper cover and lock it in place (Fig. 10).
- The mid hinge should be positioned at top of the knee cap (Fig. 11). Adjust as necessary.
- Press down on the two BOA closure systems to engage them (Fig. 12 and 13), and turn both to tighten and pre-position the orthosis (Fig. 14).
- Tighten the two BOA closure systems until the desired relief is achieved, ensuring that the brace maintains the correct position on the leg.
- While standing, make sure there is no gap between the condyle pad and your knee (Fig. 15). If there is, use a roller pin to fix the best fit.

TIPS: If you need to change the pad for a better fitting, just unscrew it from the inside of the hinge (using the hook and loop system) and attach the new pad in its place (Fig. 16)

CONSELS: As you close the brace, on perform knee bends and walk around to verify no brace slippage occurs and check the relief forces:

- If the brace slips down and/or if you don't feel relief after a few steps, the relief setting is too loose (turn the BOA dial clockwise to tighten)
- If you feel the brace is digging in the setting is too tight (turn the BOA dial anti-clockwise to loosen)

Adjust the relief as necessary.

Both upper and lower BOA closure systems can be adjusted independently for a comfortable fit.

* TAKING OFF THE DONJOY OA GO

- Pull up both BOA closure systems for quick release (Fig. 17).
- Undo the upper cover and remove it from the baseplate (Fig. 18).
- Open the upper hook/loop closure system.
- Put the Donjoy OA GO down.

TIPS: Close the hook/loop closure system before storing to avoid damaging the fabric. Lock the upper cover for easy storage.

CLEANING INSTRUCTIONS

Hand wash (30°C max) with mild detergent.
Do not bleach. Rinse thoroughly. Air dry. Do not machine dry. Do not dry clean. Do not iron.

COMPOSITION

Elastane, Stainless steel, Polyester, Polyamide, Glass fiber, Polyurethane, Polyvinylmethyle, Polycarbonate, Thermoplastic Elastomer, Natural grass, Polyethylene, Silicone.

LATEX NOT MADE WITH NATURAL RUBBER LATEX

INTENDED FOR USE ON A SINGLE PATIENT

WARRANTY DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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