
ARTROMOT®-K1

SW-Version: 2.X

- DE · Gebrauchsanleitung
- GB · Operating Instructions
- FR · Mode d'emploi
- ES · Instrucciones de uso
- IT · Istruzioni per l'uso
- NL · Gebruiksaanwijzing

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<https://www.djoglobal.de/unternehmen/qualitaet-und-verantwortung>

Contents

Gerätebeschreibung der ARTROMOT®-K1	III
Abbildungen zur Einstellung der ARTROMOT®-K1	VIII
Piktogrammübersicht ARTROMOT®-K1	V
1. How to use the CPM device	35
1.1 Fields of application	35
1.2 Therapy objectives	35
1.3 Indications	35
1.4 Contraindications	35
2. Description of the ARTROMOT®-K1	36
2.1 Description of the device components	37
2.2 Description of the control pendant	38
2.3 Explanation of symbols	40
2.4 Explanation of symbols (connections and nameplate)	41
3. Safety information	43
4. Device setup	46
4.1 Connecting the device, performance check	46
4.2 Connecting the external muscle stimulator ¹	47
4.3 Adjusting the device to the femur length	47
4.4 Adapting the leg support assemblies/footplate	47
4.5 Connecting the potential equalization conductor (SN > 35,000)	48
5. Setting the treatment values	49
5.1 General information on programming the ARTROMOT®-K1	49
5.2 Programming ARTROMOT®-K1 Standard models	50
5.3 Treatment value details – Standard models	50
5.4 Programming ARTROMOT®-K1 Comfort models	52
5.5 Protocol details – Comfort model	53
6. Care, maintenance, transport	57
6.1 Care / re-use	57
6.2 Maintenance (fuse replacement)	58
6.3. Transport	59
7. Environmental protection statement	60
8. Specifications	60
9. Contact	62
10. Technical service	63
10.1 Technical hotline	63
10.2 Shipment	63
10.3 Spare parts	63
11. CE-Declaration of Conformity	64
12. IEC 60601-1-2 and IEC 60601-1-11	202
12.1 Electromagnetic emissions	202
12.2 Electromagnetic immunity	203
12.3 RF (Radio Frequency) Interference	205

1. How to use the CPM device

1.1 Fields of application

The **ARTROMOT®-K1** is a motor-operated **Continuous Passive Motion (CPM)** device used to mobilize knee and hip joints.

Suitable for use in hospitals, clinics, general practices and rental services as well as in the patient's home, it is an important supplement to medical and therapeutic treatment.

1.2 Therapy objectives

CPM therapy with the **ARTROMOT®-K1** is mainly used to prevent the negative effects of immobilization, to allow patients to regain painless mobility of joints at an early stage and to promote healing and achieve a positive functional result.

Other objectives of therapy include:

- improvement of joint metabolism
- prevention of joint stiffness
- promotion of the regeneration and healing of cartilage and damaged ligaments
- faster hematoma/fluid resorption
- improved lymph and blood circulation
- thrombosis and embolism prophylaxis

1.3 Indications

The physiotherapy unit is indicated in the treatment of most injuries and diseases of the knee and hip joints as well as in postoperative treatment after knee and hip joint surgery. Examples:

- joint distortion and contusion
- arthrotomy and arthroscopy procedures in combination with synovectomy, arthrolysis, or other intra-articular interventions
- mobilization of joints in anesthetized patients
- operative treatment of fractures, pseudoarthrosis, and osteotomy
- cruciate ligament replacement or reconstruction
- endoprosthetic implants

1.4 Contraindications

Do NOT use the **ARTROMOT®-K1** on patients with:

- acute inflammatory processes in the joints, unless on the order of a physician
- spastic paralysis
- unstable osteosynthesis

2. Description of the ARTROMOT® -K1

The motorized CPM device permits extension and flexion of the knee joint in the range of $-10^{\circ}/0^{\circ}/120^{\circ}$, and of the hip joint in the range of $0^{\circ}/7^{\circ}/115^{\circ}$. It can be used on either side and requires no configuration change.

These are some of the outstanding features of the **ARTROMOT®-K1**:

- anatomically correct adjustment
- physiological movements
- control pendant for precise adjustment of patient-specific therapy values
- symbols for easy operation of the control pendant
- programmed therapy parameters saved to chip card¹

Biocompatibility

Those parts of the **ARTROMOT®-K1** that come into contact with the patient when the device is used as intended, are designed to fulfill the biocompatibility requirements of the applicable standards.

Essential Performance Features

- The programmed angles are accurately maintained with a tolerance of $\pm 2^{\circ}$.
- The programmed speeds are accurately maintained with a tolerance of $\pm 5\%$.
- The selected mode and mechanical settings do not change during operation.

Frequent functions:

- a) Unpacking (device and accessories)
- b) Device connections mains power connection, connecting the control pendant and the motion element)
- c) Mechanical adjustments of the motion element on the rotational axes
- d) Programming the control pendant (range of motion, speed, functions)
- e) Storage

Warning!

Patient hazard -

- If the Patient operates the **ARTROMOT®-K1**, the Patient needs to know how to **START** or **STOP** the device.
- To prevent the parameters from being changed unintentionally, lock the **MENU** button on the remote control. Contact your local DJO office for assistance.

2.1 Description of the device components

Note: See page III !

1. Thigh support assembly (type B applied part)
2. Clamping lever to adjust the height of the thigh support assembly
3. Thigh length scale (femur length scale)
4. Thigh length fixation screws (femur length)
5. Knee hinge
6. Calf support assembly (type B applied part)
7. Clamping levers to adjust the height of the calf support assembly
8. Calf length fixation screws (tibia length)
9. Calf length scale (tibia length scale)
10. Strap to secure the foot to the footplate
11. Footplate (type B applied part)
12. Footplate angle fixation screw, also allows swiveling the footplate
13. Clamping lever to adjust footplate rotation and level and to remove the footplate
14. Connection for **ARTROSTIM®-FOCUS®plus** control cable²
15. Connection for control pendant
16. Connection for power cord
17. Fuse cap
18. Power switch (ON/OFF)
19. Rating plate
20. Control pendant (type B applied part)
21. Patient chip card¹
22. Compartment for storage of control pendant
23. Potential equalization terminal³


¹ Devices of the **ARTROMOT®-K1** product family with patient chip card only.

² **ARTROMOT®-K1** Comfort devices only

³ Applies only to **ARTROMOT®-K1** devices with serial numbers > 35,000


2.2 Description of the control pendant

2.2.1 Control pendant in normal mode


¹ Devices of the ARTROMOT®-K1 family with patient chip card only.

2.2.2 Control pendant in MENU selection mode

2.2.3 Control pendant in programming mode

2.3 Explanation of symbols

Also refer to symbol overview on page V/VI !

Standard protocols:

	Extension (stretching the knee)
	Flexion (bending the knee)
	Speed
	Warm-up protocol
	Extension pause
	Flexion pause
	Timer (therapy timer)
	Reverse on load (feature for patient safety)
	Transport setting
	New patient
	Total therapy time
	Service menu

Comfort protocols:

	Stretching in extension
	Stretching in flexion
	Workout protocol
	Comfort protocol
	EROM repeat extension
	EROM repeat flexion
	EMS control
	Therapy documentation

2.4 Explanation of symbols (connections and nameplate)

-  Alternating current
-  For serial numbers < 20,000:
Protection class I equipment. The medical device must be connected to a system with protective earth conductor!
-  For serial numbers > 20,000:
Protection class II equipment. The medical device has a double or reinforced insulation.
-  Type B applied part
-  Power switch OFF
-  Power switch ON
-  The year and month of manufacture are indicated next to this factory symbol.
-  The name next to this factory symbol is the manufacturer.
-  The number next to this symbol is the article reference number.
-  Device complies with Council Directive 93/42/EEC about medical devices, tested and approved by DQS Medizinprodukte GmbH.
-  The number next to this symbol is the serial number.
-  Caution! Observe warnings set forth in instructions for use!
-  "Refer to instruction manual!" ISO 7010-M002
-  "Refer to instruction manual!" ISO 7000 -1641
-  Do not dispose with unsorted municipal waste.

Protect from moisture.

IP21

The IP rating indicates the level of protection and thus the suitability of the device for use under different ambient conditions.

The rating IP21 means:

2 is the level of protection against contact and solid objects.

The digit 2 means:

- Protection against contact: protected against contact with a finger.
- Protection against solid objects: protected against solid objects (diameter of 12.5 mm and greater).

1 is the level of protection against moisture.

The digit 1 means: protection against vertically falling water drops.

Warning! Depending on the device settings, the moving parts of the device present pinch points.

Pay particular attention to small children and babies!

"Observe temperature limits" (storage)! ISO 7000-0632

The overall weight of the device is indicated next to this symbol.

For serial numbers > 35,000: Potential equalization terminal.

3. Safety information

Definitions

Be sure to read the safety information before using the CPM device. The safety statements are classified as follows:

Danger!

indicates an imminent hazard. If not avoided, this hazard will result in death or serious injury.

Warning!

indicates a hazard. If not avoided, the hazard can result in death or serious injury.

Caution!

indicates a potential hazard. If not avoided, the hazard may result in minor injury and/or product/property damage.

Warning!

Patient hazard –

- While the **ARTROMOT®-K1** is in use, do not repair and maintenance the device! There are no parts that can be repaired during use.
- There are no user-serviceable parts inside the device. If a malfunction occurs, discontinue use immediately and contact your local DJO office for assistance and service information.

Safety information

Danger!

Explosion hazard –

The **ARTROMOT®-K1** is not designed for use in areas of rooms used for medical purposes where an explosion hazard may occur. An explosion hazard may result from the use of flammable anesthetics, skin cleansing agents, and disinfectants.

Warning!

Patient hazard –

- Only authorized individuals are allowed to operate the **ARTROMOT®-K1**. Individuals are authorized after receiving training in the operation of the device and reading these instructions for use.
- Before using the device, the operator must ascertain that it is in correct working order and operating condition. The cables and connectors, in particular, must be checked for signs of damage. Damaged parts must be replaced immediately, before use.
- **Before therapy**, a **test run** consisting of several exercise cycles must be completed, first without and then with the patient. Check that all setting screws are tightened.
- Stop therapy immediately when you have doubts about the device settings and/or the therapy protocol.

Warning!

Patient hazard –

- It is important that the patient's **position** is anatomically correct. Therefore, carefully verify the following settings/positions:
 1. femur length
 2. knee joint axis
 3. tibia length and leg rotation
 4. leg support assemblies
- Movement must always be **free of pain** and **irritation**.
- The patient must be **fully conscious** during instruction and when using the device.
- The **choice** of the **therapy parameters** to program and of the **therapy protocols** to use is restricted to the responsible **physician or therapist**. The physician or therapist must decide on a case-by-case basis whether the device can be used with a particular patient.
- The patient must be familiar with the functions of the **ARTROMOT®-K1 control pendant** and the control pendant must be **within easy reach** of the patient, allowing him or her to stop therapy if needed. **Patients unable to operate the control pendant**, e.g. paralytic patients, must always be supervised by specially trained staff during therapy.
- After data storage, write the patient's name on the **patient chip card**¹. The card should only be used for this particular patient. If the patient chip card¹ is used for another patient, be sure to **delete the previous patient's data from the card first**. (See: section 5.2 Programming: "New Patient"). Use **original chip cards** only.¹
- The **ARTROMOT®-K1** must only be used with **accessories** approved by DJO.
- Modifications to the medical device described in this document without the manufacturer's written consent are prohibited.

¹ Devices of the **ARTROMOT®-K1** product family **with patient chip card only**.

- The simultaneous treatment of both legs by simultaneous use of two CPM devices is not permitted because the motion elements might interfere with each other.
- Stability of the CPM device must always be ensured while it is in use. The **ARTROMOT®-K1** must only be set up on surfaces that guarantee its stability. Very soft or instable surfaces (such as waterbeds) are **NOT** suitable.
- Do not allow **parts of the body or any objects** (such as blankets, cushions, or cables) to **get caught in the moving parts of the CPM device**.

Warning!

- Special care must be taken when small children and babies are present while the device is in use! Sufficient distance to the device is mandatory for their safety!
- Never leave the device unattended when it is switched on! Switch the device off and unplug the power cord from the wall outlet!
- After use, store the device in a safe place! Ensure device stability also during storage!

Warning!

Shock hazard –

Strictly observe the following warnings. Failure to do so endangers the lives of the patient, the user, and other persons involved.

- Allow the **ARTROMOT®-K1** to reach room temperature **before use**:
The time required for from the minimum storage temperature between uses until **ARTROMOT®-K1** is ready for its intended use when the ambient temperature is 20 °C is 5 hours.
The time required for **ARTROMOT®-K1** to cool from the maximum storage temperature between uses until In-tellect RPW 2 is ready for its Intended use when the ambient temperature is 20 °C is 5 hours.
- The **ARTROMOT®-K1** must only be operated in **dry rooms**.

- When disconnecting the device from the power line, remove the plug from the wall outlet first, before disconnecting the cable from the device.
- When connecting the device to other equipment or when creating a medical system, check that the sum of leakage currents will not cause any hazard. Please contact DJO, if you have questions in this matter.
- Do not use multiple portable socket outlets (MPSO) to connect the device to the power line. The **ARTROMOT®-K1** must be connected to a **properly installed wall outlet with a non-fused earthed wire**.
Before connecting the power cord, it must be completely unrolled and placed such that it will not get caught in the moving parts of the device.
- Before cleaning and service interventions, **disconnect the device from the power line** by removing the power cord from the wall outlet.
- **Liquids must not be allowed** to enter the CPM device or the control pendant. If liquids have entered into the devices, the **ARTROMOT®-K1** must be immediately checked by a service technician, before it can be reused.

Warning!

Equipment malfunction –

- Magnetic and electrical fields are capable of interfering with the proper performance of the unit. For this reason make sure that all external devices operated in the vicinity of the device comply with the relevant EMC requirements. X-ray equipment, MRI devices, radio systems, cell phones, etc. are possible sources of interference as they may emit higher levels of electromagnetic radiation. Keep the device away from such equipment and verify its performance before use.
- Refer repair and maintenance to authorized persons.
Persons are authorized after training by a specialist trained and commissioned by the manufacturer.

- **Route all cables below the device frame to either side**, ensuring that they cannot get caught by the moving parts during operation.
- **Check the ARTROMOT®-K1** for damage and loose connections at least **once a year**. Damaged and worn parts must immediately be replaced with original spare parts by authorized staff.

Caution!

Preventing chafing and pressure sores –
When your patient is **adipose**, very **tall** or **very short**, be sure to prevent chafing and pressure sores. Place the leg concerned in a moderate abductive position if deemed appropriate.

Caution!

Equipment damage –

- Check that the voltage and frequency ratings of your local **power line** are those indicated on the nameplate.
- The leg support element withstands a **maximum continuous load of 20 kg**
S# < 20,000 / 25 kg S# > 20,000.
- **Do not allow any objects** (such as blankets, cushions, or cables) to get caught **in the moving parts** of the CPM device.
- Do not expose the **ARTROMOT®-K1** to direct sunlight, because some of the components may warm up excessively. Furthermore, the device must be set up at a safe distance from radiators to prevent excessive temperature rises.
- The presence of children, pets and rodents does not normally impair the functioning of the device. However, avoid contamination of the device by children or animals, from dust and lint, and keep them at a safe distance from the device. The safety statements set forth apply.

4. Device setup

Note: See pages III and VIII !

4.1 Connecting the device, performance check

The equipment supplied includes these items:

base unit, control pendant (20), patient chip card (21 - **ARTROMOT®-K1** versions with chip card only), power cord (not shown), instructions for use

1. Connect the control pendant (20) to the provided socket (15) by plugging it in and closing the bayonet lock. (With patient chip card versions, insert the chip card into the slot in the control pendant.)
2. Connect the **power cord** to socket (16) of the device and connect the **power plug** to a wall outlet with a non-fused earthed wire (100–240 Volt, 50/60 Hz).
3. Turn the **power switch** (18) on.
4. Follow these steps to set the carriage to the **home position**:

ARTROMOT®-K1 without patient chip card

Press the MENU button on the control pendant until you reach program level 3 (Standard model) or program level 5 (Comfort model).

Press the New Patient parameter button → **0** ← .

Press the parameter button again to activate the function.

Press the **START** button. The carriage automatically enters the home position.

ARTROMOT®-K1 with patient chip card

Initial setup for new patients

Insert the original patient chip card (21) into the control pendant (20).

Press the MENU button on the control pendant until you reach program level 3 (Standard model) or program level 5 (Comfort model). Press the New Patient parameter button → **0** ← .

Press the parameter button again to activate the function. Press the **START** button. The carriage automatically enters the home position.

Setup with programmed chip card

Insert the original patient chip card (21) into the control pendant (20).

Press the **START** button.

The carriage automatically enters the home position.

Performance check

If the control pendant can be operated as described above and the **ARTROMOT®-K1** enters the home position (for home position values, refer to section 5.3 and 5.5), the device has passed the performance check and functions correctly.

The device also runs performance checks regularly during operation. This is what happens, if a problem is identified:

- An audio signal sounds.
- the device switches off immediately
- The message "ERR", accompanied by a code number (e.g. ERR 5), appears on the display.

In this situation, you may attempt to restart the device by turning it briefly off and on again with the power switch. With the unit switched off, check that all plugs are correctly connected.

Switch the CPM device on again: if the error message persists, have the device inspected by a Service technician, before using it again.

4.2 Connecting the external muscle stimulator¹

⚠ Caution!


Patient hazard, equipment malfunction –

Only the **ARTROSTIM®-FOCUS®plus** muscle stimulator with EN 60601-1 approval, manufactured by ORMED GmbH, may be connected to the **ARTROMOT®-K1** CPM device. Use the "muscle stimulation" cable (part no. 2.0037.024) to connect the muscle stimulator. Other cables are not approved. The cable can be ordered separately from DJO.

In addition, all such combinations – Medical Electrical Systems – shall comply with the requirements stated in clause 16 of IEC 60601-1. Any equipment not complying with the leakage current requirements in IEC 60601-1 shall be kept outside the patient environment i.e. at least 1.5 m from the patient support or shall be supplied via a Separation transformer to reduce the leakage currents.

Any person who connects external equipment to the device has formed a Medical Electrical System and is therefore responsible for the system to comply with the requirements in clause 16 of IEC 60601-1. If in doubt, contact qualified medical technician or your local DJO representative.

1. Connect the jack plug to the **ARTROSTIM®-FOCUS®plus** muscle stimulator.

¹For control of an external muscle stimulator, an **ARTROMOT®-K1 Comfort device and the corresponding connector** are required.

2. Connect the round plug to connector (14) of the **ARTROMOT®-K1** and turn it a short way clockwise to lock.
3. Turn on the **ARTROMOT®-K1** first, then turn on the muscle stimulator.

Note!

For information on connecting and programming the muscle stimulator, refer to the operation manual that comes with the **ARTROSTIM®-FOCUS®plus**.

4.3 Adjusting the device to the femur length

1. Measure the length of the patient's thigh (femur) from the greater trochanter to the lateral knee joint cavity (Fig. A).
2. Set the carriage to the home position (see 4.1).
3. Adjust the measured value at the **femur scale** (3) of the carriage (Fig. B):
 - Loosen the two fixation screws (4).
 - Extend the scale (3) to the required length.
 - Tighten the fixation screws (4) to set the scale to the new length.

⚠ Caution!

Equipment damage –

Do not attempt to extend the femur scale beyond the stop.

4.4 Adapting the leg support assemblies/footplate

1. Set the **leg support assemblies and the footplate** (1, 6 and 11) to the expected positions before accommodating the patient.
 - Loosen fixation screws (8) to adjust the **footplate** (11) to the length of the patient's lower leg (Fig. C).

Loosen clamping lever (13) and adapt the footplate's rotation and height to the patient (Fig. D).

Loosen fixation screw (12) and adapt the angle to the patient's foot (turn the screw a few revolutions until the footplate can be easily adjusted).

For **short patients** you can reverse the footplate's bracket 180° (Fig. H) to adapt the footplate to shorter calves:

- Loosen clamping lever (13) and remove the footplate (11).
- Loosen the fixation screws (12).
- Reverse the bracket 180°.
- Screw the footplate to the bracket and tighten the clamping lever.

Note!

When reversing the footplate, ensure that the pins below the clamping lever engage with the recesses in the bracket.

- To adjust the height of the **support assemblies for calf** (6) and **thigh** (1), loosen clamping levers (2) and (7) (Figs. E/F).
2. Place the patient's leg on the carriage and repeat the steps outlined under 1 above to adjust the CPM device to the patient.

⚠ Caution!

Equipment damage –

Cover the **leg support assemblies with disposable tissues**, when using the **ARTROMOT®-K1** after surgery. This helps prevent discoloration.


⚠ Caution!

Patient hazard –

Ensure that the rotational axes of the CPM device and of the knee joint coincide both in the vertical and in the horizontal plane (Fig. G).


Symbol 1:

Measurement of the patient's femur length from the greater trochanter to the lateral knee joint cavity.


Symbol 2:

Set the carriage to the home position (see 4.1) and adjust it to the measured femur length.

Symbol 3:

Adjust height of calf and thigh support assemblies. Adjust the footplate to the height and length of the lower leg.

4.5 Connecting the potential equalization conductor

To prevent interference with monitoring devices when used in treatment rooms, the CPM device must be connected to a potential equalization system according to DIN 42801.

For this purpose, when using devices with serial numbers > 35,000, connect the equalization conductor to the potential equalization terminal (23) of **ARTROMOT®-K1** and to the potential equalization connection of the treatment room.

Observe the requirements of IEC 60601-1 for potential equalization connections.

5. Setting the treatment values

Warning!

Patient hazard –

Before therapy, a test run consisting of several exercise cycles must be completed without the patient. Then, repeat the test run with the patient and check that the movements do not cause any pain.

Note!

See also 2.2 and 2.3 as well as VIII!

5.1 General information on programming the ARTROMOT®-K1

1. You activate the programming mode by briefly pressing the **MENU** button on the control pendant.
2. The various treatment parameters and functions are allocated to three (Standard model) or five (Comfort model) programming levels (4 per level).

To be able to program a parameter you will have to access the corresponding programming level. This is also done with the **MENU** button. With each press of the button you advance one level. The code M1, M2, etc. that appears in the middle of the display indicates the programming level.

3. You activate the treatment parameters and functions with the 4 parameter buttons below the display. The symbols above the **4 parameter buttons** indicate the assigned parameters and functions.

This is what happens when you press one of the parameter buttons to select a parameter:

- The corresponding symbol appears on the display in a larger format.
- The set value is displayed.
- The symbol above the parameter button is highlighted.

4. With the **+/-** buttons (plus/minus) you change the displayed value. When you press and hold the key, the value will change at a higher rate.

Some of the (special) functions can only be enabled and disabled. This is done by pressing the corresponding parameter button or with the **+/-** buttons. Active parameters are identified with a check mark in the circle next to the symbol.

5. Having programmed all parameters, press the **STOP** button to save the values.
6. Then press the **START** button to start therapy.

Note!

- Refer to sections 5.3 and 5.5 for a description of the parameters.
- To **view the set parameter values**, press the corresponding parameter button. However, this is only possible when you press the **STOP** button first.
- To prevent accidental changes of the parameter settings, you can **lock the buttons**. To do so, simultaneously press the **+ and - buttons** for **approx. 4 seconds**.

Press both buttons again for approx. 4 seconds to unlock.

- Selecting the "New Patient" function will automatically delete the data on the patient chip card. When you have finished programming the unit and press the **STOP** button, the settings will automatically also be saved to the patient chip card.
- **Emergency stop function:** The **ARTROMOT®-K1** will stop immediately, when any of the buttons is pressed during therapy. Patient treatment can be resumed by pressing the **START** button. The unit will automatically change the direction of movement.
- The manufacturer recommends an application period of up to 1 hour.

- For a uniform, smooth and gentle transition between the two directions of motion, the speed is automatically reduced before reaching the reversal point and, on reaching the reversal point, it is continuously increased until it reaches the set speed.

Patients with a programmed chip card

- Insert the chip card (the patient is not yet positioned on the unit).
- Perform the mechanical adjustments (femur length, etc.).
- Place the patient on the CPM device and press the **START** button to initiate therapy.

5.2 Programming ARTROMOT®-K1 Standard models

Different programming levels are provided to program the **ARTROMOT®-K1 Standard** models.

You change between levels by pressing the **MENU** button.

The display always indicates the currently selected level.

The following **treatment values, settings and information** can be entered/viewed on the control pendant (20):

LEVEL 1:

- Extension (stretching the knee)
- Flexion (bending the knee)
- Speed
- Warm-up protocol

LEVEL 2:

- Extension pause
- Flexion pause

- Timer (therapy time)
- Reverse on load (feature for patient safety)

MENU

LEVEL 3:

- Transport setting
- New patient
- Total therapy time
- Service menu

Note!

- **While you adjust the extension/flexion values**, the carriage will move to the set range. This allows you to easily and quickly determine the ROM where the patient does not experience pain.
- **Special functions** can be programmed and retrieved with **ARTROMOT® -K1 Comfort** devices (see sections 5.4 and 5.5).
- Only **ARTROMOT®-K1 chip card** models allow therapy protocols to be **saved** to an inserted chip card.

5.3 Treatment value details – Standard models

- You access the different programming levels by repeated depressions of the MENU button.
- You select the treatment parameters with the corresponding parameter button.
- You change the treatment values with the +/- buttons and you enable/disable functions by pressing the corresponding parameter button again.
- You save the settings by pressing the STOP button.

LEVEL 1:

- **Extension (stretching)** 
 - Maximum **knee extension**: -10°
 - Maximum **hip extension**: 7°
- **Flexion (bending)** 
 - Maximum **knee flexion**: 120°
 - Maximum **hip flexion**: 115°

Note!

The programmed value and the value measured at the patient's knee may differ slightly.

- **Speed** 

The speed can be adjusted in steps of 5% between 5% and 100%.

Default setting: 50 %

- **Warm-up protocol** 

During warm-up, the patient will slowly become used to the set maximum extension and flexion values, starting from the center position.

The warm-up protocol starts in the middle between the two maximum values set for stretching and bending. The range of motion increases with each cycle, until the programmed maximum values are reached after a total of 15 cycles.

Default setting: disabled

LEVEL 2:

- **Extension pause** 

Pauses occur at the set extension limit, just before the bending movement starts. Pauses can be set to any value between 0 and 59 seconds in steps of 1 second, and to values between 1 and 59 minutes in steps of 1 minute.

Default setting: no pause

- **Flexion pause** 

Pauses occur at the set flexion limit, just before the stretching movement starts. Pauses can be set to any value between 0 and 59 seconds in steps of 1 second, and to values between 1 and 59 minutes in steps of 1 minute.

Default setting: no pause

- **Timer (therapy time)** 

Default setting is continuous operation. A clock symbol in the upper right-hand corner of the display identifies the continuous mode of operation. The clock indicates the elapsed therapy time. In the **continuous mode**, the device must be stopped with the **STOP** button.

However, you can also select therapy durations of **1 to 59 minutes in steps of 1 minute** and of **1 to 24 hours in steps of 30 minutes**.

When the time has elapsed, the unit switches **automatically** off and stops in the position equivalent to the set extension value $+10^{\circ}$. In this case, a circle replaces the clock symbol. The circle fills as the therapy time progresses.

- **Reverse on load**  **(safety function)**

The device automatically starts moving in the opposite direction of the last movement when the patient's resistance (load) exceeds the set value. Adjustable levels for reverse on load feature: 1 - 25. At level 1, very low resistance will cause the device to reverse; at level 25, a high resistance is required to initiate the reversal.

Default setting: Level 25

⚠ Caution!

Patient hazard –

The reverse-on-load feature is only a safety measure in the event of cramps, spasms, locked joints and similar conditions. The manufacturer cannot be held liable for misuse of this feature.

LEVEL 3:

■ Transport setting →

With this function, the carriage will move to a position optimally suited for packing the CPM device. Adjust the femur length to 49 cm and the tibia length to 45 cm. Select the function and press the **START** button. The carriage moves to the transport position. (Also refer to 6.3.)

■ New patient → 0 ←

With this function, the CPM device will move to the home position, allowing the mechanical settings to be completed. Select the function and press the **START** button. The device enters the home position and existing therapy parameters will be deleted.

With **ARTROMOT®-K1** devices with patient chip card, the factory defaults will be restored. All values stored on the chip card will be deleted.

The carriage will stop in the home position.

The "New Patient" function (home position) selects the following settings:

- Extension: 25°
- Flexion: 35°
- Speed: 50 %
- Warm-up: disabled
- Extension pause: 0
- Flexion pause: 0
- Timer: continuous operation
- Load reversal: 25
- Total therapy time: 0

■ Total therapy time Σ

ARTROMOT®-K1 models without chip card

The total therapy time is the added sum of operating hours. If the device is used by only one patient, this time is equivalent to the duration of all the patient's therapy sessions.

Under menu item "total therapy time" of **ARTROMOT®-K1 models with chip**

card, you can view each patient's total therapy time (duration of all the patient's therapy sessions).

Deleting the stored therapy time

Press and hold the parameter button for 5 seconds or select the New Patient function.

■ Service MENU

For service purposes only, refer to Service Manual.

Reminder:

You save the selected parameter values by pressing the **STOP** button.

5.4 Programming ARTROMOT®-K1 Comfort models

ARTROMOT®-K1 devices of the **Comfort** series offer two more programming levels for additional functions.

The programming levels are selected in the same way as on Standard models.

Programming levels 1 and 2 are identical with programming levels 1 and 2 of the Standard models.

With the default settings all **special functions** are disabled.

The following treatment values, settings and information can be entered/viewed on the control pendant (20):

LEVEL 1:

- Extension (stretching the knee) 
- Flexion (bending the knee) 
- Speed 
- Warm-up protocol 

LEVEL 2:

- Extension pause 
- Flexion pause 

- Timer (therapy time) 
 - Reverse on load (feature for patient safety) 
- LEVEL 3:** 
- Stretching in extension 
 - Stretching in flexion 
 - Workout protocol 
 - Comfort protocol 
- LEVEL 4:** 
- EROM repeat extension 
 - EROM repeat flexion 
 - EMS control 
 - Total therapy time 
- LEVEL 5:** 
- Transport setting 
 - New patient 
 - Therapy documentation 
 - Service menu 

5.5 Protocol details – Comfort model

- You access the different programming levels by repeated depressions of the **MENU** button.
- You select the treatment parameters with the corresponding **parameter button**.
- You change the treatment values with the **+/- buttons** and you enable/disable functions by pressing the corresponding parameter button again.

- You save the settings by pressing the **STOP** button.

With the **default settings** all **special functions** are disabled.

Note!

LEVEL 1: equivalent to level 1 of the Standard model (see: 5.3)

LEVEL 2: equivalent to level 2 of the Standard model (see: 5.3)

■ Stretching in extension

With the special "stretching in extension" function the joint will be gently stretched beyond the extension limit.

Starting at the middle position the carriage will first move to the programmed flexion limit and then to the programmed extension limit.

Subsequently the carriage reverses 5° toward the flexion angle and then moves very slowly back again to the programmed extension limit (display <=). After that it attempts to stretch the joint another 5°, moving even slower than before (display <<).

If a high resistance toward the additional 5° is sensed, the reverse-on-load function is automatically activated and the device moves in the opposite direction.

This stretch cycle is repeated 10 times. After that the carriage moves to the programmed flexion limit and restarts the stretching in extension cycle.

It is not possible to activate the special "stretching in extension" and "stretching in flexion" functions at the same time.

Note!

If an extension pause is programmed when you enable the function, the carriage will stop for the duration of the pause each time the maximum stretching value is attained.

■ Stretching in flexion

With the special "stretching in flexion" function the joint will be gently stretched beyond the flexion limit.

Starting at the middle position the car-

riage will first move to the programmed extension limit and then to the programmed flexion limit.

Subsequently the carriage reverses 5° toward the extension angle and then moves very slowly back again to the programmed flexion limit (display =>). After that it attempts to stretch the joint another 5°, moving even slower than before (display >>).

If a high resistance toward the additional 5° is sensed, the reverse-on-load function is automatically activated and the device moves in the opposite direction.

This stretch cycle is repeated 10 times. After that the carriage moves to the programmed extension limit and re-starts the stretching in flexion cycle.

It is not possible to activate the special "stretching in flexion" and "stretch in extension" functions at the same time.

Note!

If a flexion pause is programmed when you enable the function, the carriage will stop for the duration of the pause each time the maximum stretching value is attained.

■ Workout protocol

With the special "workout" function, a series of special, programmed protocols can be completed in one session.

The program includes the following protocols in a given sequence: warm-up, stretching in extension, EROM repeat extension, stretching in flexion, EROM repeat flexion, and cool-down.

The **entire workout protocol** takes **approx. 38 – 40 minutes** to complete.

- 5-minute protocol: **Warm-up**
Starting from the middle position, the range of motion is gradually increased toward extension and flexion in steps of 1°.
- 5-minute exercise according to **programmed settings**
- 5-minute protocol: **Stretching in flexion**
- 5-minute protocol: **EROM repeat flexion**

- 5-minute protocol: **Stretching in extension**

- 5-minute protocol: **EROM repeat extension**

- 5-minute exercise according to **programmed settings**

- 3-minute protocol: **Cool-down**
The cool-down protocol is the warm-up protocol of the workout mode reversed. Starting from the maximum values, the carriage reduces the range of motion by 1° until the middle position is reached.

The device switches off, when the protocol has been completed.

The indicated minutes are approximate values. Depending on the programmed maximum range of motion, the times may vary.

■ Comfort protocol

With the special "Comfort" function, the range of motion is gradually extended until the patient attains the maximum programmed extension and flexion values.

For this protocol, the maximum values are programmed first, then the special function is activated and, eventually, the treatment is started.

The **ARTROMOT®-K1 Comfort** will now complete 5 cycles in both directions with the maximum programmed values minus 5°. Then the range of motion is increased by 1° per cycle in both directions until the programmed limit values are reached. Once the limit values have been attained, the CPM device continues in the programmed range of motion until the end of the therapy session.

LEVEL 4:

■ EROM repeat extension

The special "EROM repeat extension" function allows a more efficient exercise in the last 10° before the set maximum extension value.

For this protocol, the CPM device starts midway between the set extension and flexion values. It will first move to the programmed flexion value and then to

the programmed extension value. When the extension value has been reached, the carriage reverses 10° toward the flexion angle and then moves back again to the maximum extension value. The movement through the final 10° is repeated 5 times at a slow speed.

At the end of the cycle, the carriage will again move to the maximum flexion value and then starts another cycle with 5 repetitions through the last 10° of the extension angle.

■ EROM repeat flexion

The special "EROM repeat flexion" function allows a more efficient exercise in the last 10° before the set maximum extension value.

For this protocol, the CPM device starts midway between the set extension and flexion values. It will first move to the programmed extension value and then to the programmed flexion value.

When the flexion value has been reached, the carriage reverses 10° toward the extension angle and then moves back again to the maximum flexion value. The movement through the final 10° is repeated 5 times at a slow speed.

At the end of the cycle, the carriage will again move to the maximum extension value and then starts another cycle with 5 repetitions through the last 10° of the flexion angle.

■ EMS control

This special function allows the **ARTROMOT®-K1** Comfort device to be operated in conjunction with the electrical muscle stimulator **ARTROSTIM®-FOCUS®plus**, for which it is ideally suited.

For this modality, first plug the "muscle stimulation" connection cable (part no. 2.0037.024) into the **ARTROSTIM®-FOCUS®plus** device (see **ARTROSTIM®-FOCUS®plus** instructions for use) and then into the **ARTROMOT®-K1** Comfort device (see 4.2).

Then activate the EMS control program at the **ARTROMOT®-K1** Comfort.

This function will select the following settings:

- Extension: as per your protocol

- Flexion: as per your protocol
- Extension pause: **20 seconds**
- Flexion pause: 0 seconds

Note!

All other protocols will be disabled or reset to the default values.

Start therapy by pressing the **START** button. Starting from the middle position, the carriage will move to the programmed flexion angle, then to the programmed extension setting. At this point, the carriage will stop and simultaneously send a trigger pulse to **ARTROSTIM®-FOCUS®plus**.

For the duration of the pause (**default: 20 seconds**), the **muscle stimulation protocol that you programmed at the ARTROSTIM®-FOCUS®plus** will be completed with the selected intensity (for programming of the EMS device, see: **ARTROSTIM®-FOCUS®plus instructions for use**).

After the pause, the **ARTROMOT®-K1** will turn off the EMS device and the carriage moves to the programmed maximum flexion angle.

The electrical muscle stimulation will be repeated in the next extension pause.

Note!

- **After** activation of the special "EMS control" function, the pause settings can be changed.
- Electrical muscle stimulation can **only** take place during **extension pauses**.

Caution!

Patient hazard –

The only approved device for electrical muscle stimulation is the **ARTROSTIM®-FOCUS®plus** muscle stimulator with EN 60601-1 approval, available from **ORMED GmbH**.

■ Total therapy time Σ ⌚

ARTROMOT®-K1 models without chip card

The total therapy time is the added sum of operating hours.

If the device is used by only one patient, this time is equivalent to the duration of all the patient's therapy sessions.

Under menu item "total therapy time" of **ARTROMOT®-K1 models with chip card**, you can view each patient's total therapy time (duration of all the patient's therapy sessions).

Deleting the stored therapy time

Press and hold the parameter button for 5 seconds or select the New Patient function.

LEVEL 5:

■ Transport setting →

With this function, the carriage will move to a position optimally suited for packing the CPM device. Adjust the femur length to 49 cm and the tibia length to 45 cm.

Select the function and press the **START** button. The carriage moves to the transport position. (Also refer to 6.3.)

■ New patient → 0 ←

With this function, the CPM device will move to the home position, allowing the mechanical settings to be completed.

Select the function and press the **START** button. The carriage moves to the home position.

With **ARTROMOT®-K1** devices with **patient chip card**, the factory defaults will be restored. All stored values will be deleted. The carriage will stop in the home position.

The New Patient (home position) function selects the following settings:

- Extension: 25°
- Flexion: 35°
- Speed: 50 %
- Warm-up: disabled
- Extension pause: 0
- Flexion pause: 0
- Timer: continuous operation

- Load reversal: 25
- Total therapy time: 0
- Stretching in extension: disabled
- Stretching in flexion: disabled
- EROM repeat extension: disabled
- EROM repeat flexion: disabled
- Comfort- protocol: disabled
- Workout protocol: disabled
- EMS control: disabled
- Therapy documentation: reset

■ Therapy documentation

ARTROMOT®-K1 devices of the **Comfort series with patient chip card** have a special documentation function which provides a log of all therapy sessions.

The device run times as well as the range of motion of the sessions are recorded. The collected data are presented graphically in the form of a coordinate system (X-axis = range of motion / Y-axis = time) where the upper curve illustrates the trend of the flexion movement and the lower curve the trend of the extension movement.

■ Service MENU

For service purposes only, refer to Service Manual.

Reminder:

You save the selected parameter values by pressing the **STOP button**.

6. Care, maintenance, transport

The **ARTROMOT®-K1** is suitable for re-use if the following points are observed.

6.1 Care / re-use

Warning!

Shock hazard –

Unplug the device from the power line before cleaning.

Shock hazard, equipment damage –
Liquids must not enter the device or the control pendant.

- The **ARTROMOT®-K1** must be **disinfected by wiping down** with a disinfectant. Thus, it complies with the special hygiene standards for medical technical equipment.
- **The enclosure and removable support assemblies can be cleaned with commonly used disinfectants and mild household detergents.**
- Use only a **damp cloth** to wipe the CPM device down.

Warning!

Patient hazard – patient contamination

- Before using the device on another patient, be sure to clean and disinfect it according to the instructions given here.

Note!

The manufacturer recommends using only special disinfectants approved for medical devices and with the characteristics specified below under "Caution".

Disinfection with the following disinfectants and in compliance with the instructions for use has been evaluated and approved by the manufacturer in the risk assessment:

- DESCOSEPT AF lemon
(Part no.: 00-311L-xxx)
BAuA Reg.No. N-55153, CE-0482
Manufacturer: Dr. Schumacher GmbH
- Contact time: 2 minutes minimum

Caution!

Equipment damage –

- The plastic material used is not resistant to mineral acids, formic acid, phenols, cresols, oxidants and strong organic or inorganic acids with a pH value below 4.
- Use only clear disinfectants to prevent discoloration of the material.
- Do not expose the CPM device to strong ultraviolet radiation (sunlight) or fire.

6.2 Maintenance (fuse replacement)

Check before each use

Visually inspect the device for signs of mechanical damage before each use.

If you detect damage or malfunctions that may impair the safety of the patient or of the operator, have the device repaired before using it.

Technical inspections

Due to their design, the bearings and joints of the **ARTROMOT®-K1** do not require maintenance and all materials are protected from corrosion, BUT: For safety, the devices require regular maintenance. To maintain the functional and operational safety, check all components for damage and loose connections at least **once a year**.

These checks should be performed by persons with adequate training and experience. Damaged and worn parts must immediately be replaced with original spare parts by authorized staff.

Warning!

Patient hazard –
equipment malfunction and damage

- Refer repair and maintenance to authorized persons. Persons are authorized after training by a specialist trained and commissioned by the manufacturer.

ORMED GmbH will make all documents required for servicing, such as circuit diagrams, parts lists, descriptions or calibration instructions, available to authorized experts.

The inspections can be carried out by DJO Technical Service within the framework of a service agreement. Contact DJO for details.

Other than that, the manufacturer does not require any other regular maintenance.

Note!

Regarding further technical or other inspections or inspection intervals, observe applicable local requirements, such as IEC 62353, DGUV3 or similar regulations to be observed by users of medical equipment or electric devices.

Fuse replacement

Warning!

Patient hazard –
equipment malfunction and damage

The replacement of fuses must be referred to specialists as defined in DIN VDE

0105 or IEC 60364 or other applicable standards (e.g. biomedical technicians, electricians, electronics installers).

Fuses used must be



S# < 20,000: T1A L250Vac

S# > 20,000: T2A H250Vac.

6.3. Transport

The following adjustments need to be carried out to prepare the **ARTROMOT®-K1** for transport:

1. Adjust the carriage to a femur length of 49 cm and a tibia length of 42 cm.
2. Activate the "transport setting" function" →  in the menu (also refer to 5.3).
3. Turn off the **ARTROMOT®-K1** power switch.
4. Disconnect the power cord and the control pendant.
5. Only use the original shipping box to transport the device. ORMED GmbH cannot be held liable for transport damage, if the device is not shipped in its original shipping box.
6. Set the ankle joint to the horizontal position.
7. Put the foam parts ① and the power cord ② in the box.
8. Then lower the **ARTROMOT®-K1** into the box.
9. Fold down the block ③ that secures the frame and insert the foam part ④ for the angle joint.
10. Then adjust the carriage to a tibia length of 50 cm so that the ankle joint will be fixed by the foam part ④. Insert the brace ⑤.
11. Put the control pendant (20) in the supplied small box.
12. Fold down the fixation element ⑥ for the control pendant and store the control pendant in the box.

7. Environmental protection statement

The product described in these instructions for use must not be disposed of with unsorted household or municipal waste. It requires separate disposal. Please contact DJO for information about the possible recycling of the product.

The expected service life of the device and all supplied parts and accessories is 6 years. The application beyond this period is the responsibility of the user.

8. Specifications

Model:	ARTROMOT®-K1
Device name:	ARTROMOT®-K1 Standard, Part no.: 80.00.040 ARTROMOT®-K1 Standard Chip, Part no.: 80.00.041 ARTROMOT®-K1 Comfort, Part no.: 80.00.042 ARTROMOT®-K1 Comfort Chip, Part no.: 80.00.043
Input ratings:	100–240 Vac / 50–60 Hz tolerance –15 % to +10 %
Current consumption:	100 V 240 V
Standby (ON):	5 VA 5 VA
Operation (maximum):	40 VA 40 VA
Fuses:	serial numbers < 20,000: 2 × T1A L250 Vac serial numbers > 20,000: 2 × T2A H250 Vac according to IEC 60127-2/3, UL 248-14, CSA C22.2 no. 248.14
Protection class:	serial numbers < 20,000: I serial numbers > 20,000: II
IP class:	serial numbers < 20,000: IPX0 serial numbers > 20,000: IP21
Applied part:	Type B
Max. load on carriage:	serial numbers < 20,000: 20 kg serial numbers > 20,000: 25 kg
Dimensions	
length:	96 cm
width:	35 cm
height:	min. 23 cm to max. 56 cm

Adjustment ranges (min./max.)	
femur range:	approx. 31 to 49 cm
lower leg range:	approx. 25 to 57 cm
Accuracy of measured values	
Goniometer in the measuring range:	from -10° to $+120^{\circ}$
Accuracy:	$\pm 2^{\circ}$
Weight:	11 kg
Materials used:	ABS, POM, PUR, PA, FR 4, aluminum, stainless steel, brass
MDD:	class IIa
Standards compliance:	93/42/EEC IEC 60601-1 IEC 60601-1-6 IEC 60601-1-11 IEC 62366 IEC 62304 EN ISO 14971 ANSI/UL 60601-1 (S# < 20,000) CAN/CSA C22.2 No. 601.1-M90 (S# < 20,000) ANSI AAMI ES 60601-1 (S# > 20,000) CAN CSA 22.2 No. 60601-1 (S# > 20,000)
EMC (electromagnetic compatibility)	IEC 60601-1-2
Manufactured in compliance with:	EN ISO 13485
Ambient conditions (storage, transport)	
Temperature:	-25°C to $+70^{\circ}\text{C}$
Relative humidity:	at 70°C up to 93%, no condensation
Atmospheric pressure:	500 hPa to 1060 hPa
Ambient conditions (operation)	
Temperature:	$+5^{\circ}\text{C}$ to $+40^{\circ}\text{C}$
Relative humidity:	15 % to 93 %
Atmospheric pressure:	700 hPa to 1060 hPa

9. Contact

We would be happy to answer any questions you may have about our products and services.

ARTROMOT® International

Please contact your local dealer, the global DJO Headquarters in the USA, the DJO International Headquarters in England, or DJO Germany.

DJO International

Headquarters: DJO UK Ltd.
1a Guildford Business Park
Guildford
Surrey, GU2 8XG, United Kingdom
Telephone: +44 (0) 1483 459659
Fax: +44 (0) 1483 459470
Email: info@DJOglobal.eu
Website: www.DJOglobal.eu

DJO Headquarters

DJO, LLC
2900 Lake Vista Drive
Lewisville, TX 75067, USA
Telephone: +1 800 321 9549
Fax: +1 800 936 6569
Email: webmaster@DJOglobal.com
Website: www.DJOglobal.com

Manufacturer / Headquarters Germany

ORMED GmbH
A DJO Company

Serial numbers < 45,000:
Merzhauser Straße 112
79100 Freiburg, Germany

Serial numbers > 45,000:
Bötzingen Straße 90
79111 Freiburg, Germany

Telephone: +49 (0) 761 4566 01
Fax: +49 (0) 761 456655-01
Email: medizintechnik@DJOglobal.com
Website: www.DJOglobal.de

Quarters Austria

DJO
ORTHOMED Medizintechnik GmbH
Annagasse 5/1/Top 4
1010 Vienna, Austria
Telephone: +43 1 5320834
Fax: +43 1 532083431
Email: orthomed@orthomed.cc
Website: www.orthomed.cc

Warranty:

2 years (mechanical parts)
2 years (electronics)

10. Technical service

10.1 Technical hotline

Do you have any technical questions?
Do you need technical service?

Telephone: +49-180-5-1 ormed de
+49-180-5-1 67 63 33
Fax : +49-180-5-3 ormed de
+49-180-5-3 67 63 33

10.2 Shipment

To prevent damage during transport, only use the original shipping box. These boxes can be obtained from DJO.

Before packing the therapy unit, set it to the transport position (see chapter 5).

10.3 Spare parts

Refer to the Service Manual for the most recent list of spare parts.

When ordering spare parts, always specify:

- Item
- Description
- Part number
- Quantity
- Serial number of the device

Note!

Refer repairs to authorized, specially trained staff.

DJO offers service training for your staff.

Note!

Country-specific power cords are available from DJO. Please contact DJO or your DJO distributor.

Surcharges may apply in certain cases to spare parts ordered in low quantities.

Item	Description	Part no.	Qty
1.	Patient chip card	0.0034.035	<input type="checkbox"/>
2.	Patient chip card (protocol)	0.0031.006	<input type="checkbox"/>
3.	Marker pen for patient chip card	2.0037.024	<input type="checkbox"/>
4.	Power cord EU version H05VV-F3G 3 × 1 mm ² length 4 m	0.0034.118	<input type="checkbox"/>
5.	Fuse T2A H250 Vac 1500A Breaking Capacity	0.0000.107	<input type="checkbox"/>

11. CE-Declaration of Conformity

The certificate can be download from:

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<https://www.djoglobal.de/unternehmen/qualitaet-und-verantwortung>

12. IEC 60601-1-2 and IEC 60601-1-11

For **ARTROMOT®-K1** or **ARTROMOT®-K1 Classic** devices with serial numbers >#35.000.

The product designation **ARTROMOT®-K1** used in the text below includes the **ARTROMOT®-K1** and **ARTROMOT®-K1 Classic** products with all their variants.

The **ARTROMOT®-K1** device is subject to particular precautions regarding electromagnetic compatibility (EMC). The device must be installed and put into service strictly in compliance with the EMC directives set forth in the accompanying documents.

Portable and mobile RF communication systems may affect the **ARTROMOT®-K1**. The **ARTROMOT®-K1** should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, **ARTROMOT®-K1** should be observed to verify normal operation in the configuration in which it will be used.

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If it is necessary to replace assemblies or cables, only the manufacturer's original parts must be used to ensure continued compliance with EMC requirements after repair. This requirement applies to the power supply unit, cables and cable lengths, drive unit consisting of the motor and the control system, the programming unit incl. the coiled cable and the connector.

With regards to the electromagnetic compatibility, there is no maintenance operation required during the expected product life time!

WARNING!

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

12.1 Electromagnetic emissions

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

ARTROMOT®-K1 is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-K1** unit is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions to CISPR 11	Group 1	ARTROMOT®-K1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions to CISPR 11	Class B	ARTROMOT®-K1 is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions to IEC 61000-3-2	complies with class A requirements	
Voltage fluctuations/flicker emissions to IEC 61 000-3-3	complies	

12.2 Electromagnetic immunity

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

ARTROMOT®-K1 is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-K1** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment – guidance
Electrostatic discharge (ESD) to IEC 61 000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/bursts to IEC 61 000-4-4	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power should be that of a typical commercial or hospital environment.
Surges to IEC 61 000-4-5	± 1 kV line-to-line voltage ± 2 kV line-to-ground voltage	± 1 kV line-to-line voltage ± 2 kV line-to-ground voltage	Mains power should be that of a typical commercial or hospital environment.
Voltage dips To IEC 61000-4-11	0% U_T ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T ; 1 cycle and 70% U_T ; 2 5/30 cycles Single phase at 0°	0% U_T ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T ; 1 cycle and 70% U_T ; 25/30 cycles Single phase at 0°	Mains power should be that of a typical commercial or hospital environment. If the user of the ARTROMOT®-K1 unit requires continued operation during power mains interruptions, it is recommended that the ARTROMOT®-K1 unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field to IEC 61 000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.
Voltage interruptions to IEC 61000-4-11	0% U_T ; 250/300 cycles	0% U_T ; 250/300 cycles	

Note: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

ARTROMOT®-K1 is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®-K1** unit is used in such an environment.

Immunity test	IEC 60601 test level Home healthcare	Compliance level Home healthcare	Electromagnetic environment - guidance
Conducted RF to IEC 61000-4-6	6 V 0,15 MHz-80 MHz 80% AM at 1 kHz	6 V 0,15 MHz-80 MHz 80% AM at 1 kHz	
Radiated RF to IEC 61000-4-3	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communications equipment to IEC 61000-4-3	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ARTROMOT®-K1** exceeds the applicable **ARTROMOT®-K1** unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **ARTROMOT®-K1** unit.

12.3 RF (Radio Frequency) Interference

EMC disturbance can affect the performance of the **ARTROMOT®-K1**. To prevent the user from unacceptable risks the device runs performance checks regularly during operation. This is what happens, if a problem is identified:

- an audio signal sounds
- the unit stops immediately
- the message ERROR and an error code appear on the display

In this case try to quit the error with the STOP key. If the error persists you may attempt to restart the unit by turning it briefly off and on with the power switch.

WARNING!

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the **ARTROMOT®-K1**, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Notice!

(for home use)

Portable and mobile RF communication systems may affect the **ARTROMOT®-K1** device.

For this reason, make sure that wireless communication equipment, such as wireless home network devices, cell phones, cordless telephones and their base stations and walkie-talkies operate at a minimum distance of 3.3 m from the device. (Calculated on the basis of the maximum power output of a typical cell phone of 2 W)..

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Symbol Overview standard models ARTROMOT®-K1

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therapy time

Lastumkehr

reverse on load
(feature for
patient safety)

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Symbol Overview comfort models ARTROMOT®-K1

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Lastumkehr
reverse on load
(feature for
patient safety)

Pause Flexion
pause flexion

EBENE 3 / LEVEL 3:

Dehnung Extension
stretching extension

Work-out Programm
workout protocol

Comfort Programm
Comfort protocol

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

AMENDMENT TO THE INSTRUCTION FOR USE (IFU)

This document outlines the amendment to the Instruction for Use (IFU) for necessary updates and changes.

List of the affected devices:

Reference Number	Article Name	GTIN
80.00.040	ARTROMOT® -K1 Standard	00888912403191
80.00.041	ARTROMOT® -K1 Standard Chip	00888912403207
80.00.042	ARTROMOT® -K1 Comfort	00888912403214
80.00.043	ARTROMOT® -K1 Comfort Chip	00888912403221
80.00.044	ARTROMOT® -K1 Standard Basic	00190446751823
80.00.045	ARTROMOT® -K1 Classic	00888912403238
80.00.047	ARTROMOT® -K1 Classic Basic	00190446751830
80.00.023	ARTROMOT® -S3 Standard	00888912403108
80.00.024	ARTROMOT® -S3 Comfort	00888912403115
80.00.033	ARTROMOT® -E2 Compact	00888912403146
80.00.031	ARTROMOT® -E2	00888912403139
80.00.035	ARTROMOT® -SP3 Standard	00888912403153
80.00.036	ARTROMOT® -SP3 Standard Chip	00888912403160
80.00.037	ARTROMOT® -SP3 Comfort	00888912403177
80.00.038	ARTROMOT® -SP3 Comfort Chip	00888912403184
80.00.070	ARTROMOT® ACTIV E-K	00888912403283
80.00.084	ARTROMOT® -S4 Comfort	00190446278764
80.00.022	OptiFlex®-S	00190446794783
80.00.046	OptiFlex®-K1 Classic	00190446796015
80.00.048	OptiFlex®-K1 Standard	00190446796039
80.00.049	OptiFlex®-K1 Comfort Chip	00190446796046

First change: Warning

New warnings have been implemented due to labeling changes that add or strengthen a warning, precaution, or include additional information.

Safety warnings:

Warning!


"The treatment of minors or persons with limited physical, sensory or mental abilities (e.g. movement disorders, epilepsy, seizures, mental disorders, behavioural disorders or neurodevelopmental disorders) may only be carried out under the constant observation of attendants (e.g. specialist staff, legal guardians or custodians)

– After data storage, write the patient's name on the patient chip card. The card should only be used for this particular patient. If the patient chip card is used for another patient, be sure to delete the previous patient data from the card first. Use original chip cards only."

Second change: Label

Symbols

Table 1: Symbols changes

Symbol change	Old version	New version	Description of symbol acc. ISO 15223 + 2014/
UDI	New requirement		Indicates a carrier that contains unique device identifier information. This symbol identifies the UDI carrier, including the AIDC and human readable information. (01) GTIN code of device (11) Date of production (21) serial number of device
Manufacturer + Date of production + country of manufacture			Indicates the medical device manufacturer Indicates the date when the medical device was manufactured DE is the two or three letter country code as defined in ISO 3166-1. Indicates the country of manufacture
Manufacturer's name	ORMED GmbH A DJO Company	ORMED GmbH a company of Enovis	Name of the manufacturer
Manufacturer's address spelling	Böszinger Straße 90 79111 Freiburg / Germany	Böszinger Strasse 90 79111 Freiburg / Germany	Address of the manufacturer
Norm/Standard ETL classification	Conforms to ANSI/ AAMI ES STD 60601-1 Certified CAN/CSA C22.2 No. 60601-1		Medical electrical equipment relevant standards according ETL classification – control number above

Other changes

Manufacturers contact was changed to: Tel. +49 761 4566-01 · Fax +49 761 4566-5501, medizintechnik@enovis.com, enovis-medtech.de

Third change: Electromagnetic emissions, Electromagnetic immunity and RF (radio frequency) interference

The EMC table provides a detailed comparison of electromagnetic compatibility as per the latest standards.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
ARTROMOT® CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the ARTROMOT® CPM device unit is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions to CISPR 11: 2015+A1:2016+A2:2019	Group 1	ARTROMOT® CPM device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions to CISPR 11: 2015+A1:2016+A2:2019	Class B	ARTROMOT® CPM device is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions to IEC 61000-3-2: 2005+A1:2008 +A2:2009	complies with class A requirements	
Voltage fluctuations/flicker emissions to IEC 61000-3-3: 2013	complies	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

ARTROMOT® CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®** CPM device unit is used in such an environment.

Immunity test	IEC 60601-1-2:2014 +A1:2020		Electromagnetic environment – guidance
	test level Home healthcare	Compliance level Home healthcare	
Electrostatic discharge (ESD) to IEC 61000-4-2:2008	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/bursts to IEC 61000-4-4:2012	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power should be that of a typical commercial or hospital environment.
Surges to IEC 61000-4-5:2014	± 1 kV line-to-line voltage	± 1 kV line-to-line voltage	Mains power should be that of a typical commercial or hospital environment.
Voltage dips To IEC 61000-4-11:2017	0% Ut ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% Ut ; 1 cycle and 70% Ut ; 25/30 cycles Single phase at 0°	0% Ut ; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% Ut ; 1 cycle and 70% Ut ; 25/30 cycles Single phase at 0°	Mains power should be that of a typical commercial or hospital environment. If the user of the ARTROMOT® CPM device unit requires continued operation during power mains interruptions, it is recommended that the ARTROMOT® CPM device unit be powered from an uninterruptible power supply or a battery.
Voltage interruptions to IEC 61000-4-11:2017	0% Ut; 250/300 cycles	0% Ut; 250/300 cycles	
Power frequency (50/60 Hz) magnetic field to IEC 61000-4-8:2009	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical commercial or hospital environment.

Note: Ut is the a.c. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

ARTROMOT® CPM device is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the **ARTROMOT®** CPM device unit is used in such an environment.

Immunity test	IEC 60601-1-2:2014 +A1:2020		Electromagnetic environment – guidance
	test level Home healthcare	Compliance level Home healthcare	

Conducted RF to IEC 61000-4-6:2013	3 V 0.15 MHz-80 MHz 80% AM at 1kHz 6V 0.15 MHz-80 MHz 80% AM at 1kHz at ISM frequency bands and Amateur radio bands	3 V 0.15 MHz-80 MHz 80% AM at 1kHz 6V 0.15 MHz-80 MHz 80% AM at 1kHz at ISM frequency bands and Amateur radio bands	
Radiated RF to IEC 61000-4-3: 2006+A1:2007 +A2:2010	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz-2,7 GHz 80% AM at 1 kHz	
Proximity fields from RF wireless communica- tions equipment to IEC 61000-4-3: 2006+A1:2007 +A2:2010	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	27 V/m 380-390 MHz PM 18 Hz 28 V/m 430-470 MHz FM +/- 5 kHz 1 kHz sine 9 V/m 704-787 MHz PM 217 Hz 28 V/m 800-960 MHz PM 18 Hz 28 V/m 1,7-1,99 GHz PM 217 Hz 28 V/m 2,4 GHz-2,57 GHz at PM 217 Hz 9 V/m 5,1 GHz-5,8 GHz at PM 217 Hz	
Radiated fields in close proximity IEC 61000-4-39:2017	8 A/m 30 kHz No modulation (CW) 65 A/m 134.2 kHz PM 2.1 kHz 7,5 A/m 13.56 MHz PM 50 kHz	8 A/m 30 kHz No modulation (CW) 65 A/m 134.2 kHz PM 2.1 kHz 7,5 A/m 13.56 MHz PM 50 kHz	

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **ARTROMOT®** CPM device exceeds the applicable **ARTROMOT®** CPM device unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **ARTROMOT®** CPM device unit.

The **ARTROMOT®** CPM movementsplint is intended for operation in the specified electromagnetic environment. The customer or user of the **ARTROMOT®** CPM device should ensure that it is used in such an environment. The required interference levels of the specified standards are met.

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